
EN 5	<b>Evolution® Colonic Stent System – Uncovered</b> Instructions for Use
BG 10	<b>Стент система за колон Evolution® – без обвивка</b> Инструкции за употреба
HR 17	<b>Sustav kolonskog stenta Evolution® – otkriveni</b> Upute za uporabu
CS 23	<b>Systém kolonického stentu Evolution® – nekrytý</b> Návod k použití
DA 29	<b>Evolution® colonstentsystem – udækket</b> Brugsanvisning
NL 35	<b>Evolution® colonstentsysteem – niet-bekleed</b> Gebruiksaanwijzing
ET 41	<b>Käärsoole stendisüsteem Evolution® – katteta</b> Kasutusjuhised
FR 46	<b>Système d'endoprothèse colique Evolution® – non couverte</b> Mode d'emploi
DE 53	<b>Evolution® Kolonstentsystem – Nicht gecovert</b> Gebrauchsanweisung
EL 59	<b>Σύστημα ενδοπρόσθεσης κόλου Evolution® – ακάλυπτο</b> Οδηγίες χρήσης
HU 66	<b>Evolution® colonsztentrendszer – fedetlen</b> Használati utasítás
IT 72	<b>Sistema con stent colico Evolution® – Non rivestito</b> Istruzioni per l'uso
LV 79	<b>Evolution® lokzarnas stenta sistēma — bez pārkājuma</b> Lietošanas instrukcija
LT 85	<b>„Evolution®“ gaubtinės žarnos stento sistema – nedengta</b> Naudojimo nurodymai
NO 91	<b>Evolution® stentsystem for kolon – utildekket</b> Bruksanvisning
PL 96	<b>System Evolution® do stentowania okrężnicy – nieosłonięty</b> Instrukcja użycia
PT 103	<b>Sistema de stent de cólon Evolution® – Não revestido</b> Instruções de utilização
RO 109	<b>Sistem pentru stentare colonică Evolution® – fără strat de acoperire</b> Instrucțiuni de utilizare
SK 115	<b>Systém kolonického stentu Evolution® – nekrytý</b> Návod na použitie
ES 121	<b>Sistema de stent colónico Evolution® – sin recubrimiento</b> Instrucciones de uso
SV 127	<b>Evolution®-stentsystem för kolon – utan täckning</b> Bruksanvisning


I F U 0 0 5 2 - 1 2


1.


2.

3.

4.

5.

## **EVOLUTION® COLONIC STENT SYSTEM - UNCOVERED**

**CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### **IMPORTANT INFORMATION**

Please review prior to use.

### **INTENDED USE**

This device is used for palliative treatment of colonic obstruction or colonic strictures caused by malignant neoplasms, and to relieve large bowel obstruction prior to colectomy in patients with malignant strictures.

### **STENT DESCRIPTION**

This flexible, self-expanding stent is constructed of a single, woven, nitinol wire. The stent has a flange at both ends to aid in preventing migration. The stent foreshortens due to its design. The total length of the stent is indicated by radiopaque markers on the inner catheter, indicating the actual length of the stent at nominal stent diameter. To aid in visibility under fluoroscopy, there are four main tantalum marker bands at either end of the stent. The proximal and distal ends of the stent are designed with numerous rounded "crowns" to potentially reduce tissue trauma.

### **DELIVERY SYSTEM DESCRIPTION**

The stent is mounted on an inner catheter, which accepts a .035-inch wire guide, and is constrained by an outer catheter. An endoscopically and fluoroscopically visible yellow marker defines the proximal end of the stent when constrained in the catheter. A pistol-grip delivery handle allows controlled stent deployment or recapture. The position of the red marker on top of the handle in relation to the point-of-no-return position on the handle label indicates the point after which it is not possible to recapture the deployed stent.

## **CONTRAINDICATIONS**

Those specific to GI endoscopy and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to: enteral ischemia • suspected or impending perforation • intra-abdominal abscess/perforation • inability to pass wire guide or stent through obstructed area • patients for whom endoscopic procedures are contraindicated • significant coagulopathy • benign disease.

## **POTENTIAL ADVERSE EVENTS**

Potential adverse events associated with GI endoscopy include, but are not limited to: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • fever • hemorrhage • hypotension • infection • perforation • respiratory depression or arrest.

Additional adverse events include, but are not limited to: allergic reaction to nickel • bowel impaction • constipation • death (other than due to normal disease progression) • diarrhea • erosion of the luminal mucosa • foreign body sensation • inadequate stent expansion • intestinal perforation • nausea/vomiting • pain/discomfort • peritonitis • pressure necrosis • septicemia • stent misplacement and/or migration • stent occlusion • symptoms of tenesmus or urgency/incontinence • tumor ingrowth or overgrowth • ulcerations.

## **PRECAUTIONS**

Refer to product package label for the minimum channel size required for this device.

A complete diagnostic evaluation must be performed prior to use to determine proper stent size.

If wire guide or stent cannot advance through obstructed area, do not attempt to place stent.

Stent should be placed endoscopically with fluoroscopic monitoring.

The stent should only be placed with the Cook delivery system, which is provided with each stent.

This device is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

Long-term patency of this device has not been established. Periodic evaluation is advised.

Use of this device is restricted to a trained healthcare professional.

Do not use this device for any purpose other than stated intended use.

Stent cannot be retrieved after the deployment threshold has been passed. The position of the red marker on the top of the handle in relation to the point-of-no-return position on the handle label indicates when threshold has been passed.

## **WARNINGS**

The stent is not intended to be removed or repositioned after stent placement and is intended to remain in the body permanently. Attempts to remove or reposition stent after placement may cause damage to surrounding tissue or mucosa.

This stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

To minimize pain and tenesmus, the stent end nearest to the anal canal/anus should be placed 2 cm above the anal canal or 6 cm from the anus.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

After stent placement, alternative methods of treatment such as chemotherapy and radiation should not be administered as this may increase risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

The device should be used with caution and only after careful consideration in patients with:

- Patients with radiation colitis or proctitis.
- Patients with elevated bleeding times, coagulopathies.

## MR INFORMATION

This symbol means the device is MR Conditional.

Non-clinical testing has demonstrated that the Evolution® Colonic **Stent** is MR Conditional according to ASTM F2503. A patient with this **stent** can be scanned safely anytime after placement under the following conditions:

- Static magnetic field of 3.0 Tesla or less.
- Spatial magnetic gradient field of 720 Gauss/cm. Non-clinical evaluation was conducted in an MR system (Excite, General Electric Healthcare) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

### MRI-Related Heating

- 1.5 and 3.0 Tesla Systems: It is recommended to scan in normal operation mode (whole body averaged specific absorption rate (SAR) = 2.0 W/kg) ("Normal Operating Mode" is defined as the mode of operation of the MR system in which none of the outputs have a value that cause physiological stress to the patient) for 15 minutes of scanning (i.e., per scanning sequence).

#### 1.5 Tesla Temperature Rise

In non-clinical testing, the Evolution® Colonic **Stent** produced maximum temperature rises of 1.9 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a 1.5 Tesla MR System (Magnetom, Siemens Medical Solutions, Malvern, PA, Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.1 W/kg).

#### 3.0 Tesla Temperature Rise

In non-clinical testing, the Evolution® Colonic **Stent** produced maximum temperature rises of 2.7 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a 3 Tesla MR System (Excite, HDx, General Electric Healthcare, Milwaukee, WI, Software 14X.M5) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.7 W/kg).

## **Image Artifacts**

MR image quality may be compromised if the area of interest is within the lumen or within approximately 3 mm of the position of the Evolution® Colonic **Stent**, as found during non-clinical testing using the sequence: T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI, Software 14X.M5). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic **stent**.

**Additional Information:** The safety of performing an MRI procedure in a patient with overlapping colonic stents or other MRI-conditional device(s) in direct contact with this device has not been determined. Performing MRI in such situations is not recommended.

It is recommended that patients register the conditions under which the implant can safely be scanned with the MedicAlert Foundation ([www.medicalert.org](http://www.medicalert.org)) or an equivalent organization.

## **EQUIPMENT REQUIRED**

Stent and delivery system

.035-inch wire guide

## **SELECTION OF STENT**

A stent length of at least 4 cm longer than the stricture is recommended (e.g. approximately 2 cm longer at either end of the stricture)

## **PATIENT PREPARATION**

With colonoscope in place, introduce a wire guide until it is fluoroscopically visualized in position through stricture.

## **INSTRUCTIONS FOR USE**

## **Illustrations**

1. Remove stent delivery system from the package.
2. Introduce delivery system in short increments over wire guide, into accessory channel, until it is endoscopically visualized exiting the scope.
3. To ensure that stent will bridge stricture after deployment, fluoroscopically position radiopaque markers on inner catheter beyond extremities of stricture to be crossed. **Note:** Yellow marker on the delivery system located at proximal end of stent can also serve as an endoscopic/fluoroscopic reference for positioning of the proximal (nearest user) end of the stent relative to the stricture.

- 4.** Confirm desired stent position fluoroscopically and deploy stent by removing red safety guard from the handle. (**Fig. 1**)
- 5.** Continue deploying stent by squeezing trigger. **Note:** Each trigger squeeze will deploy stent by an equal amount.
- 6.** If stent repositioning is required during deployment, it is possible to recapture stent. **Note:** It is **not** possible to recapture stent after passing point-of-no-return, indicated when red marker on top of the handle has passed the point-of-no-return position on handle label. (**Fig. 2**) The stent can be recaptured a maximum of five times.
- 7.** To recapture stent, push direction button on side of delivery system to opposite side. (**Fig. 3**) **Note:** Hold thumb on button when squeezing trigger for first time to recapture.
- 8.** Continue squeezing trigger as required to recapture stent by desired amount.
- 9.** To resume deployment, push button to opposite side again and hold thumb on button for first stroke while squeezing trigger. (**Fig. 4**)
- 10.** When stent point-of-no-return has been passed, pull safety wire out of delivery handle near wire guide port. (**Fig. 5**)
- 11.** Continue deploying stent by squeezing trigger.
- 12.** After deployment, fluoroscopically confirm full stent expansion. Once full expansion is confirmed, delivery system can be safely removed.

**Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.**

#### **HOW SUPPLIED**

Store in a dry location, away from temperature extremes.

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© COPYRIGHT COOK 2020 ©
**COOK IRELAND LTD.**