
<b>EN</b> 3	<b>Pancreatic Stent</b> Instructions for Use
<b>CS</b> 5	<b>Pankreatický stent</b> Návod na použití
<b>DA</b> 7	<b>Pancreasstent</b> Brugsanvisning
<b>NL</b> 9	<b>Pancreasstent</b> Gebruiksaanwijzing
<b>FR</b> 11	<b>Endoprothèse pancréatique</b> Mode d'emploi
<b>DE</b> 14	<b>Pankreasstent</b> Gebrauchsanweisung
<b>EL</b> 16	<b>Παγκρεατική ενδοπρόσθεση</b> Οδηγίες χρήσης
<b>HU</b> 18	<b>Hasnyálmirigysztent</b> Használati utasítás
<b>IT</b> 21	<b>Stent pancreatico</b> Istruzioni per l'uso
<b>NO</b> 23	<b>Bukspyttkjertelstent</b> Bruksanvisning
<b>PL</b> 25	<b>Stent do przewodów trzustkowych</b> Instrukcja stosowania
<b>PT</b> 28	<b>Stent pancreático</b> Instruções de utilização
<b>RO</b> 30	<b>Stent pancreatic</b> Instrucțiuni de utilizare
<b>ES</b> 32	<b>Stent pancreático</b> Instrucciones de uso
<b>SV</b> 34	<b>Pankreasstent</b> Bruksanvisning

**INTENDED USE**

This device is used to drain obstructed pancreatic ducts.

**NOTES**

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

Retain the positioning sleeve (if any)\* for use when introducing the stent flaps into the accessory channel.

**CONTRAINDICATIONS**

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Inability to pass the wire guide or stent through the obstructed area.

**POTENTIAL COMPLICATIONS**

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Those associated with pancreatic stent placement include, but are not limited to: trauma to the pancreatic tract or duodenum, obstruction of the common bile duct, stent migration.

**PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

A complete diagnostic evaluation must be performed, prior to use, to determine proper stent size.

The positioning sleeve\* or the pigtail straightener\* is not intended for use in the accessory channel of the endoscope.

Care must be exercised when straightening the pigtail curls\* in order to avoid kinking or breaking the stent.

The tapered tip end of the stent\* must be positioned in the pancreatic duct while the other end remains in the duodenum.

This device should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

Select the Cook stent introducer system (if not included) in the appropriate French size.

## **WARNINGS**

MRI compatibility of pancreatic stents has not been established.

## **INSTRUCTIONS FOR USE**

## **Illustrations**

### **I. FOR NON-PIGTAIL STENTS:**

1. Gently ensure full extension of all side flaps\*.
2. Load positioning sleeve onto duodenal flap end of stent.

#### **For 7 FR and smaller stents:**

- 3a. Introduce stent, tapered tip first, and positioning sleeve onto a pre-positioned wire guide.

#### **For 8.5 FR and larger stents:**

- 3b. Remove Tuohy-Borst adapter from end of guiding catheter\*, then introduce guiding catheter\* into accessory channel over a pre-positioned wire guide.
- 3c. Introduce stent, tapered tip first, and positioning sleeve onto guiding catheter\* and pre-positioned wire guide.

## **THEN REFER TO STEPS 4-7 BELOW**

### **II. FOR PIGTAIL STENTS:**

#### **For 7 FR and smaller stents:**

1. Advance pigtail straightener\* along shaft of stent from the tapered tip end in order to straighten pigtail curl. (*See fig. 1*)
2. Introduce stent, tapered tip first, and pigtail straightener\* onto pre-positioned wire guide.
3. Advance pushing catheter over wire guide to advance pigtail stent into accessory channel.

#### **For 8.5 FR and larger stents:**

- 3a. Advance pigtail straightener\* along shaft of stent from the tapered tip end in order to straighten pigtail curl. (*See fig. 1*)
- 3b. Introduce stent, pigtail first, and pigtail straightener\* onto guiding catheter\*. Advance guiding catheter over pre-positioned wire guide.

## **THEN REFER TO STEPS 4-7 BELOW**

4. Advance guiding catheter\* and/or pushing catheter in 1-2 cm increments until stent is in desired position.
5. Fluoroscopically and endoscopically confirm desired stent position.
6. After confirming stent position, gently remove wire guide, then the guiding catheter\* from endoscope while maintaining position of stent with pushing catheter.
7. Gently remove pushing catheter from accessory channel.

## **Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.**

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© COPYRIGHT COOK 2019©
**COOK IRELAND LTD.**