EN 8	Pancreatic Stents and Sets Instructions for Use
BG 13	Панкреатични стентове и комплекти Инструкции за употреба
CS 19	Pankreatické stenty a soupravy Návod k použití
DA 25	Pancreasstenter og stentsæt Brugsanvisning
DE 31	Pankreasstents und Pankreasstent-Sets Gebrauchsanweisung
EL 37	Παγκρεατικές ενδοπροθέσεις και σέτ Οδηγίες χρήσης
ES 43	Stents y equipos de stents pancreáticos Instrucciones de uso
ET 49	Pankrease stendid ja stendikomplektid Kasutusjuhend
FR 55	Endoprothèses pancréatiques et sets d'endoprothèses pancréatiques Mode d'emploi
HR 61	Komplet stenta za gušteraču Upute za uporabu
HU 66	Hasnyálmirigysztent és készletek Használati utasítás
IT 72	Stent e set pancreatici Istruzioni per l'uso
LT 78	Kasos stentai ir rinkiniai Naudojimo instrukcija
LV 84	Aizkuņģa dziedzerā stentu komplekti Lietošanas instrukcija
NL 90	Pancreasstents en -sets Gebruiksaanwijzing
NO 96	Bukspyttkjertelstenter og -sett Bruksanvisning
PL 101	Stenty do przewodów trzustkowych i zestawy Instrukcja użycia
PT 108	Conjuntos e stents pancreáticos Instruções de utilização
RO 113	Stenturi pancreatice și seturi de stenturi pancreatice Instrucțiuni de utilizare
SK 119	Pankreatické stenty a súpravy Návod na použitie
SL 125	Stenti za trebušno slinavko in sklopi stentov Navodila za uporabo
SR 131	Pankreasni stentovi i kompleti Uputstvo za upotrebu

SV
136

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
Bruksanvisning

TR
142

Pankreatik Stentler ve Stent Setleri

Kullanma Talimatı

EN Fig. 1: Geenen Pancreatic Stent

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GEENEN®, GEENEN® SOF-FLEX®, AND ZIMMON® PANCREATIC STENTS AND STENT SETS

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Please read all instructions before using this device.

DEVICE DESCRIPTION

This IFU covers Geenen and Zimmon Pancreatic Stents represented by Figures 1 and 2.

These devices are available in two device configurations:

- Stent only (with or without flap protector)
- Stent set, which includes a combination of a stent (with or without a flap protector) and an introducer (pushing catheter and, where applicable, a guiding catheter).

Performance Characteristics

The function and key features of these components are described below.

A variety of pancreatic stents that provide drainage are available across the device range to accommodate various patient anatomies. They are offered in French sizes of between 3 Fr and 11.5 Fr, in various shape configurations and in labelled lengths of between 2 cm and 25 cm. Pancreatic stents have the following features:

- Atraumatic tip – The ends of the stent are designed to be atraumatic to the anatomy
- Sideports – Aid drainage
- Anti-migration features – Aim to help maintain the position of the stent once it is placed by providing mechanical interference with the anatomy
- Radiopacity – The stent materials are radiopaque to facilitate visibility under fluoroscopy. Where present, marker bands further support the fluoroscopic visibility of the stent.

Flap protector – This is a tube used to collapse the flap(s) as it is introduced into the endoscope.

Guiding catheter – The stent is loaded onto the guiding catheter. The function of the guiding catheter is to guide the stent to its intended location.

Guiding catheters contain one or more of the following features:

- Atraumatic tip – The end of the guiding catheter is designed to be atraumatic to the anatomy
- Radiopacity – The guiding catheter tubing material is radiopaque. The guiding catheter also contains marker bands to further enhance its fluoroscopic visibility.
- Hub – The guiding catheter contains a hub with a connection to facilitate the attachment of a standard Luer syringe.

Pushing catheter – The function of the pushing catheter is to advance the stent over a pre-positioned wire guide and, where applicable, a guiding catheter to its intended location, and to maintain the position of the stent as it's being deployed. It contains the following feature:

- Proximal identifier – An identifier at its proximal end. It is used to identify the orientation in which the device must be used.

See product label for all dimensions.

Device Compatibility

The Geenen and Zimmon Pancreatic Stents are compatible with the following:

- Endoscope with 3.2 mm, 3.7 mm, or 4.2 mm accessory channel (refer to label for selection of the appropriate accessory channel)
- 0.018", 0.021", or 0.035" wire guide (refer to label for selection of the appropriate wire guide size)
- Endoscopic cap or wire guide locking device
- Stent retriever or forceps
- Cook Medical pushing catheter
- Cook Medical guiding catheter
- Contrast media
- Sterile water or saline
- Standard Luer syringe

Table 1 lists the compatible introducers for pancreatic stents supplied as stent only.

Table 1: Compatible introducers for pancreatic stents only.

Name	Device Code Prefix	Pushing Catheter (PC): Compatible Stent Size (Fr)	Guiding Catheter (GC): Compatible Stent Size (Fr)
Geenen Pancreatic Stent	GPSO, GPSOS	PC-3: 3 Fr PC-4: 4 Fr PC-5: 5 Fr PC-6: 6 Fr PC-7: 7 Fr PC-8.5: 8.5 Fr PC-10: 10 Fr PC-11.5: 11.5 Fr	GC-5: 8.5 Fr GC-6: 10 Fr, 11.5 Fr
	GPSO-SF, GPSOS-SF		N/A
Zimmon Pancreatic Stent	SPSOS, SPSOF, ZPSOS, ZPSOF		N/A

Note: Stent sets ZEPDF, ZEPDS, GEPD, and GPDS contain an introducer within the sets.

Qualitative and Quantitative Information

The materials for the stent implant are outlined in Table 2.

Table 2: Stent implant materials.

Product	Implant Material Qualitative Information		Weight (g)
3 Fr and 4 Fr Zimmon® Pancreatic Stent (SPSOF, SPSOS)	Polyethylene	Stent	Up to 0.26
5 Fr Zimmon® Pancreatic Stent (ZEPDF, ZEPDS, ZPSOF, ZPSOS)	Polyethylene	Stent	Up to 0.30
5 Fr Zimmon® Pancreatic Stent (SPSOF, SPSOS)	Polyethylene Platinum with 10% iridium	Stent Marker band	Up to 0.28 Up to 0.004
7 Fr Zimmon® Pancreatic Stent (ZEPDF, ZEPDS, ZPSOF, ZPSOS)	Ethylene-vinyl acetate (EVA) copolymer	Stent	Up to 0.85
7 Fr Zimmon® Pancreatic Stent (SPSOF, SPSOS)	EVA copolymer Tantalum	Stent Marker band	Up to 0.79 Up to 0.10
6 Fr, 8.5 Fr, and 10 Fr Zimmon® Pancreatic Stent (SPSOF)	Polyethylene Tantalum	Stent Marker band	Up to 1.46 Up to 0.10
3 Fr, 4 Fr, 5 Fr, 6 Fr, and 8.5 Fr Geenen® Pancreatic Stent (GEPD, GPDS, GPSO, GPSOS)	Polyethylene	Stent	Up to 0.52
7 Fr, 10 Fr, and 11.5 Fr Geenen® Pancreatic Stent (GEPD, GPDS, GPSO, GPSOS)	EVA copolymer	Stent	Up to 1.00
5 Fr Geenen® Sof-Flex® Pancreatic Stent (GPSO-SF, GPSOS-SF)	Polyurethane	Stent	Up to 0.20

Patient Population

These devices can be used in all patient populations (≥ 10 kg) requiring pancreatic stenting. The nature of the underlying pathology is prevalent in varying patients, therefore the devices are indicated for patients with pancreatitis, stricture, stones, cancer, anatomic anomalies, disrupted duct, fistula/leak, pancreatic fluid collection and patients in whom prophylactic stenting for prevention of post-ERCP (endoscopic retrograde cholangiopancreatography) pancreatitis is required.

Intended User

The product is intended for use by physicians trained and experienced in ERCP techniques. Standard techniques for the endoscopic placement of stents should be employed.

Contact with Body Tissue

The device is tissue contacting in line with intended use.

Operating Principles

The pancreatic stent operates by providing a lumen through which pancreatic fluid can drain. The stent has an anti-migration feature which provides mechanical interference with the anatomy to help

maintain the position of the stent. The stent and marker band materials are radiopaque to facilitate visibility under fluoroscopy. Sideports on the stent aid drainage by providing additional openings for fluid flow. The pancreatic stent is placed over a pre-positioned wire guide using an introducer. An introducer consists of a pushing catheter only or a pushing catheter and a guiding catheter combined. The guiding catheter contains marker bands which are radiopaque to facilitate visibility under fluoroscopy and guides stent placement. The pushing catheter operates by pushing the stent into place along a pre-positioned wire guide. To protect the flaps on some of the stents a flap protector is provided. The flap protector operates by providing a lumen through which the stent can be inserted. This collapses the flaps on the stent before inserting into the endoscope.

INTENDED USE

This device is used to drain pancreatic ducts.

INDICATIONS FOR USE

Endoscopic pancreatic stent placement is used for pancreatic drainage that could be caused by pancreatitis, stricture, pancreatic cancer, anatomic anomalies of the pancreas, pancreatic fluid collection, pancreatic stones, disrupted duct, fistula/pancreatic leak. Pancreatic stents are also used prophylactically for prevention of post-ERCP pancreatitis.

CLINICAL BENEFITS

Drainage of the pancreatic duct.

CONTRAINDICATIONS

Those specific to ERCP and any procedure to be performed in conjunction with stent placement. Inability to pass the wire guide, guiding catheter, and/or stent through the obstructed area.

WARNINGS

- This single-use device is not designed for reuse. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or migration and/or mechanical integrity failure of the device.
- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect the device with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.
- These pancreatic stents are intended up to a maximum of 3 months indwell.
- For post-ERCP pancreatitis (PEP) prevention, these pancreatic stents are intended for a minimum of 12–24 hours and up to a maximum of 10 days indwell.

PRECAUTIONS

- Refer to the label for the minimum channel size required for this device.
- Refer to the label for selection of the appropriate wire guide size.
- A complete diagnostic evaluation must be performed, prior to use, to determine proper stent size.
- Select the Cook stent introducer (if not included) in the appropriate size. See Table 1 for guidance.
- Fluoroscopic monitoring must be utilized during stent placement.
- Refer to Figures 1 to 2 for the intended placed orientation of the stent in the pancreatic duct.
- The positioning sleeve (where supplied) is not intended for use in the accessory channel of the endoscope.
- Care must be exercised when straightening the pigtail curls (where applicable) in order to avoid kinking or breaking the stent.
- Periodic evaluation of the device is recommended during indwell period.
- Do not use excessive force to advance the stent.
- Dislodgement of a placed stent is possible when attempting additional procedures.
- Sphincterotomy is not necessary for device placement.
- Do not use this device for any purpose other than stated intended use.
- Store the device in a dry location.

- Use of this device is restricted to a trained healthcare professional.
- The product is intended for use by physicians trained and experienced in ERCP techniques.
- Standard techniques for the endoscopic placement of stents should be employed.

MRI SAFETY INFORMATION

This symbol means the stent is MR Safe.

The Geenen Pancreatic Stents are MR Safe. The Zimmon Pancreatic Stents which do not contain radiopaque marker bands (ZEPDF, ZEPDS, ZPSOF, ZPSOS, 4 Fr SPSOF, and 3 Fr SPSOS) are MR Safe. In non-clinical testing, the image artifact caused by the device extends approximately 7 mm from the stent when imaged with a gradient echo pulse sequence and a 3 T MR system.

This symbol means the stent is MR Conditional.

Non-clinical testing has demonstrated that the 5 Fr, 6 Fr, 7 Fr, 8.5 Fr, and 10 Fr Zimmon Pancreatic Stents (SPSOF, SPSOS) which contain radiopaque marker bands are MR conditional. A patient with these devices may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.

Device name	Zimmon Pancreatic Stents which contain radiopaque marker bands (5, 6, 7, 8.5, and 10 Fr SPSOF and SPSOS)
Static magnetic field	1.5 T or 3 T
Maximum spatial field gradient	3000 gauss/cm (30 T/m)
RF excitation	Circularly polarized (CP)
RF transmit coil type	Whole body transmit coil, head RF transmit-receive coil
Maximum whole body specific absorption rate (SAR)	4.0 W/kg (first level operating mode)
Limits on scan duration	4.0 W/kg whole body average SAR for 60 minutes of continuous RF (a sequence or back-to-back series/scan without breaks)
MR image artifact	The presence of this implant may produce an image artifact of 7 mm.

Note: If information about a specific parameter is not included, there are no conditions associated with that parameter.

It is recommended that patients register the conditions under which the implant can safely be scanned with the MedAlert Foundation (www.medicalert.org) or an equivalent organization.

POTENTIAL ADVERSE EVENTS

Those associated with ERCP: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with pancreatic stent placement: fever • infection • obstruction of the common bile duct • pain/discomfort • pancreatitis • perforation • stent migration • stent occlusion • trauma to the pancreatic duct or duodenum.

HOW SUPPLIED

These devices are supplied ethylene oxide (EtO) sterilized in a peel-open pouch.

Following sterilization, the pouch is further packaged into a box (Sof-Flex devices only – GPSO-SF, GPSOS-SF – as these are UV light sensitive).

This device is accompanied by an implant card that should be given to the patient after it has been completed by the healthcare professional.

DEVICE PREPARATION

1. Prior to introducing any devices ensure working channel of endoscope is lubricated with water or a water based lubricant.
2. Referring to the device package label, ensure that the endoscope has a channel size greater than or equal to the minimum channel size required to operate the device.
3. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
4. Prior to using the device(s) visually inspect it/them for abnormalities that could result in improper function.
5. Attach endoscopic cap or wire guide locking device to endoscope.
6. With endoscope in place, introduce a wire guide, floppy tip first and advance until it is fluoroscopically visualized in position through stricture.

INSTRUCTIONS FOR USE

Illustrations

Note: Retain the positioning sleeve for use when introducing the stent duodenal flaps into the accessory channel.

*where applicable

I. FOR NON-PIGTAIL STENTS:

1. Gently ensure the full extension of all flaps.
2. Load positioning sleeve* onto the duodenal end of the stent.

For 7 Fr and smaller stents:

- 3a. Introduce the stent, ductal end first, and the positioning sleeve* onto a pre-positioned wire guide.
- 3b. Advance the pushing catheter over the wire guide to advance the stent into the accessory channel. As the front flap* of the stent enters the accessory channel, keep the positioning sleeve over the back flap*, until the stent is completely in the accessory channel.

For 8.5 Fr and larger stents:

- 3c. Remove the hub from the end of guiding catheter, then introduce the guiding catheter into the accessory channel over a pre-positioned wire guide.
- 3d. With the elevator open, advance the device until a sufficient length of the guiding catheter is above the obstruction.
- 3e. Introduce the stent, ductal end first, and the positioning sleeve onto the guiding catheter and the pre-positioned wire guide.
- 3f. Advance the pushing catheter over the wire guide and the guiding catheter to advance the stent into the accessory channel. As the front flap of stent enters the accessory channel, keep the positioning sleeve over the back flap until the stent is completely in the accessory channel.

THEN REFER TO STEPS 4-9 BELOW

II. FOR PIGTAIL STENTS:

1. Introduce the stent, ductal end first, onto the pre-positioned wire guide.
2. Advance the pushing catheter over the wire guide.
3. Advance the pigtail stent into the accessory channel using the pushing catheter.

THEN REFER TO STEPS 4-9 BELOW

4. As the pushing catheter advances the stent completely into the accessory channel, slide the positioning sleeve* back over the pushing catheter until it reaches the end of the catheter, keeping it clear of the accessory channel.
5. Advance the pushing catheter in 1-2 cm increments until the stent is in the desired position.
6. Fluoroscopically and endoscopically confirm the desired stent position. **Note:** For non pigtail stents of 8.5 Fr or greater only: If desired replace the hub on the guiding catheter before flushing with contrast media.
7. After confirming stent position, gently remove wire guide and the guiding catheter* from endoscope while maintaining position of stent with pushing catheter.
8. Gently remove the pushing catheter from the accessory channel.
9. These stents may be removed using standard endoscopic techniques.

DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be disposed in accordance with institutional guidelines. The introducer (if supplied) should be coiled for disposal.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of.

When available, the EUDAMED website (<https://ec.europa.eu/tools/eudamed/>), along with the BUDI for this product (0827002CIRL202007020003AX), can be used to locate the Summary of Safety and Clinical Performance (SSCP) for this product.

Patient-facing information may be accessed at cookmedical.eu/patient-implant-information.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2026
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.