

- EN 6 **Evolution® Biliary Stent System – Uncovered**
 Instructions for Use
- BG 13 **Система с билиарен стент Evolution® – без покритие**
 Инструкции за употреба
- CS 21 **Systém biliárního stentu Evolution® – nekrytý**
 Návod k použití
- DA 28 **Evolution® galdestentsystem - udækket**
 Brugsanvisning
- NL 35 **Evolution® biliair stentsysteem – niet-bekleed**
 Gebruiksaanwijzing
- ET 42 **Biliaarne stendisüsteem Evolution® – ilma katteta**
 Kasutusjuhend
- FR 49 **Système d'endoprothèse biliaire Evolution® – non couvert**
 Mode d'emploi
- DE 57 **Evolution® Gallengangstentsystem – Nicht gecouvert**
 Gebrauchsanweisung
- EL 64 **Σύστημα ενδοπρόσθεσης χοληφόρων Evolution® - ακάλυπτο**
 Οδηγίες χρήσης
- HU 72 **Evolution® epevezetéksztent-rendszer – fedetlen**
 Használati utasítás
- IT 79 **Sistema con stent biliare Evolution® – non rivestito**
 Istruzioni per l'uso
- LV 86 **Evolution® žultsvadu stenta sistēma – bez pārklājuma**
 Lietošanas norādījumi
- LT 93 **„Evolution®“ tulžies lataku stento sistema – nedengta**
 Naudojimo nurodymai
- NO 100 **Evolution® gallestentsystem – ikke tildekket**
 Bruksanvisning
- PL 107 **Stent do dróg żółciowych Evolution® – niepokryty**
 Instrukcja użycia
- PT 114 **Sistema de stent biliair Evolution® – Não coberto**
 Instruções de utilização
- RO 121 **Sistem pentru stent biliair Evolution® – neacoperit**
 Instrucțiuni de utilizare
- SK 128 **Systém biliárneho stentu Evolution® – nepokrytý**
 Návod na použitie
- ES 135 **Sistema de stent biliair Evolution® – sin recubrimiento**
 Instrucciones de uso
- SV 142 **Evolution® biliärt stentsystem – utan täckning**
 Bruksanvisning





A1


A2


E2

F

EVOLUTION® BILIARY STENT SYSTEM – UNCOVERED

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

IMPORTANT INFORMATION

Please review prior to use.

Warning: The safety and effectiveness of this device for use in the vascular system has not been established.

INTENDED USE

This device is used in palliation of malignant neoplasms in the biliary tree.

STENT DESCRIPTION

This flexible, self-expanding biliary metal stent is constructed of a single, woven, platinum core nitinol radiopaque wire. The stent has a flange at both ends to aid in preventing migration after the uncovered stent has been placed in the bile duct. The stent foreshortens due to its design. The total length of the stent is indicated by radiopaque markers on the inner catheter (*Fig. D, 2 & 4*), indicating the actual length of the stent in a collapsed state. This stent is supplied sterile. This stent is designed for single use only. Attempts to reprocess, resterilize and/or reuse may lead to device failure and/or transmission of disease.

DELIVERY SYSTEM DESCRIPTION

The stent is mounted on an inner catheter, which accepts a .035 inch wire guide, and is constrained by an outer catheter. A pistol-grip delivery handle allows stent deployment or recapture. There are four radiopaque markers to aid in the deployment of the stent while using fluoroscopy. This delivery system is supplied sterile. This delivery system is designed for single use only. Attempts to reprocess, resterilize and/or reuse may lead to device failure and/or transmission of disease.

NOTES

Do not use this device for any other purpose than the stated intended use.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Use of this device is restricted to a trained healthcare professional.

Store in a dry location, away from temperature extremes. Not compatible with THSF wire guide.

CONTRAINDICATIONS

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to: inability to pass the wire guide or stent through the obstructed area, biliary duct strictures of benign etiology, biliary obstruction preventing endoscopic cholangiography, concurrent perforated bile duct, those patients for whom endoscopic procedures are contraindicated, patients with coagulopathy, concurrent bile duct stones, very small intrahepatic ducts and any use other than those specifically outlined under Intended Use.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, cholecystitis, cholestasis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that can occur in conjunction with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum; perforation; obstruction of the pancreatic duct; stent migration; stent occlusion; ingrowth due to tumor or excessive hyperplastic tissue; tumor overgrowth; stent misplacement, pain, fever, nausea, vomiting, inflammation, recurrent obstructive jaundice, bile duct ulceration, death (other than due to normal disease progression).

PRECAUTIONS

Refer to the package label for the minimum channel size required for this device.

A complete diagnostic evaluation should be performed prior to placement to measure the stricture length and determine the proper stent length. The stent length chosen should allow for additional length on either side of the stricture.

Note: In the event a single stent will not adequately cover the stricture, a second stent of the same diameter should be placed providing adequate overlapping (minimum 1cm) of the initially placed stent to ensure a bridging of the stricture between the stents.

If wire guide or stent cannot advance through obstructed area, do not attempt to place stent.

Stent should be placed using fluoroscopic and endoscopic monitoring.

The stent should only be placed with the Cook delivery system, which is provided with each stent.

This stent is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

After stent placement, additional methods of treatment such as chemotherapy and irradiation may increase the risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

Long-term patency with this stent has not been established. Periodic evaluation of the stent is advised. Assessment must be made to determine the necessity of sphincterotomy or balloon dilation prior to stent placement. In the event sphincterotomy or balloon dilation is required, all appropriate cautions, warnings and contraindications must be observed.

WARNINGS

MRI INFORMATION

This symbol means the device is MR conditional

Non-clinical testing has demonstrated that the Evolution® Biliary Stent is MR Conditional according to ASTM F2503. A patient with this stent can be scanned safely anytime after placement under the following conditions.

Static Magnetic Field

- Static magnetic field of 3 Tesla or less
- Highest spatial magnetic gradient field of 1600 Gauss/cm or less

Non-clinical evaluation was conducted in a MR system (Excite, General Electric Healthcare) with a maximum spatial magnetic gradient field of 1600 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

- 1.5 and 3.0 Tesla Systems: It is recommended to scan in normal operation mode (whole body averaged specific absorption rate (SAR) = 2.0W/kg) ("Normal Operating Mode" is defined as the mode of operation of the MR system in which none of the outputs have a value that cause physiological stress to the patient) for 15 minutes of scanning (i.e., per scanning sequence).

1.5 Tesla Temperature Rise

In non-clinical testing, the Evolution® Biliary Stent produced maximum temperature rises of 2.1 and 2.6 °C (for a single stent and pair of overlapping stents) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 1.5 Tesla System (Siemens Medical Solutions, Malvern, PA, Software Numaris/4) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.1 W/kg).

3.0 Tesla Temperature Rise

In non-clinical testing, the Evolution® Biliary Stent produced maximum temperature rises of 3.1 and 4.2 °C (for a single stent and pair of overlapping stents) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3 Tesla System (Excite, GE Electric Healthcare, Milwaukee, WI, Software 14X.M5) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.7 W/kg).

Image Artifacts

MR image quality may be compromised if the area of interest is within the lumen or within approximately 7 mm of the position of the Evolution Biliary Stent, as found during non-clinical testing using the sequence: T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic stent.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail:	MedicAlert Foundation International 2323 Colorado Avenue Turlock, CA 95382
Phone:	888-633-4298 (toll free) 209-668-3333 from outside the US
Fax:	209-669-2450
Web:	www.medicalert.org

This stent is **not intended to be removed** and is considered a permanent implant. Attempts to remove the stent after placement may cause damage to the surrounding mucosa.

Stent cannot be recaptured after the deployment threshold has been passed. Corresponding marks on outer catheter and delivery handle indicate when threshold has been passed.

This device is not intended to be deployed through the wall of a previously placed or existing metal stent. Doing so could result in difficulty or inability to remove introducer.

The safety and effectiveness of this device for use in the vascular system has not been established.

The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

EQUIPMENT REQUIRED

Stent and Delivery System

.035 inch wire guide

Fusion Wire Lock

IMPORTANT: This stent system is intended for single use only. Do not attempt to reload a deployed stent.

PROCEDURE PREPARATION

With endoscope in place, introduce a wire guide, floppy tip first and advance until it is fluoroscopically visualized in position through stricture.

INSTRUCTIONS FOR USE

Illustrations

1. Remove stent delivery system from package and backload device over a prepositioned wire guide, ensuring that the wire guide exits catheter at Zip port (*Fig. A1, A2*).

2. Unlock wire guide from Fusion Wire Guide Locking Device.

3. Introduce device in short increments into accessory channel of endoscope until Zip port is inside of the accessory channel, then relock wire guide. Continue advancing device in short increments.

4. With elevator open, advance device until endoscopically visualized exiting endoscope.

5. For transpapillary placement reference the yellow marker for accurate stent positioning. **Note:** Stents bridging papilla should extend beyond papilla and into duodenum approximately 0.5 cm after deployment.

6. Under fluoroscopic guidance, with elevator open, continue to advance the device in short increments until the stent is fluoroscopically visualized through the stricture. Fluoroscopically visualize radiopaque markers on inner catheter at either end of stent and position the inner radiopaque markers (*Fig. D, 2 and D, 4*) a minimum of 1 cm beyond the stricture so that it bridges stricture completely.

7. Confirm desired stent position fluoroscopically. To deploy stent, remove red safety guard from handle, (*Fig. B*) then squeeze the trigger. **Note:** Each full trigger squeeze will deploy stent by an equal amount.

8. If stent repositioning is required during deployment, it is possible to recapture stent. **Note:** It is **not** possible to recapture stent after passing point-of-no-return, indicated when red marker on top of handle has passed the point-of-no-return indicator on handle label (*Fig. C*).

To fluoroscopically monitor the point-of-no-return (*Fig. D, 5*), the stent can be recaptured until the point where the outer catheter radiopaque marker (*Fig. D, 1*) is aligned with the second inner catheter radiopaque marker (*Fig D, 3*).

Note: It is **not** possible to recapture stent after the outer catheter radiopaque marker (*Fig. D, 1*) passes the second inner catheter radiopaque marker (*Fig D, 3*).

Note: The stent can be recaptured a maximum of three times.

9. To reposition the stent, the stent must first be recaptured and the elevator must be open. **Note:** Do not push forward on the delivery system with stent partially deployed. Push direction button on side of delivery handle to opposite side (*Fig. E1*). **Note:** Hold thumb on button when squeezing trigger for first time to recapture. Continue squeezing trigger as required to recapture stent by desired amount.

10. To resume deployment, push button to opposite side again and hold button for first stroke while squeezing trigger (*Fig. E2*).

11. When stent point-of-no-return has been passed, disconnect luer lock fitting and remove safety wire completely from delivery handle (*Fig. F*).

12. Continue deploying stent by squeezing trigger.

13. After deployment, fluoroscopically confirm full stent expansion. While maintaining wire guide position, push direction button on side of delivery handle to opposite side (*Fig. E1*). Squeeze the trigger to completely recapture the introduction system. **Unlock the wire guide from Fusion Wire Guide**

locking device. The device can be safely removed with the elevator of the endoscope fully down.

If placement of a second stent is required, it can be placed by following steps 1 to 13. **Caution;** passing a delivery system through a just deployed stent may cause a previously placed stent to dislodge. To reduce the risk of stent dislodgement ensure that the introduction system is completely recaptured following stent deployment.

Warning: This metal biliary stent is **not intended to be repositioned or removed after stent placement** and is considered a permanent implant. Attempts to remove the stent after placement may cause damage to the surrounding mucosa. In case of accidental deployment or improper placement (immediately following deployment), stent should be left in place and placement of a second stent of the same diameter should be attempted to achieve desired result providing adequate overlapping (minimum 1cm).

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2016 ©
Cook Ireland Ltd.