

- EN 7 **Evolution® Biliary Stent System – Partially Covered**
 Instructions for Use
- BG 13 **Система с билиарен стент Evolution® – с частично покритие**
 Инструкции за употреба
- HR 20 **Evolution® sustav bilijarnog stenta – djelomično pokriveni**
 Upute za uporabu
- CS 26 **Systém biliárního stentu Evolution® – částečně krytý**
 Návod k použití
- DA 33 **Evolution® galdestentsystem – delvist dækket**
 Brugsanvisning
- NL 40 **Evolution® biliaire stentsysteem – gedeeltelijk bekleed**
 Gebruiksaanwijzing
- ET 46 **Biliaarne stendisüsteem Evolution® – osaliselt kaetud**
 Kasutusjuhised
- FR 53 **Système d'endoprothèse biliaire Evolution® – partiellement couverte**
 Mode d'emploi
- DE 60 **Evolution® Gallengangstentsystem – Teilweise gecouvert**
 Gebrauchsanweisung
- EL 67 **Σύστημα ενδοπρόσθεσης χοληφόρων Evolution® - μερικώς καλυμμένο**
 Οδηγίες χρήσης
- HU 75 **Evolution® epeúti sztentrendszer – részlegesen fedett**
 Használati utasítás
- IT 81 **Sistema con stent biliare Evolution® – parzialmente rivestito**
 Istruzioni per l'uso
- LV 88 **Evolution® žultsvadu stenta sistēma — ar daļēju pārklājumu**
 Lietošanas instrukcija
- LT 95 **„Evolution®“ tulžies lataų stento sistema – iš dalies dengta**
 Naudojimo nurodymai
- NO 101 **Evolution® gallestentsystem – delvis tildekket**
 Bruksanvisning
- PL 107 **Stent do dróg żółciowych Evolution® — częściowo pokryty**
 Instrukcja użycia
- PT 114 **Sistema de stent biliar Evolution® — parcialmente coberto**
 Instruções de utilização
- RO 121 **Sistem pentru stent biliar Evolution® – acoperit parțial**
 Instrucțiuni de utilizare
- SK 128 **Systém biliárneho stentu Evolution® – čiastočne pokrytý**
 Návod na použitie
- ES 134 **Sistema de stent biliar Evolution® – parcialmente recubierto**
 Instrucciones de uso
- SV 141 **Evolution® biliärt stentsystem – delvis täckning**
 Bruksanvisning




A1


A2



E2

F

G

EVOLUTION® BILIARY STENT SYSTEM – PARTIALLY COVERED

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

IMPORTANT INFORMATION

Please review prior to use.

Warning: The safety and effectiveness of this device for use in the vascular system has not been established.

INTENDED USE

This device is used in palliation of malignant neoplasms in the biliary tree.

STENT DESCRIPTION

This flexible, self-expanding biliary metal stent is constructed of a single, woven, platinum core nitinol radiopaque wire. The partially covered stent has the potential to reduce tumor ingrowth through the stent. The stent has a flange at both ends to aid in preventing migration after the stent has been placed in the bile duct. There is a yellow grasping loop at the end of the stent, the purpose is to reposition distally/remove the stent during initial stent placement procedure in the event of incorrect placement. The stent foreshortens due to its design. The total length of the stent is indicated by radiopaque markers on the inner catheter (**Fig. D, 2 & 4**), indicating the actual length of the stent in a collapsed state. The proximal and distal ends of the stent are designed with numerous rounded “crowns” to potentially reduce tissue trauma.

DELIVERY SYSTEM DESCRIPTION

The stent is mounted on an inner catheter, which accepts a .035-inch wire guide, and is constrained by an outer catheter. A pistol-grip delivery handle allows controlled stent deployment or recapture. The position of the red marker on top of the handle in relation to the “point-of-no-return” position on the handle label indicates the point after which it is not possible to recapture the deployed stent. There are four radiopaque markers to aid in the deployment of the stent while using fluoroscopy, and a yellow endoscopic marker for transpapillary placement.

CONTRAINDICATIONS

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to: inability to pass the wire guide or stent through the obstructed area • biliary duct strictures of benign etiology • biliary obstruction preventing endoscopic cholangiography • concurrent perforated bile duct • those patients for whom endoscopic

procedures are contraindicated • patients with coagulopathy, concurrent bile duct stones, very small intrahepatic ducts and any use other than those specifically outlined under Intended Use.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with ERCP include, but are not limited to: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • cholecystitis • cholestasis • hemorrhage • hypotension • infection • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Additional adverse events that can occur in conjunction with biliary stent placement include, but are not limited to: allergic reaction to nickel • bile duct ulceration • death (other than due to normal disease progression) • fever • inflammation • ingrowth due to tumor or excessive hyperplastic tissue • nausea • obstruction of the pancreatic duct • pain/discomfort • perforation • recurrent obstructive jaundice • stent migration • stent misplacement • stent occlusion • trauma to the biliary tract or duodenum • tumor overgrowth • vomiting.

PRECAUTIONS

Refer to the package label for the minimum channel size required for this device.

A complete diagnostic evaluation should be performed prior to placement to measure the stricture length and determine the proper stent length. The stent length chosen should allow for additional length on either side of the stricture.

If wire guide or stent cannot advance through obstructed area, do not attempt to place stent.

Stent should be placed using fluoroscopic and endoscopic monitoring.

The stent should only be placed with the Cook delivery system, which is provided with each stent.

This stent is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

Long-term patency with this stent has not been established. Periodic evaluation of the stent is advised. Assessment must be made to determine the necessity of sphincterotomy or balloon dilation prior to stent placement. In the event sphincterotomy or balloon dilation is required, all appropriate cautions, warnings and contraindications must be observed.

The stent cannot be recaptured after the deployment threshold has been passed. The position of the red marker on the top of the handle in relation to point-of-no-return position on the handle label indicates when threshold has been passed.

Do not attempt to reload a deployed stent.

Do not use this device for any other purpose than the stated intended use.

Not compatible with THSF wire guide.

Use of this device is restricted to a trained healthcare professional.

Passing a second stent delivery system through a just deployed stent is not recommended and could cause the stent to dislodge. To reduce the risk of stent dislodgement ensure that the delivery system is completely recaptured following stent deployment.

WARNINGS

The safety and effectiveness of this device for use in the vascular system has not been established.

Placement of a partially covered stent across a branch duct or major bifurcation may result in obstruction.

This stent is not intended to be removed after initial stent placement and is considered a permanent implant. Attempts to remove the stent after placement may cause damage to the surrounding mucosa.

This device is not intended to be deployed through the wall of a previously placed or existing metal stent. Doing so could result in difficulty or inability to remove introducer.

The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

This stent system is designed for single use only. Attempts to reprocess, resterilize and/or reuse may lead to device failure and/or transmission of disease.

After stent placement, additional methods of treatment such as chemotherapy and irradiation may increase the risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Testing of overlapping stents has not been completed and is not recommended.

MRI INFORMATION

This symbol means the device is MR conditional

Non-clinical testing has demonstrated that the Evolution® Biliary Stent is MR Conditional according to ASTM F2503. A patient with this stent can be scanned safely anytime after placement under the following conditions.

Static Magnetic Field

- Static magnetic field of 3 Tesla or less
- Highest spatial magnetic gradient field of 1600 Gauss/cm or less

Non-clinical evaluation was conducted in a MR system (Excite, General Electric Healthcare) with a maximum spatial magnetic gradient field of 1600 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 and 3.0 Tesla Systems: It is recommended to scan in normal operation mode (whole body averaged specific absorption rate (SAR) = 2.0W/kg) ("Normal Operating Mode" is defined as the mode of operation of the MR system in which none of the outputs have a value that cause physiological stress to the patient) for 15 minutes of scanning (i.e., per scanning sequence).

1.5 Tesla Temperature Rise

In non-clinical testing, the Evolution® Biliary Stent produced a maximum temperature rise of 2.1°C (for a single stent) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 1.5 Tesla System (Siemens Medical Solutions, Malvern, PA, Software Numaris/4) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.1 W/kg).

3.0 Tesla Temperature Rise

In non-clinical testing, the Evolution® Biliary Stent produced a maximum temperature rise of 3.1°C (for a single stent) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3 Tesla System (Excite, GE Electric Healthcare, Milwaukee, WI, Software 14X.M5) at

an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.7 W/kg).

Image Artifacts

MR image quality may be compromised if the area of interest is within the lumen or within approximately 7 mm of the position of the Evolution Biliary Stent, as found during non-clinical testing using the sequence: T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic stent.

It is recommended that patients register the conditions under which the implant can safely be scanned with the MedicalAlert Foundation (www.medicalert.org) or an equivalent organization.

EQUIPMENT REQUIRED

Stent and Delivery System
.035 inch wire guide
Fusion Wire Lock

PROCEDURE PREPARATION

With duodenoscope in place, introduce a wire guide, floppy tip first and advance until it is fluoroscopically visualized in position through stricture.

INSTRUCTIONS FOR USE

Illustrations

1. Remove stent delivery system from package and backload device over a prepositioned wire guide, ensuring that the wire guide exits catheter at Zip port (**Fig. A1, A2**).
2. Unlock wire guide from Fusion Wire Guide Locking Device.
3. Introduce device in short increments into accessory channel of duodenoscope until Zip port is inside of the accessory channel, then relock wire guide. Continue advancing device in short increments.
4. With elevator open, advance device until endoscopically visualized exiting duodenoscope.
5. For transpapillary placement reference the yellow marker for accurate stent positioning. **Note:** Stents bridging papilla should extend beyond papilla and into duodenum approximately 0.5 cm after deployment.
6. Under fluoroscopic guidance, with elevator open, continue to advance the device in short increments until the stent is fluoroscopically visualized through the stricture. Fluoroscopically visualize radiopaque markers on inner catheter at either end of stent and position the inner radiopaque markers (**Fig. D,2 and D,4**) a minimum of 1 cm beyond the stricture so that it bridges stricture completely.

7. Confirm desired stent position fluoroscopically. To deploy stent, remove red safety guard from handle, **(Fig. B)** then squeeze the trigger. **Note:** Each full trigger squeeze will deploy stent by an equal amount.

8. If stent repositioning is required during deployment, it is possible to recapture stent. **Note:** It is **not** possible to recapture stent after passing point-of-no-return, indicated when red marker on top of handle has passed the point-of-no-return position on the handle label **(Fig. C)**.

To fluoroscopically monitor the point-of-no-return **(Fig. D, 5)**, the stent can be recaptured until the point where the outer catheter radiopaque marker **(Fig. D, 1)** is aligned with the second inner catheter radiopaque marker **(Fig. D, 3)**.

Note: It is **not** possible to recapture stent after the outer radiopaque marker **(Fig. D, 1)** passes the second inner radiopaque marker **(Fig. D, 3)**.

Note: The stent can be recaptured a maximum of three times.

9. To reposition the stent, the stent must first be recaptured and the elevator must be open. **Note:** Do not push forward on the delivery system with stent partially deployed. Push direction button on side of delivery handle to opposite side **(Fig. E1)**. **Note:** Hold thumb on button when squeezing trigger for first time to recapture. Continue squeezing trigger as required to recapture stent by desired amount.

10. To resume deployment, push button to opposite side again and hold button for first stroke while squeezing trigger **(Fig. E2)**.

11. When stent point-of-no-return has been passed, disconnect luer lock fitting to remove safety wire completely from delivery handle **(Fig. F)**.

12. Continue deploying stent by squeezing trigger.

13. After deployment, fluoroscopically confirm full stent expansion. While maintaining wire guide position, push direction button on side of delivery handle to opposite side **(Fig. E1)**. Squeeze the trigger to completely recapture the delivery system. **Unlock the wire guide from Fusion Wire Guide locking device.** The device can be safely removed with the elevator of the duodenoscope fully down.

To reposition distally/remove stent during initial stent placement procedure, grasp the yellow grasping loop at end of stent with forceps **(Fig. G)** and reposition distally/remove as required.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Store in a dry location, away from temperature extremes.

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© COPYRIGHT COOK 2019 ©
Cook Ireland Ltd