
EN 9	Zilver Flex™35 Vascular Stent Instructions for Use
BG 14	Васкуларен стент Zilver Flex™35 Инструкции за употреба
CS 20	Cévní stent Zilver Flex™35 Návod k použití
DA 26	Zilver Flex™35 vaskulær stent Brugsanvisning
NL 31	Zilver Flex™35 vasculaire stent Gebruiksaanwijzing
ET 37	Vaskulaarne stent Zilver Flex™35 Kasutusjuhised
FR 42	Endoprothèse vasculaire Zilver Flex™35 Mode d'emploi
DE 48	Zilver Flex™35 Gefäßstent Gebrauchsanweisung
EL 54	Αγγειακή ενδοπρόσθεση Zilver Flex™35 Οδηγίες χρήσης
HU 61	Zilver Flex™35 vaszkuláris sztent Használati utasítás
IT 67	Stent vascolare Zilver Flex™35 Istruzioni per l'uso
LV 73	Zilver Flex™35 vaskulārais stents Lietošanas instrukcija
LT 79	„Zilver Flex™35“ kraujagyslių stentas Naudojimo nurodymai
NO 85	Zilver Flex™35 vaskulær stent Bruksanvisning
PL 90	Stent naczyniowy Zilver Flex™35 Instrukcja użycia
PT 96	Stent vascular Zilver Flex™35 Instruções de utilização
RO 102	Stent vascular Zilver Flex™35 Instrucțiuni de utilizare
SK 108	Vaskulárny stent Zilver Flex™35 Návod na použitie
ES 114	Stent vascular Zilver Flex™35 Instrucciones de uso
SV 120	Zilver Flex™35 vaskulär stent Bruksanvisning

a. Handle
b. Hub
c. Safety Lock
d. Delivery System: Outer Sheath
e. Tip of Delivery System Inner Catheter
f. Side-arm Flushing Port
g. Metal Cannula
h. Radiopaque Marker on the Delivery System
i. Gold Radiopaque Markers

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

Fig. 5

ZILVER FLEX™35 VASCULAR STENT

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Do not re-sterilize.

DEVICE DESCRIPTION

The ZILVER FLEX™35 Vascular Stent (ZFV) is a self-expandable stent made of nitinol. The special laser-cut pattern of the nitinol tube provides a construction with a strong radial force and high flexibility.

The stent has radiopaque markers at both ends (i).

The stent is supplied preloaded in a 6.0 French (2.0mm) introducer catheter. A radiopaque marker (h) on the distal tip of the sheath is used to visualize deployment of the stent. Stent deployment is controlled by means of a hand-held device.

Refer to Stent Size Selection Table below.

Stent Length (mm)		20		30		40		60		80	
Delivery System (cm)		80	125	80	125	80	125	80	125	80	125
Stent Outer Diameter (mm)	5	x	x	x	x	x	x	x	x	x	x
	6	x	x	x	x	x	x	x	x	x	x
	7	x	x	x	x	x	x	x	x	x	x
	8	x	x	x	x	x	x	x	x	x	x
	9	x	x	x	x	x	x	x	x	x	x
	10	x	x	x	x	x	x	x	x	x	x

Stent Length (mm)		100		120		140		170		200	
Delivery System (cm)		80	125	80	125	80	125	80	125	80	125
Stent Outer Diameter (mm)	5	x	x	x	x	x	x	x	x	x	x
	6	x	x	x	x	x	x	x	x	x	x
	7	x	x	x	x	x	x	x	x	x	x
	8	x	x	x	x	x	x	x	x	x	x
	9	x	x	x	x	x	x				
	10	x	x	x	x	x	x				

INTENDED USE

The product is intended for use in the iliac, superficial femoral artery (SFA) and above-the-knee popliteal artery for the following treatments:

- Arteriosclerotic stenosis
- Total occlusions that have been recanalized

The product provides mechanical support to maintain constant blood flow of the vessel.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques.

Placement of the stent requires advanced skills in interventional procedures. The following instructions will give technical guidance but do not obviate formal training in the use of the device.

Standard techniques for placement of an arterial access sheath and wire guide should be employed during use of the vascular stent.

CONTRAINDICATIONS

Patients with bleeding diathesis or other disorder, e.g. peptic ulceration or recent cerebrovascular accident, limiting the use of antiplatelet and anticoagulant therapy.

Patients with total arterial occlusion.

Coronary use is not recommended.

WARNINGS

- Possible allergic reactions should be avoided.

PRECAUTIONS

- Manipulation of the product requires high-resolution fluoroscopic control.
- The lesion and vessel can be predilated before placement of the stent.
- Do not try to remove the stent from the introducer system before use.
- Ensure that the safety lock is not inadvertently removed prior to stent release.

MRI INFORMATION

Non-clinical testing has demonstrated that the ZILVER FLEX™35 Vascular Stent is MR conditional. A patient with this stent can be scanned safely after placement under the following conditions.

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Highest spatial magnetic gradient of 720-gauss/cm

Non-clinical evaluation was conducted in a 3-Tesla MR system with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gauss meter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual) and resulted in a deflection angle of 3 degrees. An implant that displays a deflection angle of ≤ 45 degrees is considered to pose no harm to a patient or individual.

MRI-Related Heating

1.5 and 3.0 Tesla Systems

- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.9 W/kg for 15 minutes of scanning (i.e., per pulse sequence)

Non-clinical evaluation, including analysis of the effects of blood flow, was conducted for single and two overlapped versions of the ZILVER FLEX™35 Vascular Stents (lengths = 140 and 200 mm) producing a maximum temperature change of 3.1°C during 15 minutes of MRI performed at 1.5-Tesla (Magnetom, Siemens Medical Solutions, Malvern, PA) and 3-Tesla (Excite, General Electric Healthcare, Milwaukee, WI) under the following minimum conditions:

140 mm and 200 mm lengths:

	<u>1.5-Tesla</u>	<u>3-Tesla</u>
MR system reported, whole body averaged SAR	2.9-W/kg	3.0-W/kg
Calorimetry measured values, whole body averaged SAR	2.1-W/kg	2.8-W/kg

Image Artifacts

MR image quality may be compromised if the area of interest is in the same area or within approximately 5 mm of the position of the ZILVER FLEX™35 Vascular Stent, as found during non-clinical testing using the sequence: T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant. Heating in the MRI environment for stents with fractured struts is not known.

POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur include, but are not limited to:

- Abrupt stent closure
- Allergic reaction to nitinol
- Amputation
- Angina/coronary ischemia
- Arrhythmia
- Arterial aneurysm
- Arterial rupture
- Arteriovenous fistula
- Atheroembolization (Blue Toe Syndrome)
- Death
- Embolism
- Fever
- Hematoma/hemorrhage
- Hypersensitivity reactions
- Hypotension/hypertension
- Infection/abscess formation at access site
- Intimal injury/dissection
- Ischemia requiring intervention (bypass or amputation of toe, foot or leg)
- Myocardial infarction
- Pseudoaneurysm formation
- Pulmonary embolism
- Renal failure
- Restenosis of the stented artery
- Septicemia/bacteraemia
- Stent malapposition
- Stent migration
- Stent strut fracture
- Stroke

- Spasm
- Tissue necrosis
- Worsened claudication/rest pain

PRODUCT RECOMMENDATIONS

Arterial access

For arterial access, we recommend the use of an access set that accepts a 6.0 French (2.0mm) introducer catheter.

Wire Guide Selection

The use of a .035-inch (0.89mm) wire guide is recommended.

PTA Balloon Selection

For pre- and post-dilation, an appropriately sized balloon catheter is recommended.

Stent Selection

The chosen stent diameter should be at least 1 mm larger than the diameter of the reference vessel.

PROCEDURE

Illustrations

Access

1. Gain access at the appropriate site utilizing an introducer sheath or guiding catheter. Refer to product labeling for required introducer sheath or guiding catheter size.
2. Introduce the extra or ultra stiff .035" (0.89mm) wire guide through the introducer sheath or guiding catheter across the distal segment of the target lesion.
3. Predilate if necessary. Remove the balloon catheter, leaving the wire guide in place.
4. Determine the proper stent size after complete diagnostic evaluation. The stent diameter should be at least 1 mm larger than the diameter of the reference vessel. Measure the length of the target lesion to determine the length of the stent required. Allow for the proximal and distal aspects of the stent to cover the entire target area. **NOTE:** If multiple stents are required to cover the length of the lesion, please refer to the Multiple Stent Placement section of these instructions for use for further recommendations.
5. Use the 1 ml syringe to flush the wire guide lumen with saline through the hub.
6. Immediately before placing the delivery system into the body, use a 1 ml syringe to flush the delivery system with saline through the side-arm flushing port. Flush until a few drops of saline exit the distal tip, between the delivery system outer sheath (d) and inner catheter (e).
7. Insert the delivery system over the wire guide.
8. Under fluoroscopy, advance the delivery system past the lesion.
9. Pull back on the stent delivery system under fluoroscopy until the radiopaque markers on the stent (i) are at the desired position. The stent is now ready to be deployed. **(Fig. 1)**

Deployment of the Stent

1. Before deployment, it is important to straighten the proximal part of the delivery system as much as possible and to keep the handle in a stable position.
2. The stent expansion must be performed under fluoroscopic control.
3. Hold the hub (b) on the metal cannula (g) steady. To deploy the stent, remove the red safety lock (c). **(Fig. 2)**

4. Hold the hub end stationary. The stent will be deployed as you pull the handle (a) toward the hub (b). (**Fig. 3**) The radiopaque marker on the delivery system (h) indicates the progress of the deployment. **NOTE:** Full deployment of the stent length will occur when the distal end of the sheath has been retracted past the proximal part of the stent.
5. As deployment occurs, continue sliding the handle (a) toward the hub (b) in a slow, smooth and consistent fashion. (**Fig. 4**) **NOTE: Once stent deployment has begun, the stent must be fully deployed. Repositioning of the ZILVER FLEX™35 Vascular Stent is not possible since the delivery system's outer sheath cannot be readvanced over the stent once deployment begins. Refer to the Multiple Stent Placement section of these instructions for use for information on missed lesions.**
6. The stent is fully deployed when the handle (a) reaches the hub (b). (**Fig. 5**)
7. Perform an arterial angiogram to ensure proper placement of the device. **NOTE:** If incomplete expansion exists within the stent at any point along the lesion, post-deployment balloon dilatation (standard PTA) can be performed at the discretion of the physician. Select an appropriate size PTA balloon catheter and dilate the lesion with conventional technique. The inflation diameter of the PTA balloon used for post dilatation should approximate the diameter of the reference vessel. Remove the PTA balloon from the patient.
8. Delivery System Removal - Do not advance sheath after the stent has been deployed. **Delivery System can be removed without recapturing tip.** Check for delivery system integrity post removal from the patient.
9. Introducer sheath and wire guide may be removed at this point.

Multiple Stent Placement

If placement of more than one stent is required, the following recommendations should be considered:

- When more than one stent is required, resulting in stent-to-stent contact, stent materials should be of similar composition to avoid the possibility of dissimilar metal corrosion.
- In relation to the lesion site, the distal area of narrowing should be stented first, followed by the proximal locations (i.e., a second stent should be placed proximal to the previously placed stent).
- Stents placed in tandem should slightly overlap.

HOW SUPPLIED

Supplied sterile by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

Upon completion of procedure, dispose of components per institutional guidelines for biohazardous medical waste.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COPYRIGHT COOK 2017 ©
COOK IRELAND LTD.