EchoTip Ultra Endobronchial HD Ultrasound Needle for Pentax Scopes

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English

INTENDED USE

This device is used to sample targeted submucosal and extramural lesions within or adjacent to the tracheobronchial tree through the accessory channel of an ultrasound endoscope for Fine Needle Aspiration (FNA).

NOTES

Do not use this device for any purpose other than stated intended use.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify COOK Medical for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device is restricted to a trained healthcare professional.

This device is supplied sterile and intended for single use only.

This device is intended for use with an Pentax EBUS scope.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications include, but are not limited to: coagulopathy.

POTENTIAL COMPLICATIONS

Potential complications associated with bronchial endoscopy include, but are not limited to: Perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest and tumor seeding.

WARNING

Not for use in the heart or vascular system. The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the endoscope.

SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of the syringe has a Luer slip fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (See Fig. 1)
2. Prepare syringe as follows:
 - a) With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b) Turn stopcock to "closed" position.
 - c) Pull plunger back until it is locked into place at desired setting, creating a suction.
 - d) Set prepared syringe aside until aspiration biopsy is desired.
3. To adjust device sheath length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the sheath length is at "0". **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window. (See Fig. 2) Tighten thumbscrew on the sliding sheath adjuster to maintain the sheath length.
4. Lock thumbscrew on safety ring until "0" appears in safety ring window.

INSTRUCTIONS FOR USE

Illustrations

1. Identify the desired targeted site by endoscopic ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked to hold the needle in place, introduce the ultrasound needle into the accessory channel of the endoscope. **Note:** If resistance is encountered on needle introduction, reduce angulation of scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
4. Attach the device to the endoscopes accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on

the endoscopic view, confirming that the sheath has emerged from the working channel of the scope.

6. With the endoscope and device straight, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for the needle advancement appears in the window of the safety ring. (See Fig. 3) Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the endoscope. Failure to attach the device prior to needle adjustment or extension may result in damage to the endoscope.

7. While maintaining the position of the endoscope, extend the needle by advancing the needle handle of the device to the pre-positioned safety ring.

8. Advance the needle into the lesion.

9. Remove the stylet of the ultrasound needle by gently pulling back on the plastic hub seated in the metal fitting of the needle handle. Preserve the stylet for use if additional cell collection is desired later.

10. Attach the Luer lock fitting of the previously prepared syringe securely into the fitting on the needle handle.

11. Turn the stopcock to the "open" position aligned with the syringe, allowing the negative pressure in the syringe to aspirate cells. Gently move the needle handle, in small increments, back and forth within the biopsy site. **Note:** Do not remove the needle from the biopsy site during cell collection.

12. Upon completion of cell collection, retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring to hold the needle in place.

13. Unlock the thumbscrew of sheath adjuster and remove needle while keeping sheath adjuster attached to endoscope. (See Fig. 4)

14. Exchange the locking syringe for an empty syringe and attach the new syringe to the Luer lock fitting on the needle handle. Unlock the thumbscrew on the safety ring and extend the needle.

15. Push in the syringe plunger to expel aspirate specimen, then prepare per institutional guidelines.

16. Remaining aspirate may be recovered for examination by reinserting stylet or flushing device.

17. For additional cell collection from the same lesion, gently reinsert the stylet into the metal fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water and ensure that the needle is fully retracted into sheath and locking ring is locked in place.

While supporting the stylet at the Luer lock fitting, advance it in small increments, until the hub is engaged in the fitting.

18. Readvance the sheath through the attached sheath adjuster, ensuring that both thumbscrews are aligned on same side.

19. Additional samples may be obtained by repeating Step 5 through 16 of "Instructions for Use".

20. Upon completion of procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counter clockwise, and withdraw the entire device from the endoscope.

Upon completion of the procedure, dispose of the device, per institutional guidelines for biohazardous medical waste.

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Cook is a registered trademark of Cook Incorporated.
Cook Ireland Ltd.