
EN 7	EchoTip® Ultra Endobronchial HD Ultrasound Needle for Pentax Scopes Instructions for Use
BG 10	Ендобронхиална ултразвукова игла с висока разделителна способност EchoTip® Ultra за ендоскопи Pentax Инструкции за употреба
CS 14	Endobronchiální jehla EchoTip® Ultra pro ultrazvuk s vysokým rozlišením, pro endoskopy Pentax Návod k použití
DA 17	EchoTip® Ultra endobronkial HD ultralydsnål til Pentax skoper Brugsanvisning
DE 20	Endobronchiaale kõrglahutusega ultrahelinoel EchoTip® Ultra endoskoopidele Pentax Kasutusjuhised
EL 23	Ενδοβρογχική βελόνα υπερήχων υψηλής ανάλυσης EchoTip® Ultra για ενδοσκόπια Pentax Οδηγίες χρήσης
ES 27	Aguja de ecografía endobronquial de alta definición EchoTip® Ultra para endoscopios Pentax Instrucciones de uso
ET 30	Aiguille ultrasonore HD endobronchique EchoTip® Ultra pour les endoscopes Pentax Mode d'emploi
FR 33	EchoTip® Ultra endobronchiale HD-Ultraschallnadel für Pentax-Endoskope Gebrauchsanweisung
HR 36	Endobronhijalna HD ultrazvučna igla EchoTip® Ultra za endoskope Pentax Upute za uporabu
HU 39	EchoTip® Ultra endobronchialis, nagy felbontású, echogén tű Pentax endoszkópokhoz Használati utasítás
IT 42	Ago ecografico endobronchiale ad alta definizione EchoTip® Ultra per endoscopi Pentax Istruzioni per l'uso
LT 45	„EchoTip® Ultra“ endobronchinio didelės skiriamosios gebos ultragarso adata, skirta naudoti su „Pentax“ endoskopais Naudojimo nurodymai
LV 48	EchoTip® Ultra endobronhiālā augstas precizitātes ultraskaņas adata Pentax endoskopiem Lietošanas instrukcija
NL 51	EchoTip® Ultra endobronchiale HD ultrasone naald voor Pentax endoscopen Gebruiksaanwijzing
NO 54	EchoTip® Ultra endobronkial HD-ultralýdnål for Pentax-skop Bruksanvisning
PL 57	Igła EchoTip® Ultra do ultrasonografii wewnątrzoskrzelowej o wysokiej rozdzielczości do endoskopów Pentax Instrukcja użycia
PT 60	Aguļa de ultrassons endobrônquica EchoTip® Ultra de alta definição para endoscópios Pentax Instruções de utilização
RO 63	Ac ecografic HD endobronşic EchoTip® Ultra pentru endoscoape Pentax Instrucţiuni de utilizare
SK 66	Endobronchiálna ultrazvuková ihla s vysokým rozlišením EchoTip® Ultra pre endoskopy Pentax Návod na použitie
SV 69	EchoTip® Ultra endobronkial HD ultraljudsnål för Pentax endoskop Bruksanvisning
TR 72	Pentax Skopları için EchoTip® Ultra Endobronşiyal HD Ultrason İğnesi Kullanma Talimatı


EN Plunger

BG Бутало
 CS Pist
 DA Stempel
 DE Kolben
 EL Έμβολο
 ES Émbolo
 ET Kolb
 FR Piston
 HR Klip
 HU Dugattyú
 IT Stantuffo
 LT Stūmoklis
 LV Virzulis
 NL Plunjer
 NO Stempel
 PL Tłoczek
 PT Êmbolo
 RO Piston
 SK Piest
 SV Kolv
 TR Piston

EN Stopcock "open"

BG Кранче "отворено"
 CS „Otevřený“ uzavírací kohout
 DA Stophane "åben"
 DE Absperrhahn geöffnet
 EL Στρόφιγγα ανοικτή
 ES Llave de paso «abierta»
 ET Sulgemisventiil avatud
 FR Robinet ouvert
 HR Sigurnosni ventil otvoren
 HU Elzárócsap nyitott állásban
 IT Rubinetto aperto
 LT Kramelis „atviras“
 LV Slēgkrāns atvērts
 NL Afsluitkraan 'open'
 NO Stoppekran åpen
 PL Kranik w pozycji otwartej
 PT Torneira de passagem "aberta"
 RO Robinet „deschis“
 SK Otvorený kohútik
 SV Kran "öppen"
 TR Stopkok "açık"

Fig. 2**EN Sliding Sheath Adjuster**

BG	Плъзгач се регулатор на защитната обвивка
CS	Posuvné nastavovací zařízení sheathu
DA	Glidende sheath-justering
DE	Verschiebbarer Schleuseneinsteller
EL	Συρόμενος προσαρμογέας θηκariού
ES	Ajustador deslizante de la vaina
ET	Liughülsi regulaator
FR	Ajusteur de gaine coulissant
HR	Klizni podešivac
HU	Elcsúsztható hüvelyhossz-beállító
IT	Regolatore scorrevole della guaina
LT	Slankūsis vamzdelio reguliatorius
LV	Pārbidāmais ievadapvalka regulētājs
NL	Schuiwregeling van de sheath
NO	Glidende hylsejusteringsmekanisme
PL	Przesuwany regulator koszulki
PT	Regulador deslizante da bainha
RO	Ajustor glisant al tecii
SK	Posuvný nastavovač puzdra
SV	Glidande hylsjusterare
TR	Kayar Kilif Ayarlayıcı

EN Needle Handle

BG	Дръжка на иглата
CS	Rukojeť jehly
DA	Nålegreb
DE	Nadelgriff
EL	Λαβή βελόνας
ES	Mango de la aguja
ET	Nõela käepide
FR	Poignée de l'aiguille
HR	Drška igle
HU	Tű fogantyúja
IT	Impugnatura dell'ago
LT	Adatos kotelis
LV	Adatas rokturis
NL	Veiligheidsring
NO	Nålehåndtak
PL	Uchwyt igły
PT	Punho da agulha
RO	Mânerul acului
SK	Rukoväť ihly
SV	Nålhandtag
TR	İğne Sapı

EN Safety Ring

BG	Предпазен пръстен
CS	Pojistný prstenec
DA	Sikkerhedsring
DE	Sicherheitsring
EL	Δακτύλιος ασφαλείας
ES	Anillo de seguridad
ET	Kaitserõngas
FR	Bague de sécurité
HR	Sigurnosni prsten
HU	Biztonsági gyűrű
IT	Anello di sicurezza
LT	Apsauginis žiedas
LV	Drošības gredzens
NL	Veiligheidsring
NO	Sikkerhetsring
PL	Pierścień bezpieczeństwa
PT	Anel de segurança
RO	Inel de siguranță
SK	Bezpečnostný krúžok
SV	Säkerhetslåsring
TR	Güvenlik Halkası

EN Stylet Hub

BG	Втулка на стилето
CS	Nástavec stiletu
DA	Stiletmuffe
DE	Mandrinansatz
EL	Ομφαλός του στείλεου
ES	Conector del estilete
ET	Stiletikinnitus
FR	Embase du stylet
HR	Cvorište sonde
HU	Mandrin végidoma
IT	Pomello del mandrino
LT	Zondo kaištis
LV	Stiletas galviņa
NL	Aanzetstuk van het stilet
NO	Stillettffeste
PL	Głowica mandrynu
PT	Conector do estilete
RO	Acoperitoarea stiletului
SK	Hrdlo stiletu
SV	Mandringfattning
TR	Stile Göbeği

EN Sheath Reference Mark

BG	Означение за дължината на защитната обвивка
CS	Referenční značka sheathu
DA	Referencemærke på sheath
DE	Schleusenbezugsmarke
EL	Σήμανση αναφοράς θηκariού
ES	Marca de referencia de la vaina
ET	Hülsi tähis
FR	Repère de référence sur la gaine
HR	Referentna oznaka uvodnice
HU	Hüvelyhossz referencijele
IT	Indicatore di riferimento della guaina
LT	Vamzdelio atskaitos žymė
LV	Ievadapvalka atskaites atzīme
NL	Referentiemarkering van de sheath
NO	Referansemerke på hylsen
PL	Znacznik referencyjny wysunięcia koszulki
PT	Marca de referência da bainha
RO	Marcaj de referință al tecii
SK	Referenčná značka puzdra
SV	Hylsans referensmarkering
TR	Kilif Referans İşareti

Fig. 3




Fig. 4

ECHOTIP® ULTRA ENDOBRONCHIAL HD ULTRASOUND NEEDLE FOR PENTAX SCOPES

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

INTENDED USE

This device is used to sample targeted submucosal and extramural lesions within or adjacent to the tracheobronchial tree through the accessory channel of an ultrasound endoscope for Fine Needle Aspiration (FNA).

DEVICE DESCRIPTION

The device consists of a handle with adjustable components to allow the user to adjust the extension of the needle and sheath. The device has a 4.1 Fr outer sheath for protection of the needle and the scope, and it is attached to the accessory channel of a Pentax endobronchial ultrasound (EBUS) scope. The device has a removable sheath adjuster to allow the user to insert the needle into the scope in the same orientation repeatedly. The device has a maximum needle extension length of 5 cm and sheath adjustment length of 3 cm. The device is supplied with a 10 ml syringe to utilize the standard vacuum technique for sample collection if desired. The device is available in two different needle sizes: 22 and 25 ga.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications include, but are not limited to: coagulopathy.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with bronchial endoscopy include, but are not limited to: allergic reaction to medication, allergic reaction to nickel, aspiration, cardiac arrhythmia or arrest, damage to blood vessels, death, discomfort, fever, hemorrhage, hypotension, infection, inflammation, nerve damage, pain, perforation, pneumoperitoneum, respiratory depression or arrest, sepsis, septicemia/bacteremia and tumor seeding.

WARNINGS

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the scope.

Ensure the stylet is fully inserted when advancing the needle into the target site.

Do not use this device for any purpose other than stated intended use.

Use of this device is restricted to a trained healthcare professional (Rx only).

When targeting multiple sites, adequate care should be taken to minimize the risk of tumor seeding.

This device is intended for use with a Pentax EBUS scope.

SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when the stopcock is in "open" position, aligned with syringe (**Fig. 1**).
2. Prepare syringe as follows:
 - a) With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b) Turn stopcock to "closed" position.
 - c) Pull plunger back until it is locked into place at desired setting, creating a vacuum.
 - d) Set prepared syringe aside until aspiration is required.

INSTRUCTIONS FOR USE

Illustrations

1. Identify the desired target site by endobronchial ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope. **Caution:** If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred sheath length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (**Fig. 2**). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining the position of the scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring (**Fig. 3**). Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
8. Standard vacuum syringe techniques may be applied for sample collection (Steps 9-11) or, if desired, other techniques that may or may not incorporate the use of the stylet may be used.
9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional sample collection is required later.
10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer fitting on the needle handle.
11. Turn the stopcock to the "open" position aligned with the syringe, allowing the negative pressure in the syringe to aspirate the sample. Move the needle handle back and forth within the target site. **Note:** Do not remove the needle from the target site during sample collection.
12. Upon completion of sample collection, if using the syringe, turn the stopcock to the "closed" position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
13. To remove the needle from the scope: unlock the thumbscrew on the sliding sheath adjuster while maintaining attachment of the sliding sheath adjuster to the scope and remove needle

(Fig. 4). Or, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.

Note: When reinserting device into the scope for additional passes and before advancing it into the target site, you must position the handle onto the scope such that the thumbscrews on the handle are aligned in the same direction as they were prior to device removal.

14. Unlock the thumbscrew on the safety ring and extend the needle.
15. Use a new syringe or the stylet to expel the specimen, then prepare per institutional guidelines.
16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.
17. For additional sample collection from the same target site, gently reinsert the stylet into the Luer fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.
18. If the needle was removed from the scope while maintaining attachment of the sheath adjuster to the scope in Step 13: readvance the sheath through the attached sheath adjuster, ensuring that both thumbscrews are aligned on same side. Tighten the thumbscrew on the sliding sheath adjuster at the preferred sheath length.
19. Additional samples may be obtained by repeating Step 2 through 16 of "Instructions for Use".
20. Upon completion of the procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.

Upon completion of the procedure, dispose of the device, per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Store in a dry location, away from temperature extremes.

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© COOK 2021
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.