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


Fig. 4

ECHOTIP ULTRA® ENDOBRONCHIAL HD ULTRASOUND NEEDLE FOR PENTAX SCOPES

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The EchoTip Ultra® Endobronchial HD Ultrasound Needle for Pentax Scopes is available in two different needle sizes (22 and 25 gage) and contains the following components:

- Needle
- Syringe

Performance Characteristics

The function and key features of these components are described below:

Endobronchial needle

- Handle – The handle contains adjustable components to allow the user to adjust the extension of the needle and the sheath. The handle safety ring slides and locks the needle at the desired needle length extension (0–5 cm). The '0 cm' reference mark ensures complete needle retraction within the sheath. The sliding sheath adjuster allows the user to set the required sheath length extension (0–3 cm) for their Pentax endobronchial ultrasound (EBUS) scope. The '0 cm' reference mark denotes the minimum sheath length available. The base of the sliding sheath adjuster has a Luer lock fitting, which facilitates the attachment of the device to the scope accessory channel.
- 4.1 Fr Sheath – The sheath covers the needle during introduction, advancement, and withdrawal of the device from the scope accessory channel. It protects the scope accessory channel from damage by the needle tip.
- Needle – The stainless steel needle has a dimpling pattern to provide visualization of the needle tip under endobronchial ultrasound. The needle has a beveled tip to puncture and facilitate sample collection from the target site.
- Stylet – The stylet is inserted through the length of the needle. The distal end of the stylet is beveled to match the bevel of the needle tip. The stylet provides support to the needle during advancement into the target site, prevents sample collection outside of the target site, and may also be used for sample retrieval following collection if desired.

Syringe – The device is supplied with a 10 mL syringe to utilize the standard vacuum technique for sample collection, if desired. The Luer lock fitting at the tip of the syringe can be attached to the Luer fitting on the handle. The syringe has a plunger and a stopcock to create a vacuum for sample collection. It is supplied with 5 mL and 10 mL clips attached to the plunger to allow the user to select the preferred volume for sample aspiration.

Device Compatibility

The EchoTip Ultra Endobronchial HD Ultrasound Needle for Pentax Scopes is compatible with:

- Pentax EBUS scope with 2 mm minimum accessory channel
- 10 mL Standard Luer syringe

Patient Population

Adult patients requiring fine needle aspiration (FNA) of submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree.

Intended User

Use of this device is restricted to a trained healthcare professional.

Contact with Body Tissue

This device is tissue contacting in line with intended use.

Operating Principles

Pentax EBUS scopes have a Luer fitting to allow attachment of the needle to the accessory channel port. The device is inserted into the scope accessory channel port and attached when the fittings meet.

The sheath length can be adjusted depending on the length of the EBUS scope using the sliding sheath adjuster. The needle length extension can be adjusted as required. The device is used to puncture the target site to obtain a sample. If using a standard vacuum syringe technique, the stylet is then removed, and the syringe is attached to the Luer fitting on the needle handle. Turning the syringe stopcock to the "open" position allows negative pressure to aspirate the sample. The handle is gently moved in small increments back and forth within the target site to obtain a sample. If using the syringe, the stopcock is returned to the "closed" position. The needle is retracted into the sheath and the device is disconnected from the accessory channel port.

INTENDED USE

This device is used to sample targeted submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree through the accessory channel of an ultrasound bronchoscope for fine needle aspiration (FNA).

INDICATIONS FOR USE

This device is indicated for endobronchial ultrasound guided tissue sampling of submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree, which is used for subsequent pathological examination to aid in disease diagnosis and patient management.

CLINICAL BENEFITS

This device facilitates successful sampling of submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree to aid in disease diagnosis and patient management.

CONTRAINDICATIONS

Those specific to the primary endobronchial procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications: coagulopathy.

WARNINGS

- Not for use in the heart or vascular system.
- The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.
- This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.
- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect the device with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.
- This device contains nickel, which may cause allergic reaction in individuals with nickel sensitivity.

PRECAUTIONS

- Ensure the stylet is fully inserted when advancing the needle into the target site.
- The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement, or withdrawal of the device. Failure to retract the needle may result in damage to the scope.
- Do not use this device for any purpose other than stated intended use.
- Use of this device is restricted to a trained healthcare professional.
- The product is intended for use by physicians trained and experienced in endobronchial ultrasound techniques.

- When targeting multiple sites, adequate care should be taken to minimize the risk of tumor seeding.
- This device is intended to be used with a Pentax EBUS scope.
- Store the device in a dry location.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with an endobronchial ultrasound procedure: allergic reaction to medication • aspiration • cardiac arrhythmia or arrest • fever • hemorrhage • hypotension • infection • pain/discomfort • perforation • respiratory depression or arrest.

Those associated with the device: allergic reaction to nickel • damage to blood vessels • death (other than due to normal disease progression) • fever • hemorrhage • infection • inflammation • nerve damage • pain/discomfort • perforation • pneumoperitoneum • sepsis • septicemia/bacteremia • tumor seeding.

HOW SUPPLIED

These devices are supplied ethylene oxide (EO) sterilized in a tray with a Tyvek lid.

SYSTEM PREPARATION

Illustrations

1. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
2. Visually inspect the device components confirming they are free from damage prior to use.
3. Examine the syringe. It has two plunger locks that must be depressed to advance the plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when the stopcock is in the "open" position, aligned with the syringe (**Fig. 1**).
4. Prepare the syringe as follows:
 - a. With the stopcock in the "open" position, depress the plunger locks and fully advance the plunger into the syringe.
 - b. Turn the stopcock to the "closed" position.
 - c. Pull the plunger back until it is locked into place at the desired setting, creating a vacuum.
 - d. Set the prepared syringe aside until aspiration is required.

INSTRUCTIONS FOR USE

1. Identify the desired target site by endobronchial ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope. **Caution:** If resistance is encountered on needle introduction, reduce angulation of scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the accessory channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred sheath length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (**Fig. 2**). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining the position of the scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring (**Fig. 3**). Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement,

retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.

8. Standard vacuum syringe techniques may be applied for sample collection (Steps 9–11) or, if desired, other techniques that may or may not incorporate the use of the stylet may be used.
9. Remove the stylet of the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional sample collection is required later.
10. Attach the Luer lock fitting of the previously prepared syringe securely into the fitting on the needle handle.
11. Turn the stopcock to the “open” position aligned with the syringe, allowing the negative pressure in the syringe to aspirate the sample. Move the needle handle back and forth within the target site. **Note:** Do not remove the needle from the target site during sample collection.
12. Upon completion of sample collection, if using the syringe, turn the stopcock to the “closed” position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
13. To remove the needle from the scope: unlock the thumbscrew on the sliding sheath adjuster while maintaining attachment of the sliding sheath adjuster to the scope and remove needle (**Fig. 4**). Or, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counter clockwise, and withdraw the entire device from the scope. **Note:** When reinserting device into scope for additional passes and before advancing it into the target site, you must position the handle onto the scope such that the thumbscrews on the handle are aligned in the same direction as they were prior to device removal.
14. Unlock the thumbscrew on the safety ring and extend the needle.
15. Use a new syringe or stylet to expel specimen, then prepare per institutional guidelines.
16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.
17. For additional sample collection from the same target site, gently reinsert the stylet into the metal fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.
18. If needle was removed from the scope while maintaining attachment of the sheath adjuster to the scope in Step 13: readvance the sheath through the attached sheath adjuster, ensuring that both thumbscrews are aligned on same side. Tighten the thumbscrew on the sliding sheath adjuster at the preferred sheath length.
19. Additional samples may be obtained by repeating Steps 2 through 16 of “Instructions for Use”.
20. Upon completion of procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counter clockwise, and withdraw the entire device from the scope.

DEVICE DISPOSAL

This device may be contaminated with potentially infectious substances of human origin and should be coiled for disposal in accordance with institutional guidelines.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2025
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.