EN
6

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





3.

7.

Evolution® Esophageal Stent System – Partially Covered / Fully Covered

IMPORTANT INFORMATION

Please review prior to use.

INTENDED USE

This device is used to maintain patency of malignant esophageal strictures and/or to seal tracheoesophageal fistulas.

STENT DESCRIPTION

This flexible, self-expanding stent is constructed of a single, woven, nitinol wire. The stent foreshortens due to its design. The total length of the stent is indicated by radiopaque markers on the inner catheter, indicating the actual length of the stent at nominal stent diameter. There is a lasso at the proximal end of the stent (partially covered) or the proximal and distal ends of the stent (fully covered), purpose is to reposition the stent in the event of incorrect placement.

DELIVERY SYSTEM DESCRIPTION

Stent is mounted on an inner catheter, which accepts a .035 inch wire guide, and is constrained by an outer catheter. A pistol-grip delivery handle allows stent deployment or recapture.

NOTES

Do not use this device for any purpose other than stated intended use.

Part numbers with the prefix EVO-FC are Fully Covered.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to upper GI endoscopy and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to: total esophageal obstruction, strictures that cannot be dilated a minimum size as outlined in the precautions section, placement requiring positioning of stent within 2 cm of the cricopharynx, surgical resection candidates, hiatal hernia and gastric prolapse in the esophageal area, patients with a perforated esophagus, placement in actively bleeding tumors, benign diseases.

Relative contraindications include, but are not limited to: uncooperative patient, coagulopathy, tracheal compression, recent myocardial infarction, cervical arthritis with fixed

cervical spine, large tumor mass occupying the mediastinum, nonobstructive tumor, gastric outlet obstruction, necrotic esophageal mucosa, acutely angled stenosis.

POTENTIAL COMPLICATIONS

Those associated with upper GI endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, reflux, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications include, but are not limited to: stent misplacement and/or migration; tumor ingrowth or overgrowth; esophageal ulceration and erosion; nausea; chest or retrosternal pain; foreign body sensation; food bolus impaction; gasbloat; sensitivity to metal components; fistula involving trachea, bronchi or pleural space; intestinal obstruction secondary to migration; mediastinitis or peritonitis; airway compression; tracheal obstruction; aorto and arterioesophageal fistula; erosion or perforation of stent into adjacent vascular structures.

PRECAUTIONS

A complete diagnostic evaluation must be performed prior to use to determine proper stent size.

Stent should be placed using fluoroscopic monitoring.

Stent should only be placed with the Cook delivery system, which is provided with each stent.

Note: Prior to advancing system, area to be stented should be dilated to:

For (18mm x 23mm) stent – a minimum 9 mm and a maximum of 11 mm. If area is dilated greater than 11 mm, stent may migrate.

For (20mm x 25mm) stent – a minimum of 10 mm and a maximum of 14 mm. If area is dilated greater than 14 mm, stent may migrate.

This device is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

After stent placement, alternative methods of treatment such as chemotherapy and radiation should not be administered as this may increase risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

Long-term patency of this device has not been established. Periodic evaluation is advised.

After stent placement, patients should be instructed to chew food well or eat soft or pureed food.

This device shortens upon deployment. With proximal strictures near the upper esophageal sphincter, deployment should be performed under fluoroscopic visualization as this may enhance placement accuracy.

WARNINGS

The **stent is not intended to be removed** and is intended to remain in the body permanently. Attempts to remove stent after placement may cause damage to esophageal mucosa.

Stent cannot be retrieved after the deployment threshold has been passed. Corresponding marks on outer catheter and delivery handle indicate when threshold has been passed.

The risk of perforation and erosion into adjacent vascular structures or aorto-esophageal and arterio-esophageal fistulas may be increased with pre- or post-operative chemotherapy and radiation, longer implantation times, aberrant anatomy, and/or mediastinal contamination or inflammation.

MRI Safety Information

This symbol means the device is MR Conditional

Nonclinical testing has demonstrated the Evolution® Esophageal Stents are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3 T or 1.5 T
- Maximum spatial field gradient of 1900 gauss/cm (19 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, nonclinical testing results indicate the Evolution® Esophageal Stent is expected to produce a maximum temperature rise of less than 4.4°C after 15 minutes of continuous scanning.

In nonclinical testing, the image artifact caused by the device extends approximately 7 mm from the Evolution® Esophageal Stent when imaged with a spin echo pulse sequence and a 3 T MR system.

The safety of MR imaging using a local transmit/receive RF coil in areas immediately adjacent or directly involving the stent has not been established and is therefore not recommended.

The safety of performing MRI imaging in patients with overlapping stents or other device(s) in the esophagus has not been established and is therefore not recommended.

For US Patients only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail:	MediAlert Foundation International 2323 Colorado Avenue Turlock, CA 95382
Phone:	888-633-4298 (toll free) 209-668-3333 from outside the US
Fax:	209-669-2450
Web:	www.medicalert.org

EQUIPMENT REQUIRED

Stent and delivery system

Esophageal Balloon Dilators or Bougie Dilators

.035 inch wire guide

PATIENT PREPARATION

1. Dilate stricture to a minimum size as outlined in the precautions section using Balloon Dilators or Bougie Dilators.
2. Endoscopically locate and mark upper and lower margins of lesion with radiopaque marker. Caution: Accurate marking of lesion borders is essential for proper stent placement.
3. Place a .035 inch wire guide, floppy tip first, through accessory channel, through stricture, until it is fluoroscopically visualized in stomach. Leave wire guide in place and remove endoscope.

INSTRUCTIONS FOR USE

Illustrations

1. Remove the protective tubing from stent.
2. Using fluoroscopic guidance, introduce and advance the introduction system in short increments over the pre-positioned wire guide into the esophagus.
3. To ensure that stent will bridge stricture after deployment, fluoroscopically position radiopaque markers on inner catheter beyond extremities of stricture to be crossed.
4. Confirm desired stent position fluoroscopically and deploy stent by removing red safety guard from the handle. (See Fig. 1)
5. Continue deploying stent by squeezing trigger. Note: Each trigger squeeze will deploy stent by an equal amount.
6. If stent repositioning is required during deployment, it is possible to recapture stent. Note: It is not possible to recapture stent after passing point-of-no-return, indicated when red marker on top of introducer has passed the point-of-no-return indicator on handle. (See Fig. 2)

7. To recapture stent, push direction button on side of delivery system to opposite side. (See Fig. 3) **Note:** Hold thumb on button when squeezing trigger for first time to recapture.

8. Continue squeezing trigger as required to recapture stent by desired amount.

9. To resume deployment, push button to opposite side again and hold thumb on button for first stroke while squeezing trigger. (See Fig. 4)

10. When stent point-of-no-return has been passed, pull safety wire out of delivery handle near wire guide port. (See Fig. 5)

11. Continue deploying stent by squeezing trigger.

12. After deployment, fluoroscopically confirm full stent expansion. Once full expansion is confirmed, introduction system can be safely removed. **Note:** You may introduce endoscope and advance to the upper margin of the stent to endoscopically confirm position and patency of the stent. **Warning:** Do not introduce the endoscope into the stent as displacement may occur.

For partially covered stent:

13. To reposition stent directly after placement, grip the white lasso at proximal-end of stent with forceps and reposition stent to desired area. (See Fig. 6)

For fully covered stent:

14. To reposition stent directly after placement, grip either the proximal-end lasso or the distal-end lasso with forceps and reposition stent to desired area. (See Fig. 7)

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2017 ©
COOK IRELAND LTD.