






- EN 6 **Evolution® Esophageal Stent System – Partially Covered / Fully Covered**  
 Instructions for Use
- BG 11 **Стент система за хранопровод Evolution® – частично покрита / напълно покрита**  
 Инструкции за употреба
- HR 17 **Sustav ezofagealnog stenta Evolution® – djelomično pokriven / potpuno pokriven**  
 Upute za uporabu
- CS 23 **Systém ezofageálního stentu Evolution® – částečně pokrytý / plně pokrytý**  
 Návod k použití
- DA 28 **Evolution® oesophagus-stentsystem – delvist dækket / helt dækket**  
 Brugsanvisning
- NL 34 **Evolution® oesofagusstentsysteem – gedeeltelijk bekleed/volledig bekleed**  
 Gebruiksaanwijzing
- ET 40 **Söögitorustendi süsteem Evolution® – osaliselt/täielikult kaetud**  
 Kasutusjuhised
- FR 45 **Système d'endoprothèse œsophagienne Evolution® – partiellement couverte / complètement couverte**  
 Mode d'emploi
- DE 51 **Evolution® Ösophagusstentsystem – teilweise gecouvert / vollständig gecouvert**  
 Anleitung zum Gebrauch
- EL 57 **Σύστημα ενδοπρόσθεσης οισοφάγου Evolution® – Μερικώς καλυμμένο / Πλήρως καλυμμένο**  
 Οδηγίες χρήσης
- HU 64 **Evolution® nyelökszisztentrendszer – részlegesen bevont / teljesen bevont**  
 Használati utasítás
- IT 69 **Sistema con stent esofageo Evolution® – Parzialmente rivestito / Completamente rivestito**  
 Istruzioni per l'uso
- LV 75 **Evolution® ezofageālā stenta sistēma — daļēji pārklāta/pilnībā pārklāta**  
 Lietošanas instrukcija
- LT 81 **„Evolution®“ stemplės stento sistema – iš dalies dengta / visiškai dengta**  
 Naudojimo nurodymai
- NO 87 **Evolution® stentsystem for øsofagus – delvis belagt / helt belagt**  
 Bruksanvisning

- PL**  
92 **System do stentowania przełyku Evolution® — częściowo pokrywany /  
całkowicie pokrywany**  
Instrukcja użycia
- PT**  
98 **Sistema de stent esofágico Evolution® – Parcialmente revestido / totalmente  
revestido**  
Instruções de utilização
- RO**  
104 **Sistem pentru stent esofagian Evolution® – cu acoperire parțială/cu  
acoperire integrală**  
Instrucțiuni de utilizare
- SK**  
110 **Systém ezofagového stentu Evolution® – čiastočne krytého/plne krytého**  
Návod na použitie
- ES**  
115 **Sistema de stent esofágico Evolution® – parcial o totalmente recubierto**  
Instrucciones de uso
- SV**  
121 **Evolution®-stentsystem för esofagus – delvis täckt/helt täckt**  
Bruksanvisning

## **EVOLUTION® ESOPHAGEAL STENT SYSTEM – PARTIALLY COVERED / FULLY COVERED**

**CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### **IMPORTANT INFORMATION**

Please review prior to use.

### **INTENDED USE**

This device is used to maintain patency of malignant esophageal strictures and/or to seal tracheoesophageal fistulas.

### **STENT DESCRIPTION**

This flexible, self-expanding stent is constructed of a single, woven, nitinol wire. The stent foreshortens due to its design. The stent has a flange at both ends to aid in preventing migration. A silicone membrane applied to the stent covers the complete stent in the fully covered option and the complete body of the stent along with part of both flanges in the partially covered stent. The covering is intended to reduce the risk of tissue ingrowth and to provide a seal for tracheoesophageal fistulas. The total length of the stent is indicated by radiopaque markers on the inner catheter, indicating the actual length of the stent at nominal stent diameter. To aid in visibility under fluoroscopy there are four tantalum bands at either end of the stent. There is a lasso at the proximal end of the stent (partially covered) or the proximal and distal ends of the stent (fully covered), the purpose of which is to reposition the stent in the event of incorrect placement. The proximal and distal ends of the stent are designed with numerous rounded “crowns” to potentially reduce tissue trauma.

### **DELIVERY SYSTEM DESCRIPTION**

The stent is mounted on an inner catheter, which accepts a .035-inch wire guide, and is constrained by an outer catheter. A pistol-grip delivery handle allows controlled stent deployment or recapture. The position of the red marker on top of the handle in relation to the “point-of-no-return” position on the handle label indicates the point after which it is not possible to recapture the deployed stent.

### **CONTRAINDICATIONS**

Those specific to upper GI endoscopy and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to: total esophageal obstruction • strictures that cannot be dilated to a minimum size as outlined in the precautions section • placement requiring positioning of stent within 2 cm of the cricopharynx • surgical resection candidates • hiatal hernia and gastric prolapse in the esophageal area • patients with a perforated esophagus • placement in actively bleeding tumors • benign diseases.

Relative contraindications include, but are not limited to: uncooperative patient • coagulopathy • tracheal compression • recent myocardial infarction • cervical arthritis with

fixed cervical spine • large tumor mass occupying the mediastinum • nonobstructive tumor • gastric outlet obstruction • necrotic esophageal mucosa • acutely angled stenosis.

## POTENTIAL ADVERSE EVENTS

Those associated with upper GI endoscopy include, but are not limited to: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • fever • hemorrhage • hypotension • infection • perforation • reflux • respiratory depression or arrest • vomiting.

Additional adverse events include, but are not limited to: airway compression • allergic reaction to nickel • aorto-esophageal fistula and arterio-esophageal fistula • chest or retrosternal pain • death (other than due to normal disease progression) • dysphagia • edema • erosion or perforation of stent into adjacent vascular structures • esophageal ulceration and erosion • esophagitis • fistula involving trachea, bronchi or pleural space • food bolus impaction • foreign body sensation or reaction • gas bloat • inadequate stent expansion • intestinal obstruction secondary to migration • mediastinitis or peritonitis • nausea • pain/discomfort • reocclusion • sensitivity to metal components • sepsis • stent misplacement and/or migration • tracheal obstruction • tumor ingrowth or overgrowth • wire entrapment.

## PRECAUTIONS

A complete diagnostic evaluation must be performed prior to use to determine proper stent size.

Part numbers with the prefix EVO-FC are Fully Covered.

Stent should be placed using fluoroscopic monitoring.

Stent should only be placed with the Cook delivery system, which is provided with each stent.

**Note:** Prior to advancing system, area to be stented should be dilated to:

**For (18mm x 23mm) stent** – a minimum of 9 mm and a maximum of 11 mm. If area is dilated greater than 11 mm, stent may migrate.

**For (20mm x 25mm) stent** – a minimum of 10 mm and a maximum of 14 mm. If area is dilated greater than 14 mm, stent may migrate.

This device is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

Long-term patency of this device has not been established. Periodic evaluation is advised.

This device shortens upon deployment. With proximal strictures near the upper esophageal sphincter, deployment should be performed under fluoroscopic visualization as this may enhance placement accuracy.

Do not use this device for any purpose other than stated intended use.

Stent cannot be retrieved after the deployment threshold has been passed. The position of the red marker on top of the handle in relation to the point-of-no-return position on the handle label indicates when the threshold has been passed.

Use of this device is restricted to a trained healthcare professional.

## WARNINGS

The **stent is not intended to be removed** and is considered a permanent implant. Attempts to remove stent after placement may cause damage to esophageal mucosa.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect product with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working conditions, do not use. Please notify Cook for return authorization.

This stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

Do not introduce the endoscope into the stent as displacement may occur.

After stent placement, alternative methods of treatment such as chemotherapy and radiation should not be administered as this may increase risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

The risk of perforation and erosion into adjacent vascular structures or aorto-esophageal and arterio-esophageal fistulas may be increased with pre- or post-operative chemotherapy and radiation, longer implantation times, aberrant anatomy, and/or mediastinal contamination or inflammation.

## MR INFORMATION

**This symbol means the device is MR Conditional**

Non-clinical testing has demonstrated that Evolution® Esophageal Stents are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

### Static Magnetic Field

- Static magnetic field of 3.0 Tesla or less
- Maximum spatial field gradient of 1900 gauss/cm (19 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

## **MRI-Related Heating**

### **1.5 and 3.0 Tesla Systems**

Non-clinical testing was conducted under the following conditions and produced a maximum temperature rise of 4.4°C:

- A maximum whole body average specific absorption rate (SAR) of 3.8 W/kg (corresponding to a calorimetry measured value of 3.1 W/kg) for 15 minutes of MR scanning in a 1.5 Tesla Magnetom (Siemens Medical Solutions, Malvern, PA, Software Numaris/4) MR scanner.
- A maximum whole body average specific absorption rate (SAR) of 3.0 W/kg (corresponding to a calorimetry measured value of 2.8 W/kg) for 15 minutes of MR scanning in a 3.0 Tesla Excite (GE Electric Healthcare, Milwaukee, WI, Software G3.0-052B) MR scanner.

## **Image Artifacts**

In nonclinical testing, the image artifact caused by the device extends approximately 7 mm from the Evolution® Esophageal Stent when imaged with a spin echo pulse sequence and a 3 T MR system.

The safety of MR imaging using a local transmit/receive RF coil in areas immediately adjacent or directly involving the implant has not been established and is therefore not recommended.

The safety of performing MR imaging in patients with overlapping stents or other device(s) in the esophagus has not been established and is therefore not recommended.

## **For US Patients only**

It is recommended that patients register the conditions under which the implant can safely be scanned with the MedicAlert Foundation ([www.medicalert.org](http://www.medicalert.org)) or an equivalent organization.

## **EQUIPMENT REQUIRED**

Stent and delivery system

Esophageal Balloon Dilators or Bougie Dilators

.035 inch wire guide

## **PATIENT PREPARATION**

1. Intubate the patient using an endoscope per standard procedure. Assess the stricture location and dimensions upon direct visualization.
2. If required, dilate stricture to a minimum size as outlined in the precautions section using Balloon Dilators or Bougie Dilators.
3. Endoscopically locate and mark upper and lower margins of lesion with radiopaque marker. **Caution:** Accurate marking of lesion borders is essential for proper stent placement.
4. The esophageal stent should be at least 4 cm longer than the stricture, allowing for approximately 2 cm longer at either end of the stricture.

5. Place a .035-inch wire guide, floppy tip first, through accessory channel, through stricture, until it is fluoroscopically visualized in stomach. Leave wire guide in place and remove endoscope.

## INSTRUCTIONS FOR USE

## Illustrations

1. Remove stent delivery system from the package and protective tubing.
  2. Using fluoroscopic guidance, introduce and advance the delivery system in short increments over the pre-positioned wire guide into the esophagus.
  3. To ensure that stent will bridge stricture after deployment, fluoroscopically position radiopaque markers on inner catheter beyond extremities of stricture to be crossed.
  4. Confirm desired stent position fluoroscopically and deploy stent by removing red safety guard from the handle. **(Fig. 1)**
  5. Continue deploying stent by squeezing trigger. **Note:** Each trigger squeeze will deploy stent by an equal amount.
  6. If stent repositioning is required during deployment, it is possible to recapture stent. **Note:** It is not possible to recapture stent after passing point-of-no-return, indicated when red marker on top of the handle has passed the point-of-no-return position on the handle label. **(Fig. 2)**
  7. To recapture stent, push direction button on side of delivery system to opposite side. **(Fig. 3) Note:** Hold thumb on button when squeezing trigger for first time to recapture.
  8. Continue squeezing trigger as required to recapture stent by desired amount.
  9. To resume deployment, push button to opposite side again and hold thumb on button for first stroke while squeezing trigger. **(Fig. 4)**
  10. When stent point-of-no-return has been passed, pull safety wire out of delivery handle near wire guide port. **(Fig. 5)**
  11. Continue deploying stent by squeezing trigger.
  12. After deployment, fluoroscopically confirm full stent expansion. Once full expansion is confirmed, delivery system can be safely removed. **Note:** You may introduce endoscope and advance to the upper margin of the stent to endoscopically confirm position and patency of the stent.
- For partially covered stent:**
13. To reposition stent directly after placement, grip the white lasso at proximal-end of stent with forceps and reposition stent to desired area.
- For fully covered stent:**
14. To reposition stent directly after placement, grip either the proximal-end lasso or the distal-end lasso with forceps and reposition stent to desired area. **(Fig. 7)**

### **Post Procedure**

Observe the patient for development of any complications of endoscopy, esophageal dilation or stent placement. During the first 24 hours, vital signs should be monitored, and clear liquid given in an upright position. Patients being treated for fistula should have no fluids or food by mouth until successful seal of fistula has been confirmed.

Subsequent follow-ups (such as endoscopy and/or chest x-rays) with patient should be performed at the physician's direction.

### **PATIENT COUNSELLING INFORMATION**

After stent placement, post 24 hours, the patient should be instructed to chew food well or eat soft or pureed food and to drink fluids before and after meals.

Patients with stents placed in the distal esophagus or across the GE junction may be instructed to elevate the head of the bed and may be prescribed acid suppression therapy to minimize gastric reflux into the stent.

### **HOW SUPPLIED**

Store in a dry location, away from temperature extremes.

**Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.**

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© COPYRIGHT COOK 2020 ©
**COOK IRELAND LTD.**