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





7

8

EVOLUTION® ESOPHAGEAL STENT SYSTEM – FULLY COVERED

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

IMPORTANT INFORMATION

Please review prior to use.

STENT DESCRIPTION

This flexible, self-expanding stent is constructed of a single, woven, nitinol wire. The stent foreshortens due to its design. The stent has a flange at both ends to aid in preventing migration. The stent has a silicone membrane covering. The covering is intended to reduce the risk of tissue ingrowth and to provide a seal for tracheoesophageal fistulas. The total length of the stent is indicated by radiopaque markers on the inner catheter, indicating the actual length of the stent at nominal stent diameter. To aid in visibility under fluoroscopy there are four tantalum bands at either end of the stent. There is a lasso at the proximal and distal ends of the stent that allows repositioning of the stent directly after placement or is used for stent removal in benign indications. There are radiopaque markers on the proximal and distal ends of the stent. The proximal and distal ends of the stent are designed with numerous rounded “crowns” to potentially reduce tissue trauma.

DELIVERY SYSTEM DESCRIPTION

The stent is mounted on an inner catheter, which accepts a .035-inch wire guide, and is constrained by an outer catheter. A pistol-grip delivery handle allows controlled stent deployment or recapture. The position of the red marker on top of the handle in relation to the “point-of-no-return” position on the handle label indicates the point after which it is not possible to recapture the deployed stent.

INTENDED USE

This device is used to maintain luminal patency of the esophagus in cases of: obstruction caused by intrinsic or extrinsic malignancies, refractory benign strictures, or to seal tracheoesophageal fistulas.

MALIGNANT AND BENIGN

The above symbol means the following: Malignant and Benign.

CONTRAINDICATIONS

Contraindications include those specific to upper gastrointestinal (GI) endoscopy and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to:

- total esophageal obstruction
- strictures that cannot be dilated to a minimum size as outlined in the precautions section
- placement requiring positioning of stent within 2 cm of the cricopharynx
- surgical resection patients
- hiatal hernia and gastric prolapse in the esophageal area
- patients with a perforated esophagus
- placement in actively bleeding tumors
- placement in polypoid lesions

Relative contraindications include, but are not limited to:

- uncooperative patient
- coagulopathy
- tracheal compression
- recent myocardial infarction
- cervical arthritis with fixed cervical spine
- large tumor mass occupying the mediastinum
- nonobstructive tumor
- gastric outlet obstruction
- necrotic esophageal mucosa
- acutely angled stenosis
- placement in an esophago-jejunostomy (following gastrectomy)

POTENTIAL ADVERSE EVENTS

Those associated with upper GI endoscopy include, but are not limited to: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • fever • hemorrhage • hypotension • infection • perforation • reflux • respiratory depression or arrest • vomiting.

Additional adverse events include, but are not limited to: airway compression • allergic reaction to nickel • aortoesophageal fistula and arterioesophageal fistula • chest or retrosternal pain • death (other than due to normal disease progression) • dysphagia • edema • erosion or perforation of stent into adjacent vascular structures • esophageal ulceration and erosion • esophagitis • fistula involving trachea, bronchi or pleural space • food bolus impaction • foreign body sensation or reaction • gas bloat • inadequate stent expansion • intestinal obstruction secondary to migration • mediastinitis or peritonitis • nausea • pain/discomfort • reocclusion • sensitivity to metal components • sepsis • stent misplacement and/or migration • tracheal obstruction • tumor overgrowth • wire entrapment.

PRECAUTIONS

- A complete diagnostic evaluation must be performed prior to use to determine proper stent size.
- The Evolution Esophageal Stent – Fully Covered should be used with caution and only after careful consideration in patients with:
 - Significant preexisting pulmonary or cardiac disease.
 - Conditions where esophageal stenting may cause a secondary tracheal stenosis.
 - Prior gastric pull up operations resulting in altered anatomy.
- Stent should be placed using fluoroscopic monitoring.
- Prior to advancing the system, physician discretion should be used to determine if dilation is required. If required, dilate area to be stented to a diameter that is at least equal to the diameter of the delivery system and a maximum of 11 mm.
- Long-term patency of this device has not been established. Periodic evaluation is advised.
- This device shortens upon deployment. With proximal strictures near the upper esophageal sphincter, deployment should be performed under fluoroscopic visualization as this may enhance placement accuracy.
- For Malignant use, alternate methods of therapy should be investigated prior to placement.
- The stent should only be placed with the Cook delivery system, which is provided with each stent.

- Do not use this device for any purpose other than the stated intended use.
- Use of this device is restricted to a trained healthcare professional.
- Stent cannot be retrieved after the deployment threshold has been passed. The position of the red marker on the top of the handle in relation to point-of-no-return position on the handle label indicates when threshold has been passed.

WARNINGS

- The safety and effectiveness of leaving this stent in place, or removing the Evolution Esophageal Stent – Fully Covered from a benign lesion beyond 8 weeks has not been established.
- Exercise care when removing a stent that has significant overgrowth around the stent ends.
- It is not intended to remove a stent from a malignant lesion once original placement has been completed and confirmed. Attempts to remove the stent after placement in a malignant lesion may cause damage to the esophageal mucosa.
- Stenting across the gastroesophageal junction may increase the risk of migration.
- This stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- The effectiveness of the delivery system cannot be guaranteed if the stent has been recaptured more than 5 times.
- Testing of overlapping stents has not been completed and is not recommended.
- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.
- Do not introduce the endoscope into the stent, as displacement may occur.
- After stent placement, alternative methods of treatment such as chemotherapy and radiation should not be administered as this may increase risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.
- The migration rate may be increased by absence of a tumor abutment in benign esophageal stenoses, esophageal fistula without stenoses, incorrect sizing of the stent, or in cases of malignant stenoses where radiotherapy and/or chemotherapy is employed.
- The risk of perforation and erosion into adjacent vascular structures or aorto-esophageal and arterio-esophageal fistulas may be increased with pre- or post-operative chemotherapy and radiation, longer implantation times, aberrant anatomy, and/or mediastinal contamination or inflammation.

MR INFORMATION

This symbol means the device is MR Conditional

Non-clinical testing has demonstrated the Evolution® Esophageal Stents are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

Static Magnetic Field:

- Static magnetic field of 3.0 Tesla or less
- Maximum spatial field gradient of 1900 gauss/cm (19 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/Kg (Normal Operating Mode).

MRI-Related Heating

1.5 and 3.0 Tesla Systems:

Non-clinical testing was conducted under the following conditions and produced a maximum temperature rise of 4.4 °C:

- a maximum whole body average specific absorption rate (SAR) of 3.8 W/kg (corresponding to a calorimetry measured value of 3.1 W/kg) for 15 minutes of MR scanning in a 1.5 Tesla Magnetom (Siemens Medical Solutions, Malvern, PA, Software Numaris/4) MR scanner.
- a maximum whole body average specific absorption rate (SAR) of 3.0 W/kg (corresponding to a calorimetry measured value of 2.8 W/kg) for 15 minutes of MR scanning in a 3.0 Tesla Excite (GE Electric Healthcare, Milwaukee, WI, Software G3.0-052B) MR scanner.

In non-clinical testing, the image artifact caused by the device extends approximately 7mm from the Evolution® Esophageal Stent when imaged with a spin echo pulse sequence and 3T MR system.

The safety of MR imaging using a local transmit/receive RF coil in areas immediately adjacent or directly involving the implant has not been established and is therefore not recommended.

The safety of performing MR imaging in patients with overlapping stents or other device(s) in the esophagus has not been established and is therefore not recommended.

For US Patients only

It is recommended that patients register the conditions under which the implant can safely be scanned with the MedicalAlert Foundation (www.medicalert.org) or an equivalent organization.

INSTRUCTIONS FOR USE

Equipment Required

Stent and delivery system

Esophageal Balloon Dilators or Bougie Dilators

.035 inch wire guide

Patient Preparation

1. Intubate the patient using an endoscope per standard procedure. Assess the stricture location and dimensions upon direct visualization.
2. **If required**, dilate stricture to a minimum size as outlined in the Precautions section using Balloon Dilators or Bougie Dilators.
3. Endoscopically locate and mark upper and lower margins of lesion with radiopaque marker, internal or external, or use an injectable radiopaque substance. **CAUTION:** Accurate marking of lesion borders is essential for proper stent selection and placement.
4. The esophageal stent should be at least 4 cm longer than the stricture, allowing for approximately 2 cm longer at either end of the stricture.
5. Place a .035 inch wire guide, floppy tip first, through accessory channel, through stricture, until it is fluoroscopically visualized in stomach. Leave wire guide in place and remove endoscope.

Stent Placement

1. Remove stent delivery system from the package and protective tubing.
2. Using fluoroscopic guidance, introduce and advance the delivery system in short increments over the pre-positioned wire guide into the esophagus.
3. To ensure that the stent will bridge stricture after deployment, fluoroscopically position radiopaque markers on inner catheter beyond extremities of stricture to be crossed.
4. Confirm desired stent position fluoroscopically and to enable stent deployment remove the red safety guard from the handle. **(Fig.1)**
5. Continue deploying stent by squeezing the trigger. **NOTE:** Each trigger squeeze will deploy the stent by an equal amount.
6. If stent repositioning is required during deployment, it is possible to recapture the stent. **NOTE:** It is not possible to recapture the stent after passing point-of-no-return, indicated when the red marker on top of the handle has passed the point-of-no-return position on the handle label. **(Fig.2)**
7. To recapture the stent, push the direction button on the side of the delivery system to the opposite side. **(Fig.3)**
NOTE: Hold thumb on button when squeezing the trigger for the first time to recapture.
8. Continue squeezing the trigger as required, to recapture the stent by desired amount.
9. To resume deployment, push the button to the opposite side again and hold thumb on the button for first stroke while squeezing the trigger. **(Fig.4)**
10. When stent point-of-no-return has been passed, pull and remove the safety wire completely out of the delivery handle near the wire guide port to allow full stent deployment. **(Fig.5)**
11. Continue deploying the stent by squeezing the trigger.
12. After deployment, fluoroscopically confirm full stent expansion. Once full expansion is confirmed, the delivery system and wire guide can be safely removed.
NOTE: You may introduce endoscope and advance to the upper margin of the stent to endoscopically confirm position and patency of the stent.

Stent Repositioning Technique

To reposition the stent directly after placement, grip either the proximal or distal-end green lasso with forceps and reposition the stent to the desired area. **(Fig.6)**

Stent Removal Technique (Benign Indications Only)

1. Locate the proximal and distal-end green lasso of the stent by direct endoscopic visualization.
2. With forceps, gently apply tension to each green lasso to mobilize the distal and then the proximal end of the stent from the esophageal wall. **(Fig.7)**
3. Once the stent has been mobilized, grasp the green lasso at proximal end of stent securely and pull gently to collapse proximal end of stent and gently remove the stent. **(Fig.8)**

Post Procedure

- Observe the patient for development of any complications of endoscopy, esophageal dilation or stent placement. During the first 24 hours, vital signs should be monitored and clear liquid given in an upright position. Patients being treated for fistula should have no fluids or foods by mouth until successful sealing of fistula has been confirmed.
- Subsequent follow-ups (such as endoscopy and/or chest x-rays) with patient should be performed at the physician's discretion.

Patient Counselling Information

- After stent placement, post 24 hours, the patient should be instructed to chew food well or eat soft or pureed food and to drink fluids before and after meals.
- Patients with stents placed in the distal esophagus or across the GE junction may be instructed to elevate the head of the bed, and may be prescribed acid suppression therapy to minimize gastric reflux into the stent.

HOW SUPPLIED

Store in a dry location, away from temperature extremes.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2020 ©
COOK IRELAND LTD.