**EN**  
**5**
**EchoTip ProCore® HD Ultrasound Biopsy Needle**

Instructions for Use

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Please review prior to use.

### **INTENDED USE**

This device is used with an ultrasound endoscope for fine needle biopsy, (FNB), of submucosal lesions, mediastinal masses, lymph nodes and intraperitoneal masses within or adjacent to the gastrointestinal tract.

### **NOTES**

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

### **CONTRAINDICATIONS**

Those specific to primary endoscopic procedure to be performed in gaining access to desired site.

Coagulopathy.

### **POTENTIAL COMPLICATIONS**

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, damage to blood vessels, nerve damage, and acute pancreatitis. Those associated with EUS needle biopsy include but are not limited to: pain, death, peritonitis, portal vein gas and thrombosis, pneumoperitoneum and tumor seeding of the needle tract.

### **WARNINGS**

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.

## PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Ensure the stylet is fully inserted when advancing the needle into the biopsy site.

When targeting multiple sites, replace device for each site.

Needle must be retracted into sheath and thumbscrew on safety ring must be locked to hold needle in place prior to introduction, advancement or withdrawal of device. Failure to retract needle may result in damage to endoscope.

## SYSTEM PREPARATION

### Illustrations

**1.** Examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when stopcock is in “open” position, aligned with syringe. (See fig. 1)

**2.** Prepare syringe as follows:

- a. With stopcock in “open” position, depress plunger locks and fully advance plunger into syringe.
- b. Turn stopcock 90° to “closed” position.
- c. Pull plunger back until it is locked into place at desired setting, creating suction.
- d. Set prepared syringe aside until biopsy is desired.

**3.** Advance device into ultrasound endoscope to determine preferred sheath length. To adjust length, loosen thumbscrew lock on sliding sheath adjuster and slide until preferred length is attained. **Note:** Reference mark for sheath length will appear in sliding sheath adjuster window. (See fig. 2) Tighten thumbscrew on sliding sheath adjuster to maintain preferred sheath length.

## INSTRUCTIONS FOR USE

**1.** Identify desired biopsy site by endoscopic ultrasound.

**2.** Introduce ultrasound needle into accessory channel and advance device in small increments until Luer lock fitting at base of sliding sheath adjuster meets fitting on accessory channel. (See fig. 3)

**3.** Attach device to accessory channel port by rotating device handle until fittings are connected.

**4.** With ultrasound endoscope and device straight, adjust needle to desired length by loosening thumbscrew on safety ring, and advancing

it until desired reference mark for needle advancement appears in the window of safety ring. (*See fig. 4*) Tighten thumbscrew to lock safety ring in place. **Note:** Number in safety lock ring window indicates extension of needle in centimeters. **Caution:** During needle adjustment or extension, ensure device is attached to accessory channel. Failure to attach device prior to needle adjustment or extension may result in damage to endoscope.

**5.** While maintaining position of ultrasound endoscope, extend needle by advancing needle handle to pre-positioned safety ring.

**6.** Advance needle into lesion.

**7.** Remove stylet from needle by gently pulling back on plastic hub seated in metal fitting of needle handle.

**8.** Attach Luer lock fitting of previously prepared syringe securely onto fitting on handle.

**9.** Turn stopcock to “open” position aligned with syringe, allowing negative pressure in syringe to facilitate biopsy. Gently move handle in small increments back and forth within biopsy site. **Note:** Do not remove needle from biopsy site during biopsy.

**10.** Upon completion of biopsy, turn stopcock to “close” position.

**11.** Retract needle into sheath by pulling back on handle shaft. Lock thumbscrew on safety ring at 0 cm mark.

**12.** Disconnect Luer lock fitting from accessory channel and withdraw entire device from ultrasound endoscope.

**13.** Unlock thumbscrew on safety ring and extend needle. Exchange locking syringe for an empty syringe; attach to Luer lock fitting on handle.

**14.** Push in plunger to expel specimen, then prepare per institutional guidelines.

**15.** Remaining material may be recovered for examination by reinserting stylet or flushing device.

**Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.**

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© COPYRIGHT COOK 2014 ©
**COOK IRELAND LTD.**