
EN 5	EchoTip ProCore® HD Ultrasound Biopsy Needle Instructions for Use
BG 9	Игла за биопсия под ехографски контрол с висока разграничителна способност EchoTip ProCore® Инструкции за употреба
CS 13	Bioptická jehla EchoTip ProCore® používající ultrazvuk s vysokým rozlišením Návod k použití
DA 17	EchoTip ProCore® ultralydbiopsinål med høj opløsning Brugsanvisning
DE 21	EchoTip ProCore® HD Ultraschall-Biopsienadel Gebrauchsanweisung
EL 25	Βελόνα βιοψίας υπερήχων υψηλής ευκρίνειας EchoTip ProCore® Οδηγίες χρήσης
ES 29	Aguja de biopsia por ecografía de alta definición EchoTip ProCore® Instrucciones de uso
ET 33	HD-ultraheli-biopsianõel EchoTip ProCore® Kasutusjuhised
FR 37	Aiguille de biopsie ultrasonore à haute définition EchoTip ProCore® Mode d'emploi
HR 41	EchoTip ProCore® igla za HD ultrazvučnu biopsiju Upute za uporabu
HU 45	EchoTip ProCore® nagy felbontású ultrahangos biopsziás tű Használati utasítás
IT 49	Ago di biopsia per ecografia ad alta definizione EchoTip ProCore® Istruzioni per l'uso
LT 53	EchoTip ProCore® didelės skiriamosios gebos ultragarsinė biopsijos adata Naudojimo nurodymai
LV 57	EchoTip ProCore® augstas izšķirtspējas ultraskaņas biopsijas adata Lietošanas instrukcija
NL 78	EchoTip ProCore® HD ultrasone biopsienaald Gebruiksaanwijzing
NO 63	EchoTip ProCore® HD ultralydbiopsinål Bruksanvisning
PL 66	Igła biopsyjna do USG o wysokiej rozdzielczości EchoTip ProCore® Instrukcja użycia
PT 69	Agulha de biopsia de ultra-sons HD EchoTip ProCore® Instruções de utilização
RO 72	Ac ecografic HD pentru biopsie EchoTip ProCore® Instrucțiuni de utilizare
SK 75	Ultrazvuková bioptická ihla s vysokým rozlišením EchoTip ProCore® Návod na použitie
SV 78	EchoTip ProCore® ultraljudsbiopsinål med hög upplösning Bruksanvisning
TR 81	EchoTip ProCore® HD Ultrason Biyopsi İğnesi Kullanma Talimatı


EN Plunger

BG	Буталo
CS	Pist
DA	Stempel
DE	Kolben
EL	Έμβολο
ES	Émbolo
ET	Kolb
FR	Piston
HR	Klip
HU	Dugattyú
IT	Stantuffo
LT	Stūmoklis
LV	Virzulis
NL	Plunjer
NO	Stempel
PL	Tłoczek
PT	Émbolo
RO	Piston
SK	Piest
SV	Kolv
TR	Piston

EN Stopcock "open"

BG	Кранче "отворено"
CS	„Otevřený“ uzavírací kohout
DA	Stophane "åben"
DE	Absperrhahn geöffnet
EL	Στρόφιγγα ανοικτή
ES	Llave de paso «abierta»
ET	Sulgemisventiil avatud
FR	Robinet ouvert
HR	Sigurnosni ventil otvoren
HU	Elzárócsap nyitott állásban
IT	Rubinetto aperto
LT	Kranelis „atviras“
LV	Slēgkrāns atvērts
NL	Afsluitkraan 'open'
NO	Stoppekran åpen
PL	Kranik w pozycji otwartej
PT	Torneira de passagem "aberta"
RO	Robinet „deschis“
SK	Otvorený kohútik
SV	Kran "öppen"
TR	Stopkok "açık"

EN Sliding Sheath Adjuster

BG	Плъзгащ се регулатор на защитната обвивка
CS	Posuvné nastavovací zařízení sheathu
DA	Glidende sheath-justering
DE	Verschiebbarer Schleuseneinsteller
EL	Συρόμενος προσαρμογέας θηκαρίου
ES	Ajustador deslizante de la vaina
ET	Liughülsi regulaator
FR	Ajusteur de gaine coulissant
HR	Klizni podešivac uvodnice
HU	Elcsúsztható hüvelyhossz-beállító
IT	Regolatore scorrevole della guaina
LT	Slankusis vamzdelio reguliatorius
LV	Pārbidāmais ievadpalka regulētājs
NL	Schuifregeling van de sheath
NO	Glidende hylsejusteringsmekanisme
PL	Przesuwany regulator koszulki
PT	Regulador deslizante da bainha
RO	Ajustor glisant al tecii
SK	Posuvný nastavovač puzdra
SV	Glidande hylsjusterare
TR	Kayar Kilif Ayarlayıcı

EN Sheath Reference Mark

BG	Означение за дължината на защитната обвивка
CS	Referenční značka sheathu
DA	Referencemærke på sheath
DE	Schleusenbezugsmarke
EL	Σήμανση αναφοράς θηκαρίου
ES	Marca de referencia de la vaina
ET	Hülsi tähis
FR	Repère de référence sur la gaine
HR	Referentna oznaka uvodnice
HU	Hüvelyhossz referencijele
IT	Indicatore di riferimento della guaina
LT	Vamzdelio atskaitos žymė
LV	Ievadpalka atskaites atzīme
NL	Referentiemarkering van de sheath
NO	Referansemerke på hylsen
PL	Znacznik referencyjny wysunięcia koszulki
PT	Marca de referência da bainha
RO	Marcaj de referință al tecii
SK	Referenčná značka puzdra
SV	Hylsans referensmarkering
TR	Kilif Referans İşareti

EN Needle Handle

BG	Дръжка на иглата
CS	Rukojeť jehly
DA	Nålegreb
DE	Nadelgriff
EL	Λαβή βελόνας
ES	Mango de la aguja
ET	Nõela käepide
FR	Poignée de l'aiguille
HR	Drška igle
HU	Tű fogantyúja
IT	Impugnatura dell'ago
LT	Adatos kotelis
LV	Adatas rokturis
NL	Veiligheidsring
NO	Nålehåndtak
PL	Uchwył igły
PT	Punho da agulha
RO	Mânerul acului
SK	Rukoväť ihly
SV	Nålhandtag
TR	İğne Sapi


EN Safety Ring

BG	Предпазен пръстен
CS	Pojistný prstenec
DA	Sikkerhedsring
DE	Sicherheitsring
EL	Δακτύλιος ασφαλείας
ES	Anillo de seguridad
ET	Kaitserõngas
FR	Bague de sécurité
HR	Sigurnosni prsten
HU	Biztonsági gyűrű
IT	Anello di sicurezza
LT	Apsauginis žiedas
LV	Drošības gredzens
NL	Veiligheidsring
NO	Sikkerhetsring
PL	Pierścień bezpieczeństwa
PT	Anel de segurança
RO	Inel de siguranță
SK	Bezpečnostný krúžok
SV	Säkerhetslåsring
TR	Güvenlik Halkası

EN Stylet Hub

BG	Втулка на стилето
CS	Nástavec stiletu
DA	Stiletmuffe
DE	Mandrinansatz
EL	Ομφαλός του στείλεου
ES	Conector del estilete
ET	Stiletikinnitus
FR	Embase du stylet
HR	Cvorište sonde
HU	Mandrin végidoma
IT	Pomello del mandrino
LT	Zondo kaištis
LV	Stiletas galviņa
NL	Aanzetstuk van het stilet
NO	Stilletfeste
PL	Głowica mandryny
PT	Conector do estilete
RO	Acoperitoarea stiletului
SK	Hrdlo styletu
SV	Mandringfattning
TR	Stile Göbeği

3.

4.

ECHOTIP PROCORE® HD ULTRASOUND BIOPSY NEEDLE

IMPORTANT INFORMATION

Please review prior to use.

INTENDED USE

This device is used with an ultrasound endoscope for fine needle biopsy, (FNB), of submucosal lesions, mediastinal masses, lymph nodes and intraperitoneal masses within or adjacent to the gastrointestinal tract.

NOTES

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired site.

Coagulopathy.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, damage to blood vessels, nerve damage, and acute pancreatitis. Those associated with EUS needle biopsy include but are not limited to: pain, death, peritonitis, portal vein gas and thrombosis, pneumoperitoneum and tumor seeding of the needle tract.

WARNINGS

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Ensure the stylet is fully inserted when advancing the needle into the biopsy site.

When targeting multiple sites, replace device for each site.

Needle must be retracted into sheath and thumbscrew on safety ring must be locked to hold needle in place prior to introduction, advancement or withdrawal of device. Failure to retract needle may result in damage to endoscope.

SYSTEM PREPARATION

Illustrations

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (See fig. 1)
2. Prepare syringe as follows:
 - a. With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b. Turn stopcock 90° to "closed" position.
 - c. Pull plunger back until it is locked into place at desired setting, creating suction.
 - d. Set prepared syringe aside until biopsy is desired.
3. Advance device into ultrasound endoscope to determine preferred sheath length. To adjust length, loosen thumbscrew lock on sliding sheath adjuster and slide until preferred length is attained. **Note:** Reference mark for sheath length will appear in sliding sheath adjuster window. (See fig. 2) Tighten thumbscrew on sliding sheath adjuster to maintain preferred sheath length.

INSTRUCTIONS FOR USE

1. Identify desired biopsy site by endoscopic ultrasound.
2. Introduce ultrasound needle into accessory channel and advance device in small increments until Luer lock fitting at base of sliding sheath adjuster meets fitting on accessory channel. (See fig. 3)
3. Attach device to accessory channel port by rotating device handle until fittings are connected.
4. With ultrasound endoscope and device straight, adjust needle to desired length by loosening thumbscrew on safety ring, and advancing it until desired reference mark for needle advancement appears in the window of safety ring. (See fig. 4) Tighten thumbscrew to lock safety ring in place. **Note:** Number in safety lock ring window indicates extension of needle in centimeters. **Caution:** During needle adjustment or extension, ensure device is attached to accessory channel. Failure to attach device prior to needle adjustment or extension may result in damage to endoscope.
5. While maintaining position of ultrasound endoscope, extend needle by advancing needle handle to pre-positioned safety ring.
6. Advance needle into lesion.
7. Remove stylet from needle by gently pulling back on plastic hub seated in metal fitting of needle handle.

8. Attach Luer lock fitting of previously prepared syringe securely onto fitting on handle.
9. Turn stopcock to "open" position aligned with syringe, allowing negative pressure in syringe to facilitate biopsy. Gently move handle in small increments back and forth within biopsy site. **Note:** Do not remove needle from biopsy site during biopsy.
10. Upon completion of biopsy, turn stopcock to "close" position.
11. Retract needle into sheath by pulling back on handle shaft. Lock thumbscrew on safety ring at 0 cm mark.
12. Disconnect Luer lock fitting from accessory channel and withdraw entire device from ultrasound endoscope.
13. Unlock thumbscrew on safety ring and extend needle. Exchange locking syringe for an empty syringe; attach to Luer lock fitting on handle.
14. Push in plunger to expel specimen, then prepare per institutional guidelines.
15. Remaining material may be recovered for examination by reinserting stylet or flushing device.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COOK 2021
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.