
EN 5	EchoTip ProCore® HD Ultrasound Biotip Needle Instructions for Use
CS 9	Bioptická jehla EchoTip ProCore® používající ultrazvuk s vysokým rozlišením Návod k použití
DA 13	EchoTip ProCore® HD ultralydsbiopsinål Brugsanvisning
DE 17	EchoTip ProCore® HD Ultraschall-Biopsienadel Gebrauchsanweisung
EL 21	Βελόνα βιοψίας υπερήχων υψηλής ευκρίνειας EchoTip ProCore® Οδηγίες χρήσης
ES 26	Aguja de biopsia por ecografía de alta definición EchoTip ProCore® Instrucciones de uso
ET 31	HD-ultraheli-biopsianõel EchoTip ProCore® Kasutusjuhend
FI 35	EchoTip ProCore® -HD-ultraäänibiopsianeula Käyttöohjeet
FR 39	Aiguille de biopsie ultrasonore EchoTip ProCore® HD Mode d'emploi
HR 43	HD igla za ultrazvučnu biopsiju EchoTip ProCore® Upute za uporabu
HU 47	EchoTip ProCore® nagy felbontású echogén biopsziás tű Használati utasítás
IT 52	Ago ecogenico per biopsia ad alta definizione EchoTip ProCore® Istruzioni per l'uso
LT 56	„EchoTip ProCore®“ didelės skiriamosios gebos ultragarso biopsijos adata Naudojimo instrukcija
LV 60	EchoTip ProCore® augstas precizitātes ultraskaņas biopsijas adata Lietošanas instrukcija
NL 64	EchoTip ProCore® HD ultrasone biopsienaald Gebruiksaanwijzing
NO 69	EchoTip ProCore® HD-ultralydbiopsinål Bruksanvisning
PL 73	Ultraźwiękowa igła biopsyjna o wysokiej rozdzielczości EchoTip ProCore® Instrukcja uźycia
PT 77	Agulha de biopsia de ultrassons HD EchoTip ProCore® Instruções de Utilização
RO 82	Ac ecografic HD pentru biopsie EchoTip ProCore® Instrucțiuni de utilizare
SK 86	Ultrazvuková bioptická ihla s vysokým rozlišením EchoTip ProCore® Návod na použitie
SL 90	Biopsijska igla za HD-ultrazvočne preiskave EchoTip ProCore® Navodila za uporabo
SV 95	EchoTip ProCore® ultraljudsbiopsinål med hög upplösning Bruksanvisning
TR 99	EchoTip ProCore® HD Ultrason Biyopsi İğnesi Kullanma Talimatı


EN Plunger

CS	Píst
DA	Stempel
DE	Kolben
EL	Έμβολο
ES	Émbolo
ET	Kolb
FI	Mäntä
FR	Poussoir
HR	Klip
HU	Dugattyú
IT	Stantuffo
LT	Stūmoklis
LV	Virzulis
NL	Plunjer
NO	Stempel
PL	Tłoczek
PT	Émbolo
RO	Piston
SK	Piest
SL	Bat
SV	Kolv
TR	Piston

EN Stopcock "open"

CS	„Otevřený“ uzavírací kohout
DA	Stophane "åben"
DE	Absperrhahn „geöffnet“
EL	Στρόφιγγα «ανοικτή»
ES	Llave de paso «abierta»
ET	Sulgemisventiil avatud
FI	Sulkuhana "auki"
FR	Robinet ouvert
HR	Sigurnosni ventil otvoren
HU	Elzárócsap nyitott állásban
IT	Rubinetto "aperto"
LT	„Atsuktas“ čiaupas
LV	Noslēgkrāns atvērts
NL	Afsluitkraan 'open'
NO	Stoppekran «åpen»
PL	Kranik w pozycji „otwartej”
PT	Torneira de passagem "aberta"
RO	Robinet „deschis”
SK	Otvorený kohútik
SL	Zapiralni ventil „odprt”
SV	Kran "öppen"
TR	Stopkok "açık"

EN Sliding Sheath Adjuster

CS	Posuvné nastavovací zařízení sheathu
DA	Glidende sheathregulator
DE	Verschiebbarer Schleuseneinsteller
EL	Συρόμενος προσαρμογέας θηκαρίου
ES	Ajustador deslizante de la vaina
ET	Liughülsi regulaator
FI	Holkin liukusäädin
FR	Ajusteur de gaine coulissant
HR	Klizni podešivač uvodnice
HU	Elcsúsztható hüvelyhossz-beállító
IT	Regolatore scorrevole della guaina
LT	Slankusis vamzdelio regulatorius
LV	Pārbidāmais ievadapvalka regulētājs
NL	Schuifregeling van de sheath
NO	Glidende hylsejusteringsmekanisme
PL	Przesuwany regulator koszulki
PT	Regulador deslizante da bainha
RO	Ajustor glisant al tecii
SK	Posuvný nastavovač puzdra
SL	Drčni nastavjalnik tulca
SV	Glidande hylsjusterare
TR	Kayar Kilif Ayarlayıcı

EN Sheath Reference Mark

CS	Referenční značka sheathu
DA	Referencemærke på sheath
DE	Schleusenbezugsmarke
EL	Σήμανση αναφοράς θηκαρίου
ES	Marca de referencia de la vaina
ET	Hülsi tähis
FI	Holkin viitemerkintä
FR	Repère de référence sur la gaine
HR	Referentna oznaka uvodnice
HU	Hüvelyhossz referencijale
IT	Indicatore di riferimento della guaina
LT	Vamzdelio atskaitos žymė
LV	Ievadapvalka atskaites atzīme
NL	Referentiemarkering van de sheath
NO	Referansemerke på hylsen
PL	Znacznik referencyjny wysunięcia koszulki
PT	Marca de referência da bainha
RO	Marcaj de referință al tecii
SK	Referenčná značka puzdra
SL	Referenčna oznaka tulca
SV	Hylsans referensmarkering
TR	Kilif Referans İşareti

EN Needle Handle



CS	Rukojeť jehly
DA	Nålegreb
DE	Nadelgriff
EL	Λαβή βελόνας
ES	Mango de la aguja
ET	Nõela käepide
FI	Neulan kahva
FR	Poignée de l'aiguille
HR	Drška igle
HU	Tű fogantyúja
IT	Impugnatura dell'ago
LT	Adatos kotelis
LV	Adatas rokrturis
NL	Naaldhandvat
NO	Nålhandtak
PL	Uchwyt igły
PT	Punho da agulha
RO	Mâner al acului
SK	Rukoväť ihly
SL	Ročaj igle
SV	Nålhandtag
TR	İğne Sapi

EN Safety Ring

CS	Pojistný prstenec
DA	Sikkerhedsring
DE	Sicherheitsring
EL	Διακότιλος ασφαλείας
ES	Anillo de seguridad
ET	Kaitserõngas
FI	Turvarengas
FR	Bague de sécurité
HR	Sigurnosni prsten
HU	Biztonsági gyűrű
IT	Anello di sicurezza
LT	Apsauginis žiedas
LV	Drošības gredzens
NL	Veiligheidsring
NO	Sikkerhetsring
PL	Pierścień bezpieczeństwa
PT	Anel de segurança
RO	Inel de siguranță
SK	Bezpečnostný krúžok
SL	Varnostni obroč
SV	Säkerhetslåsring
TR	Güvenlik Halkası

EN Stylet Hub

CS	Ústí styletu
DA	Stiletmuffe
DE	Mandrinansatz
EL	Ομφαλός του στείλεου
ES	Conector del estilete
ET	Stiletikinnitus
FI	Mandriinin kanta
FR	Embase du stylet
HR	Čvoršte sonde
HU	Mandrin végidoma
IT	Pomello del mandrino
LT	Zondo kaištis
LV	Stiletas galviņa
NL	Aanzetstuk van het stilet
NO	Stilettfeste
PL	Złączka mandrynu
PT	Conector do estilete
RO	Acoperitoarea stiletului
SK	Hrdlo styletu
SL	Pesto stileta
SV	Mandrängfattning
TR	Stile Göbeği

ECHOTIP PROCORE® HD ULTRASOUND BIOPSY NEEDLE

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The EchoTip ProCore® HD Ultrasound Biopsy Needle is available in four different needle sizes (19, 20, 22, and 25 gage) and contains the following components:

- Needle
- Syringe

Performance Characteristics

The function and key features of these components are described below:

Ultrasound Biopsy Needle

- Handle – The handle contains adjustable components to allow the user to adjust the extension of the needle and the sheath. The handle safety ring slides and locks the needle at the desired needle length extension (0–8 cm). The '0 cm' reference mark ensures complete needle retraction within the sheath. The sliding sheath adjuster allows the user to set the required sheath length extension (0–5 cm) for their ultrasound endoscope (EUS). The '0 cm' reference mark denotes the minimum sheath length available. The base of the sliding sheath adjuster has a Luer lock fitting, which facilitates the attachment of the device to the scope accessory channel.
- Sheath – The sheath covers the needle during introduction, advancement, and withdrawal of the device from the scope accessory channel. It protects the needle and the scope accessory channel from damage by the needle tip. ECHO-HD-3-20-C has a flexible coiled sheath (PTFE-coated stainless steel coil spring) that includes a nylon covering. Table 1 outlines the sheath sizes (Fr) for each reference part number.

Table 1: Sheath sizes (Fr) for reference part numbers.

Reference Part Number (RPN)	Sheath Size (Fr)
ECHO-HD-19-C	4.8
ECHO-HD-3-20-C	7.0–7.95
ECHO-HD-22-C	5.2
ECHO-HD-25-C	5.2

- Needle – The stainless steel needle has a dimpling pattern to provide visualisation of the needle tip under endoscopic ultrasound. The cut-off section has been added to the distal tip of the needle cannula for biopsy purposes. The cutting edge of the notch is reverse facing for ECHO-HD-19-C, ECHO-HD-22-C, and ECHO-HD-25-C and forward facing for ECHO-HD-3-20-C. The needle cannula has a bevelled tip design to puncture and facilitate sample collection from the target site.
- Stylet – The stylet is inserted through the length of the needle. The distal end of the stylet is bevelled (ECHO-HD-25-C) to match the bevel of the needle tip or has a ball tip design (ECHO-HD-19-C, ECHO-HD-3-20-C, ECHO-HD-22-C). The stylet provides support to the needle during advancement into the target site, prevents sample collection outside of the target site, and may also be used for sample retrieval following collection if desired. The stylet is straight (ECHO-HD-19-C, ECHO-HD-22-C, ECHO-HD-25-C) or coiled (ECHO-HD-3-20-C) upon removal. The stylet is removed when aspiration is required at the target site.

Syringe – The device is supplied with a 10 mL syringe to utilize the standard vacuum technique for sample collection, if desired. The Luer lock fitting at the tip of the syringe can be attached to the Luer fitting on the handle. The syringe has a plunger and a stopcock to create a vacuum for sample collection. It is supplied with 5 mL and 10 mL clips attached to the plunger to allow the user to select the preferred volume for sample aspiration.

Device Compatibility

EchoTip ProCore HD Ultrasound Biopsy Needles are compatible with:

- Ultrasound endoscope with 2.8 mm minimum accessory channel (ECHO-HD-19-C, ECHO-HD-22-C, ECHO-HD-25-C)
- Ultrasound endoscope with 3.7 mm minimum accessory channel (ECHO-HD-3-20-C)
- 10 mL standard Luer syringe

Patient Population

Adult patients requiring fine needle biopsy (FNB) of submucosal lesions or tissue, mediastinal masses, lymph nodes, and intraperitoneal masses within or adjacent to the gastrointestinal tract.

Intended User

Use of this device is restricted to a trained healthcare professional.

Contact with body tissue

This device is tissue contacting in line with intended use.

Operating Principles

The device is inserted into the ultrasound endoscope accessory channel and attached when the fittings meet. The sheath length can be adjusted depending on the length of the EUS scope using the sliding sheath adjuster. The needle length extension can be adjusted as required. The device is used to puncture the target site to obtain a sample. If using a standard vacuum syringe technique, the stylet is then removed, and the syringe is attached to the Luer fitting on the needle handle. Turning the syringe stopcock to the "open" position allows negative pressure to aspirate the sample. The handle is gently moved in small increments back and forth within the target site to obtain a sample. If using the syringe, the stopcock is returned to the "closed" position. Upon completion of the procedure the needle is retracted into the sheath and the device is disconnected from the scope accessory channel.

INTENDED USE

This device is used with an ultrasound endoscope for fine needle biopsy (FNB), of submucosal lesions or tissue, mediastinal masses, lymph nodes, and intraperitoneal masses within or adjacent to the gastrointestinal tract.

INDICATIONS FOR USE

This device is indicated for ultrasound guided tissue sampling of submucosal lesions or tissue, mediastinal masses, lymph nodes, and intraperitoneal masses within or adjacent to the GI tract, which is used for subsequent pathological examination to aid in disease diagnosis and patient management.

CLINICAL BENEFITS

This device facilitates successful sampling of submucosal lesions or tissue, mediastinal masses, lymph nodes, and intraperitoneal masses within or adjacent to the gastrointestinal tract to aid in disease diagnosis and patient management.

CONTRAINDICATIONS

- Those specific to primary endoscopic procedure to be performed to gain access to the target site.
- Coagulopathy.

WARNINGS

- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect the device with particular attention to kinks, bends, and breaks. If an abnormality is detected that would prohibit proper working condition do not use. Please notify Cook Medical for return authorization.
- Not for use in the heart or vascular system.
- The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.
- This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.

- This device contains nickel, which may cause allergic reaction in individuals with nickel sensitivity.

PRECAUTIONS

- Refer to package label for minimum channel size required for this device.
- Ensure the stylet is fully inserted when advancing the needle into the target site.
- When targeting multiple sites, replace the device for each site.
- Needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold needle in place prior to introduction, advancement, or withdrawal of the device. Failure to retract needle may result in damage to the endoscope.
- Device not indicated for use in pregnant/nursing women.
- Do not use this device for any purpose other than the stated intended use.
- Use of this device is restricted to a trained healthcare professional.
- The product is intended for use by physicians trained and experienced in ultrasound endoscopy techniques.
- Store the device in a dry location.

POTENTIAL ADVERSE EVENTS

Those associated with gastrointestinal endoscopic ultrasound procedure: allergic reaction to medication • aspiration • cardiac arrhythmia or arrest • damage to blood vessels • fever • hemorrhage • hypotension • infection • pain/discomfort • perforation • respiratory depression or arrest.

Those associated with the device: acute pancreatitis • allergic reaction to nickel • fever • hemorrhage • infection • pain/discomfort • perforation • peritonitis • portal vein gas and thrombosis • pneumoperitoneum • tumor seeding (through the needle tract).

HOW SUPPLIED

These devices are supplied ethylene oxide (EO) sterilized in a tray with a Tyvek lid.

SYSTEM PREPARATION

Illustrations

1. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
2. Visually inspect the device components confirming they are free from damage prior to use.
3. Examine syringe. It has two plunger locks that must be depressed to advance the plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when the stopcock is in the “open” position, aligned with the syringe. (**Fig. 1**)
4. Prepare the syringe as follows:
 - a. With the stopcock in the “open” position, depress the plunger locks and fully advance the plunger into syringe.
 - b. Turn the stopcock 90° to the “closed” position.
 - c. Pull the plunger back until it is locked into place at the desired setting, creating vacuum.
 - d. Set the prepared syringe aside until biopsy is required.

INSTRUCTIONS FOR USE

1. Identify desired target site by endoscopic ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope. **Caution:** If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device in small increments until Luer lock fitting at the base of sliding sheath adjuster meets the Luer fitting on the accessory channel port. (**Fig. 3**)
4. Attach the device to the scope accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference

- mark for the sheath length will appear in the sliding sheath adjuster window. Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining position of the scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring. (**Fig. 4**) Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
 7. Extend the needle by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
 8. Standard vacuum syringe techniques may be applied for biopsy (steps 9–11) or, if desired, other techniques that may or may not incorporate the use of the stylet may be used.
 9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of needle handle. Preserve the stylet for use if additional biopsy is required later.
 10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer fitting on the needle handle.
 11. Turn the stopcock to the “open” position aligned with syringe, allowing negative pressure in the syringe to facilitate biopsy. Move the needle handle back and forth in small increments within the target site. **Note:** Do not remove the needle from the target site during biopsy.
 12. Upon completion of biopsy, if using a syringe, turn the stopcock to the “closed” position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
 13. Disconnect the Luer lock fitting of the device from accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
 14. Unlock the thumbscrew on the safety ring and extend the needle.
 15. Use a new syringe or the stylet to expel the specimen, then prepare per institutional guidelines.
 16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.
 17. For additional biopsy from the same target site, gently reinsert the stylet into the Luer fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.
 18. Additional samples may be obtained by repeating steps 2 through 16 of the “Instructions for Use”.
 19. Upon completion of procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.

DEVICE DISPOSAL

This device may be contaminated with potentially infectious substances of human origin and should be coiled for disposal in accordance with institutional guidelines.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2023
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.