EN
4
OASIS® One Action Stent Introduction System

Instructions for Use

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

INTENDED USE

This device (with preloaded stent, if applicable) is intended for endoscopic biliary stent placement to drain obstructed bile ducts.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications associated with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Wire guide diameter and inner lumen of wire-guided device must be compatible.

Not compatible with pigtail stents.

Do not use this device if stent cannot advance through strictured area.

Preloaded stents should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

SYSTEM PREPARATION

1. If applicable, pre-load desired stent and positioning sleeve onto tip of introducer.

INSTRUCTIONS FOR USE

Illustrations

I. If using Intra Ductal Exchange (IDE) port and a pre-positioned short wire. (See fig. 1)

1. Unlock short wire from Wire Guide Locking Device and advance introducer with preloaded stent over pre-positioned wire guide ensuring that wire guide exits guiding catheter at IDE port.
2. Introduce device into endoscope accessory channel and relock wire guide.
3. With elevator open, advance device in short increments until endoscopically visualized exiting scope. **Note:** As front flap of stent enters accessory channel, keep positioning sleeve over back flap until completely in accessory channel. Slide positioning sleeve over pushing catheter, keeping it clear of accessory channel.
4. Advance device into appropriate duct.
5. Fluoroscopically monitor radiopaque band(s) on guiding catheter. When a sufficient length of guiding catheter is above majority of stricture, disconnect Luer lock fitting on introducer. **Important:** A sufficient portion of guiding catheter must be above stricture prior to stent positioning.
6. Maintain position of guiding catheter and advance stent into desired position using pushing catheter.
7. Fluoroscopically and endoscopically confirm desired stent placement.
8. Unlock wire guide from Wire Guide Locking Device.
9. Retract wire guide until it exits guiding catheter at IDE port.

Note: Fluoroscopically visualize radiopaque band at IDE port. When radiopaque distal tip of wire guide passes band, a disengagement from wire guide lumen will occur.

10. Advance disengaged wire guide to maintain ductal access.
11. Relock wire guide into Wire Guide Locking Device.
12. Retract guiding catheter into endoscope while maintaining position of pushing catheter.
13. Confirm position of stent once guiding catheter is removed.
14. Remove pushing catheter.

II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide. (See fig. 2)

Note: For best results, wire guide should be kept wet.

1. Remove stylet, if applicable.
2. Advance introducer with preloaded stent over pre-positioned wire guide and into accessory channel of endoscope until proximal end of wire guide exits PWP. **Note:** If applicable, ensure that wire guide does not exit at IDE port.

REFER TO STEPS 3-7 IN "SECTION I", THEN RESUME WITH STEP 3 BELOW:

3. Retract and remove both guiding catheter and wire guide into endoscope while maintaining position of pushing catheter.
4. Confirm position of stent once guiding catheter and wire guide are removed.
5. Remove pushing catheter.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2013
COOK IRELAND LTD.