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
EN Figure 1: OASIS® One Action Stent Introduction System with Preloaded Cotton-Leung® Biliary Stent (OACL)

- BG Фигура 1: Система OASIS® за едностъпково въвеждане на стент с предварително зареден билиарен стент Cotton-Leung® (OACL)
- CS Obrázek 1: Jednofázový zaváděcí systém stentu OASIS® s předinstalovaným biliárním stentem Cotton-Leung® (OACL)
- DA Figur 1: OASIS® stentindføringsystem med én funktion med forudladet Cotton-Leung® galdestent (OACL)
- DE Abbildung 1: OASIS® One Action Stent-Einführsystem mit vorgeladenem Cotton-Leung® Gallengangstent (OACL)
- EL Εικόνα 1: Σύστημα εισαγωγής ενδοπρόσθεσης μίας ενέργειας OASIS® με προτοποθετημένη ενδοπρόσθεση χοληφόρων Cotton-Leung® (OACL)
- ES Figura 1: Sistema de introducción de stents en un paso OASIS® con stent biliar Cotton-Leung® precargado (OACL)
- ET Joonis 1. Stendi ühesammuline sisestussüsteem OASIS® eelnevalt laaditud biliaarse stendiga Cotton-Leung® (OACL)
- FR Figure 1 : Système d'introduction d'endoprothèse en un temps OASIS® avec endoprothèse biliaire Cotton-Leung® préchargée (OACL)
- HR Slika 1: Sustav za uvođenje stenta OASIS® One Action s unaprijed umetnutim bilijarnim stentom Cotton-Leung® (OACL)
- HU 1. ábra: OASIS® egyműveletű sztentfelvezető rendszer előre betöltött Cotton-Leung® epeúti sztenttel (OACL)
- IT Figura 1. Sistema di introduzione di stent in un'unica azione OASIS® con stent biliare Cotton-Leung® (OACL) precaricati
- LT 1 paveikslas. „OASIS® One Action” stento įvedimo sistema su iš anksto įstatytu „Cotton-Leung®” tulžies latakų stentu (OACL)
- LV 1. attēls. OASIS® One Action stenta ievadišanas sistēma ar iepriekš ievietotu Cotton-Leung® žultsvadu stentu (OACL)
- NL Afbeelding 1: OASIS® One Action-stentintroductiesysteem met voorgeladen Cotton-Leung® biliaire stent (OACL)
- NO Figur 1: OASIS® One Action stentinnføringsystem med forhåndspåsett Cotton-Leung® gallestent (OACL)
- PL Rysunek 1: OASIS® jednostopniowy system wprowadzania stentu z wstępnie załadowanym stentem do dróg żółciowych Cotton-Leung® (OACL)
- PT Figura 1: Sistema introdutor de stent com um só movimento OASIS® com stent biliar Cotton-Leung® pré-carregado (OACL)
- RO Figura 1: Sistem de introducerea a stentului OASIS® One Action cu stent biliar Cotton-Leung® preîncărcat (OACL)
- SK Obrázok 1: Zavádzací systém stentov OASIS® One Action s vopred založeným biliárnym stentom Cotton-Leung® (OACL)
- SL Slika 1: Uvajalni sistem stenta z enim gibom OASIS® s prednastavljenim stentom za žolčevod Cotton-Leung® (OACL)
- SR Slika 1: OASIS® Jednopotezni sistem za uvođenje stenta sa unapred umetnutim „Cotton-Leung®” bilijarnim stentom (OACL)
- SV Figur 1: OASIS® One Action stentinföringsystem med förinstallerad Cotton-Leung® biliär stent (OACL)
- TR Şekil 1: Önceden Yüklemiş Cotton-Leung® Biliyer Stentli OASIS® Tek İşlemli Stent Yerleştirme Sistemi (OACL)

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
SR Slika 3: Zaštita krilaca

SV Figur 3: Flikskydd

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OASIS® ONE ACTION STENT INTRODUCTION SYSTEMS (PRELOADED)

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

This IFU covers the OASIS® One Action Stent Introduction Systems preloaded with a Cotton-Leung Biliary Stent (OACL) or preloaded with an ST-2 Soehendra Tannenbaum Biliary Stent (OATS) as represented in Figures 1 and 2. Only stents supplied preloaded and their associated delivery systems are covered by this IFU.

OACL – This device is sold as a set consisting of a stent and a dedicated introduction system (a guiding catheter and a pushing catheter). The stents are available in 7–11.5 Fr diameters and lengths of 5–15 cm. The stent has a center bend with one flap on each end of the stent and a sideport on the ductal end. The stent is provided in a preloaded configuration on the OASIS introduction system.

OATS – This device is sold as a set consisting of a stent and a dedicated introduction system (a guiding catheter and a pushing catheter). The stents are available in 8.5–11.5 Fr diameters and lengths of 5–15 cm. The stent has a center bend with four flaps on each end of the stent. The stent is provided in a clamshell tray with the OASIS introduction system.

The introduction system consists of a guiding catheter component of length 205 cm and a pushing catheter component of length 170 cm. The guiding catheter has four radiopaque marker bands.

Performance Characteristics

The function and key features of the device are described below.

A variety of biliary stents that provide drainage are available across the device range to accommodate various patient anatomies. They are offered in French sizes of between 7 and 11.5 Fr, in various shape configurations, and in labelled lengths of between 5 and 15 cm. Biliary stents have the following features:

- Atraumatic tip – The ends of the stent are designed to be atraumatic to the anatomy.
- Sideports – Aid drainage (not applicable to ST-2 Soehendra Tannenbaum).
- Anti-migration features – Stent flaps aim to help maintain the position of the stent once it is placed by providing mechanical interference with the anatomy.
- Radiopacity – The stent materials are radiopaque to facilitate visibility under fluoroscopy.

Flap protector – This is a tube used to collapse the flap(s) as it is introduced into the endoscope.

Guiding catheter – The function of the guiding catheter is to guide the stent to its intended location.

The guiding catheter contains the following features:

- Atraumatic tip – The end of the guiding catheter is designed to be atraumatic to the anatomy.
- Radiopacity – The guiding catheter tubing material is radiopaque. The guiding catheter also contains marker bands to further enhance its fluoroscopic visibility.
- Proximal wire port (hub) – Allows the wire guide to exit in the long wire configuration.

Pushing catheter – The function of the pushing catheter is to advance the stent over a pre-positioned wire guide and a guiding catheter to its intended location and to maintain the position of the stent as it is being deployed. It contains the following feature:

- Proximal identifier – An identifier at its proximal end. It is used to identify the orientation in which the device must be used.

Device Compatibility

These devices are compatible with the following:

- Endoscope with 3.2, 3.7, or 4.2 mm accessory channel (refer to label for selection of the appropriate accessory channel)
- 0.035" wire guide

- Endoscopic cap or wire guide locking device
- Stent retriever or forceps
- Contrast media
- Sterile water or saline
- Standard Luer syringe
- Water-soluble lubricant

NOTE: Preloaded stents are only compatible with the introduction systems with which they are supplied.

Qualitative and Quantitative information

The materials for the stent implant are outlined in Table 1.

Table 1: Stent implant materials

Product	Qualitative Information Device Material		Quantitative Information Weight (g)
7 Fr Cotton Leung®	Polyethylene	Stent	Up to 0.51
8.5 Fr Cotton Leung®	Polyethylene	Stent	Up to 0.76
8.5 Fr ST-2 Soehendra Tannenbaum®	Polytetrafluoroethylene	Stent	Up to 1.39
10 Fr ST-2 Soehendra Tannenbaum®	Polytetrafluoroethylene	Stent	Up to 1.42
10 Fr Cotton Leung®	Ethylene-vinyl acetate (EVA) copolymer	Stent	Up to 1.19
11.5 Fr Cotton Leung®	EVA copolymer	Stent	Up to 1.59
11.5 Fr ST-2 Soehendra Tannenbaum®	Polytetrafluoroethylene	Stent	Up to 1.92

Patient Population

Adult patients requiring biliary stenting for obstruction. The nature of the underlying pathology is prevalent in varying patients, therefore the devices are indicated for patients with biliary duct obstruction caused by common bile duct stones, malignant biliary obstruction, and benign or malignant strictures.

Intended User

This device is intended for use by physicians trained and experienced in endoscopic retrograde cholangiopancreatography (ERCP) techniques. Standard techniques for the endoscopic placement of stents should be employed.

Contact with Body Tissue

This device is tissue contacting in line with intended use.

Operating Principles

The biliary stent operates by providing a lumen through which biliary fluid can drain. The stent has an anti-migration feature (stent flaps) which provides mechanical interference with the anatomy to help maintain the position of the stent. The stent is radiopaque to facilitate visibility under fluoroscopy. A sideport (Cotton-Leung only) on the stent aids drainage by providing an additional opening for fluid flow. The biliary stent is placed over a pre-positioned wire guide using an introducer. An introducer consists of a pushing catheter and a guiding catheter combined. The guiding catheter contains marker bands which are radiopaque to facilitate visibility under fluoroscopy and guides stent placement. The pushing catheter operates by pushing the stent into place along a pre-positioned wire guide. To protect the flaps on the stents a flap protector is provided. The flap protector operates by providing a lumen through which the stent can be inserted and it collapses the flaps on the stent before insertion into the endoscope.

INTENDED USE

This device is used to drain obstructed biliary ducts.

INDICATIONS FOR USE

Endoscopic biliary stent placement for biliary drainage of obstructed ducts that could be caused by common bile duct stones, malignant biliary obstruction, benign or malignant strictures or other obstructed biliary conditions requiring drainage.

CLINICAL BENEFITS

Drainage of the biliary duct.

CONTRAINDICATIONS

Those specific to ERCP.

Inability to pass a wire guide, guiding catheter, and/or stent through the obstructed area.

WARNINGS

- This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or migration and/or mechanical integrity failure of device.
- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.
- This device is intended up to a maximum of 3 months indwell.

PRECAUTIONS

- A complete diagnostic evaluation must be performed prior to use to determine the proper stent size.
- Refer to the label for the minimum channel size required for this device.
- This device must be used under fluoroscopic monitoring.
- Do not use excessive force to advance the stent.
- The flap protector (supplied with the stent) is not intended for use in the accessory channel of the endoscope.
- Periodic evaluation of the device is recommended during indwell period.
- Sphincterotomy is not necessary for device placement.
- Dislodgement of a placed stent is possible when attempting additional procedures.
- Do not use this device for any purpose other than stated intended use.
- Store the device in a dry location.
- Use of this device restricted to a trained healthcare professional.
- The product is intended for use by physicians trained and experienced in ERCP techniques.
- Standard techniques for the endoscopic placement of stents should be employed.

MRI SAFETY INFORMATION

This symbol means the stent is MR Safe.

The ST-2 Soehendra Tannenbaum® Biliary Stent and Cotton-Leung® Biliary Stent are MR Safe.

In non-clinical testing, the image artifact caused by the device extends approximately 7 mm from the stent when imaged with a gradient echo pulse sequence and a 3 T MR system.

POTENTIAL ADVERSE EVENTS

Those associated with ERCP: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • liver abscess • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with biliary stent placement: fever • obstruction of the pancreatic duct • pain/discomfort • stent migration • stent occlusion • trauma to the biliary duct or duodenum.

HOW SUPPLIED

These devices are supplied ethylene oxide (EO) sterilized in a peel-open pouch.

This device is accompanied by an implant card that should be given to the patient after it has been completed by the healthcare professional.

DEVICE PREPARATION

1. Prior to introducing any devices ensure working channel of endoscope is lubricated with water or a water-based lubricant.
2. Referring to the device package label, ensure that the endoscope has a channel size greater than or equal to the minimum channel size required to operate the device.
3. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
4. Prior to using the device(s) visually inspect it/them for abnormalities that could result in improper function.
5. Attach endoscopic cap or wire guide locking device to endoscope.
6. Load flap protector onto duodenal end of the stent in order to collapse the flaps. The stent should fit into the flap protector with the entirety of the duodenal flaps collapsing and being covered.
7. Lubricate the wire guide locking device/endoscopic cap with a water based lubricant.

INSTRUCTIONS FOR USE

Illustrations

1. Advance the introducer with the preloaded stent over a pre-positioned wire guide and into the accessory channel of the endoscope until the proximal end of the wire guide exits PWP.
2. With the elevator open, advance the device in short increments until it is endoscopically visualized exiting the scope. **Note:** As the front flap of stent enters the accessory channel, keep the flap protector over the back flap until it is completely in the accessory channel. Slide the flap protector over the pushing catheter, keeping it clear of the accessory channel.
3. Advance the device into the appropriate duct.
4. Fluoroscopically monitor the radiopaque band(s) on the guiding catheter. When a sufficient length of the guiding catheter is above the obstruction, disconnect the guiding and pushing catheter hubs on the introducer. **Note:** A sufficient length of the guiding catheter must be above the obstruction prior to the stent positioning.
5. Maintain the position of the guiding catheter and advance the stent into the desired position using the pushing catheter.
6. Fluoroscopically and endoscopically confirm the desired stent placement.
7. Retract and remove both the guiding catheter and the wire guide while maintaining the position of the stent with the pushing catheter.
8. Confirm the position of the stent once the guiding catheter and the wire guide are removed.
9. Remove the pushing catheter.
10. These stents may be removed using standard endoscopic techniques.

DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be disposed in accordance with institutional guidelines. The introduction system should be coiled for disposal.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of.

When available, the EUDAMED website (<https://ec.europa.eu/tools/eudamed>), along with the BUDI for this product (0827002CIRL202007013010B6), can be used to locate the Summary of Safety and Clinical Performance (SSCP) for this product.

Patient-facing information may be accessed at cookmedical.eu/patient-implant-information.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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William A. Cook Australia Pty Ltd
COOK IRELAND LTD.