EN  
 3

**Fusion® Pushing Catheter**

Instructions for Use

 BG  
 6

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
## FUSION® PUSHING CATHETER

**CAUTION: US Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

This IFU covers Fusion Pushing Catheter devices as represented in **Figure 1**. Refer to the stent instructions for use for stent information. The Fusion Pushing Catheter is an introducer and consists of a pushing catheter and positioning sleeve. It is available in 5 Fr and 6 Fr in a length of 170 cm.

### Performance Characteristics

The function and key features of the device are described below.

- IDE port – The IDE port located 6 cm from the distal tip allows for use with shorter wire guide configurations.
- Pigtail straightener/flap protector/positioning sleeve - To facilitate the loading of the stent onto the wire guide.
- Stylet – Enhances the pushability of the device when used in a 'short-wire' configuration.
- Endoscopic visibility – The pushing catheter is visible endoscopically.
- Proximal wire port (hub) – Used to facilitate the attachment and detachment of the stylet and allows the wire guide to exit in the long wire configuration.

### Device Compatibility

The Fusion Pushing Catheter is compatible with the following:

- Endoscope with 3.2 mm accessory channel
- 0.035" wire guide
- Endoscopic cap/wire guide locking device
- Cook stents (**see Table 1**)

**Table 1. Compatible Cook Biliary Stents for Fusion® Pushing Catheter**

| Device Family                           | Cook Biliary Stents |              |               |                         |                        |         |
|-----------------------------------------|---------------------|--------------|---------------|-------------------------|------------------------|---------|
|                                         | Cotton-Huibregtse®  | Zimmon®      | Cotton-Leung® | Cotton-Leung® Sof-Flex® | Plastic Biliary Stent* | Compass |
|                                         | CHBSO               | ZSO          | CLSO          | CLSO-SF                 | PBS*                   | CBBSO   |
| <b>Fusion® Pushing Catheter (FS-PC)</b> | 7 Fr                | 5 Fr<br>7 Fr | 5 Fr<br>7 Fr  | 7 Fr                    | 7 Fr                   | 7 Fr    |

**\*Note:** PBS devices are not available in all markets, including the US.

Convention is as follows:

- Fr refers to the stent diameter in French size.

### Patient Population

Adult patients requiring placement of a biliary stent. The nature of the underlying pathology is prevalent in varying patients, therefore the devices are indicated for patients with biliary duct obstruction caused by common bile duct stones, malignant biliary obstruction, and benign or malignant strictures.

## Intended User

The product is intended for use by physicians trained and experienced in ERCP techniques. Standard techniques for the endoscopic placement of stents should be employed.

## Contact with Body Tissue

This device is tissue contacting in line with intended use.

## Operating Principles

The function of the pushing catheter is to advance the stent over a pre-positioned wire guide to its intended location and to maintain the position of the stent as it is being deployed.

The FS-PC may be used for either 'long-wire' or 'short-wire' ERCP techniques.

Long-wire: Upon removal of the stylet the FS-PC may be used in an over-the-wire configuration, whereby the full device is loaded onto the wire guide until the wire guide exits at the proximal wire port (hub).

Short-wire: The FS-PC contains an IDE port, located 6 cm from its distal tip. By allowing the wire guide to exit its intended lumen through this port, only the distal section of the device is loaded onto the wire guide during use, thereby facilitating the 'short-wire' ERCP technique. The stylet enhances pushability when placing the stent.

## INTENDED USE

This device is used for endoscopic biliary stent placement.

## INDICATIONS FOR USE

For endoscopic biliary stent placement for biliary drainage of obstructed ducts that could be caused by common bile duct stones, malignant biliary obstruction, benign or malignant strictures or other obstructed biliary conditions requiring drainage.

## CLINICAL BENEFITS

The biliary stent introducers have an indirect therapeutic impact as they function to facilitate the successful placement of a biliary stent. The introducers are essential for stent placement.

## CONTRAINDICATIONS

Contraindications include those specific to ERCP.

Inability to pass wire guide/stent through the obstructed area.

## WARNINGS

- Introducer contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.
- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

## PRECAUTIONS

-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a cobalt alloy (or stainless-steel alloy containing cobalt), which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.
- Only Cook biliary stents are recommended for use with this device (see Table 1 for guidance).
- The positioning sleeve\* is not intended for use in the accessory channel of the endoscope.
- Care must be exercised when straightening the pigtail curls\* in order to avoid kinking or breaking the stent.

- Do not use excessive force to advance the device.
- Dislodgement of a placed stent is possible when attempting additional procedures.
- Do not use this device for any purpose other than the stated intended use.
- Store the device in a dry location.
- Use of this device is restricted to a trained healthcare professional.
- The product is intended for use by physicians trained and experienced in ERCP techniques.
- Standard techniques for the endoscopic placement of stents should be employed.

(\*) indicates "if any"

## POTENTIAL ADVERSE EVENTS

Those associated with ERCP: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • liver abscess • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with biliary stent placement: allergic reaction to nickel • fever • obstruction of the pancreatic duct • pain/discomfort • stent migration • trauma to the biliary duct or duodenum.

## HOW SUPPLIED

These devices are supplied Ethylene Oxide (EO) sterilized in a peel open pouch.

The FS-PC is further packaged into a holder to maintain integrity of the device.

## DEVICE PREPARATION

1. Prior to introducing any devices ensure working channel of endoscope is lubricated with water or a water based lubricant.
2. Referring to the device package label, ensure that the endoscope has a channel size greater than or equal to the minimum channel size required to operate the device.
3. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
4. Prior to using the device(s) visually inspect it/them for abnormalities that could result in improper function.

## INSTRUCTIONS FOR USE (SEE FIGURE 1)

## Illustrations

### I. IF USING INTRA DUCTAL EXCHANGE (IDE) PORT AND A PRE-POSITIONED SHORT WIRE:

1. Unlock the short wire from the Wire Guide Locking Device and load the stent and positioning sleeve over a pre-positioned wire guide. Note: Use the positioning sleeve included with this device.
2. Advance the pushing catheter over a pre-positioned wire guide, ensuring that the wire guide exits the pushing catheter at the IDE port.
3. As the front flap of stent enters accessory channel, keep the positioning sleeve over the back flap until it is completely in the accessory channel. Slide the positioning sleeve over the pushing catheter and wire guide, keeping it clear of the accessory channel.
4. Relock the wire guide.
5. With the elevator open, advance the device in small increments until it is visualized exiting the endoscope.
6. Advance the device until the stent is in the desired position.
7. Fluoroscopically and endoscopically confirm the desired stent placement.
8. Unlock the wire guide and gently remove the wire guide from the endoscope while maintaining the position of the stent with the pushing catheter.
9. Re-lock the wire guide and gently remove the pushing catheter from the accessory channel.
10. When resistance is felt at the Wire Guide Locking Device when removing the pushing catheter, unlock the wire lock and perform a standard exchange over the remaining section of the pushing catheter. Once the distal tip exits from the Wire Guide Locking Device, relock the wire guide.

### II. IF USING PROXIMAL WIRE PORT (PWP) AND A PRE-POSITIONED LONG WIRE GUIDE:

1. Remove the stylet.

2. Load the stent and the positioning sleeve onto a pre-positioned wire guide.
3. Advance the pushing catheter over the pre-positioned wire guide, ensuring that the wire guide exits the pushing catheter at the PWP.
4. As the front flap of the stent enters the accessory channel, keep the positioning sleeve over the back flap until it is completely in the accessory channel. Slide the positioning sleeve over the pushing catheter, keeping it clear of the accessory channel.

**REFER TO STEPS 5-9 IN "SECTION I", THEN RESUME WITH STEP 5 BELOW:**

5. Remove the pushing catheter.

## **DISPOSAL OF DEVICES**

This device may be contaminated with potentially infectious substances of human origin and should be coiled for disposal in accordance with institutional guidelines.

## **PATIENT COUNSELLING INFORMATION**

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

## **SERIOUS INCIDENT REPORTING**

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2024
William A. Cook Australia Pty Ltd
**COOK IRELAND LTD.**