EN
 3

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INTENDED USE

This device is used to drain obstructed pancreatic ducts.

NOTES

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to: inability to dilate stricture, inability to pass wire guide or stent through obstructed area.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Other potential complications associated with pancreatic stent placement include, but are not limited to: trauma to pancreatic duct or duodenum, obstruction of common bile duct, stent migration.

WARNINGS

MRI compatibility of pancreatic stents has not been established.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

A complete diagnostic evaluation must be performed prior to use to determine proper stent size.

This stent must be placed under fluoroscopic monitoring.

This stent must only be placed using introduction system provided with set.

This stent should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

INSTRUCTIONS FOR USE**Illustrations**

1. Flush the injection port on the wedge stent introducer with sterile water or saline.

2. Advance preloaded introduction system over pre-positioned wire guide until it reaches accessory channel. **Note:** If stent is not preloaded, back load stent onto guiding catheter.
3. Advance preloaded introduction system into accessory channel until stent is endoscopically visualized exiting scope. **Note:** Endoscope elevator must be open to allow stent to exit.
4. Under fluoroscopic monitoring, slowly advance stent into desired position. Stent must be positioned with wedged end and length of stent in duct and proximal end exiting into duodenum. **Note:** Stents bridging papilla should extend beyond papilla and into duodenum approximately 0.5 cm after deployment. **Note:** Stent is constructed of radiopaque material, allowing fluoroscopic visualization of stent position.
5. Fluoroscopically and endoscopically confirm desired stent position.
6. Gently remove wire guide from introduction system while maintaining position of stent with introduction system.
7. Re-confirm position of stent once wire guide is removed.
8. Loosen Tuohy-Borst adapter on end of introduction system and gently withdraw guiding catheter portion of introduction system. (*See fig. 1*)
9. Lumen may be flushed, at flushport, to observe stent drainage.
10. Remove guiding catheter completely from endoscope while maintaining position of stent with pushing catheter.
11. Re-confirm position of stent fluoroscopically and endoscopically once guiding catheter and wire guide are removed.
12. Gently remove pushing catheter from accessory channel.
13. Re-confirm proper positioning of stent fluoroscopically and endoscopically before removing endoscope.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2015 ©
COOK IRELAND LTD.