EN
6
Johlin™ Pancreatic Wedge Stent and Introducer

Instructions for Use

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CAUTION: U.S Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Please read all instructions before using this device.

INTENDED USE

This device is used to drain obstructed pancreatic ducts.

DEVICE DESCRIPTION

This IFU covers the Johlin Pancreatic Wedge Stent and Introducer as represented by Figures 1 and 2.

CONTRAINDICATIONS

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to: inability to dilate stricture • inability to pass wire guide or stent through obstructed area.

POTENTIAL ADVERSE EVENTS

Those associated with ERCP include, but are not limited to: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with pancreatic stent placement include, but are not limited to: fever • obstruction of the biliary/pancreatic duct • obstruction of common bile duct • pain/discomfort • stent migration • stent occlusion • trauma to the biliary tract or duodenum • trauma to pancreatic duct or duodenum.

WARNINGS

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

MR SAFETY INFORMATION

This symbol means the stent is MR Safe.

The Johlin™ Pancreatic Wedge Stent is MR Safe.

In non-clinical testing, the image artifact caused by the device extends approximately 7 mm from the stent when imaged with a gradient echo pulse sequence and a 3 T MR system.

PRECAUTIONS

Refer to the label for the minimum channel size required for this device.

Refer to the label for the appropriate wire guide size required for this device.

A complete diagnostic evaluation must be performed prior to use to determine proper stent size.

This stent must be placed under fluoroscopic monitoring.

This stent must only be placed using introduction system provided with set.

This stent should not be left indwelling for more than three months. Periodic evaluation is recommended.

Sphincterotomy is not necessary for device placement.

Dislodgement of a placed stent is possible when attempting additional procedures.

Do not use excessive force to advance device.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

SYSTEM PREPARATION

1. If applicable, load the stent onto the tip of the introduction system and position the pushing catheter such that it is located adjacent to the stent.
2. Tighten the Tuohy-Borst adapter to secure the position of the pushing catheter.
3. Flush the injection port on the wedge stent introducer with sterile water or saline.

INSTRUCTIONS FOR USE

Illustrations

1. Advance the device over the pre-positioned wire guide until it reaches the accessory channel.
2. Advance the device into the accessory channel in small increments until the stent is endoscopically visualized exiting the scope. **Note:** Endoscope elevator must be open to allow the stent to exit.
3. Under fluoroscopic monitoring, slowly advance the stent into desired position. **Note:** Stents bridging papilla should extend beyond papilla and into duodenum approximately 0.5 cm after deployment.
4. Gently remove the wire guide from the introduction system while maintaining the position of the stent with the introduction system.
5. Lumen may be flushed, at the wire guide port to observe stent drainage.
6. Loosen the Tuohy-Borst adapter on the end of introduction system and gently withdraw the guiding catheter portion of the introduction system while maintaining the position of the stent with the pushing catheter.

7. Fluoroscopically and endoscopically confirm the desired stent position once the guiding catheter and wire guide are removed.
8. Gently remove the pushing catheter from the accessory channel.
9. Fluoroscopically and endoscopically confirm the desired stent position before removing the endoscope.
10. These stents may be removed using standard endoscopic techniques.

Upon completion of procedure, remove all components of the device from the working channel of the endoscope and dispose of device per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

These devices are supplied Ethylene oxide (EtO) sterilized in a peel open pouch.

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© COPYRIGHT COOK 2019 ©
COOK IRELAND LTD.