
EN 6	Johlin™ Pancreatic Wedge Stent and Introducer Instructions for Use
BG 9	Панкреасен клиновиден стент и въвеждащо устройство Johlin™ Инструкции за употреба
CS 13	Pankreatický klínový stent a zaváděč Johlin™ Návod k použití
DA 17	Johlin™ pancreaswedgestent og indfører Brugsanvisning
DE 21	Keilförmiger Johlin™ Pankreasstent mit Einführsystem Gebrauchsanweisung
EL 25	Παγκρεατική σφηνοειδής ενδοπρόσθεση και εισαγωγέας Johlin™ Οδηγίες χρήσης
ES 30	Stent pancreático angulado e introductor Johlin™ Instrucciones de uso
ET 34	Johlin™ pankrease kiilstent ja sisestusvahend Kasutusjuhend
FR 38	Endoprothèse pancréatique cunéiforme et introducteur Johlin™ Mode d'emploi
HR 42	Johlin™ klinasti stent za gušteraču i uvodnica Upute za uporabu
HU 46	Johlin™ hasnyálmirigyékszent és bevezetőeszköz Használati utasítás
IT 50	Stent pancreatico a cuneo Johlin™ con introduttore Istruzioni per l'uso
LT 54	„Johlin™“ kasos pleištinis stentas ir intubatorius Naudojimo instrukcija
LV 58	Johlin™ aizkuņģa dziedzera ķīļveida stents un ievadītājs Lietošanas instrukcija
NL 62	Johlin™-pancreaswigstent en introducer Gebruiksaanwijzing
NO 66	Johlin™ kilestent og innførsingsenhet for bukspyttkjertel Bruksanvisning
PL 70	Stent klinowy i koszulka wprowadzająca Johlin™ do przewodu trzustkowego Instrukcja użycia
PT 74	Stent de encravamento pancreático e introdutor Johlin™ Instruções de utilização
RO 78	Stent pancreatic cuneiform Johlin™ și introductor Instrucțiuni de utilizare
SK 82	Pankreatický klínový stent a zaváděč Johlin™ Návod na použití
SL 86	Zagozdni stent za trebušno slinavko in uvajalo Johlin™ Navodila za uporabo
SR 90	„Johlin™“ pankreasni klinasti stent i uvodnik Uputstvo za upotrebu
SV 93	Johlin™ pankreatisk kilstent och införare Bruksanvisning
TR 97	Johlin™ Pankreatik Kama Stent ve İntrodüser Kullanma Talimatı

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JOHLIN™ PANCREATIC WEDGE STENT AND INTRODUCER

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Please read all instructions before using this device.

DEVICE DESCRIPTION

This IFU covers the Johlin Pancreatic Wedge Stent and Introducer as represented by Figures 1 and 2.

This device is sold as a set consisting of a stent and dedicated introducer (a guiding catheter and a pushing catheter). The stents are available in 8.5 and 10 Fr diameter and lengths of 8–22 cm.

The JPWS 8.5 Fr stent is straight and is offered in a preloaded configuration and the JPWS 10 Fr stents have a bend at the duodenal end and are supplied in a set where the stent is not preloaded onto the introducer. A Tuohy-Borst adapter allows for the guiding catheter and pushing catheter to be secured in position relative to each other. It also provides a connection for injection via the injection port.

Performance Characteristics

The function and key features of these components are described below.

The stents that provide drainage are available in 8.5 and 10 Fr diameter and lengths of 8–22 cm to accommodate various patient anatomies, the size and location of the obstruction. The shape of the 8.5 Fr stent is straight and the 10 Fr stent has a bend (anti-migration feature) at the duodenal end.

They have the following features:

- Atraumatic tip – The ends of the stent are designed to be atraumatic to the anatomy.
- Sideports – Aid drainage
- Anti-migration features – Aim to help maintain the position of the stent once it is placed by providing mechanical interference with the anatomy.
- Radiopacity – The stent materials are radiopaque to facilitate visibility under fluoroscopy.

Guiding catheter – The function of the guiding catheter is to guide the stent to its intended location.

Guiding catheters contain one or more of the following features:

- Atraumatic tip – The end of the guiding catheter is designed to be atraumatic to the anatomy.
- Radiopacity – The guiding catheter tubing material is radiopaque. The guiding catheter also contains marker bands to further enhance its fluoroscopic visibility.
- Tuohy-Borst – The guiding catheter contains a Tuohy-Borst with a connection to facilitate the attachment of a standard Luer syringe.

Pushing catheter – The function of the pushing catheter is to advance the stent over a pre-positioned wire guide and a guiding catheter to its intended location, and to maintain the position of the stent as it is being deployed. It contains the following feature:

- Proximal identifier – An identifier at its proximal end. It is used to identify the orientation in which the device must be used.

Stent Size	Pushing Catheter Diameter	Guiding Catheter Diameter
8.5 Fr	8.5 Fr	5 Fr
10 Fr	10 Fr	6 Fr

See product label for all dimensions.

Device Compatibility

The Johlin pancreatic stents are compatible with the following:

- Endoscope with 3.2 mm or 3.7 mm accessory channel (refer to label for selection of the appropriate accessory channel)
- 0.035" wire guide (refer to label for selection of the appropriate wire guide size)
- Endoscopic cap or wire guide locking device
- Contrast media
- Sterile water or saline
- Standard Luer syringe
- Stent retriever or forceps

Qualitative and Quantitative Information

The materials for the stent implant are outlined in the table below:

Product	Implant Material Qualitative Information		Weight (g)
Johlin Pancreatic Wedge Stent	Polyurethane	Wedge stent	Up to 1.28

Patient Population

These devices can be used in all patient populations (≥ 10 kg) requiring pancreatic stenting for obstruction. The nature of the underlying pathology is prevalent in varying patients, therefore, the devices are indicated for patients with pancreatic duct obstruction caused by pancreatitis, stricture, stones, cancer, anatomic anomalies, disrupted duct, and pancreatic fluid collection.

Intended User

The product is intended for use by physicians trained and experienced in endoscopic retrograde cholangiopancreatography (ERCP) techniques.

Standard techniques for the endoscopic placement of stents should be employed.

Contact with Body Tissue

The device is tissue contacting in line with intended use.

Operating Principles

The Johlin stent operates by providing a lumen through which pancreatic fluid can drain, thereby relieving obstructed pancreatic ducts. The 10 Fr stent has an anti-migration feature which provides mechanical interference with the anatomy to help maintain the position of the stent. The Johlin stent is placed using an introducer which is supplied as part of the set. The introducer consists of a pushing catheter, a guiding catheter and a Tuohy-Borst. The guiding catheter contains marker bands which are radiopaque to facilitate visibility under fluoroscopy and guide stent placement. The pushing catheter operates by pushing the stent into place. Tuohy-Borst allows for the guiding catheter and pushing catheter to be secured in position relative to each other. An injection port allows the flushing of sterile or saline water through the device.

INTENDED USE

This device is used to drain pancreatic ducts.

INDICATIONS FOR USE

Endoscopic pancreatic stent placement for pancreatic drainage of obstructed ducts that could be caused by pancreatitis, stricture, pancreatic cancer, anatomic anomalies of the pancreas, pancreatic fluid collection, pancreatic stones, and disrupted duct.

CLINICAL BENEFITS

Drainage of the pancreatic duct.

CONTRAINDICATIONS

Those specific to ERCP and any procedure to be performed in conjunction with stent placement. Inability to dilate stricture, inability to pass wire guide or stent through obstructed area.

WARNINGS

This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or migration and/or mechanical integrity failure of device.

The Johlin Pancreatic Wedge Stents are intended for a maximum of 3 months indwell.

Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.

Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

The safety and effectiveness of this device for prophylactic use has not been established.

PRECAUTIONS

- Refer to the label for the minimum channel size required for this device.
- Refer to the label for the appropriate wire guide size required for this device.

- A complete diagnostic evaluation must be performed prior to use to determine proper stent size.
- Fluoroscopic monitoring must be utilized during stent placement.
- Refer to Fig. 2 for the intended placed orientation of the stent.
- This stent must only be placed using introducer provided with set.
- Periodic evaluation of the device is recommended during indwell period.
- Sphincterotomy is not necessary for device placement.
- Dislodgement of a placed stent is possible when attempting additional procedures.
- Do not use excessive force to advance device.
- Do not use this device for any purpose other than the stated intended use.
- Store the device in a dry location.
- Use of this device restricted to a trained healthcare professional.
- The product is intended for use by physicians trained and experienced in ERCP techniques.
- Standard techniques for the endoscopic placement of stents should be employed.

MRI SAFETY INFORMATION

This symbol means the stent is MR Safe.

The Johlin™ Pancreatic Wedge Stent is MR Safe.

In non-clinical testing, the image artifact caused by the device extends approximately 7 mm from the stent when imaged with a gradient echo pulse sequence and a 3 T MR system.

POTENTIAL ADVERSE EVENTS

Those associated with ERCP: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with pancreatic stent placement: fever • infection • obstruction of the common bile duct • pain/discomfort • pancreatitis • perforation • stent migration • stent occlusion • trauma to the pancreatic duct or duodenum.

HOW SUPPLIED

These devices are supplied ethylene oxide (EtO) sterilized in a peel-open pouch. Following sterilization, the pouch is further packaged into an outer protective box, which protects the UV-sensitive device from light.

This device is accompanied by an implant card that should be given to the patient after it has been completed by the healthcare professional.

DEVICE PREPARATION

1. Prior to introducing any devices ensure working channel of endoscope is lubricated with water or a water-based lubricant.
2. Referring to the device package label, ensure that the endoscope has a channel size greater than or equal to the minimum channel size required to operate the device.
3. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
4. Prior to using the device(s) visually inspect it/them for abnormalities that could result in improper function.
5. Attach an endoscopic cap to the endoscope working channel, apply a small amount of water-based lubricant to the opening leading to the accessory channel on the endoscopic cap.
6. If applicable, load the stent onto the tip of the introduction system and position the pushing catheter such that it is located adjacent to the stent.
7. Tighten the Tuohy-Borst adapter to secure the position of the pushing catheter.
8. Flush the injection port on the wedge stent introducer with sterile water or saline.

INSTRUCTIONS FOR USE

Illustrations

1. Advance the device over the pre-positioned wire guide until it reaches the accessory channel.

2. Advance the device into the accessory channel in small increments until the stent is endoscopically visualized exiting the scope. **Note:** Endoscope elevator must be open to allow the stent to exit.
3. Under fluoroscopic monitoring, slowly advance the stent into desired position. **Note:** Stents bridging papilla should extend beyond papilla and into duodenum approximately 0.5 cm after deployment.
4. Gently remove the wire guide from the introduction system while maintaining the position of the stent with the introduction system.
5. Lumen may be flushed, at the wire guide port to observe stent drainage.
6. Loosen the Tuohy-Borst adapter on the end of introduction system and gently withdraw the guiding catheter portion of the introduction system while maintaining the position of the stent with the pushing catheter.
7. Fluoroscopically and endoscopically confirm the desired stent position once the guiding catheter and wire guide are removed.
8. Gently remove the pushing catheter from the accessory channel.
9. Fluoroscopically and endoscopically confirm the desired stent position before removing the endoscope.
10. These stents may be removed using standard endoscopic techniques.

DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be disposed in accordance with institutional guidelines. The introducer should be coiled for disposal.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

When available, the EUDAMED website (<https://ec.europa.eu/tools/eudamed>), along with the BUDI for this product (0827002CIRL202007020003AX), can be used to locate the Summary of Safety and Clinical Performance (SSCP) for this product.

Patient-facing information may be accessed at cookmedical.eu/patient-implant-information.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2026
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.