EN
3**Nasal Biliary Drainage Set**

Instructions for Use

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

Bruksanvisning

INTENDED USE

This device is used for temporary endoscopic drainage of the biliary duct through the nasal passage by use of an indwelling catheter.

NOTES

Do not use this device for any purpose other than stated intended use.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

Each set includes a drainage catheter, nasal transfer tube and a drainage connecting tube to facilitate irrigation and drainage.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

Contraindications of nasal biliary drainage include, but are not limited to: coagulopathy, sepsis, active pancreatitis or other infection, recent food ingestion, inability to pass a wire guide or drainage catheter through the obstructed area.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Complications associated with nasal biliary drainage include, but are not limited to: trauma to the common bile duct or duodenum, obstruction of the pancreatic duct, blockage of the drainage catheter, nasal irritation and/or sore throat.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Pigtail straightener, if any, is not intended for use in accessory channel of endoscope.

SYSTEM PREPARATION

If tip of drainage catheter is configured with a pigtail and is packaged with a pigtail straightener, advance pigtail straightener along shaft and to tip of catheter to straighten pigtail curl. Drainage catheter is now ready for placement. **Note:** Care must be exercised when straightening pigtail curl in order to avoid kinking or breaking catheter.

1. Flush drainage catheter using a syringe filled with sterile water or saline attached to Luer lock fitting on end of catheter.
2. Carefully remove Tuohy-Borst connector from end of drainage catheter.
(See fig. 1)
3. Introduce and advance drainage catheter and pigtail straightener, if any, in small increments into accessory channel of endoscope over pre-positioned wire guide until catheter is endoscopically visualized exiting endoscope. **Note:** Elevator of endoscope must be open to allow drainage catheter to exit.
4. Under fluoroscopic monitoring, slowly advance drainage catheter to desired position in biliary duct above obstruction.
5. When drainage catheter is in desired position, slowly remove endoscope while holding drainage catheter and wire guide in place at patient's mouth.
6. Fluoroscopically confirm that drainage catheter and wire guide remain positioned in biliary duct, then slowly remove wire guide in small increments while maintaining position of drainage catheter.
7. Lubricate and advance nasal transfer tube through nostril until it exits in back of throat.
8. Grasp connection end of nasal transfer tube in throat and pull it out through mouth. (See fig. 2)
9. If placing an 8.5 French or smaller tube, thread the tip of the tube through the oral end of the nasal transfer tube and advance until it exits through the nasal end. If placing a 10 French or larger tube, securely connect the tip of the tube to the oral end of the nasal transfer tube. (See fig. 3)
10. Hold drainage catheter, to prevent kinking and removal, and slowly pull nasal transfer tube until drainage catheter emerges through nostril. (See fig. 4)
11. Replace Tuohy-Borst connector onto end of drainage catheter. Contrast may be injected to check placement of catheter within duct.
12. Secure drainage catheter to patient, then connect end of catheter to Luer lock fitting on drainage connecting tube. Attach universal connector to appropriate drainage bag to begin drainage, per institutional guidelines.

Upon completion of procedure, dispose of components per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2013 ©
COOK IRELAND LTD.