
EN 8	Solus® Double Pigtail Stent with Introducer Instructions for Use
BG 12	Двустранно завит като свинска опашка стент с въвеждащо устройство Solus® Инструкции за употреба
CS 16	Double pigtailový stent Solus® se zavaděčem Návod k použití
DA 20	Solus® stent med dobbelt grisehale og indfører Brugsanvisning
DE 24	Solus® Doppelpigtail-Stent mit Einführschleuse Gebrauchsanweisung
EL 29	Ενδοπρόσθεση διπλού σπειροειδούς άκρου Solus® με εισαγωγή Οδηγίες χρήσης
ES 34	Stent de doble pigtail con introductor Solus® Instrucciones de uso
ET 38	Kahe rõngasse keeratud otsaga stent koos sisestajaga Solus® Kasutusjuhend
FR 42	Endoprothèse double pigtail Solus® avec introducteur Mode d'emploi
HR 47	Stent s dvostrukim „Pigtailom“ s uvodnicom Solus® Upute za uporabu
HU 51	Solus® kettős pigtail sztent bevezetőeszközzel Használati utasítás
IT 55	Stent a doppio pigtail con introduttore Solus® Istruzioni per l'uso
LT 60	„Solus®“ stentas dviem riestais galiukais su intubatoriumi Naudojimo instrukcija
LV 64	Solus® stents ar abusēji saritinātiem galiem ar ievadītāju Lietošanas instrukcija
NL 68	Solus®-stent met dubbele pigtail met introducer Gebruiksaanwijzing
NO 72	Solus® stent med dobbel pigtail og innføringsenhet Bruksanvisning
PL 76	Stent z dwiema końcówkami pigtail Solus® z introduktorem Instrukcja użycia
PT 81	Stent de duplo pigtail Solus® com introdutor Instruções de utilização
RO 85	Stent dublu-pigtail Solus® cu introductor Instrucţiuni de utilizare
SK 90	Dvojité špirálový stent Solus® so zavádzačom Návod na použitie
SL 94	Stent z dvojnim zavitim koncem in z uvajalom Solus® Navodila za uporabo

SR 98	„Solus®” stent sa dvostrukim „pigtail” sa uvodnikom Uputstvo za upotrebu
SV 102	Solus® stent med dubbla pigtails och introducer Bruksanvisning
TR 106	Solus® İntrodüserli Çift Pigtail’li Stent Kullanma Talimatı

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SR Slika 2: „Pigtail” ispravljač

SV Figur 2: Pigtailuträtare

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EN The indicated stent length is the distance between the pigtails

- BG Посочваната дължина на стента е разстоянието между извитите краища с форма на свинска опашка (pigtail)
- CS Uvedená délka stentu představuje vzdálenost mezi pigtaily
- DA Den angivne stentlængde er afstanden mellem grisehalerne
- DE Die angegebene Stentlänge bezieht sich auf den Abstand zwischen den Pigtails.
- EL Το μήκος ενδοπρόθεσης που αναφέρεται είναι η απόσταση μεταξύ των σπειροειδών άκρων
- ES La longitud indicada del stent es la distancia entre los pigtails
- ET Stendi näidatud pikkus on rõngassabade vahekaugus
- FR La longueur indiquée de l'endoprothèse est la distance entres les pigtails
- HR Indicirana duljina stenta predstavlja udaljenost između zavojnica
- HU A megadott sztenthossz a pigtailek közötti távolsággal egyenlő
- IT La lunghezza dello stent indicata corrisponde alla distanza tra pigtail
- LT Nurodytas stento ilgis yra atstumas tarp riestų galiukų
- LV Indicētais stenta garums ir attālums starp saritinātajiem galiem
- NL De aangegeven stentlengte is de afstand tussen de pigtails
- NO Den angitte stentlengden er avstanden mellom pigtailene
- PL Wskazana długość stentu to odległość między końcówkami pigtail
- PT O comprimento de stent indicado é a distância entre os pigtails
- RO Lungimea indicată a stentului reprezintă distanța dintre regiunile pigtail
- SK Uvádzaná dĺžka stentu je vzdialenosť medzi špirálami
- SL Navedena dolžina stenta je razdalja med zavitima koncema
- SR Namenjena dužina stenta je udaljenost između „pigtails“
- SV Den indikerade stentlängden är avståndet mellan pigtails
- TR Belirtilen stent uzunluğu, pigtail'ler arasındaki mesafedir

SOLUS® DOUBLE PIGTAIL STENT WITH INTRODUCER

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

This IFU covers Solus Double Pigtail Stents with Introduction System (ZSS) as represented in **Figure 1**. Only stents supplied with this introduction system are compatible for use with the Solus Stent Introduction System.

Solus Double Pigtail Stent with Introducer is a Biliary Stent Set which includes a combination of a double-pigtail Biliary stent (with a pigtail straightener except for the 1 cm stent which does not include a pigtail straightener) and an introduction system.

The stents are available in 10 Fr diameter and lengths of 1–15 cm.

The introduction system consists of a guiding catheter component of length 210 cm and a pushing catheter component of length 165 cm. The guiding catheter has 4 radiopaque marker bands.

Performance Characteristics

The function and key features of the device are described below.

A variety of Solus biliary stents that provide drainage are available across the device range to accommodate various patient anatomies. They are offered in 10 Fr and in labelled lengths of various sizes between 1 cm and 15 cm. Biliary stents have the following features:

- Atraumatic Ends – The ends of the stent are designed to be atraumatic to the anatomy.
- Side ports – Aid drainage.
- Anti-migration features – the double pigtails aim to help maintain the position of the stent once it is placed by providing mechanical interference with the anatomy.
- Radiopacity – The stent materials are radiopaque to facilitate visibility under fluoroscopy. Marker bands further support the fluoroscopic visibility of the stent.

Pigtail Straightener – This is a tube used to straighten the pigtails.

Guiding Catheter – The function of the guiding catheter is to guide the stent to its intended location. The guiding catheter contains the following features:

- Atraumatic tip – The end of the guiding catheter is designed to be atraumatic to the anatomy.
- Radiopacity – The guiding catheter tubing material is radiopaque. The guiding catheter also contains marker bands to further enhance its fluoroscopic visibility.
- Proximal Wire Port (Hub) – allows the wire guide to exit in the long wire configuration.

Pushing Catheter – The function of the Pushing Catheter is to advance the Stent over a pre-positioned wire guide and a Guiding Catheter, to its intended location, and to maintain the position of the Stent as it is being deployed. It contains the following feature:

- Proximal Identifier – An identifier at its proximal end. It is used to identify the orientation in which the device must be used.

Device Compatibility

These devices are compatible with the following:

- Endoscope with 3.7 mm accessory channel
- 0.035" Wire Guide
- Endoscopic Cap
- Stent retriever/forceps
- Contrast media
- Sterile water/Saline
- Standard Luer syringe
- Water Soluble Lubricant

Qualitative and Quantitative Information

The materials for the stent implant are outlined in **Table 1**.

Table 1: Stent Implant Materials

Product	Qualitative Information Device Material		Quantitative Information Weight (g)
10 Fr Solus Double Pigtail Biliary Stent (ZSS)	Tantalum	Stent Marker Band	Up to 0.10
	Polyurethane	Stent	Up to 2.20

Patient Population

Adult patients requiring biliary stenting for obstruction. The nature of the underlying pathology is prevalent in varying patients, therefore the devices are indicated for patients with biliary duct obstruction caused by common bile duct stones, malignant biliary obstruction and benign or malignant strictures.

Intended User

This product is intended for use by physicians trained and experienced in endoscopic retrograde cholangiopancreatography (ERCP) techniques.

Standard techniques for the endoscopic placement of stents should be employed.

Contact with Body Tissue

This device is tissue contacting in line with intended use.

Operating Principles

The Solus biliary stent operates by providing a lumen through which biliary fluid can drain. The stent pigtails act as an anti-migration feature which provides mechanical interference with the anatomy to help maintain the position of the stent. The stent is radiopaque and contains marker bands to facilitate visibility under fluoroscopy. Side-ports on the stent aid drainage by providing additional openings for fluid flow. The biliary stent is placed over a prepositioned wire guide using an introducer. The introducer consists of a pushing catheter and a guiding catheter combined. The guiding catheter contains marker bands which are radiopaque to facilitate visibility under fluoroscopy and guides stent placement. The pushing catheter operates by pushing the stent into place along a pre-positioned wire guide. The pigtail straightener operates by providing a lumen through which the stent can be inserted thus straightening the pigtail.

INTENDED USE

This device is used to drain obstructed biliary ducts.

INDICATIONS FOR USE

Endoscopic biliary stent placement for biliary drainage of obstructed ducts that could be caused by common bile duct stones, malignant biliary obstruction, benign or malignant strictures or other obstructed biliary conditions requiring drainage.

CLINICAL BENEFITS

Drainage of the biliary duct.

CONTRAINDICATIONS

Those specific to ERCP.

Inability to pass a wire guide, guiding catheter and/or stent through obstructed area.

WARNINGS

This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or migration and/or mechanical integrity failure of device.

Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.

Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working conditions, do not use. Please notify Cook for return authorization. The Solus Double Pigtail Stent is intended up to a maximum of 3 months indwell.

PRECAUTIONS

A complete diagnostic evaluation must be performed prior to use to determine the proper stent size. This stent must be placed under fluoroscopic monitoring.

This stent must only be placed using the Solus Double Pigtail Stent Introduction System, which is included in this device.

Periodic evaluation of the device is recommended during indwell period.

The pigtail straightener is used to straighten the pigtail curls on the stent. A pigtail straightener is not intended for use in the accessory channel.

An assessment must be made to determine the necessity of sphincterotomy prior to stent placement. If sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed. Dislodgement of a placed stent is possible when attempting additional procedures.

Do not use excessive force to advance the device.

Do not use this device for any purpose other than stated intended use.

Store the device in a dry location.

Use of this device restricted to a trained healthcare professional.

The product is intended for use by physicians trained and experienced in ERCP techniques. Standard techniques for the endoscopic placement of stents should be employed.

MRI SAFETY INFORMATION

This symbol means the stent is MR Conditional.

Non-clinical testing has demonstrated that the Solus Double Pigtail Stents are MR Conditional. A patient with the Solus® Double Pigtail Stent which contains radiopaque marker bands may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.

Device name	Solus Double Pigtail Stent with introducer
Static magnetic field	1.5 T or 3 T
Maximum spatial field gradient	3000 gauss/cm (30 T/m)
RF excitation	Circularly polarized (CP)
RF transmit coil type	Whole body transmit coil, head RF transmit-receive coil
Maximum whole body specific absorption rate (SAR)	4.0 W/kg (first level operating mode)
Limits on scan duration	4.0 W/kg whole body average SAR for 60 minutes of continuous RF (a sequence or back-to-back series/scan without breaks)
MR image artifact	The presence of this implant may produce an image artifact of 7 mm.

Note: If information about a specific parameter is not included, there are no conditions associated with that parameter.

For US Patients Only:

It is recommended that patients register the conditions under which the implant can safely be scanned with the MedicAlert Foundation (www.medicalert.org) or an equivalent organization.

POTENTIAL ADVERSE EVENTS

Those associated with ERCP: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • liver abscess • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with biliary stent placement: fever • obstruction of the pancreatic duct • pain/discomfort • stent migration • stent occlusion • trauma to the biliary duct or duodenum

HOW SUPPLIED

These devices are supplied ethylene oxide (EO) sterilized in a peel-open pouch.

Following sterilization, the pouch is further packaged into a box as these are UV sensitive.

This device is accompanied by an implant card that should be given to the patient after it has been completed by the healthcare professional.

DEVICE PREPARATION

1. Prior to introducing any devices, ensure working channel of endoscope is lubricated with water or a water-based lubricant.
2. Referring to the device package label, ensure that the endoscope has a channel size greater than or equal to the minimum channel size required to operate the device.
3. Prior to using the device(s) visually inspect it/them for abnormalities that could result in improper function.
4. Attach endoscopic cap to endoscope.
5. Lubricate the visible portion of the guiding catheter with a water-soluble lubricant.

INSTRUCTIONS FOR USE

Illustrations

1. Advance the pigtail straightener, if applicable, along the shaft of the stent to the end.
2. Using the pigtail straightener as an aid, load the stent onto the introduction system until it reaches the pushing catheter.
3. Advance the device over a prepositioned wire guide into the accessory channel. Slide the pigtail straightener over the pushing catheter, keeping it clear of the accessory channel.
4. Advance the device into the accessory channel in small increments until endoscopically visualized exiting the scope. **Note:** The endoscope elevator must be open to allow the stent to exit.
5. Advance the device into the duct. Fluoroscopically monitor the radiopaque band(s) on the guiding catheter. Advance the device into the duct until a sufficient length of the guiding catheter is above the obstruction.
6. Disconnect the guiding and pushing catheter hubs on to the introducer. Maintain the position of the guiding catheter and advance the stent into the desired position using the pushing catheter.
7. Fluoroscopically and endoscopically confirm the desired stent placement.
8. Retract and remove both the guiding catheter and the wire guide while maintaining the position of the stent with the pushing catheter.
9. Remove the pushing catheter.
10. These stents may be removed using standard endoscopic techniques.

DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be disposed in accordance with institutional guidelines. The introducer should be coiled for disposal.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

When available, the EUDAMED website (<https://ec.europa.eu/tools/eudamed>), along with the BUDI for this product (0827002CIRL202007013010B6), can be used to locate the Summary of Safety and Clinical Performance (SSCP) for this product.

Patient-facing information may be accessed at cookmedical.eu/patient-implant-information.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2026
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.