
- EN**
5 **EchoTip® Ultra Ultrasound Needle**
Instructions for Use
- CS**
8 **Ultrazvuková jehla EchoTip® Ultra**
Návod k použití
- DA**
11 **EchoTip® Ultra ultralydnål**
Brugsanvisning
- NL**
14 **EchoTip® Ultra ultrasone naald**
Gebruiksaanwijzing
- FR**
17 **Aiguille ultrasonore EchoTip® Ultra**
Mode d'emploi
- DE**
20 **EchoTip®-Ultra Ultraschallnadel**
Gebrauchsanweisung
- EL**
23 **Βελόνα υπερήχων EchoTip® Ultra**
Οδηγίες χρήσης
- HU**
26 **EchoTip® Ultra echogén tű**
Használati utasítás
- IT**
29 **Ago EchoTip® Ultra ecogenico**
Istruzioni per l'uso
- NO**
32 **EchoTip® Ultra ultralydnål**
Bruksanvisning
- PL**
35 **Igła ultrasonograficzna EchoTip® Ultra**
Instrukcja użycia
- PT**
38 **Agulha de ultra-sons EchoTip® Ultra**
Instruções de utilização
- ES**
41 **Aguja de ecografía EchoTip® Ultra**
Instrucciones de uso
- SV**
44 **EchoTip® Ultra ultraljudsnål**
Bruksanvisning

1.


Plunger

Píst
Stempel
Plunjer
Piston
Kolben
Έμβολο
Dugattyú
Stantuffo
Stempel
Tłoczek
Έμβολο
Émbolo
Kolb

Stopcock "open"

Uzavírací kohout v otevřené poloze
Hane "åben"
Afsluitkraan 'open'
Robinet ouvert
Absperrhahn geöffnet
Στρόφιγγα ανοικτή
Elzárócsap „nyitott” állásban
Rubinetto aperto
Stoppekran åpen
Kranik w pozycji otwartej
Torneira de passagem "aberta"
Llave de paso «abierta»
Kran "öppen"

Sliding Sheath Adjuster

Posuvné nastavovací zařízení sheathu
 Glidende sheath-justering
 Schuifafstelling van de sheath
 Ajusteur de gaine coulissant
 Verschiebbarer Schleuseneinsteller
 Συρόμενος προσαρμογέας θηκαριού
 Elcsúsztható hüvelyhossz-beállító
 Regolatore della guaina scorrevole
 Glidende hylsejusteringsmekanisme
 Przesuwany regulator koszulki
 Regulador deslizante da bainha
 Ajustador deslizante de la vaina
 Flyttbar justerare till höljet

Sheath Reference Mark

Referenční značka sheathu
 Referencemærke på sheath
 Referentiemarkering van de sheath
 Repère de référence sur la gaine
 Schleusenbezugsmarke
 Σήμανση αναφοράς θηκαριού
 Hüvelyhossz referencijele
 Indicatore di riferimento della guaina
 Referansemerke på hylsen
 Znacznik referencyjny wysunięcia koszulki
 Marca de referència da bainha
 Marca de referencia de la vaina
 Höljets referensmarkering

Needle Handle

Rukojeť jehly
 Nålegreb
 Naaldhandvat
 Poignée de l'aiguille
 Nadelgriff
 Λαβή βελόνας
 Tű fogantyúja
 Impugnatura dell'ago
 Nålehåndtak
 Uchwyt igły
 Punho da agulha
 Mango de la aguja
 Nålhandtag

Safety Ring

Pojistný prstenec
 Sikkerhedsring
 Veiligheidsring
 Bague de sécurité
 Sicherheitsring
 Δακτύλιος ασφάλειας
 Biztonsági gyűrű
 Anello di sicurezza
 Sikkerhetsring
 Pierścień bezpieczeństwa
 Anel de sigurança
 Anillo de seguridad
 Säkerhetslåsring

Stylet Hub

Nástavec stiletu
 Stiletmuffe
 Aanzetstuk van de stilet
 Embase du stylet
 Mandrinansatz
 Ομφαλός στείλεου
 Mandrin végidoma
 Pomello del mandrino
 Stilettfeste
 Głowica mandrynu
 Conector do estilete
 Conector del estilete
 Stilettfattning

3.

4.

INTENDED USE

This device is used to sample targeted submucosal gastrointestinal lesions through the accessory channel of an ultrasound endoscope.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired site.

Coagulopathy.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Needle must be retracted into sheath and thumbscrew on safety ring must be locked to hold needle in place prior to introduction, advancement or withdrawal of device. Failure to retract needle may result in damage to endoscope.

SYSTEM PREPARATION

Illustrations

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of syringe has a Luer slip fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (**See fig. 1**)
2. Prepare syringe as follows:
 - a. With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b. Turn stopcock 90° to "closed" position.

- c. Pull plunger back until it is locked into place at desired setting, creating a suction.
 - d. Set prepared syringe aside until aspiration biopsy is desired.
3. Advance device into ultrasound endoscope to determine preferred sheath length. To adjust length, loosen thumbscrew lock on sliding sheath adjuster and slide until preferred length is attained. **Note:** Reference mark for sheath length will appear in sliding sheath adjuster window. (See fig. 2) Tighten thumbscrew on sliding sheath adjuster to maintain preferred sheath length.

INSTRUCTIONS FOR USE

Illustrations

1. Identify desired biopsy site by endoscopic ultrasound.
2. Introduce ultrasound needle into accessory channel and advance device in small increments until Luer lock fitting at base of sliding sheath adjuster meets fitting on accessory channel. (See fig. 3)
3. Attach device to accessory channel port by rotating device handle until fittings are connected.
4. With ultrasound endoscope and device straight, adjust needle to desired length by loosening thumbscrew on safety ring, and advancing it until desired reference mark for needle advancement appears in the window of safety ring. (See fig. 4) Tighten thumbscrew to lock safety ring in place. **Note:** Number in safety lock ring window indicates extension of needle in centimeters.
Caution: During needle adjustment or extension, ensure device has been attached to accessory channel. Failure to attach device prior to needle adjustment or extension may result in damage to endoscope.
5. While maintaining position of ultrasound endoscope, extend needle by advancing needle handle to pre-positioned safety ring.
6. Advance needle into lesion.
7. Remove stylet from needle by gently pulling back on plastic hub seated in metal fitting of needle handle. Preserve stylet for use if additional cell collection is desired.
8. Attach Luer lock fitting of previously prepared syringe securely onto fitting on handle.
9. Turn stopcock to "open" position aligned with syringe, allowing negative pressure in syringe to aspirate cells. Gently move handle in small increments back and forth within biopsy site. **Note:** Do not remove needle from biopsy site during cell collection.
10. Upon completion of cell collection, retract needle into sheath by pulling back on handle shaft. Lock thumbscrew on safety ring at 0 cm mark.
11. Disconnect Luer lock fitting from accessory channel and withdraw entire device from ultrasound endoscope.

12. Unlock thumbscrew on safety ring and extend needle. Exchange locking syringe for an empty syringe and attach new syringe to Luer lock fitting on handle.
13. Push in plunger to expel aspirate specimen, then prepare per institutional guidelines.
14. Aspirate remaining can be recovered for examination by flushing device.
15. For additional cell collection from same lesion, gently reinsert stylet into metal fitting on handle. **Note:** Prior to reinserting stylet, wipe with saline or sterile water. While supporting sheath at Luer lock fitting, advance stylet in small increments until stylet hub is engaged in fitting.
16. Additional samples may be obtained by repeating Step 2 of "System Preparation" through Step 14 of the "Instructions for Use."

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2013 ©
COOK IRELAND LTD.