
EN 5	EchoTip® Ultra Endoscopic Ultrasound Needle Instructions for Use
CS 8	Endoskopická ultrazvuková jehla EchoTip® Ultra Návod k použití
DA 11	EchoTip® Ultra endoskopisk ultralydsnål Brugsanvisning
DE 14	EchoTip® Ultra endoskopische Ultraschallnadel Gebrauchsanweisung
EL 17	Ενδοσκοπική βελόνα υπερήχων EchoTip® Ultra Οδηγίες χρήσης
ES 20	Aguja ecoendoscópica EchoTip® Ultra Instrucciones de uso
FR 23	Aiguille ultrasonore endoscopique EchoTip® Ultra Mode d'emploi
HR 26	Endoskopska ultrazvučna igla EchoTip® Ultra Upute za uporabu
HU 29	EchoTip® Ultra endoszkópiás ultrahangos tű Használati utasítás
IT 32	Ago ecoendoscopico EchoTip® Ultra Istruzioni per l'uso
NL 35	EchoTip® Ultra endoscopische ultrasone naald Gebruiksaanwijzing
NO 38	EchoTip® Ultra endoskopisk ultralydnål Bruksanvisning
PL 41	Endoskopowa igła ultrasonograficzna EchoTip® Ultra Instrukcja użycia
PT 44	Agulha de ultrassons endoscópica EchoTip® Ultra Instruções de utilização
RO 47	Ac pentru ecografie endoscopică EchoTip® Ultra Instrucțiuni de utilizare
SK 50	Endoskopická ultrazvuková ihla EchoTip® Ultra Návod na použitie
SV 53	EchoTip® Ultra endoskopisk ultraljudsnål Bruksanvisning
TR 56	EchoTip® Ultra Endoskopik Ultrason İğnesi Kullanma Talimatı

1.


Plunger

Pist
Stempel
Kolben
Έμβολο
Émbolo
Poussoir
Klip
Dugattyú
Stantuffo
Plunjer
Stempel
Tłoczek
Έμβολο
Piston
Piest
Kolv
Piston

Stopcock "open"

„Otvorený“ uzavírací kohout
Stophane „åben“
Absperrhahn geöffnert
Στρόφιγγα ανοικτή
Llave de paso «abierta»
Robinet ouvert
Sigumosni ventil otvoren
Elzárócsap nyitott állásban
Rubinetto aperto
Afsluitkraan 'open'
Stoppekran åpen
Kranik w pozycji otwartej
Torneira de passagem "aberta"
Robinet „deschis“
Otvorený kohútik
Kran "öppen"
Stopkok "açık"

2.

Sliding Sheath Adjuster

Posuvné nastavovací zařízení sheathu
 Glidende sheath-justering
 Verschiebbarer Schleuseneinsteller
 Συρόμενος προσαρμογέας θηκαιοῦ
 Ajustador deslizante de la vaina
 Ajusteur de gaine coulissant
 Klizni podešivač uvodnice
 Elcsúsztatható hüvelyhossz-beállító
 Regolatore scorrevole della guaina
 Schuifregeling van de sheath
 Glidende hylsejusteringsmekanisme
 Przesuwany regulator koszulki
 Regulador deslizante da bainha
 Ajustor glisant al tecii
 Posuvný nastavovač puzdra
 Glidande hylsjusterare
 Kayar Kilif Ayarlayıcı

Sheath Reference Mark

Referenční značka sheathu
 Referencemærke på sheath
 Schleusenbezugsmarke
 Σήμανση αναφοράς θηκαιοῦ
 Marca de referencia de la vaina
 Repère de référence sur la gaine
 Referentna oznaka uvodnice
 Hüvelyhossz referencijele
 Indicatore di riferimento della guaina
 Referentiemarkering van de sheath
 Referansemerke på hylsen
 Znacznik referencyjny wysunięcia koszulki
 Marca de referència da bainha
 Marcaj de referință al tecii
 Referenčná značka puzdra
 Hylsans referensmarkering
 Kilif Referans İşareti

Needle Handle

Rukojet' jehly
 Nålegreb
 Nadelgriff
 Λαβή βελόνας
 Mango de la aguja
 Poignée de l'aiguille
 Drška igle
 Tű fogantyúja
 Impugnatura dell'ago
 Naaldhandvat
 Nähåndtak
 Uchwyt igły
 Punho da agulha
 Måner al acului
 Rukováť ihly
 Nähåndtag
 İğne Sapi

Safety Ring

Pojistný prsteneč
 Sikkerhedsring
 Sicherheitsring
 Δακτύλιος ασφαλείας
 Anillo de seguridad
 Bague de sécurité
 Sigurnosni prsten
 Biztonsági gyűrű
 Anello di sicurezza
 Veiligheidsring
 Sikkerhetsring
 Pierścień bezpieczeństwa
 Anel de sigurança
 Inel de siguranță
 Bezpečnostný krúžok
 Säkerhetslåsring
 Güvenlik Halkası

Stylet Hub

Nástavec styletu
 Stiletmuffe
 Mandrinansatz
 Ομφαλός του σπειλεού
 Conector del estilete
 Embase du stylet
 Čvoršte sonde
 Mandrin végidoma
 Pomello del mandrino
 Aanzetstuk van het stilet
 Stilettfeste
 Złącza mandryny
 Conector do estilete
 Acoperitoarea stiletului
 Hrdlo styletu
 Mandrängfattning
 Stile Göbeği

3.

ECHOTIP® ULTRA ENDOSCOPIC ULTRASOUND NEEDLE

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

INTENDED USE

This device is used to sample targeted submucosal gastrointestinal lesions through the accessory channel of an ultrasound endoscope.

DEVICE DESCRIPTION

The device consists of a handle with adjustable components to allow the user to adjust the extension of the needle and sheath. The device has a 5.2 Fr outer sheath for protection of the needle and the scope, and it is attached to the accessory channel of an ultrasound endoscope. The device has a maximum needle extension length of 8cm and sheath adjustment length of 5cm. The device is supplied with a 10ml syringe to utilize standard vacuum technique for sample collection if desired. The device is available in three different needle sizes: 19, 22 and 25 ga.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired site. Coagulopathy.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with gastrointestinal endoscopy: • acute pancreatitis • allergic reaction to medication • allergic reaction to nickel • aspiration • cardiac arrhythmia or arrest • damage to blood vessels • fever • hemorrhage • hypotension • infection • pain / discomfort • perforation • peritonitis • respiratory depression or arrest • tumor seeding (through the needle tract).

WARNINGS

- Not for use in the heart or vascular system.
- The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.
- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical and request return authorization.

PRECAUTIONS

- Refer to package label for minimum channel size required for this device.
- The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the endoscope.
- Do not use this device for any purpose other than the stated intended use.
- Use of this device is restricted to a trained healthcare professional (Rx only).
- When targeting multiple sites, adequate care should be taken to minimize the risk of tumor seeding.
- The stylet must be retracted (approximately 1 cm) from the Luer fitting on the needle handle prior to puncture.

SYSTEM PREPARATION

Illustrations

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when the stopcock is in "open" position, aligned with syringe. (**Fig. 1**)
2. Prepare syringe as follows:
 - a. With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b. Turn stopcock 90° to "closed" position.
 - c. Pull plunger back until it is locked into place at desired setting, creating a vacuum.
 - d. Set prepared syringe aside until aspiration is required.

INSTRUCTIONS FOR USE

Illustrations

1. Identify the desired target site by endoscopic ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope. **Caution:** If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device in small increments until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting on the accessory channel port. (**Fig. 3**)
4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window. (**Fig. 2**) Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining position of the scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring. (**Fig. 3**) Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters.
Caution: During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
7. Prior to puncture, retract the stylet slightly (approximately 1 cm) from the Luer fitting on the needle handle, to ensure the sharpened needle tip is exposed. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
8. Standard vacuum syringe techniques may be applied for sample collection (Steps 9 – 11) or, if desired, other techniques that may or may not incorporate the use of the stylet may be used.
9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional sample collection is required later.
10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer fitting on the handle.
11. Turn the stopcock to the "open" position aligned with the syringe, allowing negative pressure in the syringe to aspirate the sample. Move the needle handle in small increments back and forth within the target site.
Note: Do not remove the needle from the target site during sample collection.

12. Upon completion of sample collection, if using the syringe, turn the stopcock to the “closed” position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
13. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
14. Unlock the thumbscrew on the safety ring and extend the needle.
15. Use a new syringe or the stylet to expel the specimen, then prepare per institutional guidelines.
16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.
17. For additional sample collection from the same lesion, gently reinsert the stylet into the Luer fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer lock fitting, advance the stylet in small increments until the hub is engaged in the fitting.
18. Additional samples may be obtained by repeating Step 2 through Step 16 of the “Instructions for Use.”
19. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Store in a dry location, away from temperature extremes.

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© COOK 2021
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.