EN 2	14Fr Colon Decompression Set Instructions for Use
BG 5	Комплект за декомпресия на дебелото черво 14 Fr Инструкции за употреба
HR 9	Komplet za dekompresiju debelog crijeva od 14 Fr Upute za uporabu
CS 12	Souprava pro kolonoskopickou dekompresi velikosti 14 Fr Návod k použití
DA 15	14 Fr Colondekompressionsæt Brugsanvisning
NL 18	14Fr Colondecompressieset Gebruiksaanwijzing
FR 21	Set pour décompression du côlon 14 Fr Mode d'emploi
DE 24	Kolon-Dekompressions-Set (14 Fr) Gebrauchsanweisung
EL 27	Σετ αποσυμπίεσης παχέος εντέρου 14 Fr Οδηγίες χρήσης
IT 31	Set di decompressione del colon da 14 Fr Istruzioni per l'uso
PL 34	Zestaw 14 Fr do odbarczania okrężnicy Instrukcja użycia
PT 37	Conjunto de descompressão do cólon 14 Fr Instruções de utilização
SK 40	Súprava na dekompresiu hrubého čreva veľkosti 14 Fr Návod na použitie
ES 43	Equipo de descompresión del colon de 14 Fr Instrucciones de uso
SV 46	14 Fr Kolondekompressionsset Bruksanvisning

14FR COLON DECOMPRESSION SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

INTENDED USE

This device is intended for treatment of acute non-toxic megacolon, pseudo-obstruction (Ogilvie's syndrome) and colonic strictures.

DEVICE DESCRIPTION

Set includes a colon decompression tube, guiding catheter and .035 inch (0.89mm) wire guide. The colon decompression tube, guiding catheter and wire guide are radiopaque to aid with fluoroscopic placement.

CONTRAINDICATIONS

Those for colonoscopy include, but are not limited to: fulminant colitis • possible existing perforation • acute severe diverticulitis • presence of barium or inadequate bowel preparation.

Those for colon decompression include, but are not limited to: ischemic bowel • recent rectal surgery • organic colon obstruction • recent surgical anastomosis • acute myocardial infarction.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with colonoscopy include, but are not limited to: allergic reaction to medication • cardiac arrhythmia or arrest • fever • hemorrhage • hypotension • infection • perforation • respiratory depression or arrest.

Those associated with colon decompression include, but are not limited to: clogging of the tube • dislodgment of the colon decompression tube / wire guide • fever and discomfort • hemorrhage • hypotension • infection • migration • pain • perforation • sepsis.

PRECAUTIONS

The product is intended for use by physicians trained and experienced in colonoscopy techniques.

Standard technique for the treatment of acute non-toxic megacolon, pseudo-obstruction (Ogilvie's syndrome) and colonic strictures should be employed.

Refer to package label for minimum channel size required for this device.

Careful history and physical examination must be performed prior to colon decompression tube insertion in order to determine underlying cause.

Do not use this device for any purpose other than stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device is restricted to a trained healthcare professional.

It is recommended that drainage collection bag be available for use during procedure.

Colon should not be prepped unless there is an excessive amount of stool in rectum.

This device should not be left indwelling for more than 72 hours.

WARNINGS

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Visually inspect the packaging and device. If the packaging is opened or damaged when received, do not use. If a device abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

MR SAFETY INFORMATION

This symbol means the device is MR Safe.

Non-clinical testing performed to establish the MR safety was completed using a static magnet field of 3 T.

In non-clinical testing, the image artifact caused by the device extends approximately 4 mm from the 14 Fr Colon Decompression Tube when imaged with a gradient echo pulse sequence and a 3 T MR system.

INSTRUCTIONS FOR USE

1. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
2. Visually inspect the device confirming it is free from damage prior to use.
3. Advance colonoscope, insufflating minimal air or CO₂, until it has reached the desired location. Use suction intermittently throughout procedure to help decompress colon. **Caution: Failure to perform suction may result in further colonic distention.**
4. Fluoroscopically observe location of colonoscope.
5. Insert wire guide into accessory channel until it is endoscopically visualized exiting colonoscope and advance in short increments until it is fluoroscopically visualized in position at the desired location.
6. Carefully withdraw colonoscope, making certain wire guide remains stationary within colon by periodic fluoroscopic monitoring. **Caution: Failure to monitor the wire guide may result in dislodgement of the wire guide.**
7. Disassemble colon decompression tube and pre-loaded guiding catheter.

8. Flush colon decompression tube and guiding catheter with sterile water or lubricate with water-soluble lubricant, then reassemble and advance into colon over pre-positioned wire guide under periodic fluoroscopic monitoring.
Caution: Failure to monitor the position of the wire guide fluoroscopically may result in dislodgement of the wire guide.
9. Once colon decompression tube is advanced to target area, under periodic fluoroscopic monitoring withdraw both guiding catheter and wire guide, making certain the colon decompression tube remains stationary within colon.
Caution: Failure to monitor the colon decompression tube fluoroscopically may result in dislodgement of the colon decompression tube.
10. Secure external portion of colon decompression tube to patient. Attach drainage connector to appropriate drainage collection bag to begin drainage, per institutional guidelines.

Removal

11. Carefully remove the colon decompression tube by gently pulling it through the colon until it exits.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages.

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This symbol on the label indicates that this device contains phthalates. Specific phthalates contained in the device are identified beside or below the symbol by the following acronyms:

- BBP: Benzyl butyl phthalate
- DBP: Di-n-butyl phthalate
- DEHP: Di(2-ethylhexyl) phthalate
- DIDP: Diisodecyl phthalate
- DINP: Diisononyl phthalate
- DIPP: Diisopentyl phthalate
- DMEP: Di(methoxyethyl) phthalate
- DNOP: Di-n-Octyl phthalate
- DNPP: Di-n-pentyl phthalate

The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

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© COPYRIGHT COOK 2020 ©
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.