EN
3
Marcon Colon Decompression Set

Instructions for Use

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
Bruksanvisning

MARCON COLON DECOMPRESSION SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

INTENDED USE

This device is intended for treatment of acute non-toxic megacolon, pseudo-obstruction (Ogilvie's syndrome) and colonic strictures.

DEVICE DESCRIPTION

Set includes a colon decompression tube and a .035 inch (0.89 mm) wire guide. The colon decompression tube and wire guide are radiopaque to aid with fluoroscopic placement.

CONTRAINDICATIONS

Those for colonoscopy include, but are not limited to: fulminant colitis • possible existing perforation • acute severe diverticulitis • presence of barium or inadequate bowel preparation.

Those for colon decompression include, but are not limited to: ischemic bowel • recent rectal surgery • organic colon obstruction • recent surgical anastomosis • acute myocardial infarction.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with colonoscopy include, but are not limited to: allergic reaction to medication • cardiac arrhythmia or arrest • fever • hemorrhage • hypotension • infection • perforation • respiratory depression or arrest.

Those associated with colon decompression include, but are not limited to: clogging of the tube • dislodgment of the colon decompression tube / wire guide • fever and discomfort • hemorrhage • hypotension • infection • migration • pain • perforation • sepsis.

PRECAUTIONS

The product is intended for use by physicians trained and experienced in colonoscopy techniques.

Standard technique for the treatment of acute non-toxic megacolon, pseudo-obstruction (Ogilvie's syndrome) and colonic strictures should be employed.

Refer to package label for minimum channel size required for this device.

Careful history and physical examination must be performed prior to colon decompression tube insertion in order to determine underlying cause.

Do not use this device for any purpose other than stated intended use.

Use of this device is restricted to a trained healthcare professional.

It is recommended that drainage collection bag be available for use during procedure.

Colon should not be prepped unless there is an excessive amount of stool in rectum.

The device should not be left indwelling for more than 72 hours.

Store in a dry location, away from temperature extremes.

WARNINGS

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Visually inspect the packaging and device. If packaging is opened or damaged when received, do not use. If a device abnormality is detected what would prohibit proper working condition, do not use. Please notify Cook for return authorization.

MR SAFETY INFORMATION

This symbol means the device is MR Safe.

Non-clinical testing performed to establish the MR safety was completed using a static magnet field of 3 T.

In non-clinical testing, the image artifact caused by the device extends approximately 4 mm from the Marcon Colon Decompression Tubes when imaged with a gradient echo pulse sequence and a 3 T MR system.

INSTRUCTIONS FOR USE

Illustrations

1. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
2. Visually inspect the device confirming it is free from damage prior to use.
3. Advance colonoscope, insufflating minimal air or CO₂, until it has reached the desired location. **Caution: Failure to perform suction may result in further colonic distention.**
4. Fluoroscopically observe location of colonoscope.
5. Flush the colon decompression tube if required with saline or sterile water.
6. Preload the colon decompression tube with wire guide until wire guide exits tip of tube.
7. Advance preloaded colon decompression tube through colonoscope to target area, carefully observing the wire guide tip and colon decompression tube fluoroscopically.

8. Remove drainage connector. **(See fig. 1)**
9. Under periodic fluoroscopic monitoring, carefully withdraw colonoscope, making certain the colon decompression tube and wire guide remain stationary within colon. **Caution: Failure to periodically monitor the colon decompression tube fluoroscopically may result in dislodgement of the colon decompression tube/wire guide.**
10. Under periodic fluoroscopic monitoring, remove wire guide making certain the colon decompression tube remains stationary within colon. **Caution: Failure to periodically monitor the colon decompression tube fluoroscopically may result in dislodgement of the colon decompression tube.**
11. Reattach drainage connector to the colon decompression tube.
12. Secure external portion of decompression tube to patient. Attach drainage connector to appropriate drainage collection bag to begin drainage, per institutional guidelines.

Removal

13. Carefully remove the colon decompression tube by gently pulling it through the colon until it exits.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages.

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© COOK 2021
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.