

EN 3	Nasal Pancreatic Drainage Set Instructions for Use
CS 5	Souprava pro nasální pankreatickou drenáž Návod k použití
DA 7	Nasal pancreasdrænaesæt Brugsanvisning
NL 9	Nasale pancreasdrainageset Gebruiksaanwijzing
FR 11	Set de drainage pancréatique par voie nasale Mode d'emploi
DE 13	Nasales Pankreas-Drainage-Set Gebrauchsanweisung
EL 15	Σετ ρινικής παγκρεατικής παροχέτευσης Οδηγίες χρήσης
HU 17	Nazális hasnyáldrenázskészlet Használati utasítás
IT 19	Set per il drenaggio naso-pancreatico Istruzioni per l'uso
PL 21	Zestaw do przeznosowego drenażu trzustki Instrukcja użycia
PT 23	Conjunto de drenagem naso-pancreática Instruções de utilização
ES 25	Equipo de drenaje pancreático nasal Instrucciones de uso
SV 27	Nasalt pankreasdränageset Bruksanvisning

4.

3.

INTENDED USE

This device is used for temporary endoscopic drainage of the pancreatic duct through the nasal passage by use of an indwelling catheter.

NOTES

Do not use this device for any purpose other than stated intended use.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

Each set includes a drainage catheter, nasal transfer tube and a drainage connecting tube to facilitate irrigation and drainage.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

Contraindications of nasal pancreatic drainage include, but are not limited to: coagulopathy, sepsis, active pancreatitis or other infection, recent food ingestion, inability to pass the drainage catheter through the obstructed area.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Complications associated with nasal pancreatic drainage include, but are not limited to: trauma to the pancreatic duct or duodenum, obstruction of the common bile duct, blockage of the drainage catheter, nasal irritation and/or sore throat.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

INSTRUCTIONS FOR USE

Illustrations

1. Flush drainage catheter using a syringe filled with sterile water or saline attached to Luer lock fitting on end of catheter.
2. Carefully remove Tuohy-Borst connector from end of drainage catheter. **(See fig. 1)**
3. Introduce and advance drainage catheter in small increments into accessory channel of endoscope over pre-positioned wire guide until it is endoscopically visualized exiting endoscope. **Note:** Elevator of endoscope must be open to allow drainage catheter to exit.
4. Under fluoroscopic monitoring, slowly advance drainage catheter to desired position in pancreatic duct above obstruction.
5. When drainage catheter is in desired position, slowly remove endoscope while holding drainage catheter and wire guide in place at patient's mouth.

6. Fluoroscopically confirm that drainage catheter and wire guide remain positioned in pancreatic duct, then slowly remove wire guide in small increments while maintaining position of drainage catheter.
7. Lubricate and advance nasal transfer tube through nostril until it exits in back of throat.
8. Grasp connection end of nasal transfer tube in throat and pull it out through mouth.
(See fig. 2)
9. Thread proximal end of drainage catheter through oral end of nasal transfer tube and advance until it exits through nasal end. (See fig. 3)
10. Hold drainage catheter, to prevent kinking and removal, and slowly pull nasal transfer tube until drainage catheter emerges through nostril. (See fig. 4)
11. Replace Tuohy-Borst connector onto end of drainage catheter. Contrast may be injected to check placement of catheter within duct.
12. Secure drainage catheter to patient, then connect end of catheter to Luer lock fitting on drainage connecting tube. Attach universal connector to appropriate drainage bag to begin drainage, per institutional guidelines.

Upon completion of procedure, dispose of components per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2013 ©
COOK IRELAND LTD.