EN
8
EchoTip ProCore® Endobronchial HD Ultrasound Biopsy Needle for Olympus Scopes

Instructions for Use

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



Fig. 7

INTENDED USE

This device is used with an ultrasound endoscope for fine needle biopsy, (FNB), of submucosal and extramural lesions within or adjacent to the tracheobronchial tree or gastrointestinal tract.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If the package is opened or damaged when received, do not use.

Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify COOK Medical for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device is restricted to a trained healthcare professional.

This device is intended for use with an Olympus EBUS scope.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications include, but are not limited to: coagulopathy.

POTENTIAL COMPLICATIONS

Potential complications associated with bronchial endoscopy and gastrointestinal endoscopy include, but are not limited to: Perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, bleeding, damage to blood vessels, nerve damage, pain, pneumoeritoneum, tumor seeding, respiratory depression or arrest, cardiac arrhythmia or arrest and death.

WARNING

Not for use in the heart or vascular system. The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Ensure the stylet is fully inserted when advancing the needle into the biopsy site.

When targeting multiple sites, replace device for each site.

The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the endoscope.

PLACEMENT OF ACCESSORY CHANNEL ADAPTER

Illustrations

1. Examine the adapter and note features. (See Fig. 1)
2. Place the adapter, in the open position, in line with the endoscope accessory channel's metal hub.
3. Slide the adapter on to the accessory channel ensuring that the base of the adapter slides under the accessory channel's metal hub. (See Fig. 2)
4. Close the adapter ensuring it is locked. (See Fig. 3)
5. The adapter will now function as a Luer lock to which the needle can be connected. **Note:** If irrigation of accessory channel is desired/necessary, irrigate through closed adapter prior to insertion or after removal of needle.

SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of the syringe has a Luer slip fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (See Fig. 4)
2. Prepare syringe as follows:
 - a) With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b) Turn stopcock to "closed" position.
 - c) Pull plunger back until it is locked into place at desired setting, creating a suction.
 - d) Set prepared syringe aside until biopsy is desired.
3. Advance the device into the endoscope to determine preferred sheath length. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window. (See Fig. 5) Tighten thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.

INSTRUCTIONS FOR USE

1. Identify the desired site by endoscopic ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the endoscope. (See Fig. 6)
Note: If resistance is encountered on needle introduction, reduce angulation of scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
4. Attach the device to the endoscopes accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is endoscopically visible, confirming that the sheath has emerged from the working channel of the scope.
6. With the endoscope and device straight, set the needle extension to the desired length by loosening the thumbscrew on the safety ring and advancing it until the desired reference mark for the needle advancement appears in the window of the safety ring. (See Fig. 7) Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the endoscope. Failure to attach the device prior to needle adjustment or extension may result in damage to the endoscope.
7. While maintaining the position of the endoscope, extend the needle by advancing the needle handle of the device to the pre-positioned safety ring.
8. Advance the needle into the lesion.
9. Standard vacuum syringe techniques may be applied for biopsy (see steps 10 – 12) or, if desired, other techniques that may or may not incorporate the use of the stylet may be used.
10. Remove the stylet of the ultrasound needle by gently pulling back on the plastic hub seated in the metal fitting of the needle handle. Preserve the stylet for use if additional biopsy is desired later.
11. Attach the Luer lock fitting of the previously prepared syringe securely into the fitting on the needle handle.
12. Turn the stopcock to the “open” position aligned with the syringe, allowing the negative pressure in the syringe to facilitate biopsy. Move the needle handle back and forth within the biopsy site. **Note:** Do not remove the needle from the biopsy site during biopsy.
13. Upon completion of biopsy, if using a syringe, turn the stopcock to the “closed” position, retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.

14. Disconnect the Luer lock fitting of the device from the accessory channel port by rotating the device handle counter clockwise, and withdraw the entire device from the endoscope.

15. Exchange the locking syringe for an empty syringe and attach the new syringe to the Luer lock fitting on the needle handle. Unlock the thumbscrew on the safety ring and extend the needle.

16. Use syringe or stylet to expel specimen, then prepare per institutional guidelines.

17. For additional biopsy from the same lesion, gently reinsert the stylet into the metal fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer lock fitting, advance the stylet in small increments, until the hub is engaged in the fitting.

18. Additional samples may be obtained by repeating Step 3 of "System Preparation" through step 16 of the "Instructions for Use."

19. To remove adapter from endoscope, pull back on the locking clip, lift the top and slide off endoscope.

Upon completion of the procedure, dispose of the device, including the adapter, per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2015 ©
COOK IRELAND LTD.