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


Fig. 7

ECHOTIP PROCORE® ENDOBRONCHIAL HD ULTRASOUND BIOPSY NEEDLE FOR OLYMPUS SCOPES

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

INTENDED USE

This device is used with an ultrasound endoscope for fine needle biopsy, (FNB), of submucosal and extramural lesions within or adjacent to the tracheobronchial tree or gastrointestinal tract.

DEVICE DESCRIPTION

The device consists of a handle with adjustable components to allow the user to adjust the extension of the needle and sheath. The device has a 4.1 Fr outer sheath for protection of the needle and the scope, and it is attached to the accessory channel of an Olympus endobronchial ultrasound (EBUS) scope using an adapter. The device has a maximum needle extension length of 5 cm and sheath adjustment length of 3 cm. The device is supplied with a 10 mL syringe to utilize the standard vacuum technique for biopsy, if desired, and an accessory channel adapter. The device is available in two different needle sizes: 22 and 25 ga.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications include, but are not limited to: coagulopathy.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with bronchial endoscopy and gastrointestinal endoscopy include, but are not limited to: allergic reaction to medication, allergic reaction to nickel, aspiration, cardiac arrhythmia or arrest, damage to blood vessels, death, discomfort, fever, hemorrhage, hypotension, infection, inflammation, nerve damage, pain, perforation, pneumoperitoneum, respiratory depression or arrest, sepsis, septicemia/bacteremia and tumor seeding.

WARNINGS

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.

PRECAUTIONS

Refer to package label for minimum channel size required for this device. Ensure the stylet is fully inserted when advancing the needle into the target site.

The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the scope.

Do not use this device for any purpose other than stated intended use.

Use of this device is restricted to a trained healthcare professional (Rx only).

When targeting multiple sites, adequate care should be taken to minimize the risk of tumor seeding.

This device is intended for use with an Olympus EBUS scope.

PLACEMENT OF ACCESSORY CHANNEL ADAPTER

Illustrations

1. Examine the adapter and note features (**Fig. 1**).
2. Place the adapter, in the open position, in line with the scope accessory channel metal hub.
3. Slide the adapter on to the accessory channel ensuring that the base of the adapter slides under the accessory channel metal hub (**Fig. 2**).
4. Close the adapter ensuring it is locked (**Fig. 3**).
5. The adapter will now function as a Luer lock to which the needle can be connected. **Note:** If irrigation of accessory channel is desired/necessary, irrigate through closed adapter prior to insertion or after removal of needle.

SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when the stopcock is in "open" position, aligned with syringe (**Fig. 4**).
2. Prepare syringe as follows:
 - a. With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b. Turn stopcock to "closed" position.
 - c. Pull plunger back until it is locked into place at desired setting, creating a vacuum.
 - d. Set prepared syringe aside until biopsy is required.

INSTRUCTIONS FOR USE

1. Identify the desired target site by endobronchial ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope (**Fig. 5**).
Caution: If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (**Fig. 6**). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining position of the scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring (**Fig. 7**). Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters.
Caution: During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
8. Standard vacuum syringe techniques may be applied for biopsy (Steps 9 – 11) or, if desired, other techniques that may or may not incorporate the use of the stylet may be used.
9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if an additional biopsy is required later.
10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer fitting on the needle handle.

11. Turn the stopcock to the “open” position aligned with the syringe, allowing the negative pressure in the syringe to facilitate biopsy. Move the needle handle back and forth within the biopsy site. **Note:** Do not remove the needle from the biopsy site during biopsy.
12. Upon completion of biopsy, if using the syringe, turn the stopcock to the “closed” position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
13. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
14. Unlock the thumbscrew on the safety ring and extend the needle.
15. Use a new syringe or the stylet to expel the specimen, then prepare per institutional guidelines.
16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.
17. For additional biopsy from the same target site, gently reinsert the stylet into the Luer fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.
18. Additional samples may be obtained by repeating Step 2 through 16 of the “Instructions for Use.”
19. Upon completion of the procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
20. To remove the adapter from the scope, pull back on the locking clip, lift the top and slide off the scope.

Upon completion of the procedure, dispose of the device, including the adapter, per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Store in a dry location, away from temperature extremes.

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© COOK 2021
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.