EN 9	EchoTip ProCore® Endobronchial HD Ultrasound Biopsy Needle for Olympus Scopes Instructions for Use
BG 13	Ендобронхиална ултразвукова игла за биопсия с висока разделителна способност EchoTip ProCore® за ендоскопи Olympus Инструкции за употреба
CS 18	Endobronchiální biptická jehla EchoTip ProCore® pro ultrazvuk s vysokým rozlišením, pro endoskopy Olympus Návod k použití
DA 22	EchoTip ProCore® endobronkial HD ultralydsbiopsinål til Olympus skoper Brugsanvisning
DE 26	EchoTip ProCore® endobronchiale HD-Ultraschallbiopsienadel für Olympus-Endoskope Gebrauchsanweisung
EL 31	Ενδοβρογχική βελόνα βιοψίας υπερήχων υψηλής ανάλυσης EchoTip ProCore® για ενδοσκόπια Olympus Οδηγίες χρήσης
ES 36	Aguja de biopsia por ecografía endobronquial de alta definición EchoTip ProCore® para endoscopios Olympus Instrucciones de uso
ET 41	Endobronhiaalne kõrglahutusega ultrahelinõel EchoTip ProCore® endoskoopidele Olympus Kasutusjuhend
FI 45	Endobronkiaalinen EchoTip ProCore® -HD-ultraäänibiopsianeula Olympus-skooppeja varten Käyttöohjeet
FR 49	Aiguille de biopsie ultrasonore endobronchique EchoTip ProCore® HD pour les endoscopes Olympus Mode d'emploi
HR 54	Endobronhijalna HD igla za ultrazvučnu biopsiju EchoTip ProCore® za endoskope Olympus Upute za uporabu
HU 58	EchoTip ProCore® endobronchialis nagy felbontású ultrahangos biopsziás tű Olympus endoszkópokhoz Használati utasítás
IT 63	Ago ecogenico per biopsia endobronchiale ad alta definizione EchoTip ProCore® per endoscopi Olympus Istruzioni per l'uso
LT 68	„EchoTip ProCore®“ endobronchinio didelės skiriamosios gebos ultragarso biopsijos adata, skirta naudoti su „Olympus“ endoskopais Naudojimo instrukcija
LV 72	EchoTip ProCore® endobronhiālā augstas precizitātes ultraskaņas biopsijas adata Olympus endoskopiem Lietošanas instrukcija
NL 76	EchoTip ProCore® endobronchiale HD ultrasone biopsienaald voor Olympus-scopeen Gebruiksaanwijzing
NO 81	EchoTip ProCore® endobronkial HD-ultralydbiopsinål for Olympus-skop Bruksanvisning
PL 85	Igła biopsyjna EchoTip ProCore® do ultrasonografii wewnątrzoskrzelowej o wysokiej rozdzielczości do endoskopów Olympus Instrukcja użycia

PT
90

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


Fig. 7

ECHOTIP PROCORE® ENDOBRONCHIAL HD ULTRASOUND BIOPSY NEEDLE FOR OLYMPUS SCOPES

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The EchoTip ProCore® Endobronchial HD Ultrasound Biopsy Needle for Olympus Scopes is available in two different needle sizes (22 and 25 gage) and contains the following components:

- Needle
- Syringe
- Accessory channel adapter

Performance Characteristics

The function and key features of these components are described below:

Endobronchial needle

- Handle – The handle contains adjustable components to allow the user to adjust the extension of the needle and the sheath. The handle safety ring slides and locks the needle at the desired needle length extension (0–5 cm). The '0 cm' reference mark ensures complete needle retraction within the sheath. The sliding sheath adjuster allows the user to set the required sheath length extension (0–3 cm) for their Olympus endobronchial ultrasound (EBUS) scope. The '0 cm' reference mark denotes the minimum sheath length available. The base of the sliding sheath adjuster has a Luer lock fitting, which facilitates the attachment of the device to the scope accessory channel using the accessory channel adapter.
- 4.1 Fr Sheath – The sheath covers the needle during introduction, advancement, and withdrawal of the device from the scope accessory channel. It protects the scope accessory channel from damage by the needle tip.
- Needle – The stainless steel needle has a dimpling pattern to provide visualization of the needle tip under endobronchial ultrasound. The needle has a beveled tip that allows puncture and a core trap that facilitates biopsy from the target site.
- Stylet – The stylet is inserted through the length of the needle. The distal end of the stylet has a ball tip that is recessed inside the lumen of the needle (22 gage only) or it is beveled to match the bevel of the needle tip (25 gage only). The stylet provides support to the needle during advancement into the target site, prevents sample collection outside of the target site, and may also be used for sample retrieval following biopsy if desired.

Syringe – The device is supplied with a 10 mL syringe to utilize the standard vacuum technique for biopsy, if desired. The Luer lock fitting at the tip of the syringe can be attached to the Luer fitting on the handle. The syringe has a plunger and a stopcock to create a vacuum for biopsy. It is supplied with 5 mL and 10 mL clips attached to the plunger to allow the user to select the preferred volume for biopsy.

Accessory channel adapter – The accessory channel adapter consists of a Luer fitting, a biopsy port seal, and a locking clip. The accessory channel adapter attaches to the scope accessory channel port of an Olympus EBUS scope to function as a Luer fitting, to which the needle can be connected.

Device Compatibility

The EchoTip ProCore Endobronchial HD Ultrasound Biopsy Needle for Olympus Scopes is compatible with:

- Olympus EBUS scope with 2 mm minimum accessory channel
- Accessory channel adapter
- 10 mL Standard Luer syringe

Patient Population

Adult patients requiring fine needle biopsy (FNB) of submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree or the gastrointestinal tract.

Intended User

Use of this device is restricted to a trained healthcare professional.

Contact with Body Tissue

This device is tissue contacting in line with intended use.

Operating Principles

Olympus EBUS scopes do not have a Luer fitting to allow needle attachment directly to the accessory channel port. The EchoTip ProCore Endobronchial HD Ultrasound Biopsy Needle for Olympus Scopes is supplied with an accessory channel adapter to function as a Luer fitting, to which the needle can be connected. The device is inserted into the scope accessory channel port and attached when the fittings meet.

The sheath length can be adjusted depending on the length of the EBUS scope using the sliding sheath adjuster. The needle length extension can be adjusted as required. The device is used to puncture the target site to obtain a biopsy. If using a standard vacuum syringe technique, the stylet is then removed, and the syringe is attached to the Luer fitting on the needle handle. Turning the syringe stopcock to the "open" position allows negative pressure to biopsy the sample. The handle is gently moved in small increments back and forth within the target site to obtain a sample. If using the syringe, the stopcock is returned to the "closed" position. The needle is retracted into the sheath and the device is disconnected from the scope accessory channel. Upon completion of the procedure the adapter is removed from the accessory channel port.

INTENDED USE

This device is used with an ultrasound bronchoscope for fine needle biopsy (FNB) of submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree or gastrointestinal tract.

INDICATIONS FOR USE

This device is indicated for endobronchial ultrasound guided tissue sampling of submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree or the gastrointestinal tract, which is used for subsequent pathological examination to aid in disease diagnosis and patient management.

CLINICAL BENEFITS

This device facilitates successful sampling of submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree or gastrointestinal tract to aid in disease diagnosis and patient management.

CONTRAINDICATIONS

Those specific to the primary endobronchial procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications: coagulopathy.

WARNINGS

- Not for use in the heart or vascular system.
- The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.
- This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.
- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect the device with particular attention to kinks, bends, and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.
- This device contains nickel, which may cause allergic reaction in individuals with nickel sensitivity.

PRECAUTIONS

- Ensure the stylet is fully inserted when advancing the needle into the target site.
- The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement, or withdrawal of the device. Failure to retract the needle may result in damage to the scope.
- Do not use this device for any purpose other than stated intended use.
- Use of this device is restricted to a trained healthcare professional.
- The product is intended for use by physicians trained and experienced in endobronchial ultrasound techniques.
- When targeting multiple sites, adequate care should be taken to minimize the risk of tumor seeding.
- This device is intended for use with an Olympus EBUS scope.
- Store the device in a dry location.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with an endobronchial ultrasound procedure: allergic reaction to medication • aspiration • cardiac arrhythmia or arrest • fever • hemorrhage • hypotension • infection • pain/discomfort • perforation • respiratory depression or arrest.

Those associated with the device: allergic reaction to nickel • damage to blood vessels • death (other than due to normal disease progression) • fever • hemorrhage • infection • inflammation • nerve damage • pain/discomfort • perforation • pneumoperitoneum • sepsis • septicemia/bacteremia • tumor seeding.

HOW SUPPLIED

These devices are supplied ethylene oxide (EO) sterilized in a tray with a Tyvek lid.

PLACEMENT OF ACCESSORY CHANNEL ADAPTER

Illustrations

1. Examine the adapter and note features (**Fig. 1**).
2. Place the adapter, in the open position, in line with the scope accessory channel metal hub.
3. Slide the adapter on to the accessory channel ensuring that the base of the adapter slides under the accessory channel metal hub (**Fig. 2**).
4. Close the adapter ensuring it is locked (**Fig. 3**).
5. The adapter will now function as a Luer fitting to which the needle can be connected. **Note:** If irrigation of the accessory channel is desired/necessary, irrigate through the closed adapter prior to insertion or after removal of the needle.

SYSTEM PREPARATION

1. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
2. Visually inspect the device components confirming they are free from damage prior to use.
3. Examine the syringe. It has two plunger locks that must be depressed to advance the plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when the stopcock is in the "open" position, aligned with the syringe (**Fig. 4**).
4. Prepare the syringe as follows:
 - a. With the stopcock in the "open" position, depress the plunger locks and fully advance the plunger into the syringe.
 - b. Turn the stopcock to the "closed" position.
 - c. Pull the plunger back until it is locked into place at the desired setting, creating a vacuum.
 - d. Set the prepared syringe aside until biopsy is required.

INSTRUCTIONS FOR USE

1. Identify the desired target site by endobronchial ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory

- channel of the scope (**Fig. 5**). **Caution:** If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
 4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
 5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the accessory channel of the scope. To adjust the length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (**Fig. 6**). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
 6. While maintaining position of the scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring (**Fig. 7**). Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
 7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
 8. Standard vacuum syringe techniques may be applied for biopsy (Steps 9–11) or, if desired, other techniques that may or may not incorporate the use of the stylet may be used.
 9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional biopsy is required later.
 10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer fitting on the needle handle.
 11. Turn the stopcock to the “open” position aligned with the syringe, allowing the negative pressure in the syringe to facilitate biopsy. Move the needle handle back and forth within the biopsy site. **Note:** Do not remove the needle from the biopsy site during biopsy.
 12. Upon completion of biopsy, if using the syringe, turn the stopcock to the “closed” position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
 13. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
 14. Unlock the thumbscrew on the safety ring and extend the needle.
 15. Use a new syringe or the stylet to expel the specimen, then prepare per institutional guidelines.
 16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.
 17. For additional biopsy from the same target site, gently reinsert the stylet into the Luer fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.
 18. Additional samples may be obtained by repeating Steps 2 through 16 of the “Instructions for Use.”
 19. Upon completion of the procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
 20. To remove the adapter from the scope, pull back on the locking clip, lift the top, and slide off the scope.

DEVICE DISPOSAL

This device may be contaminated with potentially infectious substances of human origin and should be coiled for disposal in accordance with institutional guidelines.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2025
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.