COOK®
MEDICAL

CE
0123

EN
2
BG
8
CS
14
DA
20
NL
26
ET
32
FR
38
DE
44
EL
50
HU
56
IT
62
LV
68
LT
74
NO
80
PL
86
PT
92
RO
98
SK
104
ES
110
SV
116

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ZILVER® PTX® DRUG-ELUTING PERIPHERAL STENT

DEVICE DESCRIPTION

The Zilver PTX Drug-Eluting Peripheral Stent is a self-expanding stent made of nitinol and coated with the drug paclitaxel. The stent is intended as a permanent implant. It is a flexible, slotted tube that is designed to provide support while maintaining flexibility in the vessel upon deployment. Post-deployment, the stent is designed to impart an outward radial force upon the inner lumen of the vessel, establishing patency in the stented region.

The stent is preloaded in a 2.1mm (6 French) delivery sheath. Hand-loading of the stent is not possible. Stent deployment is controlled by rotating the thumbwheel on the handle.

Table 1 presents the stent sizes and the nominal total quantity of paclitaxel on each stent based on the established dose density of 3 µg/mm².

Table 1. Stent sizes and paclitaxel total quantity

Stent size (diameter x length, mm)	6 Fr
	Total paclitaxel (µg/stent)
5 x 40	383
6 x 40	383
7 x 40	383
8 x 40	356
5 x 60	552
6 x 60	552
7 x 60	552
8 x 60	579
5 x 80	722
6 x 80	722
7 x 80	722
8 x 80	757
5 x 100	934
6 x 100	934
7 x 100	934
8 x 100	935
5 x 120	1103
6 x 120	1103
7 x 120	1103
8 x 120	1112

MRI INFORMATION

Non-clinical testing has demonstrated that the Zilver PTX Drug-Eluting Peripheral Stent is MR conditional for overlapping lengths up to 395 mm. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla or 3.0 Tesla
- Maximum spatial field gradient of 1600 Gauss/cm (16 T/m)
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode).

Under the scan conditions defined above, the Zilver PTX Drug-Eluting Peripheral Stent is expected to produce a maximum temperature rise of less than 4.9°C after 15 minutes of continuous scanning. In non-clinical testing, the image artifact caused by the device extends approximately 5 mm from the Zilver PTX Drug-Eluting Peripheral Stent when imaged with a gradient echo pulse sequence and a 3.0 Tesla MR system.

It is recommended that patients register the conditions under which the implant can safely be scanned with the MedicAlert Foundation (www.medicalert.org) or an equivalent organization.

INTENDED USE

The Zilver PTX Drug-Eluting Peripheral Stent is intended for use in the treatment of symptomatic vascular disease of the above-the-knee femoropopliteal arteries having reference vessel diameter from 4 mm to 7 mm. To avoid involvement of the common femoral artery, the proximal end of the stent should be placed at least 1 cm below the origin of the superficial femoral artery. To avoid involvement of the below-the-knee popliteal artery, the distal end of the stent should be placed above the plane of the femoral epicondyles.

CONTRAINDICATIONS

- Stenosis that cannot be dilated to permit passage of the introducer sheath.
- Stenting of an arterial vessel where leakage from the artery could be exacerbated by placement of a stent.
- Patients with bleeding disorders.
- Women who are pregnant, breastfeeding, or plan to become pregnant in the next 5 years should not receive a Zilver PTX Drug-Eluting Peripheral Stent. It is unknown whether paclitaxel will be excreted in human milk, and there is a potential for adverse reaction in nursing infants from paclitaxel exposure.

WARNINGS

- **A signal for increased risk of late mortality has been identified following the use of paclitaxel-coated balloons and paclitaxel-eluting stents for femoropopliteal arterial disease beginning approximately 2-3 years post-treatment compared with the use of non-drug coated devices. There is uncertainty regarding the magnitude and mechanism for the increased late mortality risk, including the impact of repeat paclitaxel-coated device exposure. Physicians should discuss this late mortality signal and the benefits and risks of available treatment options with their patients.**
- Persons with allergic reactions to nickel may suffer an allergic reaction to this implant.
- Persons allergic to paclitaxel may suffer an allergic reaction to this implant.
- For adult use only.
- Use of this device is restricted to a trained healthcare professional.

- This device is designed for single use only. Attempts to reprocess, re-sterilize and/or reuse may lead to device failure and/or transmission of disease.
- Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Inspect the product to ensure no damage has occurred.
- Flow restrictions remaining after stent deployment (e.g., residual proximal or distal stenosis or dissection, or poor distal outflow) may increase the risk of stent thrombosis. Inflow and outflow should be assessed at procedure completion and additional measures considered (e.g., additional PTA, adjunctive stenting, or distal bypass) if necessary to maintain good inflow and outflow.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional vascular techniques. Standard techniques for interventional vascular procedures should be employed.
- Manipulation of the Zilver PTX Drug-Eluting Peripheral Stent requires fluoroscopic control.
- If resistance is met during advancement of the delivery system, do not force passage. Remove the delivery system and replace with a new device. Continuing to force passage may ultimately cause partial deployment of the stent. Possible outcomes of partial deployment are
 - Prolonged procedure (e.g. additional exposure to radiation / contrast)
 - Additional stent required
 - Secondary intervention (e.g. vessel incision)
- Do not try to remove the stent from the introducer system before use.
- Ensure that the red safety lock is not inadvertently depressed before stent deployment is desired.
- A 0.89 mm (0.035 inch) wire guide should be used during tracking, deployment, and removal to ensure adequate support of the system. If hydrophilic wire guides are used, they must be kept fully activated. Using a smaller diameter wire guide may ultimately result in partial deployment of the stent. Possible outcomes of partial deployment are
 - Prolonged procedure (e.g. additional exposure to radiation / contrast)
 - Additional stent required
 - Secondary intervention (e.g. vessel incision)
- Prior to disengaging the device safety lock ensure the distal end of the stability sheath is inside the introducer sheath. Failure to do so may result in stent damage and/or stent compression upon deployment.
- Do not use excessive force to deploy the stent. If excessive resistance is felt when beginning deployment, remove the delivery system without deploying the stent and replace with a new device.
- Do not expose the delivery system to organic solvents (e.g., alcohol).
- Do not use power injection systems with the delivery system.
- Do not torque the delivery system during introduction or deployment. Torquing the device may lead to difficulty or inability to deploy and/or move device over the wire guide.
- Once the stent has apposed the vessel wall repositioning of the stent is not recommended. The stent retraction sheath cannot be re-advanced nor can the stent be recaptured following the start of stent deployment.
- Following stent deployment, if resistance is met during the withdrawal of the delivery system, carefully remove the delivery system and wire guide together as a unit.
- If resistance is still encountered during removal of the delivery system and wire guide as a unit, remove the wire guide, delivery system and introducer sheath together as a unit.

POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur include the following:

- Allergic reaction to anticoagulant and/or antithrombotic therapy or contrast medium
- Allergic reaction to nickel
- Atheroembolization (Blue Toe Syndrome)
- Arterial aneurysm
- Arterial rupture
- Arterial thrombosis
- Arteriovenous fistula
- Death
- Dissection
- Embolism
- Fever
- Hematoma/hemorrhage
- Hypersensitivity reactions
- Infection
- Infection/abscess formation at access site
- Ischemia requiring intervention (bypass or amputation of toe, foot or leg)
- Occlusion
- Pain/discomfort
- Pseudoaneurysm formation
- Renal failure
- Restenosis of the stented artery
- Stent embolization
- Stent malapposition
- Stent migration
- Stent strut fracture
- Vessel perforation or rupture
- Vessel spasm
- Worsened claudication/rest pain

Although systemic effects are not anticipated, refer to the Physicians' Desk Reference for more information on the potential adverse events observed with paclitaxel. Potential adverse events, not described in the above source, may be unique to the paclitaxel drug coating:

- Allergic/immunologic reaction to the drug coating
- Alopecia
- Anemia
- Blood product transfusion
- Gastrointestinal symptoms

- Hematologic dyscrasia (including leukopenia, neutropenia, thrombocytopenia)
- Hepatic enzyme changes
- Histologic changes in vessel wall, including inflammation, cellular damage, or necrosis
- Myalgia/Arthralgia
- Myelosuppression
- Peripheral neuropathy

LATE MORTALITY SIGNAL FOR PACLITAXEL-COATED DEVICES

A meta-analysis of randomized controlled trials published in December 2018 by Katsanos et al.¹ identified an increased risk of late mortality at 2 years and beyond for paclitaxel-coated balloons and paclitaxel-eluting stents used to treat femoropopliteal arterial disease.

In response to these data, the FDA performed a patient-level meta-analysis of long-term follow-up data from the pivotal premarket randomized trials of paclitaxel-coated devices used to treat femoropopliteal disease using available clinical data through May 2019. The meta-analysis also showed a late mortality signal in study subjects treated with paclitaxel-coated devices compared to patients treated with uncoated devices. Specifically, in the 3 randomized trials with a total of 1090 patients and available 5-year data, the crude mortality rate was 19.8% (range 15.9%–23.4%) in patients treated with paclitaxel-coated devices compared to 12.7% (range 11.2%–14.0%) in subjects treated with uncoated devices. The relative risk for increased mortality at 5 years was 1.57 (95% confidence interval 1.16–2.13), which corresponds to a 57% relative increase in mortality in patients treated with paclitaxel-coated devices.

As presented at the June 2019 FDA Advisory Committee Meeting, an independent meta-analysis of similar patient-level data provided by VIVA Physicians, a vascular medicine organization, reported similar findings with a hazard ratio of 1.38 (95% confidence interval 1.06–1.80). Additional analyses have been conducted and are underway that are specifically designed to assess the relationship of mortality to paclitaxel-coated devices. The presence and magnitude of the late mortality risk should be interpreted with caution because of multiple limitations in the available data, including wide confidence intervals due to a small sample size, pooling of studies of different paclitaxel-coated devices that were not intended to be combined, substantial amounts of missing study data, no clear evidence of a paclitaxel dose effect on mortality, and no identified pathophysiologic mechanism for the late deaths.

Paclitaxel-coated balloons and stents improve blood flow to the legs and decrease the likelihood of repeat procedures to reopen blocked blood vessels compared to uncoated devices. The benefits of paclitaxel-coated devices (e.g., reduced reinterventions) should be considered in individual patients along with potential risks (e.g., late mortality).

In the Zilver PTX randomized clinical study,² Kaplan Meier mortality estimates at 2, 3 and 5 years are 5.4% (95% confidence interval 3.0%–7.9%), 10.3% (7.0%–13.6%) and 19.1% (14.2%–24.1%), respectively, for all patients treated with the Zilver PTX treatment device and 5.7% (1.8%–9.5%), 10.7% (5.6%–15.9%) and 17.1% (10.0%–24.2%), respectively, for patients treated only with the PTA/BMS control device ($P = 0.60$). Publications of large observational studies have not identified an increased risk of mortality with the Zilver PTX DES; Secemsky et al.³ reported comparable mortality rates through 4 years for DES and BMS in claudicants ($P = 0.55$) and CLI patients ($P = 0.52$); Freisinger et al.⁴ reported a hazard ratio of 1.01 (95% confidence interval 0.83–1.23) at 5 years for DES vs. BMS; Katsuki et al.⁵ reported a hazard ratio of 0.89 (95% confidence interval 0.66–1.19) at 5 years for DES vs. non-paclitaxel devices.

DRUG INTERACTIONS

Formal drug interaction studies have not been conducted with the Zilver PTX Drug-Eluting Peripheral Stent. In the absence of formal clinical drug interaction studies, caution should be exercised when administering paclitaxel concomitantly with known substrates or inhibitors of the cytochrome P450 isoenzymes CYP2C8 and CYP3A4.

PRODUCT RECOMMENDATIONS

Placement of this vascular stent requires advanced skills in interventional vascular procedures. The following instructions will give technical guidance, but do not obviate formal training in the use of the device. Delivery system components are shown in **Figure 1**.


Fig. 1

1. Ergonomic thumbwheel operated handle
2. Safety lock
3. Thumbwheel
4. Luer hub (flushing port)
5. Retraction sheath
6. Stability sheath
7. Tip of delivery system inner catheter
8. Radiopaque marker band on retraction sheath
9. Gold radiopaque markers on stent
10. Radiopaque stent stop

Multiple Stent Placement

If placements of multiple stents are required in a patient, to cover the length of the lesion, the following recommendations should be considered:

- In relation to the lesion site, the distal area of narrowing should be stented first, followed by the proximal locations (i.e., a second stent should be placed proximally to the previously placed stent).
- Stents placed in tandem must overlap to allow for complete coverage of the lesion.
- The safety and effectiveness of implanting multiple Zilver PTX stents with a total drug quantity of greater than approximately 3 mg paclitaxel per patient, which is equivalent to the amount of paclitaxel for a total stented length of 360 mm, has not been established.

1. Katsanos K, et al. Risk of death following application of paclitaxel-coated balloons and stents in the femoropopliteal artery of the leg: a systematic review and meta-analysis of randomized controlled trials. *J Am Heart Assoc.* 2018;7:e011245.
2. Dake MD, et al. Paclitaxel-coated Zilver PTX Drug-Eluting Stent treatment does not result in increased long-term all-cause mortality compared to uncoated devices. *Cardiovasc Intervent Radiol.* 2019; <https://doi.org/10.1007/s00270-019-02324-4>
3. Secemsky EM, et al. Drug-eluting stent implantation and long-term survival following peripheral artery revascularization. *J Am Coll Cardiol.* 2019;73:2636–2638.
4. Freisinger E, et al. Mortality after use of paclitaxel-based devices in peripheral arteries: a real-world safety analysis. *Eur Heart J.* 2019; doi:10.1093/eurheartj/ehz698.
5. Katsuki T, et al. Mortality risk following application of a paclitaxel-coated stent in femoropopliteal lesions. *J Endovasc Ther.* 2019; doi:10.1177/1526602819870309

Stent Sizing

1. Determine the proper stent size after complete diagnostic evaluation. The stent deployment must be performed under fluoroscopic control. Measure the length of the target lesion to determine the length of the stent required. Allow for the proximal and distal aspects of the stent to cover the entire target area.

Note: The Zilver PTX Drug-Eluting Peripheral Stent is designed not to shorten upon deployment.

The stent is recommended for use in above-the-knee femoropopliteal arteries having reference vessel diameter from 4 mm to 7 mm.


Measure the diameter of the reference vessel (proximal and distal to the lesion) and use the LARGEST reference diameter as your basis for choosing the appropriate stent size.

The stent size should be selected so that the unconstrained stent diameter is at least 1 mm larger than the reference vessel diameter and no more than 2 mm larger than the reference vessel diameter.

Introduction of the Stent

1. Gain access to the site using a 6 French sheath with a minimum I.D. of 2.1 mm.
2. To ensure adequate support of the system, introduce a 0.89 mm (0.035 inch) wire guide. If hydrophilic wire guides are used, they must be kept fully activated.
3. Pre-dilatation before stent placement is optional and at the discretion of the physician.
4. Prior to advancing the system over the wire and introducing the system into the body use a syringe to flush the device wire guide lumen and stent with saline or heparinized saline through the hub of the device. Flushing the device with contrast media is not recommended.

Flush the wire lumen until a few drops of saline exit the tip of the device (**Figure 2**). Next, gently place a finger over the device end hole and flush again until a few drops of saline can be seen exiting between the end of the sheath and device tip (**Figure 3**).

5. Under fluoroscopy, advance the delivery system over the 0.89 mm (0.035 inch) wire guide through the introducer sheath until the distal gold radiopaque markers on the stent are beyond the target lesion site (**Figure 4**).

Note: If resistance is met during advancement of the delivery system do not force passage. Remove the delivery system and replace with a new device.

Note: Ensure the distal end of the stability sheath is inside the introducer sheath.

Caution: Torquing the device may lead to difficulty or inability to deploy and/or move device over the wire guide.


Fig. 4

Alignment and Deployment of the Stent

1. Under fluoroscopy, pull the delivery system back until there is no slack in the delivery system and the radiopaque markers on the stent are aligned with the lesion (**Figure 5**).


Fig. 5

2. Disengage the device safety lock by gently depressing the red safety button in the direction indicated in **Figure 6**.


Fig. 6

3. Initiate stent deployment by slowly rotating the device thumbwheel in the direction shown in **Figure 7**.

Note: If high resistance is felt on the thumbwheel prior to stent deployment, do not force deployment. Carefully withdraw the stent system without deploying the stent.


Fig. 7

4. Rotating the thumbwheel will withdraw the retraction sheath from over the stent.

Note: Visually confirm that the radiopaque stent markers are still aligned with the desired position. Reposition if necessary.

Note: Stent deployment will begin once the distal radiopaque stent markers have advanced beyond the sheath marker band (Refer to **Figure 1** for sheath and stent marker identification). Once the stent has apposed the vessel wall repositioning of the stent is not recommended. Recapture or re-sheathing of the Zilver PTX Drug-Eluting Peripheral Stent is not possible because the stent retraction sheath cannot be re-advanced. Refer to the **Multiple Stent Placement** section of these Instructions for Use for information on missed lesions.

5. Full deployment of the stent length will occur when the distal end of the retraction sheath has been retracted past the proximal radiopaque markers of the stent (**Figure 8**).


Fig. 8

Post Stent Deployment

1. Under fluoroscopy, withdraw the entire delivery system as one unit, over the wire guide, into the introducer sheath and out of the body. Remove the delivery system from the wire guide.

Note: Following stent deployment, if resistance is met during the withdrawal of the delivery system, carefully remove the delivery system and wire guide together as a unit.

If resistance is still encountered during removal of the delivery system and wire guide as a unit, remove the wire guide, delivery system and introducer sheath together as a unit.

2. Perform an arterial angiogram to verify full deployment of the device. If incomplete expansion exists within the stent at any point along the lesion, post-deployment balloon dilatation (standard PTA) can be performed at the discretion of the physician.

Note: Flow restrictions remaining after stent deployment (e.g., residual proximal or distal stenosis or dissection, or poor distal outflow) may increase the risk of stent thrombosis. Inflow and outflow should be assessed at procedure completion and additional measures considered (e.g., additional PTA, adjunctive stenting, or distal bypass) if necessary to maintain good inflow and outflow.

3. Remove the wire guide and introducer sheath from the patient.
4. Close the entry wound as appropriate.
5. Dispose in accordance with institutional guidelines

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in an outer non-sterile foil pouch and an inner peel-open package. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Do not store above 25 °C (77°F), excursions permitted to 15°C to 30°C (59°F to 86°F). Protect from moisture. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

STERILE - DO NOT RESTERILIZE - SINGLE USE ONLY

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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1. Katsanos K, et al. Risk of death following application of paclitaxel-coated balloons and stents in the femoropopliteal artery of the leg: a systematic review and meta-analysis of randomized controlled trials. *J Am Heart Assoc.* 2018;7:e011245.
2. Dake MD, et al. Paclitaxel-coated Zilver PTX Drug-Eluting Stent treatment does not result in increased long-term all-cause mortality compared to uncoated devices. *Cardiovasc Intervent Radiol.* 2019; <https://doi.org/10.1007/s00270-019-02324-4>
3. Secemsky EM, et al. Drug-eluting stent implantation and long-term survival following peripheral artery revascularization. *J Am Coll Cardiol.* 2019;73:2636–2638.
4. Freisinger E, et al. Mortality after use of paclitaxel-based devices in peripheral arteries: a real-world safety analysis. *Eur Heart J.* 2019; [doi:10.1093/eurheartj/ehz698](https://doi.org/10.1093/eurheartj/ehz698).
5. Katsuki T, et al. Mortality risk following application of a paclitaxel-coated stent in femoropopliteal lesions. *J Endovasc Ther.* 2019; [doi:10.1177/1526602819870309](https://doi.org/10.1177/1526602819870309)

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1. Katsanos K, et al. Risk of death following application of paclitaxel-coated balloons and stents in the femoropopliteal artery of the leg: a systematic review and meta-analysis of randomized controlled trials. *J Am Heart Assoc*. 2018;7:e011245.
2. Dake MD, et al. Paclitaxel-coated Zilver PTX Drug-Eluting Stent treatment does not result in increased long-term all-cause mortality compared to uncoated devices. *Cardiovasc Intervent Radiol*. 2019; <https://doi.org/10.1007/s00270-019-02324-4>
3. Secemsky EM, et al. Drug-eluting stent implantation and long-term survival following peripheral artery revascularization. *J Am Coll Cardiol*. 2019;73:2636–2638.
4. Freisinger E, et al. Mortality after use of paclitaxel-based devices in peripheral arteries: a real-world safety analysis. *Eur Heart J*. 2019; doi:10.1093/eurheartj/ehz698.
5. Katsuki T, et al. Mortality risk following application of a paclitaxel-coated stent in femoropopliteal lesions. *J Endovasc Ther*. 2019; doi:10.1177/1526602819870309

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1. Katsanos K, et al. Risk of death following application of paclitaxel-coated balloons and stents in the femoropopliteal artery of the leg: a systematic review and meta-analysis of randomized controlled trials. *J Am Heart Assoc.* 2018;7:e011245.
2. Dake MD, et al. Paclitaxel-coated Zilver PTX Drug-Eluting Stent treatment does not result in increased long-term all-cause mortality compared to uncoated devices. *Cardiovasc Intervent Radiol.* 2019; <https://doi.org/10.1007/s00270-019-02324-4>
3. Secemsky EM, et al. Drug-eluting stent implantation and long-term survival following peripheral artery revascularization. *J Am Coll Cardiol.* 2019;73:2636-2638.
4. Freisinger E, et al. Mortality after use of paclitaxel-based devices in peripheral arteries: a real-world safety analysis. *Eur Heart J.* 2019; doi:10.1093/eurheartj/ehz698.
5. Katsuki T, et al. Mortality risk following application of a paclitaxel-coated stent in femoropopliteal lesions. *J Endovasc Ther.* 2019; doi:10.1177/1526602819870309

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COOK IRELAND LTD.