# COOK®

---

MEDICAL

# CE

0123

EN  
7

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## ZILVER® PTX® DRUG-ELUTING PERIPHERAL STENT

### DEVICE DESCRIPTION

The Zilver PTX Drug-Eluting Peripheral Stent is a self-expanding stent made of nickel titanium alloy (nitinol) and coated with the drug paclitaxel. Gold markers at each end of the stent ensure visibility under fluoroscopy. The stent is a flexible, slotted tube which is preloaded in a 2.1 mm (6 French) delivery system. Stent deployment is controlled by rotating the thumbwheel on the handle of the delivery system.

**Table 1: Zilver PTX Drug Eluting Peripheral Stent and Delivery System features overview**

Zilver PTX Drug-Eluting Peripheral Stent	
Available stent lengths	40, 60, 80, 100, 120 mm
Available stent diameters	5, 6, 7, 8 mm
Stent material	Nickel Titanium alloy (Nitinol) with Gold Markers
Drug coating	Paclitaxel
Delivery System	
Available delivery system lengths	80, 125 cm
Wire guide compatibility	0.89 mm (0.035 inch)
Delivery System Outer Diameter	2.1 mm (6 Fr)
Introducer sheath	6 Fr (minimum ID of 2.1 mm)

Zilver PTX Drug-Eluting Peripheral Stents are coated with paclitaxel API (active pharmaceutical ingredient) using a proprietary process. No excipients, polymers, carriers, binding agents, other materials, or other device modifications are involved.

Paclitaxel may be extracted from the bark, branches, needles of the yew tree or paclitaxel may be produced via plant cell fermentation using a dedicated plant cell strain. In both cases, paclitaxel is purified and concentrated by column chromatography, crystallization and recrystallization.

### Performance Characteristics

The function and key features of the device are described below.

**Delivery System** - the function of the delivery system is to advance the stent through a 6 Fr introducer sheath, over a pre-positioned 0.035 inch wire guide to its intended location. It maintains the position of the stent during deployment.

The delivery system has the following features (Fig.1):

- Delivery system lengths of 80 cm or 125 cm for approach to the target location.
- Tip of Delivery System
  - Atraumatic - the tip is designed to be atraumatic to the anatomy.
  - Flexible - the tip is designed to be flexible during introduction and withdrawal from the body.
- Radiopacity - the distal marker band enables fluoroscopic visibility of delivery system.
- Retraction Sheath - the stent is deployed by controlled withdrawal of the retraction sheath using the thumbwheel. The stent stop holds stent position as the sheath withdraws proximally to reveal the stent.
- Stability sheath - provides flexible support and protection for the retraction sheath and the inner catheter during device introduction and withdrawal.
- Thumbwheel - enables single-handed, step-by-step deployment of the stent. It does not enable recapture of the stent.
- Safety lock - prevents thumbwheel rotation prior to intended stent deployment. It must be disengaged to initiate deployment at the target location.
- Luer Hub - enables attachment of a syringe to flush the wire guide lumen and stent location inside the delivery system, prior to introduction to the body. The hub allows the wire guide to exit proximally during introduction.

**Stent** - the stent is designed to impart an outward radial force upon the inner lumen of the vessel, establishing patency in the stented region. It is coated in the drug paclitaxel. Table 3 presents the stent sizes and the nominal total quantity of paclitaxel on each stent based on the established dose density of 3 µg/mm<sup>2</sup>.

The Zilver PTX stent has the following features:

- The stent is manufactured from nickel titanium alloy (nitinol) enabling it to self-expand at body temperature upon deployment.
- Radiopacity - gold markers at each end of the stent provide visibility under fluoroscopy.
- Paclitaxel (PTX) drug coating on the stent reduces restenosis by inhibiting smooth muscle cell proliferation and neointimal hyperplasia post deployment.
- The stent is a flexible slotted tube that is designed to provide support while maintaining flexibility in the vessel upon deployment.
- The stent is intended for permanent implant.

### Device Compatibility

- 0.89 mm (0.035 inch) wire guide
- 6 French introducer sheath with a minimum ID of 2.1 mm
- Syringe
- Dilatation balloons. For pre- and post- dilation, an appropriately sized balloon catheter is recommended
- Saline or Heparinized Saline
- Contrast media\*

\* Note: The outer part of the delivery system is compatible with contrast media. The device is not to be flushed using contrast media. (refer to Precautions section)

### Qualitative and Quantitative Information

Each stent is composed of nickel titanium alloy (nitinol), gold markers and paclitaxel drug.

**Table 2: Qualitative and quantitative information**

Qualitative Information Device Material		Quantitative Information Weight (g)
Nickel Titanium alloy (Nitinol)	Stent	Up to 0.50
Gold	Gold Markers	Up to 0.0080
Paclitaxel	Drug	See below

Table 3 presents the stent sizes and the nominal total quantity of paclitaxel on each stent based on the established dose density of 3 µg/mm<sup>2</sup>.

**Table 3: Stent sizes and paclitaxel total quantity**

Stent size (diameter x length, mm)	6 Fr
	Total paclitaxel (µg/stent)
5 x 40	383
6 x 40	383
7 x 40	383
8 x 40	356
5 x 60	552
6 x 60	552
7 x 60	552
8 x 60	579
5 x 80	722
6 x 80	722
7 x 80	722
8 x 80	757
5 x 100	934
6 x 100	934
7 x 100	934
8 x 100	935
5 x 120	1103
6 x 120	1103
7 x 120	1103
8 x 120	1112

### Patient Population

Adult patients with symptomatic vascular disease of the above-the-knee femoropopliteal arteries.

### Contact with Body Tissue

This device has contact with body tissue and blood, in line with its intended use.

## Operating Principles

The Zilver PTX device is placed over a prepositioned wire guide and advanced through an introducer sheath, to the target location. Radiopaque marker bands at the distal end of the delivery system provide visibility under fluoroscopy during introduction and deployment of the stent.

Hand-loading of the stent is not possible. Once in position, stent deployment is controlled by rotating the thumbwheel on the handle of the delivery system. Post deployment, the stent is designed to impart an outward radial force upon the inner lumen of the vessel, establishing patency in the stented region. The stent is designed to provide support as well as maintain flexibility in the vessel upon deployment.

Paclitaxel (PTX) drug coating reduces restenosis by inhibiting smooth muscle cell proliferation and neointimal hyperplasia post deployment.

## INTENDED USE

The Zilver® PTX® Drug-Eluting Peripheral Stent is intended to improve luminal patency by imparting an outward radial force on the inner wall of the above-the-knee femoropopliteal arteries.

## INDICATIONS FOR USE

The Zilver® PTX® Drug-Eluting Peripheral Stent is indicated for the treatment of *de novo* or restenotic symptomatic lesions in native vascular disease of the above-the-knee femoropopliteal arteries having reference vessel diameter from 4 mm to 7 mm.

To avoid involvement of the common femoral artery, the proximal end of the stent should be placed at least 1 cm below the origin of the superficial femoral artery. To avoid involvement of the below-the-knee popliteal artery, the distal end of the stent should be placed above the plane of the femoral epicondyles.

## CLINICAL BENEFITS

- Improves luminal patency (blood flow) through the above-the-knee femoropopliteal arteries
- Reduces reinterventions
- Improves patient's quality of life

## CONTRAINDICATIONS

- Stenosis that cannot be dilated to permit passage of the introducer sheath.
- Stenting of an arterial vessel where leakage from the artery could be exacerbated by placement of a stent.
- Patients who cannot receive appropriate antiplatelet and/or anti-coagulant therapy.

## WARNINGS

- **A signal for increased risk of late mortality has been identified following the use of paclitaxel-coated balloons and paclitaxel-eluting stents for femoropopliteal arterial disease beginning approximately 2-3 years post-treatment compared with the use of non-drug coated devices. There is uncertainty regarding the magnitude and mechanism for the increased late mortality risk, including the impact of repeat paclitaxel-coated device exposure. Physicians should discuss this late mortality signal and the benefits and risks of available treatment options with their patients.**
- Persons with allergic reactions to nickel titanium alloy (nitinol) may suffer an allergic reaction to this implant.
- Persons allergic to paclitaxel or structurally-related compounds may suffer an allergic reaction to this implant.
- Do not implant stents with a total drug quantity greater than 3 mg of Paclitaxel. The safety and effectiveness of implanting multiple Zilver PTX Drug-Eluting Peripheral Stents with a total maximum drug quantity of up to 3 mg paclitaxel (i.e., a total stented length of 320 mm) has been established in clinical studies.<sup>1,2,3,4</sup> A total drug dose of 3 mg paclitaxel is also not associated with genotoxic or reproductive toxicological effects; although adverse genotoxic or reproductive toxicological effects have been reported for much higher paclitaxel doses more than 20 times the maximum allowed in treatment with the Zilver PTX Drug-Eluting Stent.<sup>5,6</sup> The safety margin for Zilver PTX Drug-Eluting Stent was further verified in animal studies which demonstrated no systemic or regional effects after placement of Zilver PTX Drug-Eluting Stents with paclitaxel dose densities up to four times the nominal level (i.e., 12 µg/mm<sup>2</sup> instead of the nominal 3 µg/mm<sup>2</sup>), or after placement of Zilver PTX Drug-Eluting stents with more than three times the total maximum clinical dose (i.e., 10.5 mg instead of the clinical maximum of 3 mg).
- This single use device is not designed for re-use. Attempts to reprocess, re-sterilize, and/or to re-use may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.
- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect the device with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for a return authorization.
- Flow restrictions remaining after stent deployment (e.g., residual proximal or distal stenosis or dissection, or poor distal outflow) may increase the risk of stent thrombosis. Inflow and outflow should be assessed at procedure completion and additional measures considered (e.g., additional PTA, adjunctive stenting, or distal bypass) if necessary to maintain good inflow and outflow.
- The stent is intended for permanent implant.
- It is not recommended that women who are pregnant or intending to become pregnant in the next 12 months, or men intending to father children in the next 9 months receive a Zilver PTX Drug-Eluting Peripheral Stent as the reproductive and genotoxic effects on human subjects are unknown. Breastfeeding should be avoided for 6 months after product implantation as there is a potential for adverse reaction in nursing infants from paclitaxel exposure.

## PRECAUTIONS

### Prior to use

- Use of this device is restricted to a trained healthcare professional.
- This product is intended for use by physicians trained and experienced in diagnostic and interventional vascular techniques. Standard techniques for interventional vascular procedures should be employed.
- Device not indicated for use in the treatment of in-stent restenosis.
- Do not store above 25°C (77°F), excursions permitted to 15°C to 30°C (59°F to 86°F). Store in a dark, dry, cool place. Avoid extended exposure to light. Device is to be protected from moisture.
- Caution should be exercised when administering paclitaxel concomitantly with known substrates or inhibitors of the cytochrome P450 isoenzymes CYP2C8 and CYP3A4.
- Antiplatelet/anticoagulant therapy should be administered during and after the procedure according to institutional standard of care.
- If multiple stents are to be placed in an overlapping fashion, they should be of similar composition (i.e. nickel titanium alloy (nitinol)).

This symbol on the label indicates that components of the delivery system contain cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the components of the delivery system contain cobalt as a cobalt alloy (or stainless-steel alloy containing cobalt), which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

## Device handling

- Manipulation of the Zilver PTX Drug-Eluting Peripheral Stent requires fluoroscopic control.
- Do not expose the delivery system to organic solvents (e.g., alcohol).
- Do not use power injection systems with the delivery system.
- Do not flush the device with contrast media.
- Do not try to remove the stent from the introducer system before use.
- A 0.89 mm (0.035 inch) wire guide should be used during tracking, deployment, and removal to ensure adequate support of the system. If hydrophilic wire guides are used, they must be kept fully activated. Using a smaller diameter wire guide may ultimately result in partial deployment of the stent. Possible outcomes of partial deployment are
  - Prolonged procedure (e.g. additional exposure to radiation/contrast)
  - Additional stent required
  - Secondary intervention (e.g. vessel incision)
- If resistance is met during advancement of the delivery system, do not force passage. Remove the delivery system and replace with a new device. Continuing to force passage may ultimately cause partial deployment of the stent. Possible outcomes of partial deployment are
  - Prolonged procedure (e.g. additional exposure to radiation/contrast)
  - Additional stent required
  - Secondary intervention (e.g. vessel incision)
- Do not torque the delivery system during introduction or deployment. Torquing the device may lead to difficulty or inability to deploy and/or move device over the wire guide.
- Ensure that the red safety lock is not inadvertently depressed before stent deployment is desired.

## Stent Placement

- Prior to disengaging the device safety lock ensure the distal end of the stability sheath is inside the introducer sheath. Failure to do so may result in stent damage and/or stent compression upon deployment.
- If high resistance is felt when beginning deployment, do not force deployment. Carefully withdraw the delivery system without deploying the stent and replace with a new device.
- Once the stent has apposed the vessel wall repositioning of the stent is not recommended. The stent retraction sheath cannot be re-advanced nor can the stent be recaptured following the start of stent deployment.

1. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Zeller T, Roubin GS, Burket MW, Khatib Y, Snyder SA, Ragheb AO, White JK, Machan LS; Zilver PTX Investigators. Paclitaxel-eluting stents show superiority to balloon angioplasty and bare metal stents in femoropopliteal disease: twelve-month Zilver PTX randomized study results. *Circ Cardiovasc Interv*. 2011;4:495-504. doi: 10.1161/CIRCINTERVENTIONS.111.962324.

2. Dake MD, Scheinert D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kavteladze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther*. 2011;18:613-623. doi: 10.1583/11-3560.1.

3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv*. 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.

4. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Machan LS, Snyder SA, O'Leary EE, Ragheb AO, Zeller T. Durable clinical effectiveness with paclitaxel-eluting stents in the femoropopliteal artery: 5-year results of the Zilver PTX randomized trial. *Circulation*. 2016;133:1472-1483. doi: 10.1161/CIRCULATIONAHA.115.016900.

5. Tinwell H and Ashby J. Genetic toxicity and potential carcinogenicity of taxol. *Carcinogenesis* 1994; 15 (8): 1499-1501.

6. Physician's Desk Reference, 52nd edition, Medical Economics Inc, 1998, pp 762-766.

## Device removal

- Following stent deployment, if resistance is met during the withdrawal of the delivery system, carefully remove the delivery system and wire guide together as a unit.
- If resistance is still encountered during removal of the delivery system and wire guide as a unit, remove the wire guide, delivery system and introducer sheath together as a unit.

## Post procedure

- Periodic evaluation of the device is recommended during the indwell period.
- Damage/dislodgement of a placed stent may occur with additional procedures.

## MRI SAFETY INFORMATION

This symbol means the device is MR conditional.

A patient with a Zilver PTX Drug-Eluting Peripheral Stent may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.

- Nominal value(s) of Static Magnetic Field: 1.5 T or 3 T
- Maximum Spatial Field Gradient: 24 T/m (2400 gauss/cm)
- RF Excitation: Circularly Polarized (CP)
- RF Transmit Coil Type: Whole body transmit coil, Head RF transmit-receive coil
- Maximum Whole Body SAR: 2.0 W/kg (Normal Operating Mode)
- Limits on Scan Duration: 2.0 W/kg whole body average SAR for 15 minutes of continuous RF (a sequence or back to back series/scan without breaks)
- MR Image Artifact: The presence of this implant may produce an image artifact of 15 mm

If information about a specific parameter is not included, there are no conditions associated with that parameter.

## POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur include the following: • allergic reaction to anticoagulant and/or antithrombotic therapy or contrast medium • allergic reaction to nickel titanium alloy (nitinol) • atheroembolization (Blue Toe Syndrome) • arterial aneurysm • arterial rupture • arterial thrombosis • arteriovenous fistula • death • dissection • embolism • fever • hematoma/hemorrhage • hypersensitivity reactions • infection • infection/abscess formation at access site • ischemia requiring intervention (bypass or amputation of toe, foot or leg) • occlusion • pain/discomfort • pseudoaneurysm formation • renal failure • restenosis of the stented artery • stent embolization • stent malapposition • stent migration • stent strut fracture • vessel perforation or rupture • vessel spasm • worsened claudication/rest pain.

Although systemic effects are not anticipated, refer to the Physicians' Desk Reference for more information on the potential adverse events observed with paclitaxel. Potential adverse events, not described in the above source, may be unique to the paclitaxel drug coating: • allergic/immunologic reaction to the drug coating • alopecia • anemia • blood product transfusion • gastrointestinal symptoms • hematologic dyscrasia (including leukopenia, neutropenia, thrombocytopenia) • hepatic enzyme changes • histologic changes in vessel wall, including inflammation, cellular damage, or necrosis • myalgia/arthralgia • myelosuppression • peripheral neuropathy.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide. The outer foil pouch is non sterile and protects the contents from light and moisture. The inner pouch is the sterile barrier and is a peel open package. Following sterilization, the pouch is further packaged into a box.

This device is accompanied by an implant card that should be given to the patient after it has been completed by the healthcare professional.

## LATE MORTALITY SIGNAL FOR PACLITAXEL-COATED DEVICES

A meta-analysis of randomized controlled trials published in December 2018 by Katsanos et al.<sup>7</sup> identified an increased risk of late mortality at 2 years and beyond for paclitaxel-coated balloons and paclitaxel-eluting stents used to treat femoropopliteal arterial disease.

In response to these data, the FDA performed a patient-level meta-analysis of long-term follow-up data from the pivotal premarket randomized trials of paclitaxel-coated devices used to treat femoropopliteal disease using available clinical data through May 2019. The meta-analysis also showed a late mortality signal in study subjects treated with paclitaxel-coated devices compared to patients treated with uncoated devices. Specifically, in the 3 randomized trials with a total of 1090 patients and available 5-year data, the crude mortality rate was 19.8% (range 15.9%–23.4%) in patients treated with paclitaxel-coated devices compared to 12.7% (range 11.2%–14.0%) in subjects treated with uncoated devices. The relative risk for increased mortality at 5 years was 1.57 (95% confidence interval 1.16–2.13), which corresponds to a 57% relative increase in mortality in patients treated with paclitaxel-coated devices.

As presented at the June 2019 FDA Advisory Committee Meeting, an independent meta-analysis of similar patient-level data provided by VIVA Physicians, a vascular medicine organization, reported similar findings with a hazard ratio of 1.38 (95% confidence interval 1.06–1.80). Additional analyses have been conducted and are underway that are specifically designed to assess the relationship of mortality to paclitaxel-coated devices. The presence and magnitude of the late mortality risk should be interpreted with caution because of multiple limitations in the available data, including wide confidence intervals due to a small sample size, pooling of studies of different paclitaxel-coated devices that were not intended to be combined, substantial amounts of missing study data, no clear evidence of a paclitaxel dose effect on mortality, and no identified pathophysiologic mechanism for the late deaths.

Paclitaxel-coated balloons and stents improve blood flow to the legs and decrease the likelihood of repeat procedures to reopen blocked blood vessels compared to uncoated devices. The benefits of paclitaxel-coated devices (e.g., reduced reinterventions) should be considered in individual patients along with potential risks (e.g., late mortality).

In the Zilver PTX randomized clinical study,<sup>8</sup> Kaplan Meier mortality estimates at 2, 3 and 5 years are 5.4% (95% confidence interval 3.0%–7.9%), 10.3% (7.0%–13.6%) and 19.1% (14.2%–24.1%), respectively, for all patients treated with the Zilver PTX treatment device and 5.7% (1.8%–9.5%), 10.7% (5.6%–15.9%) and 17.1% (10.0%–24.2%), respectively, for patients treated only with the PTA/BMS control device ( $P = 0.60$ ). Publications of large observational studies have not identified an increased risk of mortality with the Zilver PTX DES; Secemsky et al.<sup>9</sup> reported comparable mortality rates through 4 years for DES and BMS in claudicants ( $P = 0.55$ ) and CLI patients ( $P = 0.52$ ); Freisinger et al.<sup>10</sup> reported a hazard ratio of 1.01 (95% confidence interval 0.83–1.23) at 5 years for DES vs. BMS; Katsuki et al.<sup>11</sup> reported a hazard ratio of 0.89 (95% confidence interval 0.66–1.19) at 5 years for DES vs. non-paclitaxel devices.

## DRUG INTERACTIONS

Because systemic levels of paclitaxel post stent placement in clinical trials were extremely low and rapidly cleared ( $C_{max} < 10$  ng/mL;  $< 1$  ng/mL after 8 hours), possible interactions of paclitaxel with concomitantly administered medications are unlikely to be detectable.

The metabolism of paclitaxel is catalyzed by CYP2C8 and CYP3A4 which are cytochrome P450 isoenzymes. Formal drug interaction studies have not been conducted with the Zilver PTX Drug-Eluting Peripheral Stent. In the absence of formal clinical drug interaction studies, caution should be exercised when administering paclitaxel concomitantly with known substrates or inhibitors of the cytochrome P450 isoenzymes CYP2C8 and CYP3A4.

## PRODUCT RECOMMENDATIONS

### Illustrations

Placement of this vascular stent requires advanced skills in interventional vascular procedures. The following instructions will give technical guidance, but do not obviate formal training in the use of the device. Delivery system components are shown in **Figure 1**.

### Multiple Stent Placement

If placements of multiple stents are required in a patient, to cover the length of the lesion, the following recommendations should be considered:

- In relation to the lesion site, the distal area of narrowing should be stented first, followed by the proximal locations (i.e., a second stent should be placed proximally to the previously placed stent).
- Stents placed in tandem must overlap to allow for complete coverage of the lesion.
- If multiple stents are to be placed in an overlapping fashion, they should be of similar composition (i.e. nickel titanium alloy (nitinol)).
- Do not implant more than a total maximum drug quantity of 3 mg paclitaxel (i.e., a total stented length of 320 mm).

### Stent Sizing

- Determine the proper stent size after complete diagnostic evaluation. The stent deployment must be performed under fluoroscopic control. Measure the length of the target lesion to determine the length of the stent required. Allow for the proximal and distal aspects of the stent to cover the entire target area. **Note:** The Zilver PTX Drug-Eluting Peripheral Stent is designed not to shorten upon deployment.
- The stent is recommended for use in above-the-knee femoropopliteal arteries having reference vessel diameter from 4 mm to 7 mm.
- Measure the diameter of the reference vessel (proximal and distal to the lesion) and use the LARGEST reference diameter as your basis for choosing the appropriate stent size.
- The stent size should be selected so that the unconstrained stent diameter is at least 1 mm larger than the reference vessel diameter and no more than 2 mm larger than the reference vessel diameter.

### Introduction of the Stent

1. Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
2. Visually inspect the device with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for a return authorization.
3. Gain access to the site using a 6 French sheath with a minimum I.D. of 2.1 mm.

7. Katsanos K, et al. Risk of death following application of paclitaxel-coated balloons and stents in the femoropopliteal artery of the leg: a systematic review and meta-analysis of randomized controlled trials. *J Am Heart Assoc.* 2018;7:e011245.

8. Dake MD, et al. Paclitaxel-coated Zilver PTX Drug-Eluting Stent treatment does not result in increased long-term all-cause mortality compared to uncoated devices. *Cardiovasc Intervent Radiol.* 2019; <https://doi.org/10.1007/s00270-019-02324-4>

9. Secemsky EM, et al. Drug-eluting stent implantation and long-term survival following peripheral artery revascularization. *J Am Coll Cardiol.* 2019;73:2636–2638.

10. Freisinger E, et al. Mortality after use of paclitaxel-based devices in peripheral arteries: a real-world safety analysis. *Eur Heart J.* 2019; doi:10.1093/eurheartj/ehz698.

11. Katsuki T, et al. Mortality risk following application of a paclitaxel-coated stent in femoropopliteal lesions. *J Endovasc Ther.* 2019; doi:10.1177/1526602819870309

4. To ensure adequate support of the system, introduce a 0.89 mm (0.035 inch) wire guide. If hydrophilic wire guides are used, they must be kept fully activated.
5. Pre-dilatation before stent placement is optional and at the discretion of the physician.
6. Prior to advancing the system over the wire and introducing the system into the body use a syringe to flush the device wire guide lumen and stent with saline or heparinized saline through the hub of the device. Do not flush the device with contrast media.  
Flush the wire lumen until a few drops of saline exit the tip of the device (Figure 2). Next, gently place a finger over the device end hole and flush again until a few drops of saline can be seen exiting between the end of the sheath and device tip (Figure 3).
7. Under fluoroscopy, advance the delivery system over the 0.89 mm (0.035 inch) wire guide through the introducer sheath until the distal gold radiopaque markers on the stent are beyond the target lesion site (Figure 4).  
**Note:** If resistance is met during advancement of the delivery system do not force passage. Remove the delivery system and replace with a new device.  
**Note:** Ensure the distal end of the stability sheath is inside the introducer sheath.  
**Caution:** Torquing the device may lead to difficulty or inability to deploy and/or move device over the wire guide.

#### Alignment and Deployment of the Stent

1. Under fluoroscopy, pull the delivery system back until there is no slack in the delivery system and the radiopaque markers on the stent are aligned with the lesion (Figure 5).
2. Disengage the device safety lock by gently depressing the red safety button in the direction indicated in Figure 6.
3. Initiate stent deployment by slowly rotating the device thumbwheel in the direction shown in Figure 7.  
**Note:** If high resistance is felt on the thumbwheel prior to stent deployment, do not force deployment. Carefully withdraw the stent system without deploying the stent.
4. Rotating the thumbwheel will withdraw the retraction sheath from over the stent.  
**Note:** Visually confirm that the radiopaque stent markers are still aligned with the desired position. Reposition if necessary.  
**Note:** Stent deployment will begin once the distal radiopaque stent markers have advanced beyond the sheath marker band (Refer to Figure 1 for sheath and stent marker identification). Once the stent has apposed the vessel wall repositioning of the stent is not recommended. Recapture or re-sheathing of the Zilver PTX Drug-Eluting Peripheral Stent is not possible because the stent retraction sheath cannot be re-advanced. Refer to the Multiple Stent Placement section of these Instructions for Use for information on missed lesions.
5. Full deployment of the stent length will occur when the distal end of the retraction sheath has been retracted past the proximal radiopaque markers of the stent (Figure 8).

#### Post Stent Deployment

1. Under fluoroscopy, withdraw the entire delivery system as one unit, over the wire guide, into the introducer sheath and out of the body. Remove the delivery system from the wire guide.  
**Note:** Following stent deployment, if resistance is met during the withdrawal of the delivery system, carefully remove the delivery system and wire guide together as a unit.  
If resistance is still encountered during removal of the delivery system and wire guide as a unit, remove the wire guide, delivery system and introducer sheath together as a unit.
2. Perform an arterial angiogram to verify full deployment of the device. If incomplete expansion exists within the stent at any point along the lesion, post-deployment balloon dilatation (standard PTA) can be performed at the discretion of the physician.  
**Note:** Flow restrictions remaining after stent deployment (e.g., residual proximal or distal stenosis or dissection, or poor distal outflow) may increase the risk of stent thrombosis. Inflow and outflow should be assessed at procedure completion and additional measures considered (e.g., additional PTA, adjunctive stenting, or distal bypass) if necessary to maintain good inflow and outflow.
3. Remove the wire guide and introducer sheath from the patient.
4. Close the entry wound as appropriate.
5. Dispose in accordance with institutional guidelines

#### REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

#### DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be coiled for disposal in accordance with institutional guidelines.

#### PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

When available, the EUDAMED website (<https://ec.europa.eu/tools/eudamed>), along with the BUDI for this product (0827002CIRL202011054026AM), can be used to locate the Summary of Safety and Clinical Performance (SSCP) for this product.

Patient-facing information may be accessed at [cookmedical.eu/patient-implant-information](http://cookmedical.eu/patient-implant-information).

#### SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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- Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Zeller T, Roubin GS, Burket MW, Khatib V, Snyder SA, Ragheb AO, White JK, Machan LS; Zilver PTX Investigators. Paclitaxel-eluting stents show superiority to balloon angioplasty and bare metal stents in femoropopliteal disease: twelve-month Zilver PTX randomized study results. *Circ Cardiovasc Interv*. 2011;4:495-504. doi: 10.1161/CIRCINTERVENTIONS.111.962324.
- Dake MD, Scheinert D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kavteladze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther*. 2011;18:613-623. doi: 10.1583/11-3560.1.
- Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv*. 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.
- Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Machan LS, Snyder SA, O'Leary EE, Ragheb AO, Zeller T. Durable clinical effectiveness with paclitaxel-eluting stents in the femoropopliteal artery: 5-year results of the Zilver PTX randomized trial. *Circulation*. 2016;133:1472-1483. doi: 10.1161/CIRCULATIONAHA.115.016900.
- Tinwell H and Ashby J. Genetic toxicity and potential carcinogenicity of taxol. *Carcinogenesis* 1994; 15 (8): 1499-1501.
- Physician's Desk Reference, 52nd edition, Medical Economics Inc, 1998, pp 762-766.

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2. Dake MD, Scheinert D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kaveladze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther.* 2011;18:613-623. doi: 10.1583/11-3560.1.

3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv.* 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.

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2. Dake MD, Scheiner D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kavteladze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther.* 2011;18:613-623. doi: 10.1583/11-3560.1.

3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv.* 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.

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1. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Zeller T, Roubin GS, Burket MW, Khatib Y, Snyder SA, Ragheb AO, Zeller T, Machan LS; Zilver PTX Investigators. Paclitaxel-eluting stents show superiority to balloon angioplasty and bare metal stents in femoropopliteal disease: twelve-month Zilver PTX randomized study results. *Circ Cardiovasc Interv.* 2011;4:495-504. doi: 10.1161/CIRCINTERVENTIONS.111.962324.
2. Dake MD, Scheinerer D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kavteladze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther.* 2011;18:613-623. doi: 10.1583/11-3560.1.
3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv.* 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.
4. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Machan LS, Snyder SA, O'Leary EE, Ragheb AO, Zeller T. Durable clinical effectiveness with paclitaxel-eluting stents in the femoropopliteal artery: 5-year results of the Zilver PTX randomized trial. *Circulation.* 2016;133:1472-1483. doi: 10.1161/CIRCULATIONAHA.115.016900.
5. Timwell H and Ashby J. Genetic toxicity and potential carcinogenicity of taxol. *Carcinogenesis* 1994; 15 (8): 1499-1501.
6. Physician's Desk Reference, 52nd edition, Medical Economics Inc, 1998, pp 762-766.

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2. Dake MD, Scheinert D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kaveladze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther*. 2011;18:613-623. doi: 10.1583/11-3560.1.

3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv*. 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.

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1. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Zeller T, Roubin GS, Burket MW, Khatib Y, Snyder SA, Ragheb AO, White JK, Machan LS; Zilver PTX Investigators. Paclitaxel-eluting stents show superiority to balloon angioplasty and bare metal stents in femoropopliteal disease: twelve-month Zilver PTX randomized study results. *Circ Cardiovasc Interv.* 2011;4:495-504. doi: 10.1161/CIRCINTERVENTIONS.111.962324.

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2. Dake MD, Scheinert D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kavteldadze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther.* 2011;18:613-623. doi: 10.1583/11-3560.1.

3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv.* 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.

4. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Machan LS, Snyder SA, O'Leary EE, Ragheb AO, Zeller T. Durable clinical effectiveness with paclitaxel-eluting stents in the femoropopliteal artery: 3-year results of the Zilver PTX randomized trial. *Circulation.* 2016;133:1472-1483. doi: 10.1161/CIRCULATIONAHA.115.016900.

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1. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Zeller T, Roubin GS, Burket MW, Khatib Y, Snyder SA, Ragheb AO, White JK, Machan LS, Zilver PTX Investigators. Paclitaxel-eluting stents show superiority to balloon angioplasty and bare metal stents in femoropopliteal disease: twelve-month Zilver PTX randomized study results. *Circ Cardiovasc Interv.* 2011;4:495-504. doi: 10.1161/CIRCINTERVENTIONS.111.962324.
2. Dake MD, Scheinert D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kavteladze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther.* 2011;18:613-623. doi: 10.1583/11-3560.1.
3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Lottes SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv.* 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.
4. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Machan LS, Snyder SA, O'Leary EE, Ragheb AO, Zeller T. Durable clinical effectiveness with paclitaxel-eluting stents in the femoropopliteal artery: 5-year results of the Zilver PTX randomized trial. *Circulation.* 2016;133:1472-1483. doi: 10.1161/CIRCULATIONAHA.115.016900.
5. Tinwell H and Ashby J. Genetic toxicity and potential carcinogenicity of taxol. *Carcinogenesis* 1994; 15 (8): 1499-1501.
6. Physician's Desk Reference, 52nd edition, Medical Economics Inc, 1998, pp 762-766.

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1. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Zeller T, Roubin GS, Burket MW, Khatib Y, Snyder SA, Ragheb AO, White JK, Machan LS; Zilver PTX Investigators. Paclitaxel-eluting stents show superiority to balloon angioplasty and bare metal stents in femoropopliteal disease: twelve-month Zilver PTX randomized study results. *Circ Cardiovasc Interv.* 2011;4:495-504. doi: 10.1161/CIRCINTERVENTIONS.111.962324.
2. Dake MD, Scheiner D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kavteldaze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther.* 2011;18:613-623. doi: 10.1583/11-3560.1.
3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv.* 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.
4. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Machan LS, Snyder SA, O'Leary EE, Ragheb AO, Zeller T. Durable clinical effectiveness with paclitaxel-eluting stents in the femoropopliteal artery: 5-year results of the Zilver PTX randomized trial. *Circulation.* 2016;133:1472-1483. doi: 10.1161/CIRCULATIONAHA.115.016900.
5. Tinwell H and Ashby J. Genetic toxicity and potential carcinogenicity of taxol. *Carcinogenesis* 1994; 15 (8): 1499-1501.
6. Physician's Desk Reference, 52nd edition, Medical Economics Inc, 1998, pp 762-766.

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1. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Zeller T, Roubin GS, Burket MW, Khatib Y, Snyder SA, Ragheb AO, White JK, Machan LS; Zilver PTX Investigators. Paclitaxel-eluting stents show superiority to balloon angioplasty and bare metal stents in femoropopliteal disease: twelve-month Zilver PTX randomized study results. *Circ Cardiovasc Interv.* 2011;4:495-504. doi: 10.1161/CIRCINTERVENTIONS.111.962324.
2. Dake MD, Scheinert D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kavteladze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther.* 2011;18:613-623. doi: 10.1583/11-3560.1.
3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv.* 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.
4. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Machan LS, Snyder SA, O'Leary EE, Ragheb AO, Zeller T. Durable clinical effectiveness with paclitaxel-eluting stents in the femoropopliteal artery: 5-year results of the Zilver PTX randomized trial. *Circulation.* 2016;133:1472-1483. doi: 10.1161/CIRCULATIONAHA.115.016900.
5. Tinwell H and Ashby J. Genetic toxicity and potential carcinogenicity of taxol. *Carcinogenesis* 1994; 15 (8): 1499-1501.
6. Physician's Desk Reference, 52nd edition, Medical Economics Inc, 1998, pp 762-766.

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7. Katsanos K, et al. Risk of death following application of paclitaxel-coated balloons and stents in the femoropopliteal artery of the leg: a systematic review and meta-analysis of randomized controlled trials. *J Am Heart Assoc*. 2018;7:e011245.

8. Dake MD, et al. Paclitaxel-coated Zilver PTX Drug-Eluting Stent treatment does not result in increased long-term all-cause mortality compared to uncoated devices. *Cardiovasc Intervent Radiol*. 2019; <https://doi.org/10.1007/s00270-019-02324-4>

9. Secemsky EM, et al. Drug-eluting stent implantation and long-term survival following peripheral artery revascularization. *J Am Coll Cardiol*. 2019;73:2636–2638.

10. Freisinger E, et al. Mortality after use of paclitaxel-based devices in peripheral arteries: a real-world safety analysis. *Eur Heart J*. 2019; doi:10.1093/eurheartj/ehz698.

11. Katsuki T, et al. Mortality risk following application of a paclitaxel-coated stent in femoropopliteal lesions. *J Endovasc Ther*. 2019; doi:10.1177/1526602819870309

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1. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Zeller T, Roubin GS, Burket MW, Khatib Y, Snyder SA, Ragheb AO, White JK, Machan LS; Zilver PTX Investigators. Paclitaxel-eluting stents show superiority to balloon angioplasty and bare metal stents in femoropopliteal disease: twelve-month Zilver PTX randomized study results. *Circ Cardiovasc Interv.* 2011;4:495–504. doi: 10.1161/CIRCINTERVENTIONS.111.962324.

2. Dake MD, Scheinert D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kavteladze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther.* 2011;18:613–623. doi: 10.1583/11-3560.1.

3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv.* 2016;9:271–277. doi: 10.1016/j.jcin.2015.09.035.

4. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Machan LS, Snyder SA, O'Leary EE, Ragheb AO, Zeller T. Durable clinical effectiveness with paclitaxel-eluting stents in the femoropopliteal artery: 5-year results of the Zilver PTX randomized trial. *Circulation.* 2016;133:1472–1483. doi: 10.1161/CIRCULATIONAHA.115.016900.

5. Tinwell H and Ashby J. Genetic toxicity and potential carcinogenicity of taxol. *Carcinogenesis* 1994; 15 (8): 1499–1501.

6. Physician's Desk Reference, 52nd edition, Medical Economics Inc, 1998, s. 762–766.

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1. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Zeller T, Roubin GS, Burket MW, Khatib Y, Snyder SA, Ragheb AO, White JK, Machan LS; Zilver PTX Investigators. Paclitaxel-eluting stents show superiority to balloon angioplasty and bare metal stents in femoropopliteal disease: twelve-month Zilver PTX randomized study results. *Circ Cardiovasc Interv.* 2011;4:495-504. doi: 10.1161/CIRCINTERVENTIONS.111.962324.
2. Dake MD, Scheinert D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kavteladze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther.* 2011;18:613-623. doi: 10.1583/11-3560.1.
3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv.* 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.
4. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Machan LS, Snyder SA, O'Leary EE, Ragheb AO, Zeller T. Durable clinical effectiveness with paclitaxel-eluting stents in the femoropopliteal artery: 5-year results of the Zilver PTX randomized trial. *Circulation.* 2016;133:1472-1483. doi: 10.1161/CIRCULATIONAHA.115.016900.
5. Tinwell H and Ashby J. Genetic toxicity and potential carcinogenicity of taxol. *Carcinogenesis* 1994; 15 (8): 1499-1501.
6. Physician's Desk Reference, wydanie 52., Medical Economics Inc, 1998, pp 762-766.

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7. Katsanos K, et al. Risk of death following application of paclitaxel-coated balloons and stents in the femoropopliteal artery of the leg: a systematic review and meta-analysis of randomized controlled trials. *J Am Heart Assoc.* 2018;7:e011245.
8. Dake MD, et al. Paclitaxel-coated Zilver PTX Drug-Eluting Stent treatment does not result in increased long-term all-cause mortality compared to uncoated devices. *Cardiovasc Intervent Radiol.* 2019; <https://doi.org/10.1007/s00270-019-02324-4>
9. Secemsky EM, et al. Drug-eluting stent implantation and long-term survival following peripheral artery revascularization. *J Am Coll Cardiol.* 2019;73:2636–2638.
10. Freisinger E, et al. Mortality after use of paclitaxel-based devices in peripheral arteries: a real-world safety analysis. *Eur Heart J.* 2019; doi:10.1093/eurheartj/ehz698.
11. Katsuki T, et al. Mortality risk following application of a paclitaxel-coated stent in femoropopliteal lesions. *J Endovasc Ther.* 2019; doi:10.1177/1526602819870309

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1. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Zeller T, Roubin GS, Burket MW, Khatib Y, Snyder SA, Ragheb AO, White JK, Machan LS; Zilver PTX Investigators. Paclitaxel-eluting stents show superiority to balloon angioplasty and bare metal stents in femoropopliteal disease: twelve-month Zilver PTX randomized study results. *Circ Cardiovasc Interv.* 2011;4:495-504. doi: 10.1161/CIRCINTERVENTIONS.111.962324.
2. Dake MD, Scheinert D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kaveladze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther.* 2011;18:613-623. doi: 10.1583/11-3560.1.
3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv.* 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.
4. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Machan LS, Snyder SA, O'Leary EE, Ragheb AO, Zeller T. Durable clinical effectiveness with paclitaxel-eluting stents in the femoropopliteal artery: 5-year results of the Zilver PTX randomized trial. *Circulation.* 2016;133:1472-1483. doi: 10.1161/CIRCULATIONAHA.115.016900.
5. Tinwell H and Ashby J. Genetic toxicity and potential carcinogenicity of taxol. *Carcinogenesis* 1994; 15 (8): 1499-1501.
6. Physician's Desk Reference, 52nd edition, Medical Economics Inc, 1998, pp 762-766.

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1. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Zeller T, Roubin GS, Burkett MW, Khatib Y, Snyder SA, Ragheb AO, White JK, Machan LS; Zilver PTX Investigators. Paclitaxel-eluting stents show superiority to balloon angioplasty and bare metal stents in femoropopliteal disease: twelve-month Zilver PTX randomized study results. *Circ Cardiovasc Interv.* 2011;4:495-504. doi: 10.1161/CIRCINTERVENTIONS.111.962324.

2. Dake MD, Scheinert D, Tepe G, Tessarek J, Fanelli F, Bosiers M, Ruhlmann C, Kavteladze Z, Lottes AE, Ragheb AO, Zeller T; Zilver PTX Single-Arm Study Investigators. Nitinol stents with polymer-free paclitaxel coating for lesions in the superficial femoral and popliteal arteries above the knee: twelve-month safety and effectiveness results from the Zilver PTX single-arm clinical study. *J Endovasc Ther.* 2011;18:613-623. doi: 10.1583/11-3560.1.

3. Yokoi H, Ohki T, Kichikawa K, Nakamura M, Komori K, Nanto S, O'Leary EE, Lottes AE, Snyder SA, Dake MD. Zilver PTX post-market surveillance study of paclitaxel-eluting stents for treating femoropopliteal artery disease in Japan: 12-month results. *J Am Coll Cardiol Interv.* 2016;9:271-277. doi: 10.1016/j.jcin.2015.09.035.

4. Dake MD, Ansel GM, Jaff MR, Ohki T, Saxon RR, Smouse HB, Machan LS, Snyder SA, O'Leary EE, Ragheb AO, Zeller T. Durable clinical effectiveness with paclitaxel-eluting stents in the femoropopliteal artery: 5-year results of the Zilver PTX randomized trial. *Circulation.* 2016;133:1472-1483. doi: 10.1161/CIRCULATIONAHA.115.016900.

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8.Dake MD, et al. Paclitaxel-coated Zilver PTX Drug-Eluting Stent treatment does not result in increased long-term all-cause mortality compared to uncoated devices. *Cardiovasc Intervent Radiol.* 2019; <https://doi.org/10.1007/s00270-019-02324-4>

9.Secemsky EM, et al. Drug-eluting stent implantation and long-term survival following peripheral artery revascularization. *J Am Coll Cardiol.* 2019;73:2636–2638.

10.Freisinger E, et al. Mortality after use of paclitaxel-based devices in peripheral arteries: a real-world safety analysis. *Eur Heart J.* 2019; doi:10.1093/eurheartj/ehz698.

11.Katsuki T, et al. Mortality risk following application of a paclitaxel-coated stent in femoropopliteal lesions. *J Endovasc Ther.* 2019; doi:10.1177/1526602819870309

EN A symbol glossary can be found at [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • BG Речник на символите може да се намери на [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • CS Glosář značek naleznete na [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • DA Der findes en symbolforklaring på [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • DE Eine Symbollegende finden Sie auf [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • EL Μπορείτε να βρείτε ένα γλωσσάρι των συμβόλων στη διεύθυνση [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • ES En [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary), puede consultarse un glosario de símbolos • ET Sümboleite sõnastik on esitatud aadressil [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • FR Pour un glossaire des symboles, consulter le site Web [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • HR Pojmovnik simbola dostupan je na [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • HU A szimbólumok magyarázata a [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) weboldalon található • IT Un glossario dei simboli è disponibile alla pagina [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • LT Simbolių žodynelį galite rasti svetainėje [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • LV Simbolu vārdnīca ir pieejama tīmekļa vietnē [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • NL Een verklarende lijst van symbolen is te vinden op [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • NO En symbolforklaring finnes på [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • PL Słownik symboli można znaleźć pod adresem [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • PT Existe um glossário de símbolos em [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • RO Se poate accesa un glosar de simboluri la [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • SK Slovník symbolov nájdete na adrese [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • SL Slovarček simbolov najdete na [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • SV En symbolordlista finns på [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) • TR Sembol sözlüğüne [cookmedical.com/symbol-glossary](http://cookmedical.com/symbol-glossary) adresinden ulaşılabilir

EN Medical device • BG Медицинско изделие • CS Zdravotnický prostředek • DA Medicinsk udstyr • DE Medizinprodukt • EL Ιατροτεχνολογικό προϊόν • ES Producto sanitario • ET Meditsiiniseade • FR Dispositif médical • HR Medicinski proizvod • HU Orvostechnikai eszköz • IT Dispositivo medico • LT Medicinos priemonė • LV Mediciniskā ierīce • NL Medisch hulpmiddel • NO Medisinsk anordning • PL Wyrób medyczny • PT Dispositivo médico • RO Dispozitiv medical • SK Zdravotnícka pomôcka • SL Medicinski pripomoček • SV Medicinteknisk produkt • TR Tıbbi cihaz

EN Do not use if package is damaged and consult instructions for use • BG Да не се използва, ако опаковката е повредена и да се направи консултация с инструкциите за употреба • CS Nepoužívejte, pokud je obal poškozený, a prostudujte si návod k použití • DA Må ikke anvendes, hvis emballagen er beskadiget. Se brugsanvisningen • DE Nicht verwenden, wenn die Verpackung beschädigt ist, und Gebrauchsanweisung beachten • EL Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης • ES No utilizar si el envase está dañado y consultar las instrucciones de uso • ET Ärge kasutage kahjustatud pakendi korral ja tutvuge kasutusjuhendiga • FR Ne pas utiliser si l'emballage est endommagé et consulter le mode d'emploi • HR Nemojte upotrebljavati ako je ambalaža oštećena te pročítajte upute za uporabu • HU Ne használja, ha a csomagolás sérült, és olvassa el a használati utasítást • IT Non usare se la confezione è danneggiata e consultare le Istruzioni per l'uso • LT Nenaudoti, jei pakuotė pažeista, ir žiūrėti naudojimo instrukciją • LV Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju • NL Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing • NO Skal ikke brukes dersom pakningen er skadet; se bruksanvisningen • PL Nie używać, jeśli opakowanie zostało uszkodzone; sprawdź w instrukcji użycia • PT Não utilizar caso a embalagem esteja danificada, e consultar as instruções de utilização • RO A nu se utiliza dacă ambalajul este deteriorat și a se consulta instrucțiunile de utilizare • SK Ak je balenie poškodené, nepoužívajte a prečítajte si návod na použitie • SL Če je embalaža poškodovana, pripomočka ne uporabite, glejte navodila za uporabo • SV Får inte användas om förpackningen är skadad, läs bruksanvisningen • TR Ambalaj hasarlıysa kullanmayın ve kullanma talimatına başvurun

EN Single sterile barrier system • BG Единична стерилна бариерна система • CS Systém jednoduché sterilní bariéry • DA Enkelt sterilt barriersystem • DE Einfaches Sterilbarrierensystem • EL Μονό σύστημα φραγμού αποστείρωσης • ES Sistema de barrera estéril individual • ET Ühekordne steriilne kaitsemeetod • FR Système de barrière stérile simple • HR Sustav jednostruke sterilne barijere • HU Egyszeres sterilgát-rendszer • IT Sistema a barriera sterile singola • LT Vieno sterilaus barjero sistema • LV Vienas sterilās barjeras sistēma • NL Stelsysteem met enkele steriele barrière • NO System med én steril barriere • PL System pojedynczej bariery sterylnej • PT Sistema de barreira única estéril • RO Sistem de barieră sterilă unică • SK Systém jednej sterilnej bariéry • SL Enojni sterilni pregradni sistem • SV System med en steril barriär • TR Tekli steril bariyer sistemi

EN Single sterile barrier system with protective packaging outside • BG Единична стерилна бариерна система с външна защитна опаковка • CS Systém jednoduché sterilní bariéry s vnějším ochranným obalem • DA Enkelt sterilt barriersystem med beskyttelsesemballage udvendigt • DE Einfaches Sterilbarrierensystem mit äußerer Schutzverpackung • EL Μονό σύστημα φραγμού αποστείρωσης με προστατευτική εξωτερική συσκευασία • ES Sistema de barrera estéril individual con embalaje protector exterior • ET Ühekordne steriilne kaitsemeetod koos kaitsva välispakendiga • FR Système de barrière stérile simple avec emballage protecteur externe • HR Sustav jednostruke sterilne barijere s vanjskim zaštitnim pakiranjem • HU Egyszeres steril védőzáras rendszere külső védőcsomagolással • IT Sistema a barriera sterile singola con confezionamento protettivo esterno • LT Vieno sterilaus barjero sistema su apsaugine išorine pakuote • LV Vienas sterilās barjeras sistēma ar aizsargiapakoju ārpusē • NL Stelsysteem met enkele steriele barrière en beschermende buitenverpakking • NO System med én steril barriere med beskyttende emballasje på utsiden • PL System pojedynczej bariery sterylnej z opakowaniem ochronnym na zewnątrz • PT Sistema de barreira única estéril com embalagem externa protetora • RO Sistem de barieră sterilă unică cu ambalaj de protecție în exterior • SK Systém jednej sterilnej bariéry s vonkajším ochranným obalom • SL Enojni sterilni pregradni sistem z zunanjo zaščitno embalažo • SV System med en steril barriär med skyddsförpackning på utsidan • TR Dişinda koruyucu ambalaj içeren tekli steril bariyer sistemi

**Australian Sponsor**

EN Australian sponsor • BG Възложител за Австралия • CS Australský zadavatel • DA Australsk sponsor • DE Australischer Sponsor • EL Χορηγός για την Αυστραλία • ES Promotor australiano • ET Australia sponsor • FR Partenaire australien • HR Australski sponsor • HU Ausztráliai partner • IT Sponsor australiano • LT Australijos užsakovas • LV Austrālijas sponsors • NL Australische opdrachtgever • NO Australsk sponsor • PL Sponsor australijski • PT Promotor australiano • RO Sponsor australian • SK Zadávateľ z Austrálie • SL Avstralski sponsor • SV Australiensk sponsor • TR Avustralya Sponsoru

EN Contains a medicinal substance • BG Съдържа лекарствено вещество • CS Obsahuje léčivou látku • DA Indeholder et lægemiddel • DE Enthält ein Arzneimittel • EL Περιέχει φαρμακευτική ουσία • ES Contiene un fármaco • ET Sisaldab ravimainet • FR Contient une substance médicinale • HR Sadržaji ljekovitu tvar • HU Gyógyhatású anyagot tartalmaz • IT Contiene una sostanza medicinale • LT Sudėtyje yra medicininės medžiagos • LV Satur zāles • NL Bevat een geneeskrachtige substantie • NO Inneholder et legemiddel • PL Zawiera substancję leczniczą • PT Contém uma substância medicamentosa • RO Conține o substanță medicinală • SK Obsahuje liečivo • SL Vsebuje zdravilno učinkovino • SV Innehåller en läkemedelssubstans • TR Medisinal madde içerir

EN Contains the hazardous substance cobalt • BG Съдържа опасното вещество кобалт • CS Obsahuje nebezpečnou látku kobalt • DA Indeholder det skadelige stof kobolt • DE Enthält den Gefahrstoff Kobalt • EL Περιέχει την επικίνδυνη ουσία κοβάλτιο • ES Contiene la sustancia peligrosa cobalto • ET Sisaldab ohtlikku ainet koobaltit • FR Contient une substance dangereuse, le cobalt • HR Sadržaji opasnu tvar kobalt • HU Veszélyes anyagot, név szerint kobaltot tartalmaz • IT Contiene la sostanza pericolosa cobalto • LT Sudėtyje yra pavojingosios medžiagos kobalto • LV Satur bīstamu vielu kobaltu • NL Bevat de gevaarlijke stof kobalt • NO Inneholder det farlige stoffet kobolt • PL Zawiera substancję niebezpieczną: kobalt • PT Contém a substância perigosa cobalto • RO Conține substanța periculoasă cobalt • SK Obsahuje nebezpečnú látku kobalt • SL Vsebuje nevarno snov kobalt • SV Innehåller det farliga ämnet kobolt • TR Tehlikeli madde içerir

Rx ONLY

STERILE	EO
---------	----

**Australian Sponsor**

William A. Cook Australia Pty Ltd  
95 Brandl St  
Eight Mile Plains QLD 4113  
Australia

**COOK IRELAND LTD.**  
O'Halloran Road  
National Technology Park  
Limerick, Ireland

[cookmedical.com](http://cookmedical.com)

**IFU0117-6**