**EN** Zilver® PTX® Drug-Eluting Peripheral Stent

2

Instructions for Use

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Mode d'emploi

**ZILVER® PTX® DRUG-ELUTING PERIPHERAL STENT**

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

**DEVICE DESCRIPTION**

The Zilver PTX Drug-Eluting Peripheral Stent is a self-expanding stent made of nitinol and coated with the drug paclitaxel. The stent is intended for permanent implant. It is a flexible, slotted tube that is designed to provide support while maintaining flexibility in the vessel upon deployment. Post-deployment, the stent is designed to impart an outward radial force upon the inner lumen of the vessel, establishing patency in the stented region.

The stent is preloaded in a 2.1mm (6 French) delivery sheath. Hand-loading of the stent is not possible. Stent deployment is controlled by rotating the thumbwheel on the handle.

**Table 1. Zilver PTX Drug-Eluting Peripheral Stent and delivery system product description**


Zilver PTX Drug-Eluting Peripheral Stent	
Available stent lengths	40, 60, 80, 100, 120, 140 mm
Available stent diameters	5, 6, 7, 8 mm
Stent material	Nitinol
Drug Coating	Paclitaxel
Delivery System	
Available delivery system lengths	80, 125 cm
Wire guide compatibility	0.035 in
Catheter shaft outer diameter	6 Fr (2.1 mm)
Introducer sheath (minimum)	6 Fr

**Drug Component Description**

Paclitaxel may be extracted from the bark, branches, needles of the yew tree or paclitaxel may be produced via plant cell fermentation using a dedicated plant cell strain. In both cases, paclitaxel is purified and concentrated by column chromatography, crystallization and recrystallization.

Zilver PTX Drug-Eluting Peripheral Stents are coated with paclitaxel API (active pharmaceutical ingredient) using a proprietary process. No excipients, polymers, carriers, binding agents, other materials, or other device modifications are involved. The chemical description of paclitaxel is provided in **Figure 1**.

Paclitaxel
<ul style="list-style-type: none"> <li>• Synonyms: Taxol, Taxol A, Hunxol I, Paclitaxelum</li> <li>• IUPAC systematic name: <math>\beta</math>-(benzoylamino)-<math>\alpha</math>-hydroxy-6,12b-bis(acetyloxy)-12-(benzyloxy)-2a,3,4,4a,5,6,9,10,11,12,12a,12b-dodecahydro-4,11-dihydroxy-4a,8,13,13-tetramethyl-5-oxo-7,11-methano-1H-cyclodeca(3,4)benz(1,2-b)oxet-9-yl ester,(2aR-(2a-a,4-<math>\beta</math>,4a-<math>\beta</math>,6-<math>\beta</math>,9-a)-(a-R*,<math>\beta</math>-S*))11-a,12-a,12a-a,2b-a)-benzenepranoic acid</li> <li>• CAS registry number: 33069-62-4</li> <li>• Chemical formula: C<sub>47</sub>H<sub>51</sub>NO<sub>14</sub></li> <li>• Structure of paclitaxel:</li> </ul>

**Figure 1. Chemical description of paclitaxel**

**Table 2** presents the stent sizes and the nominal total quantity of paclitaxel on each stent based on the established dose density of 3  $\mu$ g/mm<sup>2</sup>.

**Table 2. Stent sizes and paclitaxel total quantity**

Model Number		Nominal Stent Size		Nominal Total Paclitaxel ( $\mu$ g/stent)
80 cm delivery system	125 cm delivery system	Diameter (mm)	Length (mm)	
ZISV6-35-80-5-40-PTX	ZISV6-35-125-5-40-PTX	5	40	383
ZISV6-35-80-6-40-PTX	ZISV6-35-125-6-40-PTX	6	40	383
ZISV6-35-80-7-40-PTX	ZISV6-35-125-7-40-PTX	7	40	383
ZISV6-35-80-8-40-PTX	ZISV6-35-125-8-40-PTX	8	40	356
ZISV6-35-80-5-60-PTX	ZISV6-35-125-5-60-PTX	5	60	552
ZISV6-35-80-6-60-PTX	ZISV6-35-125-6-60-PTX	6	60	552
ZISV6-35-80-7-60-PTX	ZISV6-35-125-7-60-PTX	7	60	552
ZISV6-35-80-8-60-PTX	ZISV6-35-125-8-60-PTX	8	60	579
ZISV6-35-80-5-80-PTX	ZISV6-35-125-5-80-PTX	5	80	722
ZISV6-35-80-6-80-PTX	ZISV6-35-125-6-80-PTX	6	80	722
ZISV6-35-80-7-80-PTX	ZISV6-35-125-7-80-PTX	7	80	722
ZISV6-35-80-8-80-PTX	ZISV6-35-125-8-80-PTX	8	80	757
ZISV6-35-80-5-100-PTX	ZISV6-35-125-5-100-PTX	5	100	934
ZISV6-35-80-6-100-PTX	ZISV6-35-125-6-100-PTX	6	100	934
ZISV6-35-80-7-100-PTX	ZISV6-35-125-7-100-PTX	7	100	934
ZISV6-35-80-8-100-PTX	ZISV6-35-125-8-100-PTX	8	100	935
ZISV6-35-80-5-120-PTX	ZISV6-35-125-5-120-PTX	5	120	1103
ZISV6-35-80-6-120-PTX	ZISV6-35-125-6-120-PTX	6	120	1103
ZISV6-35-80-7-120-PTX	ZISV6-35-125-7-120-PTX	7	120	1103
ZISV6-35-80-8-120-PTX	ZISV6-35-125-8-120-PTX	8	120	1112
ZISV6-35-80-5-140-PTX	ZISV6-35-125-5-140-PTX	5	140	1273
ZISV6-35-80-6-140-PTX	ZISV6-35-125-6-140-PTX	6	140	1273
ZISV6-35-80-7-140-PTX	ZISV6-35-125-7-140-PTX	7	140	1273

**MRI INFORMATION**

Non-clinical testing has demonstrated the Zilver PTX Drug-Eluting Peripheral Stent is MR Conditional for overlapping lengths up to 395 mm. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla or 3.0 Tesla
- Maximum spatial field gradient of 1600 Gauss/cm (16 T/m)
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the Zilver PTX Drug-Eluting Peripheral Stent is expected to produce a maximum temperature rise of less than 4.9°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 5 mm from the Zilver PTX Drug-Eluting Peripheral Stent when imaged with a gradient echo pulse sequence and a 3.0 Tesla MR system.

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International  
 2323 Colorado Avenue  
 Turlock, CA 95382  
 Phone: 888-633-4298 (toll free)  
 209-668-3333 from outside the US  
 Fax: 209-669-2450  
 Web: www.medicalert.org

## INDICATIONS FOR USE

The Zilver® PTX® Drug-Eluting Peripheral Stent is indicated for improving luminal diameter for the treatment of *de novo* or restenotic symptomatic lesions in native vascular disease of the above-the-knee femoropopliteal arteries having reference vessel diameter from 4 mm to 7 mm and total lesion lengths up to 300 mm per patient.

## CONTRAINDICATIONS

- Women who are pregnant, breastfeeding, or plan to become pregnant in the next 5 years should not receive a Zilver PTX Drug-Eluting Peripheral Stent. It is unknown whether paclitaxel will be excreted in human milk, and there is a potential for adverse reaction in nursing infants from paclitaxel exposure.
- Patients who cannot receive recommended anti-platelet and/or anti-coagulant therapy.
- Patients judged to have a lesion that prevents proper placement of the stent or stent delivery system.

## WARNINGS

- **A signal for increased risk of late mortality has been identified following the use of paclitaxel-coated balloons and paclitaxel-eluting stents for femoropopliteal arterial disease beginning approximately 2-3 years post-treatment compared with the use of non-drug coated devices. There is uncertainty regarding the magnitude and mechanism for the increased late mortality risk, including the impact of repeat paclitaxel-coated device exposure. Physicians should discuss this late mortality signal and the benefits and risks of available treatment options with their patients. See SUMMARY OF CLINICAL INVESTIGATIONS section for further information.**
- Persons with allergic reactions to nitinol, or its components, nickel and titanium, may suffer an allergic reaction to this implant.
- Persons allergic to paclitaxel or structurally-related compounds may suffer an allergic reaction to this implant.
- Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Inspect the product to ensure no damage has occurred.
- The use of this Drug-Eluting Peripheral Stent carries the risks associated with peripheral artery stenting, including vascular complications, and/or bleeding events.
- The safety and effectiveness of implanting multiple Zilver PTX Drug-Eluting Peripheral Stents with a total drug coating quantity of greater than 3 mg paclitaxel (i.e., an additive stent length greater than 160mm per limb) has not been established.
- Flow restrictions remaining after stent deployment (e.g., residual proximal or distal stenosis or dissection, or poor distal outflow) may increase the risk of stent thrombosis. Inflow and outflow should be assessed at procedure completion and additional measures considered (e.g., additional PTA, adjunctive stenting, or distal bypass) if necessary to maintain good inflow and outflow.
- The device is designed for single use only. Attempts to reprocess, resterilize and/or reuse may lead to device failure and/or transmission of disease.

## PRECAUTIONS

- To avoid involvement of the common femoral artery, the most proximal stent end should be placed at least 1 cm below the origin of the superficial femoral artery. To avoid involvement of the below-the-knee popliteal artery, the most distal stent end should be placed above the plane of the femoral epicondyles.
- This product is intended for use by physicians trained and experienced in diagnostic and interventional vascular techniques. Standard techniques for interventional vascular procedures should be employed.
- Manipulation of the Zilver PTX Drug-Eluting Peripheral Stent requires fluoroscopic control.
- Do not try to push the delivery system through stenoses that cannot be dilated to permit passage of the introducer catheter.
- If resistance is met during advancement of the delivery system, do not force passage. Remove the delivery system and replace with a new device. Continuing to force passage may ultimately cause partial deployment of the stent. Possible outcomes of partial deployment are:
  - Prolonged procedure (e.g. additional exposure to radiation/contrast)
  - Additional stent required
  - Secondary intervention (e.g. vessel incision)
- Do not try to remove the stent from the introducer system before use.
- Ensure that the red safety lock is not inadvertently depressed before stent deployment is desired.
- A 0.035 inch (0.89mm) diameter wire guide should be used during tracking, deployment, and removal in order to ensure adequate support of the system. If hydrophilic wire guides are used, they must be kept fully activated. Using a smaller diameter wire guide may ultimately result in partial deployment of the stent. Possible outcomes of partial deployment are:
  - Prolonged procedure (e.g. additional exposure to radiation/contrast)
  - Additional stent required
  - Secondary intervention (e.g. vessel incision)
- Do not use excessive force to deploy the stent. If excessive resistance is felt when beginning deployment, remove the delivery system without deploying the stent and replace with a new device.
- Do not expose the delivery system to organic solvents (e.g., alcohol).
- Do not use power injection systems with the delivery system.
- Do not torque the delivery system during introduction or deployment. Torquing the device may lead to difficulty or inability to deploy and/or move device over the wire guide.
- Appropriate antiplatelet/anticoagulant therapy should be administered pre- and post-procedure (see section entitled **PRE- and POST-PROCEDURE ANTIPLATELET REGIMEN** in these Instructions for Use). Use in patients who are unable to tolerate the appropriate antiplatelet therapy is not recommended.
- Safety and effectiveness of the Zilver PTX Drug-Eluting Peripheral Stent has not been demonstrated in patients with a history of bleeding disorders.
- Use of the Zilver PTX Drug-Eluting Peripheral Stent in an arterial vessel where leakage from the artery could be exacerbated by placement of the stent is not recommended.
- A low incidence of stent fracture has been reported (0.9% at 12 months in the randomized pivotal study). Although no clinical sequelae were associated with stent fracture in the randomized study through 12 months, the long-term clinical consequence of stent fracture is not yet established. The majority of stent fractures were associated with stent elongation  $\geq$  10% at deployment. Therefore, care should be taken when deploying the stent to minimize the risk of stent fracture due to elongation at implant.
- Once the stent has apposed the vessel wall repositioning of the stent is not recommended. Recapture or re-sheathing of the Zilver PTX Drug-Eluting Peripheral Stent is not possible because the stent retraction sheath cannot be re-advanced. Refer to the Multiple Stent Placement section of these Instructions for Use for information on missed lesions.
- If multiple stents are placed in an overlapping fashion, they should be of similar composition (i.e., nitinol).
- Do not use the stent after the "Use By" date specified on the package.
- Following stent deployment, if resistance is met during the withdrawal of the delivery system, carefully remove the delivery system and wire guide as a unit. If resistance is still encountered during removal of the delivery system and wire guide as a unit, remove the wire guide, delivery system and introducer sheath together as a unit.
- Prior to disengaging the device safety lock ensure the distal end of the stability sheath is inside the introducer sheath. Failure to do so may result in stent damage and/or stent compression upon deployment.

## POTENTIAL ADVERSE EVENTS

**Potential adverse events that may occur include, but are not limited to, the following:**

- Allergic reaction to anticoagulant and/or antithrombotic therapy or contrast medium
- Allergic reaction to nitinol
- Atheroembolization (Blue Toe Syndrome)
- Arterial aneurysm
- Arterial rupture
- Arterial thrombosis
- Arteriovenous fistula
- Death
- Dissection
- Embolism
- Hematoma/hemorrhage
- Hypersensitivity reactions
- Infection
- Infection/abscess formation at access site
- Ischemia requiring intervention (bypass or amputation of toe, foot, or leg)
- Occlusion
- Pseudoaneurysm formation
- Renal failure
- Restenosis of the stented artery
- Stent embolization
- Stent malapposition
- Stent migration
- Stent strut fracture
- Vessel perforation or rupture
- Vessel spasm
- Worsened claudication/rest pain

Although systemic effects are not anticipated, refer to the Physicians' Desk Reference for more information on the potential adverse events observed with paclitaxel. Potential adverse events, not described in the above source, may be unique to the paclitaxel drug coating:

- Allergic/immunologic reaction to the drug coating
- Alopecia
- Anemia
- Blood product transfusion
- Gastrointestinal symptoms
- Hematologic dyscrasia (including leukopenia, neutropenia, thrombocytopenia)
- Hepatic enzyme changes
- Histologic changes in vessel wall, including inflammation, cellular damage, or necrosis
- Myalgia/Arthralgia
- Myelosuppression
- Peripheral neuropathy

## PRE- AND POST-PROCEDURE ANTIPLATELET REGIMEN

In the Zilver PTX randomized clinical trial, the recommended antiplatelet therapy was clopidogrel or ticlopidine starting at least 24 hours prior to the procedure, or a loading dose during the procedure. The most common loading dose was 300 mg clopidogrel or 200 mg ticlopidine. It was recommended that the patient remain on clopidogrel or ticlopidine for 60 days post-procedure along with aspirin indefinitely. Nearly all patients were discharged taking aspirin (dose ranged from 81 to 325 mg) and clopidogrel (most common dose was 75 mg) or ticlopidine (most common dose was 200 mg twice daily) and approximately 90% of patients continued to take aspirin from 1 month through 2 years. Usage of clopidogrel or ticlopidine was lower, with approximately 90% of patients remaining on this medication at 1 month, 70% at 3 months, and 60% through 2 years.

It is recommended that patients receive aspirin as well as clopidogrel or ticlopidine daily. Aspirin should be administered indefinitely; the optimal duration of antiplatelet therapy, specifically clopidogrel or ticlopidine, is not known and the regimens recommended and used in the clinical study are described in the previous paragraph. Stent thrombosis may still occur despite continued therapy.

Prior to peripheral intervention, if a surgical or dental procedure is anticipated that requires early discontinuation of antiplatelet therapy, the treating physician and patient should carefully consider whether a DES and its associated recommended antiplatelet therapy is the appropriate treatment option. Following peripheral intervention, should a surgical or dental procedure be recommended that requires suspension of antiplatelet therapy, the risks and benefits of the procedure should be weighed against the possible risk associated with premature discontinuation of antiplatelet therapy.

## DRUG INFORMATION

### Mechanism of Action

The mechanism by which a Zilver PTX Drug-Eluting Peripheral Stent affects neointimal production has not been fully established. Paclitaxel is known to bind to microtubules and inhibit their molecular disassembly into tubulin, thus arresting mitosis. This action can prevent the smooth muscle cell proliferation and migration known to occur during the restenotic process in arteries. Several studies in animal models have shown that paclitaxel applied locally reduces restenosis by inhibiting smooth muscle cell proliferation and neointimal hyperplasia.

### Pharmacokinetics

The pharmacokinetics are described in the section entitled **SUMMARY OF CLINICAL INVESTIGATIONS** in these Instructions for Use.

### Drug Interactions

Because systemic levels of paclitaxel post-stent placement in clinical trials were extremely low and rapidly cleared ( $C_{max} < 10$  ng/mL;  $< 1$  ng/mL after 8 hours), possible interactions of paclitaxel with concomitantly administered medications are unlikely to be detectable. The metabolism of paclitaxel is catalyzed by CYP2C8 and CYP3A4 which are cytochrome P450 isoenzymes. Formal drug interaction studies have not been conducted with the Zilver PTX Drug-Eluting Peripheral Stent. In the absence of formal clinical drug interaction studies, caution should be exercised when administering paclitaxel concomitantly with known substrates or inhibitors of the cytochrome P450 isoenzymes CYP2C8 and CYP3A4.

### Carcinogenicity, Genotoxicity, and Reproductive Toxicology

Formal long-term carcinogenicity testing of paclitaxel in animals has not been performed. The mechanism of action of the paclitaxel drug is microtubule stabilization and subsequent cell growth inhibition. One consequence is the potential loss of whole chromosomes via interactions with spindle microtubules during cell division. Due to this, paclitaxel carries a risk of chromosomal disruption and therefore is classified as an aneugen. It is not known whether paclitaxel has a separate direct action on DNA in the generation of DNA strand breaks or fragments. *In vivo* and *in vitro* micronucleus genotoxicity studies which are capable of detecting DNA fragments have demonstrated a dose-dependent genotoxic potential of paclitaxel, whereas *in vitro* assays for gene mutation (bacterial reverse mutation and chromosomal aberration test) showed no evidence of genotoxicity. *In vivo*, no genotoxic effects were observed at a dose of 1 mg/kg; however, effects were noted at 50 mg/kg. A dose of 50 mg/kg is approximately 1000 times greater than the total 3 mg paclitaxel amount (adjusted for a 70 kg adult) available from the maximum recommended number of Zilver PTX Drug-Eluting Peripheral Stents for a patient. In addition, the total quantity of paclitaxel with the maximum recommended number of Zilver PTX Drug-Eluting Peripheral Stents provides a total delivered quantity that is 20 times lower than a dose reported to have no observed genotoxic effect in animals (1 mg/kg). *In vivo* reproductive toxicity of paclitaxel has previously been evaluated in rabbits and rats. Although no teratogenic effects were observed in rabbits, impairment of fertility was noted at 1 mg/kg of paclitaxel in rats<sup>2</sup>. This daily dose is approximately 55 times greater than the total 1273 µg paclitaxel quantity from the single largest Zilver PTX Drug-Eluting Peripheral Stent (adjusted for a 70 kg adult).

<sup>1</sup>Tinwell H and Ashby J. Genetic toxicity and potential carcinogenicity of taxol. *Carcinogenesis* 1994;15(8):1499-1501.

<sup>2</sup>Physician's Desk Reference, 52<sup>nd</sup> edition, *Medical Economics Company Inc.*, 1998, pp 762-766.

## USE IN SPECIAL POPULATIONS

### Pregnancy and Lactation

Pregnancy Category C: There are no adequate and well controlled studies in pregnant women of paclitaxel or Zilver PTX Drug-Eluting Peripheral Stents. The Zilver PTX Drug-Eluting Peripheral Stent is contraindicated in women who are pregnant, breastfeeding, or plan to become pregnant in the next 5 years. It is unknown whether paclitaxel will be excreted in human milk, and there is a potential for adverse reaction in nursing infants from paclitaxel exposure.

### Gender

Clinical studies of the Zilver PTX Drug-Eluting Peripheral Stent did not include formal analysis of differences in safety and effectiveness between male and female patients.

### Ethnicity

Clinical studies of the Zilver PTX Drug-Eluting Peripheral Stent did not include sufficient numbers of patients to assess for differences in safety and effectiveness due to ethnicity, either by individual category or when grouped by Caucasian and non-Caucasian.

### Pediatric Use

The safety and effectiveness of the Zilver PTX Drug-Eluting Peripheral Stent in pediatric patients have not been established.

### Geriatric Use

Clinical studies of the Zilver PTX Drug-Eluting Peripheral Stent did not have an upper age limit. In the pivotal study, there were 148 patients in the primary Zilver PTX treatment group who were age 65 or older, and 511 in the single arm study. There were 25 primary Zilver PTX Drug-Eluting Peripheral Stent patients in the pivotal study who were age 80 or older, and 71 in the single arm study. Clinical studies of the Zilver PTX Drug-Eluting Peripheral Stent did not include formal analysis of differences in safety and effectiveness between patients under 65 and over 65 years of age.

## SUMMARY OF CLINICAL INVESTIGATIONS

### LATE MORTALITY SIGNAL FOR PACLITAXEL-COATED DEVICES

A meta-analysis of randomized controlled trials published in December 2018 by Katsanos K, et al.<sup>1</sup> identified an increased risk of late mortality at 2 years and beyond for paclitaxel-coated balloons and paclitaxel-eluting stents used to treat femoropopliteal arterial disease. In response to these data, FDA performed a patient-level metaanalysis of long-term follow-up data from the pivotal premarket randomized trials of paclitaxel-coated devices used to treat femoropopliteal disease using available clinical data through May 2019. The meta-analysis also showed a late mortality signal in study subjects treated with paclitaxel-coated devices compared to patients treated with uncoated devices. Specifically, in the 3 randomized trials with a total of 1090 patients and available 5-year data, the crude mortality rate was 19.8% (range 15.9% - 23.4%) in patients treated with paclitaxel-coated devices compared to 12.7% (range 11.2% - 14.0%) in subjects treated with uncoated devices. The relative risk for increased mortality at 5 years was 1.57 (95% confidence interval 1.16 - 2.13), which corresponds to a 57% relative increase in mortality in patients treated with paclitaxel-coated devices. As presented at the June 2019 FDA Advisory Committee Meeting, an independent meta-analysis of similar patient-level data provided by VIVA Physicians, a vascular medicine organization, reported similar findings with a hazard ratio of 1.38 (95% confidence interval 1.06 - 1.80). Additional analyses have been conducted and are underway that are specifically designed to assess the relationship of mortality to paclitaxel-coated devices.

The presence and magnitude of the late mortality risk should be interpreted with caution because of multiple limitations in the available data, including wide confidence intervals due to a small sample size, pooling of studies of different paclitaxel-coated devices that were not intended to be combined, substantial amounts of missing study data, no clear evidence of a paclitaxel dose effect on mortality, and no identified pathophysiologic mechanism for the late deaths.

Paclitaxel-coated balloons and stents improve blood flow to the legs and decrease the likelihood of repeat procedures to reopen blocked blood vessels compared to uncoated devices. The benefits of paclitaxel-coated devices (e.g., reduced reinterventions) should be considered in individual patients along with potential risks (e.g., late mortality).

In the Zilver PTX randomized clinical study, Kaplan Meier mortality estimates at 2, 3 and 5 years are 5.9% (95% confidence interval 3.3% - 8.6%), 11.3% (7.7% - 14.9%) and 20.6% (15.2% - 26.0%), respectively, for patients treated at the time of the procedure with the Zilver PTX treatment device and 4.7% (1.5% - 7.8%), 8.8% (4.5% - 13.1%) and 14.7% (8.6% - 20.7%), respectively, for patients treated at the time of the procedure with the PTA/BMS control device. The Zilver PTX randomized controlled trial protocol also specified crossover to treatment with the Zilver PTX stent for patients who experienced PTA failure within the first year of the study. Kaplan Meier mortality estimates at 2, 3 and 5 years are 5.4% (95% confidence interval 3.0% - 7.9%), 10.3% (7.0% - 13.6%) and 19.1% (14.2% - 24.1%), respectively, for all patients treated with the Zilver PTX treatment device and 5.7% (1.8% - 9.5%), 10.7% (5.6% - 15.9%) and 17.1% (10.0% - 24.2%), respectively, for patients treated only with the PTA/BMS control device. Additional information regarding long-term outcomes can be found in the **RANDOMIZED CLINICAL STUDY** section.

<sup>1</sup> Katsanos K, et al. Risk of death following application of paclitaxel-coated balloons and stents in the femoropopliteal artery of the leg: A systematic review and meta-analysis of randomized controlled trials. *J Am Heart Assoc.* 2018;7(24):e011245.

## OVERVIEW OF CLINICAL STUDIES

The Zilver PTX Drug-Eluting Peripheral Stent has been the subject of a multicenter randomized clinical study (IDE #G030251), a multicenter single arm clinical study and a supporting multicenter randomized clinical study to confirm deployability and safety of 100 and 120 mm stent lengths. As detailed below, results from the randomized IDE study support the safety and effectiveness of the Zilver PTX Drug-Eluting Peripheral Stent. Specifically, the event-free survival rate for the Zilver PTX treatment group was non-inferior (i.e., equivalent or superior) to the PTA control group. Moreover, the rate of event-free survival in the Zilver PTX group was significantly higher than in the PTA group, suggesting that primary stenting with the Zilver PTX Drug-Eluting Peripheral Stent is associated with a lower MAE rate than the current standard care of PTA with provisional stenting. The primary patency rate (patency defined as duplex ultrasound measured peak systolic velocity ratio < 2.0 indicating less than 50% diameter stenosis) for the Zilver PTX Drug-Eluting Peripheral Stent was significantly higher than the primary patency rate for PTA, demonstrating that primary stenting with the Zilver PTX Drug-Eluting Peripheral Stent is significantly more effective than PTA. Additionally, the Zilver PTX Drug-Eluting Peripheral Stent demonstrated superior effectiveness to the bare Zilver stent in a randomized comparison in lesions with acute PTA failure. Results from the single arm clinical study provide additional evidence supporting the safety of the Zilver PTX Drug-Eluting Peripheral Stent in a broader patient population including more complex lesions. Results from the supporting clinical study of 100 and 120 mm stents confirmed the safety and deployability of Zilver PTX Drug-Eluting Peripheral Stent.

## RANDOMIZED CLINICAL STUDY

### Primary Objective

The primary objective of the Zilver PTX randomized study was to demonstrate the safety and effectiveness of the Zilver PTX Drug-Eluting Peripheral Stent compared to percutaneous balloon angioplasty (PTA) for the treatment of *de novo* or restenotic lesions of the above-the-knee femoropopliteal artery. The patients in the PTA control group included those with optimal PTA and suboptimal (failed) PTA that underwent a secondary randomization to stenting with either Zilver PTX or bare Zilver stents.

### Study overview

The Zilver PTX randomized study is a prospective, controlled, multi-center, multinational study enrolling patients in the United States, Japan, and Germany with *de novo* or restenotic native lesions of the above-the-knee femoropopliteal artery. Patients were randomized 1:1 to treatment with the Zilver PTX Drug-Eluting Peripheral Stent (treatment group) or with PTA (control group). Recognizing that balloon angioplasty may not be successful acutely, the trial design mandated provisional stent placement immediately after failure of balloon angioplasty in instances of acute PTA failure. Therefore, patients with suboptimal (failed) PTA underwent a secondary randomization (1:1) to stenting with either Zilver PTX or bare Zilver stents (Figure 2). This secondary randomization allows evaluation of the Zilver PTX Drug-Eluting Peripheral Stent compared to a bare metal stent.


Figure 2. Patient enrollment

The study was designed to assess primary safety and effectiveness hypotheses regarding the Zilver PTX treatment group compared to the PTA control group. Specifically, the primary safety hypothesis was non-inferior (i.e., equivalent or superior) event-free survival (defined as freedom from the clinical events committee (CEC) adjudicated major adverse events of death, target lesion revascularization, target limb ischemia requiring surgical intervention or surgical repair of the target vessel, and freedom from worsening of the Rutherford classification by 2 classes or to class 5 or 6) at 12 months. The primary effectiveness hypothesis was superior primary patency at 12 months for the Zilver PTX treatment group compared to the PTA control group. Secondary analyses included evaluation of the effectiveness of the Zilver PTX Drug-Eluting Peripheral Stent compared to a bare metal stent.

Subjects eligible to be enrolled in the study had single or bilateral stenotic or occluded atherosclerotic lesions ( $\leq 14$  cm long) of the above-the-knee femoropopliteal artery with a reference vessel diameter of 4 mm to 9 mm and Rutherford classification 2 to 6. Fifty-five (55) institutions enrolled 479 patients, including 238 in the PTA control group and 241 in the Zilver PTX treatment group. Five patients in the Zilver PTX group were enrolled as live cases (i.e., with no randomization) and are included in analyses of the as-treated population but not the intent-to-treat or per-protocol populations. Acute PTA failure was common, occurring in 120 patients in the control group, and these patients underwent a second randomization to provisional stenting with either Zilver PTX Drug-Eluting Peripheral Stents or bare Zilver stents.

Follow-up included clinical assessment and ultrasound imaging prior to discharge, at 6 and 12 months, and annually thereafter. Additionally, x-rays were required prior to discharge and at 1, 3, and 5 years to assess stent integrity. Telephone contact was scheduled for 1, 3, 9, and 18 months. Patient subsets were assigned to a pharmacokinetic substudy and to an IVUS/angiography substudy.

The study was overseen by an independent data safety monitoring board (DSMB) comprised of physicians and a biostatistician. An independent CEC adjudicated major adverse events, including all patient deaths, and independent core laboratories provided uniformly defined imaging analysis.

### Demographics

Baseline patient characteristics, angiographic data, lesion location, and lesion characteristics were generally similar between the PTA control group and Zilver PTX treatment group with a more frequent history of hypertension and more severe calcification and inflow tract stenosis in the Zilver PTX Drug-Eluting Peripheral Stent treatment group (Tables 3 – 7).

Table 3. Demographics

Demographic	Control Group	Treatment Group <sup>3</sup>	Diff. (95% CI) <sup>1</sup>	P-value <sup>2</sup>
Age (years)	67.7 $\pm$ 10.6 (238) (28 – $\geq$ 90)	67.9 $\pm$ 9.6 (236) (39 – 87)	-0.1 (-2.0, 1.7)	0.88
Gender				0.70
Male	63.9% (152/238)	65.7% (155/236)	-1.8 (-10.4, 6.8)	
Female	36.1% (86/238)	34.3% (81/236)	1.8 (-6.8, 10.4)	
Ethnicity				0.81
Asian	14.1% (29/206)	11.9% (25/210)	2.2 (-4.3, 8.6)	
Black/African American	11.2% (23/206)	11.9% (25/210)	-0.7 (-6.9, 5.4)	
Hispanic/Latino	5.3% (11/206)	7.1% (15/210)	-1.8 (-6.5, 2.8)	
White/Caucasian	69.4% (143/206)	69.0% (145/210)	0.4 (-8.5, 9.2)	
Height (in)	66.4 $\pm$ 4.4 (238)	66.7 $\pm$ 3.6 (236)	-0.2 (-1.0, 0.5)	0.55
Weight (lbs)	178.5 $\pm$ 44.3 (238)	180.4 $\pm$ 40.0 (236)	-1.9 (-9.5, 5.8)	0.62
Body mass index	28.2 $\pm$ 5.6 (238)	28.4 $\pm$ 5.3 (236)	-0.2 (-1.2, 0.8)	0.71

<sup>1</sup> Confidence interval is the difference in means for continuous variables and difference in percentages for categorical variables.

<sup>2</sup> P-values are based on t-test for continuous variables and Fisher's exact test for categorical variables.

<sup>3</sup> Five patients treated as live cases not included

Table 4. Medical history

Condition	Control Group	Treatment Group <sup>4</sup>	Diff. (95% CI) <sup>1</sup>	P-value <sup>2</sup>
Diabetes	42.0% (100/238)	49.6% (117/236)	-7.6 (-16.5, 1.4)	0.11
Diabetes type				0.56
Type I	13.0% (13/100)	16.2% (19/117)	-3.2 (-12.6, 6.2)	
Type II	87.0% (87/100)	83.8% (98/117)	3.2 (-6.2, 12.6)	
Hypercholesterolemia	69.7% (166/238)	76.3% (180/236)	-6.5 (-14.5, 1.5)	0.12
Hypertension	81.5% (194/238)	89.0% (210/236)	-7.5 (-13.8, -1.1)	0.02*
Carotid disease	20.2% (48/238)	18.2% (43/236)	2.0 (-5.1, 9.0)	0.64
Renal disease	10.5% (25/238)	10.2% (24/236)	0.3 (-5.2, 5.8)	> 0.99
Pulmonary disease	16.0% (38/238)	19.1% (45/236)	-3.1 (-9.9, 3.7)	0.39
Congestive heart failure	10.5% (25/238)	11.9% (28/236)	-1.4 (-7.0, 4.3)	0.66

Condition	Control Group	Treatment Group <sup>4</sup>	Diff. (95% CI) <sup>1</sup>	P-value <sup>2</sup>
Previous cardiac arrhythmia	13.0% (31/238)	10.6% (25/236)	2.4 (-3.4, 8.2)	0.47
Previous MI	17.2% (41/238)	21.2% (50/236)	-4.0 (-11.0, 3.1)	0.29
Smoking status				0.70
Never smoked	15.5% (37/238)	13.6% (32/236)	2.0 (-4.4, 8.3)	
Quit	51.7% (123/238)	55.5% (131/236)	-3.8 (-12.8, 5.1)	
Still smokes	32.4% (77/238)	30.9% (73/236)	1.4 (-7.0, 9.8)	
Unknown	0.4% (1/238)	0.0% (0/236)	0.4 (N/A)	
Existing tissue loss <sup>3</sup>	8.4% (20/238)	9.4% (22/235)	-1.0 (-6.1, 4.2)	0.74
Rutherford Classification				0.59
Class 1	0.8% (2/236)	0.8% (2/236)	N/A	
Class 2	46.2% (109/236)	52.5% (124/236)		
Class 3	44.5% (105/236)	37.7% (89/236)		
Class 4	4.7% (11/236)	5.9% (14/236)		
Class 5	3.4% (8/236)	3.0% (7/236)		
Class 6	0.4% (1/236)	0.0% (0/236)		
Currently taking medications	99.2% (236/238)	99.6% (235/236)	-0.4 (-1.8, 1.0)	> 0.99

<sup>1</sup> Confidence interval is the difference in means for continuous variables and difference in percentages for categorical variables.

<sup>2</sup> P values are based on t-test for continuous variables and Fisher's exact test for categorical variables.

<sup>3</sup> Tissue loss includes amputations, gangrene, and ischemic ulcers.

<sup>4</sup> Five patients treated as live cases were not included

\* Statistically significant

**Table 5. Baseline angiographic data (core lab reported)**

Baseline Angiographic Data	Control Group	Treatment Group	Diff. (95% CI) <sup>1</sup>	P-value <sup>2</sup>
Stenosed lesion length (mm) <sup>3</sup>	53.2 ± 40.3 (248) (5.8 – 220.3)	54.6 ± 40.7 (242) (5.9 – 262.3)	-1.3 (-8.5, 5.9)	0.71
Normal-to-normal Lesion length (mm) <sup>4</sup>	63.2 ± 40.5 (251)(2.0 – 140.0)	66.4 ± 38.9 (246) (5.5 – 140.0)	-3.2 (-10.2, 3.8)	0.36
Proximal RVD (mm)	5.0 ± 1.0 (249) (2.7 – 8.1)	5.1 ± 0.9 (242) (2.5 – 7.9)	-0.05 (-0.2, 0.1)	0.58
Distal RVD (mm)	5.0 ± 1.0 (249) (2.9 – 10.4)	5.0 ± 1.0 (242) (2.9 – 8.0)	-0.01 (-0.2, 0.2)	0.95
MLD in lesion (mm)	1.1 ± 0.9 (249) (0.0 – 4.1)	1.0 ± 0.9 (242) (0.0 – 3.1)	0.1 (-0.1, 0.2)	0.38
Percent diameter stenosis (%)	78.4 ± 17.1 (249) (31.5 – 100.0)	79.8 ± 17.0 (242) (47.1 – 100.0)	-1.3 (-4.4, 1.7)	0.38

<sup>1</sup> Confidence interval is the difference in means for continuous variables and difference in percentages for categorical variables.

<sup>2</sup> P values are based on t-test for continuous variables and Fisher's exact test for categorical variables.

<sup>3</sup> Region with > 20% diameter stenosis.

<sup>4</sup> Site reported diseased lesion length.

**Table 6. Lesion location**

Vessel <sup>1</sup>	Control Group	Treatment Group	Diff. (95% CI) <sup>2</sup>	P-value <sup>3</sup>
Left proximal SFA	10.8% (27/251)	8.9% (22/247)	2.9 (-4.3, 10.1)	0.63
Right proximal SFA	12.0% (30/251)	10.9% (27/247)		
Left proximal SFA/distal SFA	3.6% (9/251)	5.7% (14/247)	-2.1 (-6.7, 2.5)	
Right proximal SFA/distal SFA	2.8% (7/251)	2.8% (7/247)		
Left distal SFA	34.7% (87/251)	30.4% (75/247)	-1.0 (-9.5, 7.4)	
Right distal SFA	28.7% (72/251)	34.0% (84/247)		
Left distal SFA/popliteal artery	0.4% (1/251)	0.8% (2/247)	-1.3 (-4.3, 1.8)	
Right distal SFA/popliteal artery	2.0% (5/251)	2.8% (7/247)		
Left popliteal artery	2.0% (5/251)	0.8% (2/247)	1.5 (-2.1, 5.1)	
Right popliteal artery	3.2% (8/251)	2.8% (7/247)		

<sup>1</sup> Bilateral lesions were treated in 13 patients in the control group and 11 patients in the treatment group.

<sup>2</sup> Confidence interval is the difference in means for continuous variables and difference in percentages for categorical variables.

<sup>3</sup> P values are based on t-test for continuous variables and Fisher's exact test for categorical variables.

**Table 7. Lesion characteristics**

Characteristics	Control Group	Treatment Group	Diff. (95% CI) <sup>2</sup>	P-value <sup>3</sup>	
Lesion class (TASC) <sup>1</sup>	A	36.0% (86/239)	29.4% (69/235)	6.6 (-1.8, 15.04)	0.07
	B	25.9% (62/239)	22.6% (53/235)	3.4 (-4.3, 11.1)	
	C	31.0% (74/239)	42.6% (100/235)	-11.6 (-20.2, -3.0)	
	D	7.1% (17/239)	5.5% (13/235)	1.6 (-2.8, 6.0)	
Accessibility	Readily accessible	100% (215/215)	100% (215/215)	0 (0, 0)	N/A
	Moderate tortuosity	0.0% (0/215)	0.0% (0/215)	N/A	
	Excessive tortuosity	0.0% (0/215)	0.0% (0/215)	N/A	
Lesion angulation	Non-angulated	95.2% (237/249)	95.4% (228/241)	-0.3 (-4.0, 3.5)	> 0.99
	Moderate	4.8% (12/249)	4.6% (11/241)	0.3 (-3.5, 4.0)	

Characteristics		Control Group	Treatment Group	Diff. (95% CI) <sup>2</sup>	P-value <sup>3</sup>
Calcification	None	4.8% (12/249)	1.7% (4/241)	3.2 (0.05, 6.3)	< 0.01*
	Little	38.2% (95/249)	25.7% (62/241)	12.4 (4.3, 20.6)	
	Moderate	22.1% (55/249)	35.3% (85/241)	-13.2 (-21.1, -5.3)	
	Severe	34.9% (87/249)	37.3% (90/241)	-2.4 (-10.9, 6.1)	
Other stenosis in artery	None	51.4% (111/216)	51.7% (107/207)	-0.3 (-9.8, 9.2)	0.71
	≤ 50%	34.3% (74/216)	31.4% (65/207)	2.9 (-6.1, 11.8)	
	> 50%	14.4% (31/216)	16.9% (35/207)	-2.6 (-9.5, 4.4)	
Inflow tract stenosis	None	41.6% (96/231)	37.1% (76/205)	4.5 (-4.7, 13.7)	0.03*
	≤ 50%	45.5% (105/231)	40.5% (83/205)	5.0 (-4.3, 14.3)	
	> 50%	13.0% (30/231)	22.4% (46/205)	-9.5 (-16.6, -2.3)	
Patent runoff vessels	0	17.3% (26/150)	14.8% (22/149)	2.6 (-5.8, 10.9)	0.47
	1	52.7% (79/150)	47.7% (71/149)	5.0 (-6.3, 16.3)	
	2	21.3% (32/150)	22.8% (34/149)	-1.5 (-10.9, 7.9)	
	3	8.0% (12/150)	14.1% (21/149)	-6.1 (-13.2, 1.0)	
Ulceration		19.0% (47/248)	16.7% (40/240)	2.3 (-4.5, 9.1)	0.55
Total Occlusion		27.4% (68/248)	32.8% (79/241)	-5.4 (-13.5, 2.8)	0.20

<sup>1</sup> TASC lesion class was determined by the site and was not evaluated by the core lab.

<sup>2</sup> Confidence interval is the difference in means for continuous variables and difference in percentages for categorical variables.

<sup>3</sup> P values are based on t-test for continuous variables and Fisher's exact test for categorical variables.

\* Statistically significant.

## Results

Patient availability for study follow-up through 24 months is summarized in **Table 8**.

**Table 8. Clinical and imaging follow-up data**

Follow-up	Eligible for follow-up <sup>1</sup>	Percent of Data Available			Events Occurring Before Next Visit			
		Clinical Follow-up <sup>2</sup>	Core Laboratory X-ray Follow-up	Core Laboratory Ultrasound Follow-up <sup>3</sup>	Death	Withdrawn	Lost to Follow-up	Other Endpoint <sup>4</sup>
<b>PTA Control Group</b>								
Procedure	238	100.0% (238/238)	97.5% (117/120) <sup>5</sup>	86.6% (206/238)	2	3	0	1
6-month	232	94.0% (218/232)	76.5% (13/17) <sup>5,6</sup>	80.6% (187/232)	2	4	2	2
12-month	222	97.3% (216/222)	76.0% (114/150) <sup>5</sup>	83.3% (185/222)	4	5	0	3
2-year	210	83.3% (175/210)	n/a	66.7% (84/126)	n/a			
<b>Zilver PTX Treatment Group</b>								
Procedure	236 <sup>7</sup>	100.0% (236/236)	96.2% (227/236)	89.4% (211/236)	0	1	1	0
6-month	234	94.0% (220/234)	86.2% (25/29) <sup>6</sup>	80.3% (188/234)	9	4	4	0
12-month	217	97.7% (212/217)	84.3% (183/217)	87.1% (189/217)	9	4	3	3
2-year	198	83.8% (166/198)	n/a	51.0% (101/198)	n/a			

<sup>1</sup> Eligible for follow-up = previous eligibility for follow-up – (previous death + withdrawn + LTF).

<sup>2</sup> Includes cases with at least one of any of the following submitted: clinical form, death form, withdrawn form, or lost to follow-up form.

<sup>3</sup> Includes only ultrasound studies considered diagnostic by the core lab.

<sup>4</sup> Patients who reached an "other endpoint" include 2 patients who received a non-study stent during reintervention and 1 patient who moved but has not formally withdrawn from the trial.


<sup>5</sup> Only patients implanted with stents (i.e., acute PTA failure) were required to have x-ray follow-up.

<sup>6</sup> Only first 60 patients enrolled were required to have 6-month x-ray follow-up.

<sup>7</sup> Five patients treated as live cases not included.

## Safety

The primary safety endpoint of non-inferior (i.e., equivalent or superior) safety for the Zilver PTX Drug-Eluting Peripheral Stent compared to PTA (including provisional stenting (bare and PTX-coated) in the PTA group) was met with an event-free survival rate at 12 months of 90.4% for the Zilver PTX treatment group and 83.9% for the PTA control group ( $p < 0.01^3$ ), as illustrated in **Figure 3** and **Table 9**. The most common major adverse event was TLR (**Table 10**), which occurred approximately 70% more often in the PTA group relative to the Zilver PTX group (16.3% vs. 9.6%, respectively;  $p = 0.04$ ). No patient deaths were adjudicated by the CEC as related to the device or procedure. Non-inferiority of the Zilver PTX Drug-Eluting Peripheral Stent was maintained through 24 months.

**Figure 3. Kaplan-Meier curve for event-free survival**

<sup>3</sup> Adjusted for multiplicity.

**Table 9. Kaplan-Meier estimates for event-free survival at 0, 1, 6, 12, and 24 months**

Months Post-procedure	Event-free Survival Estimate		Standard Error		Cumulative Failed		Cumulative Censored		Number Remaining	
	Control (PTA)	Treatment (Zilver PTX)	Control (PTA)	Treatment (Zilver PTX)	Control (PTA)	Treatment (Zilver PTX)	Control (PTA)	Treatment (Zilver PTX)	Control (PTA)	Treatment (Zilver PTX)
0	100.0%	100.0%	0.0%	0.0%	0	0	0	0	236	235
1	100.0%	99.1%	0.0%	0.6%	0	2	0	0	236	233
6	94.4%	97.0%	1.5%	1.1%	13	7	6	3	217	225
12	83.9%	90.4%	2.4%	1.9%	37	22	15	16	184	197
24	77.9%	86.6%	2.8%	2.3%	50	30	22	33	164	172

**Table 10. Rates of individual major adverse events at 12 months**

Major Adverse Event	Control (PTA) <sup>1</sup>	Treatment (Zilver PTX) <sup>1</sup>	P-value	Diff. (95% CI) <sup>3</sup>
Clinically-driven TLR	16.3% (36/221) <sup>2</sup>	9.6% (21/219)	0.04	(0.5, 12.9)
Worsening of Rutherford classification by 2 classes or to a class 5 or 6	0.9% (2/221) <sup>2</sup>	0.0% (0/219)	0.49	-
Amputation	0.0% (0/221)	0.5% (1/219)	0.49	-

<sup>1</sup> Denominator is the number of patients remaining free from MAE at 12 months plus the number that have experienced a MAE or died prior to 12 months.

<sup>2</sup> One patient experienced a worsening Rutherford and a TLR and is included in both categories in this table.

<sup>3</sup> Confidence interval is the difference in percentages.

**Table 11** provides a summary of other adverse events through 24 months, not including major adverse events, occurring in the control and treatment groups. No patient deaths were adjudicated by the CEC as related to the study device or procedure.

**Table 11. Adverse events**

Event Type	Control		Treatment	
	% (n) N = 236 patients	Events	% (n) N = 235 patients <sup>1</sup>	Events
<b>Occurring within 12 months of the study procedure</b>				
Death	1.7% (4)	4	3.4% (8)	8
<b>Major Adverse Events</b>				
Clinically-driven TLR	15.3% (36)	36	8.9% (21)	21
Worsening of Rutherford classification by 2 classes or to a class 5 or 6	0.8% (2)	2	0.0% (0)	0
Amputation	0.0% (0)	0	0.4% (1)	1
<b>Cardiovascular Events</b>				
Cardiac ischemia requiring intervention	3.0% (7)	7	3.8% (9)	13
Non-Q-Wave MI	0.4% (1)	1	0.9% (2)	2
Congestive heart failure	0.8% (2)	2	3.0% (7)	7
Refractory hypertension	0.4% (1)	1	0.0% (0)	0
Arrhythmia requiring intervention or new treatment	1.7% (4)	5	3.4% (8)	9
Other cardiovascular events	12.3% (29)	41	7.7% (18)	20
<b>Pulmonary Events</b>				
Pulmonary edema requiring treatment	0.0% (0)	0	0.4% (1)	1
Ventilation greater than 24 hours in duration	0.4% (1)	1	1.3% (3)	3
Pneumonia requiring antibiotics	2.1% (5)	6	5.1% (12)	14
Supplemental oxygen at time of discharge (exclude if for high altitude)	0.0% (0)	0	0.4% (1)	1
Re-intubation	0.0% (0)	0	0.9% (1)	1
COPD	1.7% (4)	5	0.9% (2)	4
Other pulmonary events	3.0% (7)	11	8.5% (20)	28
<b>Renal Events</b>				
UTI requiring antibiotic treatment	0.8% (2)	3	2.6% (6)	7
Serum creatinine rise greater than 30% above baseline resulting in persistent value greater than 2 mg/dl	0.8% (2)	2	0.0% (0)	0
Other renal events	2.5% (6)	8	4.3% (10)	12
<b>Gastrointestinal Events</b>				
Other gastrointestinal events	3.8% (9)	12	10.2% (24)	30
<b>Wound Events</b>				
Wound infection/abscess formation	1.3% (3)	3	1.7% (4)	5
Tissue necrosis requiring debridement	1.3% (3)	3	0.4% (1)	1
Wound complication requiring return to operating room	1.7% (4)	5	0.9% (2)	3
Other wound events	2.5% (6)	6	0.9% (2)	2
<b>Vascular Events</b>				
Post-procedure percutaneous intervention (e.g., PTA and/or stent) to the study vessel	8.1% (19)	22	5.5% (13)	14
Post-procedure percutaneous intervention (e.g., PTA and/or stent) to another vessel	16.9% (40)	52	21.3% (50)	62
Ischemia requiring surgical intervention (i.e., bypass or amputation) of another vessel	2.5% (6)	6	1.7% (4)	4
Embolism distal to treated study vessel	0.4% (1)	1	0.9% (2)	2
Embolism within another vessel	0.0% (0)	0	1.3% (3)	3
Thrombosis of the study lesion	0.8% (2)	2	2.6% (6)	6
Thrombosis other than study lesion	0.0% (0)	0	0.9% (2)	2
Blue toe syndrome	0.0% (0)	0	0.4% (1)	1
Aneurysm (other)	0.0% (0)	0	0.9% (2)	2
Deep vein thrombosis	0.4% (1)	1	0.9% (2)	2
Hematoma requiring intervention at access site	0.4% (1)	1	0.4% (1)	1
Pulmonary embolism	0.4% (1)	1	0.0% (0)	0
Pseudoaneurysm or AV fistula of the study vessel	0.8% (2)	2	0.4% (1)	1
Pseudoaneurysm or AV fistula of another vessel	0.4% (1)	1	0.4% (1)	2
Study vessel spasm	0.0% (0)	0	0.4% (1)	1
Worsened claudication/rest pain	6.8% (16)	21	4.3% (10)	12
Stroke	0.8% (2)	2	0.4% (1)	1
Vascular/surgical repair of injury to another vessel (other than amputation or bypass)	0.4% (1)	1	0.9% (2)	3
Post-procedure transfusion	4.2% (10)	15	5.5% (13)	13
Other vascular events	16.5% (39)	47	9.8% (23)	29


Event Type	Control		Treatment	
	% (n) N = 236 patients	Events	% (n) N = 235 patients <sup>1</sup>	Events
<b>Miscellaneous Events</b>				
Drug reaction (including contrast reaction)	2.1% (5)	5	3.8% (9)	10
Hypersensitivity/allergic reaction	2.1% (5)	5	2.1% (5)	6
Other miscellaneous events	26.7% (63)	120	28.1% (66)	115
<b>Occurring between 12 and 24 months following the study procedure</b>				
Death	1.7% (4)	4	4.3% (10)	10
<b>Major Adverse Events</b>				
Clinically-driven TLR	4.7% (11)	11	3.4% (8)	8
Worsening of Rutherford classification by 2 classes or to a class 5 or 6	1.7% (4)	4	0.0% (0)	0
<b>Cardiovascular Events</b>				
Cardiac ischemia requiring intervention	3.0% (7)	7	3.4% (8)	8
Non-Q-Wave MI	0.8% (2)	2	0.4% (1)	1
Congestive heart failure	1.3% (3)	4	1.3% (3)	6
Refractory hypertension	0.4% (1)	1	0.0% (0)	0
Arrhythmia requiring intervention or new treatment	0.8% (2)	2	1.7% (4)	5
Other cardiovascular events	7.6% (18)	30	3.8% (9)	10
<b>Pulmonary Events</b>				
Pulmonary edema requiring treatment	0.0% (0)	0	0.4% (1)	1
Pneumonia requiring antibiotics	1.3% (3)	3	1.7% (4)	4
Supplemental oxygen at time of discharge (exclude if for high altitude)	0.0% (0)	0	0.4% (1)	1
COPD	0.8% (2)	2	0.0% (0)	0
Pleural effusion requiring treatment	0.0% (0)	0	0.4% (1)	2
Other pulmonary events	3.0% (7)	7	3.4% (8)	10
<b>Renal Events</b>				
UTI requiring antibiotic treatment	0.8% (2)	2	0.4% (1)	1
Serum creatinine rise greater than 30% above baseline resulting in persistent value greater than 2 mg/dl	0.8% (2)	2	0.0% (0)	0
Other renal events	1.7% (4)	8	2.1% (5)	5
<b>Gastrointestinal Events</b>				
Other gastrointestinal events	3.4% (8)	9	2.6% (6)	6
<b>Wound Events</b>				
Wound infection/abscess formation	0.4% (1)	2	1.3% (3)	3
Wound complication requiring return to operating room	0.4% (1)	1	0.0% (0)	0
Other wound events	0.0% (0)	0	1.3% (3)	3
<b>Vascular Events</b>				
Post-procedure percutaneous intervention (e.g., PTA and/or stent) to the study vessel	5.1% (12)	14	3.8% (9)	9
Post-procedure percutaneous intervention (e.g., PTA and/or stent) to another vessel	5.9% (14)	16	6.0% (14)	16
Ischemia requiring surgical intervention (i.e., bypass or amputation) of another vessel	0.8% (2)	2	1.7% (4)	4
Embolism distal to treated study vessel	0.0% (0)	0	0.4% (1)	1
Embolism within another vessel	0.4% (1)	1	0.0% (0)	0
Thrombosis of the study lesion	0.0% (0)	0	0.4% (1)	1
Blue toe syndrome	0.0% (0)	0	0.4% (1)	1
Deep vein thrombosis requiring surgical or lytic therapy	0.0% (0)	0	0.4% (1)	1
Worsened claudication/rest pain	2.5% (6)	6	3.4% (8)	9
Stroke	1.7% (4)	4	1.3% (3)	3
Vascular/surgical repair of injury to the study vessel (other than amputation or bypass)	0.0% (0)	0	0.4% (1)	2
Vascular/surgical repair of injury to another vessel (other than amputation or bypass)	0.4% (1)	1	0.0% (0)	0
Post-procedure transfusion	0.0% (0)	0	1.7% (4)	5
Other vascular events	4.7% (11)	12	3.4% (8)	10
<b>Miscellaneous Events</b>				
Hypersensitivity/allergic reaction	0.0% (0)	0	0.4% (1)	1
Other miscellaneous events	14.0% (33)	44	15.7% (37)	55

<sup>1</sup> One patient was not treated with a Zilver PTX Stent so is not included in this analysis

**Effectiveness**

The primary effectiveness endpoint of superior primary patency (conservatively defined as a PSV ratio < 2.0) for the Zilver PTX Drug-Eluting Peripheral Stent compared to PTA was met ( $p < 0.01^4$ ) with a primary patency rate at 12 months of 82.7% for the Zilver PTX treatment group and 32.7% for the PTA control group, as illustrated in **Figure 4** and **Table 12**. The benefit of the Zilver PTX Drug-Eluting Peripheral Stent was maintained through 24 months.

<sup>4</sup>Adjusted for multiplicity


**Figure 4. Kaplan-Meier curve for primary patency**

**Table 12. Kaplan-Meier estimates for primary patency at 0, 1, 6, 12, and 24 months**

Months Post-procedure	Primary Patency Estimate		Standard Error		Cumulative Failed		Cumulative Censored		Number Remaining	
	Control (PTA)	Treatment (Zilver PTX)	Control (PTA)	Treatment (Zilver PTX)	Control (PTA)	Treatment (Zilver PTX)	Control (PTA)	Treatment (Zilver PTX)	Control (PTA)	Treatment (Zilver PTX)
0	50.2%	99.6%	3.2%	0.4%	125	1	0	0	126	245
1	50.2%	99.6%	3.2%	0.4%	125	1	0	1	126	244
6	41.6%	95.1%	3.1%	1.4%	146	12	5	4	100	230
12	32.7%	82.7%	3.0%	2.5%	167	41	11	23	73	182
24	26.5%	74.8%	3.1%	2.9%	177	58	41	39	33	149

Comparison of the results for patients treated with the Zilver PTX Drug-Eluting Peripheral Stent to those treated with the bare metal Zilver stent when used in a similar patient population (i.e., those patients who had acute failure of PTA) provides an evaluation of the paclitaxel drug effect. Both patient populations were selected in the same way (and randomized), and both stents have the identical stent platform; therefore, this comparison provides a direct measurement of the effectiveness of the Cook PTX® drug coating on the Zilver PTX Drug-Eluting Peripheral Stent. As illustrated in **Figure 5** and **Table 13**, there was a significant difference in patency outcomes between the groups ( $p < 0.01^5$ ), with the Zilver PTX group exhibiting a higher primary patency rate at 12 months of 90.2% compared to 72.9% for the bare Zilver group. The benefit of the Zilver PTX Drug-Eluting Peripheral Stent was maintained through 24 months with a 24-month patency rate of 83.4% for the Zilver PTX group compared to 64.1% for the bare Zilver group.

<sup>5</sup> Adjusted for multiplicity.

**Figure 5. Kaplan-Meier estimate plot for primary patency for provisional bare Zilver vs. provisional Zilver PTX (ITT)**

**Table 13. Kaplan-Meier estimates for bare Zilver vs. Zilver PTX**

Months Post-procedure	Primary Patency Estimate		Standard Error		Cumulative Failed		Cumulative Censored		Number Remaining	
	Bare Zilver	Zilver PTX	Bare Zilver	Zilver PTX	Bare Zilver	Zilver PTX	Bare Zilver	Zilver PTX	Bare Zilver	Zilver PTX
0	100.0%	100.0%	0.0%	0.0%	0	0	0	0	62	63
1	100.0%	100.0%	0.0%	0.0%	0	0	0	0	62	63
6	88.4%	96.8%	4.1%	2.2%	7	2	2	1	53	60
12	72.9%	90.2%	5.8%	3.8%	16	6	5	3	41	54
24	64.1%	83.4%	6.3%	4.8%	21	10	6	7	35	46

In summary, the primary patency rate at 12 months was 82.7% in the Zilver PTX treatment group and 32.7% in the PTA control group. The effect of covariates, including diabetes, lesion length, and occluded/stenosed lesions, was not significantly different between the Zilver PTX and PTA groups.

In conclusion, the primary effectiveness hypothesis of the study was met, indicating that treatment with the Zilver PTX Drug-Eluting Peripheral Stent is significantly more effective than treatment with PTA. Additionally, results from patient groups within the secondary randomization demonstrated that stenting with the paclitaxel-coated Zilver PTX Drug-Eluting Peripheral Stent is significantly more effective in maintaining primary patency than stenting with the same bare (uncoated) stent, indicating that the PTX® coating has a significant effect and further supporting the effectiveness of the Zilver PTX Drug-Eluting Peripheral Stent.

**Secondary Endpoints**

Secondary endpoint analyses include procedural success, the patient-centered measure of freedom from worsening symptoms of ischemia (i.e., sustained clinical benefit), clinical status (i.e., clinical improvement defined as improvement by 1 Rutherford class and clinical success defined as improvement by 2 Rutherford classes), and functional status (measured by ABL and Walking Impairment Questionnaire).

Procedural success (< 30% residual stenosis) was significantly higher ( $p < 0.01$ ) in the Zilver PTX treatment group (95.0%; 229/241) compared to the PTA control group (57.7%; 143/248), demonstrating that the Zilver PTX Drug-Eluting Peripheral Stent is effective in establishing acute patency.

Sustained clinical benefit (i.e., freedom from worsening symptoms of ischemia) was evaluated in terms of freedom from the patient-centered measures of worsening claudication, worsening Rutherford class, tissue loss, and other symptoms indicating the need for reintervention (e.g., rest pain, ulcer, persistent claudication), and was considered to provide a clinically-based evaluation of patient benefit in this study, without the inclusion of surrogate endpoints. In this evaluation, patients in the Zilver PTX treatment group achieved a higher sustained clinical benefit compared to patients in the PTA control group (86.9% vs. 75.4% at 12 months). Similarly, in subgroups within the PTA control group, patients receiving provisional Zilver PTX Drug-Eluting Peripheral Stents achieved a higher sustained clinical benefit compared to patients receiving optimal PTA (87.3% vs. 69.4% at 12 months) or provisional bare Zilver stents (87.3% vs. 72.3% at 12 months).

**Substudies**

**Pharmacokinetic Substudy**

A subgroup of 60 patients from the Zilver PTX treatment group was included in the pharmacokinetic (PK) substudy to evaluate systemic paclitaxel delivery. Each patient was assigned 3 of a possible 11 time points, which included post-procedure (time 0), 20 min, 40 min, and 1, 1.5, 2, 3, 4, 6, 8, and 12 hours. Patients were divided into two groups based on the number of stents with which they were implanted. The number of patients and total quantity of paclitaxel for each group are shown in **Table 14**.

**Table 14. Pharmacokinetic substudy groups**

# of Stents	# of Patients <sup>1</sup>	# of Samples	Paclitaxel Dose Range (µg)
1	42	125	312 - 864 (mean ± SD = 694 ± 200)
2	16	48	1083 - 1728 (mean ± SD = 1398 ± 228)

<sup>1</sup> Two patients were not included in the analysis; samples from one patient were assayed beyond the known stability timeframe; one patient received three Zilver PTX Drug-Eluting Peripheral Stents.

A parametric curve was fit to the data and the maximum observed plasma paclitaxel concentration ( $C_{max}$ ), time to maximum concentration ( $T_{max}$ ), area under the plasma concentration-time curve (AUC), half-life ( $t_{1/2}$ ), and paclitaxel total clearance ( $CL_{plasma}$ ) were estimated, along with 95% confidence intervals (**Table 15**). Additionally, a curve was fit to previously reported pharmacokinetic results obtained for animals implanted with Zilver PTX Drug-Eluting Peripheral Stents with a total of 876 µg paclitaxel coating per animal, and  $C_{max}$ ,  $T_{max}$ , AUC,  $t_{1/2}$ , and  $CL_{plasma}$  were estimated.

**Table 15. Pharmacokinetic parameters (with 95% confidence intervals)**

Parameter	One Stent (n = 42)	Two Stents (n = 16)	Animal PK Study (n = 2) <sup>1</sup>
$T_{max}$ (min)	20	22	20
$C_{max}$ (ng/mL)	4.4 (4.2 - 4.6)	6.6 (6.3 - 6.9)	7.1
$AUC_{0-last}^{-2}$ (ng·h/mL)	6.5 (4.7 - 8.5)	14.0 (10.7 - 17.2)	12.8
$AUC_{0-inf}^{-3}$ (ng·h/mL)	6.5 (4.7 - 8.5)	14.9 (11.2 - 18.7)	12.8
$t_{1/2}$ (h)	2.4 (1.8 - 3.3)	7.0 (5.2 - 10.8)	1.6
$CL_{plasma}$ (L/h)	107 (81.4 - 147.3)	93.3 (74.6 - 124.7)	68.5

<sup>1</sup> Confidence intervals were not calculated for this group because of the sample size of 2.

<sup>2</sup> AUC from time zero to time of last measured concentration.

<sup>3</sup> AUC from time zero to infinity.

**Angiographic/IVUS Substudy**

IVUS was intended to confirm the lack of aneurysm or stent malapposition related to the Zilver PTX Drug-Eluting Peripheral Stent, consistent with previous results from animal studies, and angiography was intended to confirm the concordance of angiography and duplex ultrasound for assessing vessel patency at follow-up. For the first 60 patients enrolled in the study, IVUS was required post-stenting and at 6 months for those patients receiving stents. An additional eligible (80) patients, 40 each from the control and treatment groups, were assigned to an angiographic/IVUS substudy. Patients in the substudy required angiography at the 12-month follow-up and substudy patients in the Zilver PTX treatment group also required IVUS post-stenting and at the 12-month follow-up.

The angiographic analysis includes all patients with both angiography and duplex ultrasound obtained at the 12-month follow-up and analyzed by the core laboratories. The agreement of angiography and duplex ultrasound for assessing lesion patency was evaluated. At 12-month follow-up, approximately 89% of lesions assessed by both angiography and duplex ultrasound were determined to be either patent or not patent by both measures and this excellent concordance between angiography and ultrasound confirms the validity of duplex ultrasound for assessing patency in the absence of angiography.

IVUS results demonstrate that no aneurysm or stent malapposition was detected at the 6-month or 12-month follow-ups, and the rate of aneurysm and stent malapposition was low immediately after stenting (2.4%, 2/85 aneurysm, 2/84 malapposition) and was not associated with clinical sequelae. Moreover, both events of aneurysm and stent malapposition detected immediately post-stenting were associated with a pre-existing aneurysm visible on the angiogram prior to stent implantation.

**Stent Integrity**

Zilver stent integrity (including both the Zilver PTX stent and the bare metal Zilver stent used in bailout stenting procedures) was evaluated prior to discharge and at 12 months by high resolution stent x-rays intended to provide visualization of the individual stent struts. As shown in **Table 16**, no stent fractures were detected upon procedure completion (0/528) and only four stent fractures (4/457) were detected at 12 months, for a 12-month stent fracture rate of 0.9%. Two of the four stent fractures were Type I (single strut) and two were Type III (multiple strut fractures resulting in complete transection of the stent, without displacement of the stent segments). All four of the patients with a stent fracture maintained primary patency and remained free from TLR through 12-month follow-up, indicating that the stent fractures did not adversely affect patient safety or the effectiveness of the Zilver PTX Drug-Eluting Peripheral Stent.

**Table 16. Stent integrity prior to discharge and at 12 months**

Time Point	Stents Visualized by X-ray	Stents with Fracture	Stent Fracture Rate
Prior to discharge	528	0	0%
12 months	457	4	0.9%

**Stent Thrombosis**

Through 24 months, stent thrombosis (or abrupt stent closure) was reported for 3.6% (2/56) of patients receiving bare Zilver stents (after provisional stenting in the PTA control group) and for 2.3% (7/305) of patients receiving Zilver PTX Drug-Eluting Peripheral Stents (in the Zilver PTX treatment group and after provisional stenting in the PTA control group).

**Results Summary**

The primary safety and effectiveness hypotheses of the study were both met and the study was successful. The results indicate that primary treatment with the Zilver PTX Drug-Eluting Peripheral Stent is as safe or safer than treatment with PTA with provisional stenting ( $p < 0.01$ ) and that the Zilver PTX Drug-Eluting Peripheral Stent provides a significantly higher rate of primary patency compared to PTA ( $p < 0.01$ ). Additionally, stenting with the paclitaxel-coated Zilver PTX Drug-Eluting Peripheral Stent is significantly more effective ( $p < 0.01$ ) in maintaining primary patency through 12 months than stenting with the bare (uncoated) stent.

**SINGLE ARM CLINICAL STUDY**

The Zilver PTX single arm study was a prospective, non-randomized, open-label, multicenter study enrolling patients in Europe, Canada, and Korea with *de novo* or restenotic lesions of the above-the-knee femoropopliteal artery. All patients enrolled in the study were treated with the Zilver PTX Drug-Eluting Peripheral Stent.

**Study Overview**

The study was designed to quantify the safety and performance of the Zilver PTX Drug-Eluting Peripheral Stent. The primary hypothesis was that the event-free survival rate at 6 months for patients receiving the Zilver PTX Drug-Eluting Peripheral Stent was comparable to that seen in prior published endovascular studies of similar patient populations. Furthermore, this study was designed to contribute to a multi-study pool of more than 1000 Zilver PTX patients to provide broad clinical data and to establish the rate of potentially rare (device- or drug-related) adverse events with reasonable precision.

A total of 787 patients were enrolled at 30 sites. The study entry criteria were similar to the randomized study with the exception that there was no limitation on lesion length or the number of lesions treated per patient (up to 4 Zilver PTX Drug-Eluting Peripheral Stents could be used per patient) and the inclusion of lesions with in-stent restenosis was allowed. The study follow-up schedule included clinical assessment at pre-discharge and at 1, 6, 12, and 24 months, and ultrasound imaging and stent X-rays at pre-discharge and at 6 and 12 months. Telephone contact was scheduled for 3, 9, and 18 months. Clinical follow-up was available for 740 patients at 12 months and 500 patients at 24 months at the time of data analysis.

## Demographics

Patient medical history included a high incidence of diabetes (36%, 285/787), hypercholesterolemia (58%, 458/787), hypertension (80%, 627/787), and past or current smoking (80%, 632/787). A total of 1722 Zilver PTX Drug-Eluting Peripheral Stents were implanted in 900 lesions during the study procedure, with more than 60% of patients being treated with at least 2 stents. Lesions treated in the study had a mean length of 100 ± 82 mm with 25% of lesions > 14 cm in length (224/900), 38% (345/900) classified as total occlusions, and 24% (219/900) having been previously treated, including 130 lesions that had been previously stented. **NOTE:** Previously stented lesions are outside of the approved indication for use.

## Results

The 6-month event-free survival (EFS) rate for patients with a Zilver PTX Drug-Eluting Peripheral Stent was 97.4%, a rate superior to the prospectively defined objective performance criterion of 75% EFS ( $p < 0.01$ ). EFS was 89.0% at 12 months and 79.3% at 24 months. Similarly, freedom from TLR was 89.3% at 12 months and 80.5% at 24 months.

Secondary outcomes were also favorable. Specifically, primary patency (conservatively defined as a PSV ratio < 2.0) was 83.0% at 12 months, and clinical and functional status measures (i.e., Rutherford classification, ABI, walking scores, and quality of life scores) improved significantly from pre-procedure to 12 months. Stent fractures were detected in 1.5% of stents (22/1432) at 12 months.

The results from the Zilver PTX single arm study provide supporting evidence confirming the safety of the Zilver PTX Drug-Eluting Peripheral Stent for the treatment of symptomatic vascular disease of the above-the-knee femoropopliteal arteries.

## COMBINED ANALYSIS OF RANDOMIZED AND SINGLE ARM STUDIES

The combined randomized and single arm studies include more than 1000 patients treated with more than 2000 Zilver PTX Drug-Eluting Peripheral Stents.

Stent integrity and the potential for rare adverse events were evaluated in this combined group of patients receiving Zilver PTX Drug-Eluting Peripheral Stents, including patients with long lesions treated with overlapped stents. Event-free survival at 12 months was 90.4% for the primary Zilver PTX group in the randomized study and 89.0% in the single arm study. Primary patency at 12 months was also similar between the two studies at 82.7% for the primary Zilver PTX group in the randomized study and 83.0% in the single arm study. These similar results demonstrate the complementary and supportive nature of the two studies. Whereas the randomized study was the pivotal study supporting safety and effectiveness and enrolled moderate length lesions, the single arm study provides additional clinical evidence supporting safety in a broader patient population.

Regarding stent integrity, high resolution radiographs demonstrated that 98.6% of stents remained free from stent fracture at 12 months (1863/1889), for a stent fracture rate of 1.4% (Table 17).

**Table 17: Stent integrity at 12 months**

# of Stents Visualized by X-ray	# of Stents with No Fracture	# of Stents with Stent Fracture	12-month Stent Integrity Rate	12-month Stent Fracture Rate
1889	1863	26*	98.6%	1.4%

\* 3 Type I (single strut fracture), 5 Type II (multiple strut fractures resulting in complete transection of the stent, without displacement of stent segments), 4 Type III (multiple single strut fractures), 14 Type IV (multiple strut fractures resulting in displacement of segments of the stent).

Importantly, 17 of the 26 stents with fractures (65%) were associated with quantifiable stent elongation during the implant procedure, including 11 (79%) of the Type IV fractures.

Regarding potential rare adverse events, there were no reported drug or hypersensitivity reactions attributed to the paclitaxel drug coating or the nitinol stent; the majority of reactions (21/33) were attributed to the antiplatelet medications ticlopidine or clopidogrel. Neutropenia was reported for 3.6% (2/56) of patients receiving bare Zilver stents (after provisional stenting in the PTA control group) and for 0.4% (4/1092) of patients receiving Zilver PTX Drug-Eluting Peripheral Stents (all Zilver PTX Drug-Eluting Peripheral Stent patients from combined randomized and single arm studies). None of the rare cases of neutropenia were determined to be due to the paclitaxel coating on the Zilver PTX Drug-Eluting Peripheral Stent or to participation in the study. The incidence of adverse events with a potential to be related to particulate matter (i.e., embolism distal to the study vessel and blue toe syndrome) was low for patients receiving bare Zilver (1.8%, 1/56) or Zilver PTX Drug-Eluting Peripheral Stents (0.6%, 7/1092), and was not increased for patients receiving a Zilver PTX Drug-Eluting Peripheral Stent. The low rate of stent thrombosis through 24 months was not increased for patients receiving a Zilver PTX Drug-Eluting Peripheral Stent (3.4%, 37/1092) compared to patients receiving a bare Zilver stent (3.6%, 2/56), and was also consistent with previous femoropopliteal stenting outcomes reported in the literature. These results, including more than 1,000 patients and 2,000 Zilver PTX Drug-Eluting Peripheral Stents, provide additional evidence supporting the safety of the polymer-free PTX\* coating and of the Zilver PTX Drug-Eluting Peripheral Stent.

## ZILVER PTX LONGER LENGTH CLINICAL STUDY

The Zilver PTX Longer Length clinical study was a prospective, multi-center, randomized study enrolling patients in Europe with *de novo* or restenotic lesions of the above-the-knee femoropopliteal artery. All patients enrolled in the study were treated with the Zilver PTX Drug-Eluting Peripheral Stent.

### Study Overview

The study was designed to confirm the deployability and safety of the 100 and 120mm lengths of the Zilver PTX Drug-Eluting Peripheral Stent. The hypothesis was that the change in stent length (from pre- to post- deployment) would be equivalent for longer ( $\geq 100$ mm) and shorter ( $\leq 80$ mm) Zilver PTX Drug-Eluting Peripheral Stents. Additional analyses included 30-day mortality rate.

A total of 45 patients were randomized in a 2:1 fashion to treatment with either longer ( $\geq 100$ mm) Zilver PTX Drug-Eluting Peripheral Stents or overlapping shorter ( $\leq 80$ mm) Zilver PTX Drug-Eluting Peripheral Stents. Patients eligible to be enrolled in the study had symptoms of peripheral artery disease (Rutherford category  $\geq 2$ ) and a stenotic or occluded atherosclerotic lesion ( $\geq 8$ cm and  $\leq 28$ cm in length) of the above-the-knee femoropopliteal artery with a reference vessel diameter of 4mm to 9mm. **NOTE:** Treatment of arteries with reference vessel diameter >7mm are outside of the approved indication for use.

Stent length measurements were obtained pre- and post-deployment and clinical follow-up was available for 44 patients at 1 month.

### Demographics

Demographics, medical history, lesion characteristics, and baseline angiographic data were similar for the longer and shorter stent groups. Patient medical history included a high incidence of diabetes (36% and 50% for the longer and shorter stent groups), hypercholesterolemia (71% and 86%), hypertension (84% and 93%), and past or current smoking (77% and 50%). Lesions treated in the study had a mean length of 180 ± 83mm in the longer stent group and 210 ± 108mm in the shorter stent group, 61% of lesions in the longer stent group and 64% in the shorter stent group were >14cm in length, and 84% of lesions in the longer stent group and 79% in the shorter stent group were classified as total occlusions.

### Results

The primary hypothesis of equivalent change in stent length upon deployment for longer stents compared to shorter stents was met with a mean change in stent length upon deployment of 0.5% ± 4.2% for longer stents and -0.2% ± 3.8% for shorter stents ( $p < 0.01$ ; based on two one-sided t-tests with an equivalence margin of 10%), and the study is considered successful.

There were no patient deaths within 30 days of the study procedure. Therefore, the 30-day mortality rate is 0.0% (0/45).

The results from the Zilver PTX Longer Length clinical study support the safety and deployability of 100 and 120mm Zilver PTX Drug-Eluting Peripheral Stents.

## ZILVER PTX THUMBWHEEL DELIVERY SYSTEM CLINICAL STUDY

The Zilver PTX Thumbwheel Delivery System clinical study was a prospective, multi-center, single arm study enrolling patients in Europe, New Zealand, and Australia with *de novo* or restenotic lesions of the above-the-knee femoropopliteal artery. All patients enrolled in the study were treated with the Zilver PTX Drug-Eluting Peripheral Stent with the thumbwheel delivery system.

### Study Overview

The study was designed to evaluate the performance of the Zilver PTX stent thumbwheel delivery system. The hypothesis tested was the absolute mean change in stent length from pre- to post-deployment is less than 10%. The endpoint is defined as a ratio and the hypothesis is testing for the ratio of the change between pre- and post-deployment stent length relative to the pre-deployment stent length. Endpoints in this study included acute effectiveness, acute safety, placement accuracy, and technical success.

A total of 40 patients were enrolled. Patients eligible to be enrolled in the study had symptoms of peripheral artery disease (Rutherford category 2-4) and a stenotic or occluded atherosclerotic lesion of the above-the-knee femoropopliteal artery with a reference vessel diameter of 4mm to 7mm.

Stent length measurements were obtained pre- and post-deployment and clinical follow-up was available for 39 patients at 1 month.

### Demographics

Demographics, medical history, lesion characteristics, and baseline angiographic data were collected. Patient medical history included a high incidence of diabetes (35.0%), hypercholesterolemia (47.5%), hypertension (80.0%), and past or current smoking (77.5%). Lesions treated in the study had a mean length of 106 ± 76mm, 58.5% of lesions were classified as total occlusions, and 51.2% of lesions had severe calcification.

### Results

The hypothesis of mean change in stent length of less than 10% upon deployment was met with a mean change in stent length upon deployment of -1.0% ± 2.1% (median -0.4%, range -13.3% to 2.4%;  $p < 0.01$ ), and the study is considered successful. The absolute mean stent length change using the absolute values of stent length change results in a mean change in stent length upon deployment of 1.3% ± 2.0% (median 0.6%, range 0.0% to 13.3%).

Acute safety was 100%, placement accuracy was rated as optimal for all stents, and technical success was 98.4%.

The results from the Zilver PTX Thumbwheel Delivery System clinical study support the performance of the Zilver PTX Thumbwheel Delivery System.

## RETROSPECTIVE ANALYSIS OF LONG LESIONS IN THE ZILVER PTX SINGLE ARM STUDY

Additional analysis of the longer lesion length patient cohorts from the Zilver PTX single-arm study was performed to support an expanded lesion length indication. This included analysis of subgroups of patients from the single-arm study with lesion lengths up to 140 mm (493 patients, 574 lesions), >140 mm to 240 mm (110 patients, 116 lesions), and >240 mm to 300 mm (62 patients, 65 lesions) to support expanding the indication for use to include treatment of lesions up to 300 mm in length. All patients and lesions in the Zilver PTX single-arm study were included in the proposed analyses (n=665), with the exception of lesions with length > 300mm or missing and lesions enrolled for treatment of in-stent restenosis (ISR), since treatment of ISR is not included in the currently approved indication.

## Study Overview

Safety was assessed by evaluation of the combined rate of freedom from all-cause death, amputation, or TVR at 30 days post-procedure for the Zilver PTX stent. Effectiveness was assessed by evaluation of the primary patency rate (< 50% diameter stenosis by angiography or by duplex ultrasound with patency defined as PSVR < 2.4) for the Zilver PTX stent at 12 months post-procedure.

Additional outcomes included longer-term safety (freedom from all-cause death and amputation at 12 and 24 months in patients with lesion lengths up to 140 mm, > 140 mm to 240 mm, and > 240 mm to 300 mm), longer-term safety and effectiveness (freedom from TLR at 12 and 24 months in patients with lesion lengths up to 140 mm, > 140 mm to 240 mm, and > 240 mm to 300 mm) and device integrity (rate of stent fracture at 12 months for Zilver PTX stents implanted in lesions with lengths up to 140 mm, > 140 mm to 240 mm, and > 240 mm to 300 mm).

## Demographics

Demographics, medical history, and lesion characteristics are summarized in Tables 18 and 19 for each lesion length subgroup.

**Table 18: Demographics and medical history**

Demographic		Lesion Length ≤ 140 mm (n=493 patients)	Lesion Length > 140 mm to 240 mm (n=110 patients)	Lesion Length > 240 mm to 300 mm (n=62 patients)
Age (years)		66.9 ± 9.6	66.9 ± 9.5	66.4 ± 9.2
Gender	Male	75.1% (370/493)	74.5% (82/110)	82.3% (51/62)
	Female	24.9% (123/493)	25.5% (28/110)	17.7% (11/62)
Height (in)		66.5 ± 3.2	66.1 ± 3.2	67.2 ± 3.4
Weight (lb)		172.9 ± 31.5	174.3 ± 37.3	177.1 ± 32.3
Diabetes		36.9% (182/493)	39.1% (43/110)	19.4% (12/62)
Diabetes type	Type I	8.2% (15/182)	9.3% (4/43)	25.0% (3/12)
	Type II	91.8% (167/182)	90.7% (39/43)	75.0% (9/12)
Hypercholesterolemia		54.0% (266/493)	62.7% (69/110)	59.7% (37/62)
Hypertension		78.3% (386/493)	82.7% (91/110)	75.8 (47/62)
Carotid disease		10.3% (51/493)	8.2% (9/110)	11.3% (7/62)
Renal disease		10.3% (51/493)	13.6% (15/110)	3.2% (2/62)
Pulmonary disease		8.5% (42/493)	9.1% (10/110)	8.1% (5/62)
Congestive heart failure		8.8% (43/491)	14.5% (16/110)	14.5% (9/62)
Previous MI		14.3% (70/490)	20.0% (22/110)	19.4% (12/62)
Smoking status	Never smoked	17.4% (86/493)	18.2% (20/110)	9.7% (6/62)
	Quit	45.2% (223/493)	40.0% (44/110)	50.0% (31/62)
	Still smokes	34.3% (169/493)	39.1% (43/110)	35.5% (22/62)
	Unknown	3.0% (15/493)	2.7% (3/110)	4.8% (3/62)

**Table 19: Lesion characteristics**

Characteristic		Lesion Length ≤ 140 mm (n=574 lesions)	Lesion Length > 140 mm to 240 mm (n=116 lesions)	Lesion Length > 240 mm to 300 mm (n=65 lesions)
Lesion length (mm)		55.3 ± 34.4	185.7 ± 28.3	265.8 ± 18.0
Proximal RVD (mm)		5.3 ± 0.8	5.4 ± 0.9	5.8 ± 1.0
Distal RVD (mm)		5.2 ± 0.8	5.2 ± 0.8	5.4 ± 1.0
MLD in lesion (mm)		1.0 ± 0.9	0.4 ± 0.7	0.1 ± 0.2
Percent diameter stenosis (%)		80.7 ± 17.1	93.4 ± 12.9	99.2 ± 3.2
Previous treatment of study lesion		11.0% (63/574)	12.1% (14/116)	15.4% (10/65)
Calcification	None	18.5% (106/574)	13.8% (16/116)	9.2% (6/65)
	Little	32.6% (187/574)	39.7% (46/116)	23.1% (15/65)
	Moderate	30.7% (176/574)	29.3% (34/116)	36.9% (24/65)
	Severe	18.3% (105/574)	17.2% (20/116)	30.8% (20/65)
Total occlusion		26.1% (150/574)	70.7% (82/116)	92.3% (60/65)
Other stenosis in artery	None	68.6% (394/574)	61.2% (71/116)	63.1% (41/65)
	≤ 50%	18.5% (106/574)	28.4% (33/116)	30.8% (20/65)
	> 50%	11.8% (68/574)	8.6% (10/116)	3.1% (2/65)
Inflow tract stenosis	None	78.7% (452/574)	65.5% (76/116)	64.6% (42/65)
	≤ 50%	18.1% (104/574)	31.0% (36/116)	27.7% (18/65)
	> 50%	0.9% (5/574)	0% (0/116)	0% (0/65)
Patent runoff vessels	1	18.5% (106/574)	16.4% (19/116)	16.9% (11/65)
	2	32.1% (184/574)	43.1% (50/116)	43.1% (28/65)
	3	48.1% (276/574)	38.8% (45/116)	35.4% (23/65)
	> 3	0.7% (4/574)	0% (0/116)	1.5% (1/65)

## Results

At 30 days, 98.8% of patients with lesions ≤ 140 mm, 97.3% of patients with lesions > 140 to 240 mm, and 96.8% of patients with lesions > 240 to 300 mm remained free from death, amputation, or TVR (Table 20). Adverse events observed for patients treated with longer lesion lengths were captured and categorized. There were no concerning safety trends across lesion lengths or compared to pivotal cohorts.

**Table 20: Analysis of primary safety endpoint**

Measure	Lesion Length	Result	
		Kaplan-Meier Estimate*	% (n/N)
Freedom from all-cause death, amputation, or TVR at 30 days post-procedure	≤ 140 mm	98.8%	98.8% (482/488)
	> 140 mm to 240 mm	97.3%	97.3% (107/110)
	> 240 mm to 300 mm	96.8%	96.8% (60/62)

\* 11 events occurred, 4 patients were censored, and 645 patients remained at risk at 30 days.

The 12-month patency results for lesion length subgroups  $\leq 140$  mm,  $> 140$  mm to 240 mm, and  $> 240$  mm to 300 mm were 91.3%, 73.2%, and 73.5%, respectively (Table 21). These results support the effectiveness of the Zilver PTX stent for treating lesions up to 300 mm in length.

**Table 21: Analysis of primary effectiveness endpoint**

Measure	Lesion Length	Result	
		Kaplan-Meier Estimate*	% (n/N)
Primary patency at 12 months post-procedure	$\leq 140$ mm	91.3%	91.6% (488/533)
	$> 140$ mm to 240 mm	73.2%	73.6% (81/110)
	$> 240$ mm to 300 mm	73.5%	74.1% (43/58)

\* 89 events occurred, 42 lesions were censored, and 570 lesions remained at risk at 12 months.

Additional analyses of longer-term safety, longer-term safety and effectiveness, and device integrity were performed to provide further information supporting the performance of the Zilver PTX stent in lesions up to 300 mm. Longer-term safety, assessed by freedom from all-cause death and amputation, and longer-term safety and effectiveness, assessed by freedom from TLR, were both shown to have clinically acceptable rates through 24 months; the rates of freedom from TLR at 12 and 24 months were 92.9% and 88.1% for lesions up to 140 mm, 75.6% and 71.3% for lesions  $> 140$  mm to 240 mm, and 84.9% and 74.7% for lesions  $> 240$  mm to 300 mm. The Zilver stent was found to have excellent integrity in the above-the-knee femoropopliteal artery in lesions up to 300 mm in length. The results from the retrospective analysis of long lesions in the Zilver PTX Single Arm Study support the safety and effectiveness of the Zilver PTX stent for the treatment of lesions up to 300mm in length.

#### Additional Supporting Clinical Evidence from the Zilver PTX Longer Length Study

Results from the Zilver PTX Longer Length Clinical Study (further described above) provided additional evidence supporting use of the Zilver PTX stent in lesions  $> 140$  mm.

#### Additional Supporting Clinical Evidence from the Zilver PTX Japan Post-Market Surveillance Study

The Japan PMS had no inclusion/exclusion criteria and enrolled 905 patients in Japan who were treated with the Zilver PTX stent. Patient medical history included a high incidence of diabetes (58.8%), hypercholesterolemia (60.8%), hypertension (85.4%), renal disease (43.6%), and renal failure (35.5%). Critical limb ischemia (Rutherford classification 4-6) was present in 21.2% of limbs. Lesions treated in the study had a mean length of 146  $\pm$  96 mm, 41.5% of lesions were classified as total occlusions, and 18.6% of lesions entered the study with ISR. Freedom from TLR was 90.9%, 83.6%, and 79.8% at 1, 2, and 3 years, respectively.

Of the total 905 patients enrolled in the study, 717 patients with 842 lesions were included in an analysis by lesion length (ISR was excluded) including 391 patients with 494 lesions up to 140 mm in length, 183 patients with 201 lesions  $> 140$  mm to 240 mm in length, and 143 patients with 147 lesions  $> 240$  mm to 300 mm in length. Average lesion length was 68  $\pm$  34 mm, 187  $\pm$  27 mm, and 266  $\pm$  23 mm in each length group, respectively, and total occlusions were present in 26.3%, 53.2%, and 72.8% of lesions within each group. Freedom from TLR at 1, 2, and 3 years was 95.5%, 91.4%, and 87.8% for lesions up to 140 mm, 92.2%, 83.4%, and 81.4% for lesions  $> 140$  mm to 240 mm, and 88.5%, 78.8%, and 70.4% for lesions  $> 240$  mm to 300 mm. Primary patency at 1 year was 90.5%, 83.2%, and 80.0% for lesions up to 140 mm,  $> 140$  mm to 240 mm, and  $> 240$  mm to 300 mm, respectively. These results provide additional evidence supporting the real world safety and effectiveness of the Zilver PTX stent for the treatment of lesions up to 300 mm in length.

#### ANALYSIS OF 5 MM DIAMETER STENTS FROM THE ZILVER PTX CHINA CLINICAL STUDY

Analysis of the 5 mm diameter stent patient cohort from the Zilver PTX China clinical study was performed to support approval of the 5 mm diameter Zilver PTX stent.

#### Study Overview

This study was a multicenter, prospective, single-arm clinical study in China that included 68 patients treated with the 5 mm diameter Zilver PTX stent. Patients had symptomatic peripheral arterial disease (Rutherford 2-4) with a single *de novo* or restenotic atherosclerotic lesion of the above-the-knee femoropopliteal artery with  $> 50\%$  angiographically documented diameter stenosis and lesion length  $\leq 140$  mm.

#### Demographics

Baseline patient demographics and lesion characteristics are reasonably similar to a United States patient population with PAD and included an average age of 66.8  $\pm$  8.9, 67.6% male patients, and 57.4% diabetic patients. Lesions treated in the study had a mean length of 71.6  $\pm$  38.8 mm, 58.8% of lesions were classified as total occlusions, and the proximal reference vessel diameter was 4.4  $\pm$  0.6 mm.

#### Results

Procedural, safety, and effectiveness outcomes through 1 year were analyzed as described in **Table 22**. There were no deaths or amputations in this study, and the rate of adverse events was acceptable.

**Table 22: Outcomes assessed and results**

Outcome assessed	Definition	Results
<i>Procedural outcomes</i>		
Placement accuracy	No geographic miss reported by the implanting physician.	100%
Procedural success	Study lesion, including $\pm 5$ mm, with $< 30\%$ residual stenosis determined angiographically immediately after the procedure.	100% (site) 89.7% (core)
Acute effectiveness	Change in stent length from pre- to post-deployment.	1.6% $\pm$ 1.7% (absolute mean value) 0.4% $\pm$ 2.3% (mean value)
Technical success	Composite of placement accuracy and absolute mean change in stent length from pre- to post-deployment of less than 10%.	100%
<i>Safety outcomes</i>		
30-day safety	Freedom from all-cause death, amputation, or target vessel revascularization (TVR) at 30 days post-procedure.	97.1%
1-year safety	Freedom from all-cause death or amputation.	100%
1-year event-free survival	Freedom from device- or procedure-related MAE, including death, amputation, clinically-driven TLR, target limb ischemia requiring surgical intervention, or surgical repair of the target vessel; and freedom from worsening of the Rutherford Classification by two classes (or to class 5 or 6).	91.2%
<i>Effectiveness outcomes</i>		
1-year freedom from TLR	Freedom from re-intervention performed for a $\geq 50\%$ stenosis of the target lesion after documentation of recurrent clinical symptoms of PAD.	92.6%
1-year primary patency	Based on duplex ultrasound imaging reviewed by an independent core laboratory.	70.6%
1-year clinical benefit	Freedom from persistent or worsening symptoms of ischemia (i.e. claudication, rest pain, ulcer or tissue loss) after the initial treatment.	92.6%

The results from the analysis of patients treated with the 5 mm diameter stent in the Zilver PTX China clinical study demonstrate the safety and effectiveness of the 5 mm diameter Zilver PTX stent for the treatment of above-the-knee femoropopliteal artery lesions.

#### US ZILVER PTX POST-APPROVAL STUDY (PAS 1)

This study was required as a condition of approval and was a multicenter, prospective, single-arm clinical study of newly enrolled patients in the United States. All patients enrolled in the study were treated with the Zilver PTX Drug-Eluting Peripheral Stent.

#### Study Overview

Patients with symptomatic vascular disease of the native above-the-knee femoropopliteal arteries with a reference vessel diameter from 4 mm to 7 mm and total lesion lengths up to 140 mm per limb and 280 mm total were eligible for enrollment. A total of 200 patients from 17 sites were enrolled in this study. The objective of this post-approval study was to confirm that post-market results were similar to those observed in the pivotal, pre-market randomized clinical study. The primary study endpoint was the occurrence of TLR at 1 year, and the hypothesis tested the 1-year freedom from TLR rate to a performance goal of 83.1%. Follow-up was through 1 year.

#### Demographics

Baseline patient demographics are consistent with a patient population with peripheral arterial disease and included an average age of 67.9  $\pm$  9.5, 63.0% male patients, and 46.0% diabetic patients. The average lesion length was 76.1  $\pm$  44.2 mm.

#### Results

At 1-year, 175 of 194 eligible patients (90.2%) completed their clinical visit; four patients died, one was lost to follow-up, and one withdrew prior to follow-up. The 1-year freedom from TLR rate was 91.8%, and the primary hypothesis of non-inferiority to the performance goal of 83.1% was met ( $p < 0.001$ ). These results confirm device safety and performance in the post-approval setting.

## US ZILVER PTX POST-APPROVAL STUDY 2 (PAS 2)

This study was required as a condition of approval and included a combination of the patients treated with the Zilver PTX stent in the randomized pivotal study, the Japanese post-market surveillance study (Japan PMS), and the US PAS 1 study. All patients in the study were treated with the Zilver PTX Drug-Eluting Peripheral Stent.

### Study Overview

Patients with symptomatic vascular disease of the native above-the-knee femoropopliteal arteries with a reference vessel diameter from 4 mm to 7 mm and total lesion lengths up to 140 mm per limb and 280 mm total were eligible for enrollment. US PAS 2 included 904 patients: 305 Zilver PTX patients enrolled in the randomized pivotal study, 399 Zilver PTX patients enrolled in the Japan PMS, and 200 Zilver PTX patients enrolled in the US PAS 1. At least 50% of the patients enrolled in US PAS 2 were from the United States. The inclusion/exclusion criteria for US PAS 1 were the same as the criteria for the randomized pivotal study; the Japan PMS study was an all-comers study and did not have any inclusion/exclusion criteria. The objective of the US PAS 2 study was to evaluate longer-term stent integrity and other emerging safety signals in the combined patient cohort. The study endpoints were (1) TLR and stent thrombosis (defined as total occlusion of the study lesion with suspected thrombotic origin) assessed annually through 5 years post-procedure and (2) stent integrity assessed by X-ray at 1, 3, and 5 years post-procedure.

### Demographics

Baseline patient demographics and lesion characteristics were consistent with a patient population with peripheral arterial disease and included an average age of 70.1 ± 9.6 years, 68.1% male patients, and 51.7% diabetic patients. The study included 990 lesions with an average lesion length of 74.2 ± 49.1 mm. Lesions had an average baseline percent diameter stenosis of 84.9%, and 30.1% of lesions were characterized as total occlusions.

### Results

Through 5 years, 451 patients completed the study. Endpoint data through 5 years are reported in Table 23. The rates of freedom from TLR, freedom from thrombosis, and freedom from stent fracture are consistent with rates reported in other Zilver PTX studies and support the safety and performance of the device.


**Table 23. Outcomes assessed and results through 5 years**

Outcome assessed		Kaplan Meier Estimate ± Standard Error
1-Year Outcomes	1-year rate of freedom from TLR	93.2% ± 0.9%
	1-year freedom from thrombosis	97.8% ± 0.5%
	1-year freedom from stent fracture	99.7% ± 0.1%
2-Year Outcomes	2-year rate of freedom from TLR	87.0% ± 1.2%
	2-year freedom from thrombosis	97.5% ± 0.5%
3-Year Outcomes	3-year rate of freedom from TLR	83.9% ± 1.3%
	3-year freedom from thrombosis	97.2% ± 0.6%
	3-year freedom from stent fracture	98.8% ± 0.4%
4-Year Outcomes	4-year rate of freedom from TLR	82.3% ± 1.4%
	4-year freedom from thrombosis	97.0% ± 0.6%
5-Year Outcomes	5-year rate of freedom from TLR	80.8% ± 1.5%
	5-year freedom from thrombosis	96.5% ± 0.7%
	5-year freedom from stent fracture	98.5% ± 0.5%

There were no safety signals identified based on the adverse event data obtained for the US PAS 2 study. The results confirm device safety and performance in the post-approval setting.

### PRODUCT RECOMMENDATIONS

Placement of this vascular stent requires advanced skills in interventional vascular procedures. The following instructions will give technical guidance, but do not obviate formal training in the use of the device. Delivery system components are shown in **Figure 6**.

1. Ergonomic thumbwheel operated handle
2. Safety lock
3. Thumbwheel
4. Luer hub (flushing port)
5. Retraction sheath
6. Stability sheath
7. Tip of delivery system inner catheter
8. Radiopaque marker band on retraction sheath
9. Gold radiopaque markers on stent
10. Stent stop

**Figure 6**

### Multiple Stent Placement

If placements of multiple stents are required in a patient, to cover the length of the lesion, the following recommendations should be considered:

- In relation to the lesion site, the distal area of narrowing should be stented first, followed by the proximal locations (i.e., a second stent should be placed proximally to the previously placed stent).
- Stents placed in tandem must overlap to allow for complete coverage of the lesion.
- The safety and effectiveness of implanting multiple Zilver PTX Drug-Eluting Peripheral Stents with a total drug coating quantity of greater than 3mg paclitaxel (i.e., an additive length greater than 320mm) has not been established.

### INSTRUCTIONS FOR USE

#### Stent Sizing

1. Determine the proper stent size after complete diagnostic evaluation. The stent deployment must be performed under fluoroscopic control. Measure the length of the target lesion to determine the length of the stent required. Allow for the proximal and distal aspects of the stent to cover the entire target area.

**NOTE:** The Zilver PTX Drug-Eluting Peripheral Stent is designed not to shorten upon deployment. Bench testing has shown that on average, the stent length increases from pre-deployment to post-deployment by 0.7-1.3% (Table 24).

The stent is recommended for use in above-the-knee femoropopliteal arteries having reference vessel diameter from 4 mm to 7 mm.

Measure the diameter of the reference vessel (proximal and distal to the lesion) and use the LARGEST reference diameter as your basis for choosing the appropriate stent size.

The stent size should be selected so that the unconstrained stent diameter is at least 1 mm larger than the reference vessel diameter and no more than 2 mm larger than the reference vessel diameter (Table 24).

**Table 24: Stent sizing and length change**

Reference Vessel Diameter (mm)	Nominal Stent Diameter (mm)	Mean Stent Length Change upon Deployment*
4	5	0.7-1.3%
4-5	6	
5-6	7	
6-7	8	

\*Positive number indicates lengthening

**Introduction of the Stent**

1. Gain access to the site using a 6 French sheath with a minimum I.D. of 2.1 mm.
2. To ensure that the system is adequately supported, introduce a 0.035 inch (0.89 mm) diameter wire guide. If hydrophilic wire guides are used, they must be kept fully activated.
3. Pre-dilatation before stent placement is optional and at the discretion of the physician.
4. Prior to advancing the system over the wire and introducing the system into the body use a syringe to flush the device's wire guide lumen and stent with saline or heparinized saline through the hub of the device. Flushing the device with contrast media is not recommended.

Flush the wire lumen until a few drops of saline exit the tip of the device (Figure 7). Next, gently place a finger over the device's end hole and flush again until a few drops of saline can be seen exiting between the end of the sheath and the device tip (Figure 8).


Figure 7


Figure 8

5. Under fluoroscopy, advance the delivery system over the 0.035 inch (0.89mm) diameter wire guide through the introducer sheath until the distal gold radiopaque markers on the stent are beyond the target lesion site (Figure 9).

**NOTE: If resistance is met during the advancement of the delivery system do not force passage. Remove the delivery system and replace with a new device.**

**NOTE: Ensure the distal end of the stability sheath is inside the introducer sheath.**

**Caution: Torquing the device may lead to difficulty or inability to deploy and/or move device over the wire guide.**


Figure 9

**Alignment and Deployment of the Stent**

1. Under fluoroscopy, pull the delivery system back until there is no slack in the delivery system and the radiopaque markers on the stent are aligned with the lesion (Figure 10).



Figure 10


2. Disengage the device's safety lock by gently depressing the red safety button in the direction that is indicated in **Figure 11**.

**Figure 11**

3. Initiate stent deployment by rotating the device thumbwheel in the direction that is shown in **Figure 12**.

**Note:** If high resistance is felt on the thumbwheel prior to stent deployment, do not force deployment. Carefully withdraw the stent system without deploying the stent.


**Figure 12**

**Note:** Rotating the thumbwheel will withdraw the retraction sheath from over the stent.

4. Visually confirm that the radiopaque stent markers are still aligned with the desired position. Reposition if necessary.

**NOTE:** Stent deployment will begin once the distal marker band on the stent retraction sheath has been retracted past the stent's distal markers (Refer to **Figure 7** for sheath and stent marker identification). Once the stent has apposed the vessel wall repositioning of the stent is not recommended. Recapture or re-sheathing of the Zilver PTX Drug-Eluting Peripheral Stent is not possible because the stent retraction sheath cannot be re-advanced. Refer to the **Multiple Stent Placement** section of these Instructions for Use for information on missed lesions.

5. Full deployment of the stent length will occur when the distal end of the retraction sheath has been retracted past the proximal radiopaque markers (**Figure 13**).

**Figure 13**

### Post Stent Deployment

1. Under fluoroscopy, withdraw the delivery system over the wire guide, into the introducer sheath and out of the body. Remove the delivery system from the wire guide.

**Note:** Following stent deployment, if resistance is met during the withdrawal of the delivery system, carefully remove the delivery system and wire guide together as a unit.

If resistance is still encountered during removal of the delivery system and wire guide as a unit, remove the wire guide, delivery system and introducer sheath together as a unit.

2. Perform an arterial angiogram to verify full deployment of the device. If incomplete expansion exists within the stent at any point along the lesion, post-deployment balloon dilatation (standard PTA) can be performed at the discretion of the physician.

**Note:** Flow restrictions remaining after stent deployment (e.g., residual proximal or distal stenosis or dissection, or poor distal outflow) may increase the risk of stent thrombosis. Inflow and outflow should be assessed at procedure completion and additional measures considered (e.g., additional PTA, adjunctive stenting, or distal bypass) if necessary to maintain good inflow and outflow.

3. Remove the wire guide and introducer sheath from the patient.

4. Close the entry wound as appropriate.

5. Dispose in accordance with institutional guidelines.

### HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in an outer non-sterile foil pouch and inner peel-open package. Intended for one-time use. Sterile if inner package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

### REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2025
**COOK IRELAND LTD.**