
- EN**
 3 **Geenen® Sof-Flex® Pancreatic Stent**
 Instructions for Use
- BG**
 5 **Панкреасен стент Geenen® Sof-Flex®**
 Инструкции за употреба
- CS**
 7 **Pankreatický stent Geenen® Sof-Flex®**
 Návod k použití
- DA**
 9 **Geenen® Sof-Flex® pancreasstent**
 Brugsanvisning
- NL**
 11 **Geenen® Sof-Flex®-pancreasstent**
 Gebruiksaanwijzing
- ET**
 13 **Pankrease stent Geenen® Sof-Flex®**
 Kasutusjuhised
- FR**
 15 **Endoprothèse pancréatique Geenen® Sof-Flex®**
 Mode d'emploi
- DE**
 17 **Geenen® Sof-Flex® Pankreasstent**
 Gebrauchsanweisung
- EL**
 19 **Παγκρεατική ενδοπρόσθεση Geenen® Sof-Flex®**
 Οδηγίες χρήσης
- HU**
 21 **Geenen® Sof-Flex® hasnyálmirigysztent**
 Használati utasítás
- IT**
 23 **Stent pancreatico Geenen® in Sof-Flex®**
 Istruzioni per l'uso
- LV**
 25 **Geenen® Sof-Flex® aizkuņģa dziedzerā stents**
 Lietošanas instrukcija
- LT**
 27 **„Geenen® Sof-Flex®“ kasos lataų stentas**
 Naudojimo nurodymai
- NO**
 29 **Geenen® Sof-Flex® bukspyttkjertelstent**
 Bruksanvisning
- PL**
 31 **Stent do przewodów trzustkowych Geenen® Sof-Flex®**
 Instrukcja użycia
- PT**
 33 **Stent pancreático Geenen® Sof-Flex®**
 Instruções de utilização
- RO**
 35 **Stent pancreatic Geenen® Sof-Flex®**
 Instrucțiuni de utilizare
- SK**
 37 **Pankreatický stent Geenen® Sof-Flex®**
 Návod na použitie
- ES**
 39 **Stent pancreático Geenen® Sof-Flex®**
 Instrucciones de uso
- SV**
 41 **Geenen® Sof-Flex® pankreasstent**
 Bruksanvisning

Ductal end

- Дуктален край
- Duktální konec
- Duktale ende
- Ductusuiteinde
- Juhapoolne ots
- Extrémité canalaire
- Duktales Ende
- Άκρο για τους πόρους
- Ductalis vég
- Estremità duttale
- Vada gals
- Latako galas
- Ductus-ende
- Koniec prowadzący
- Extremidade com ponta cónica
- Capăt ductal
- Duktálny koniec
- Extremo del conducto
- Gångände

INTENDED USE

This device is used to drain obstructed pancreatic ducts.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

Retain the positioning sleeve for use when introducing the stent flaps into the accessory channel.

CONTRAINDICATIONS

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Inability to pass the wire guide or stent through the obstructed area.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Those associated with pancreatic stent placement include, but are not limited to: trauma to the pancreatic tract or duodenum, obstruction of the common bile duct, stent migration.

PRECAUTIONS**Illustrations**

Refer to package label for minimum channel size required for this device.

A complete diagnostic evaluation must be performed, prior to use, to determine the proper stent size.

The stent must be placed under fluoroscopic monitoring.

The positioning sleeve is not intended for use in the accessory channel of the endoscope.

The ductal end of the stent must be positioned in the pancreatic duct while the other end remains in the duodenum (see **Fig. 1**).

This device should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

Select the Cook stent introducer system in the appropriate French size.

WARNINGS

MR SAFETY INFORMATION

MR

This symbol means the device is MR Safe.

The Geenen Sof-Flex Pancreatic Stent is MR Safe.

In non-clinical testing, the image artifact caused by the device extends approximately 7 mm from the stent when imaged with a gradient echo pulse sequence and a 3 T MR system.

INSTRUCTIONS FOR USE

1. Gently ensure full extension of all side flaps.
2. Load positioning sleeve onto duodenal flap end of stent.
3. Introduce stent, ductal end first, and positioning sleeve onto a pre-positioned wire guide.
4. Advance pushing catheter in 1-2 cm increments until stent is in desired position.
5. Fluoroscopically and endoscopically confirm desired stent position.
6. After confirming stent position, gently remove wire guide from endoscope while maintaining position of stent with pushing catheter.
7. Gently remove pushing catheter from accessory channel.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2019©
COOK IRELAND LTD.