EN
4
BG
9
CS
15
DA
20
NL
25
ET
31
FR
36
DE
42
EL
48
HU
54
IT
59
LV
65
LT
70
NO
75
PL
80
PT
86
RO
92
ES
97
SK
103
SV
108

Nasal Biliary Drainage Set

Instructions for Use

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


I F U 0 1 2 9 - 1

NASAL BILIARY DRAINAGE SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

INTENDED USE

The Nasal Biliary Drainage Set is intended for temporary endoscopic drainage of the biliary duct through the nasal passage by use of an indwelling catheter.

DEVICE DESCRIPTION

Each set includes a drainage catheter, drainage connecting tube, nasal transfer tube and a pigtail straightener (if supplied). The drainage catheter is radiopaque to aid with fluoroscopic placement. Devices with a –RT suffix on their part numbers contain radiopaque marker bands.

For the drainage catheters with preloaded pigtail straightener, the pigtail straightener is intended to be used for straightening the pigtail prior to insertion of drainage catheter into the accessory channel of endoscope.

The drainage connector on the drainage connecting tube is intended to be connected to a drainage collection bag. In the middle of the drainage connecting tube is a three way stopcock (**See Fig. 5**) which allows for a flow through the drainage connecting tube during the procedure.

The nasal transfer tube enables the drainage catheter to be threaded through the oral cavity and out through the nostril. (**See Figs. 2-4**)

For Drainage Only

The above symbol means the following: For Drainage Only. (**See Fig. 1**)

OFF

The above symbol means the following: OFF. (**See Fig. 5**)

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

Contraindications of nasal biliary drainage include, but are not limited to: active pancreatitis or other infection • coagulopathy • inability to pass a wire guide or drainage catheter through the obstructed area • recent food ingestion • sepsis.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with ERCP include, but are not limited to: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • pancreatitis • perforation • respiratory depression or arrest • sepsis • trauma / edema to upper GI tract

Potential adverse events associated with nasal biliary drainage include, but are not limited to: blockage of the drainage catheter • cholangitis • fever • hemorrhage • infection • nasal irritation and/or sore throat • obstruction of the pancreatic duct • pain / discomfort • perforation • stasis • trauma / edema of duct • trauma / edema of papilla / papillary fronds • trauma to duodenum • trauma to nasal cavity.

PRECAUTIONS

The product is intended for use by physicians trained and experienced in upper GI endoscopy and ERCP techniques.

Standard techniques for the endoscopic placement of a drainage catheter in the biliary duct should be employed.

Refer to package label for minimum channel size required for this device.

A complete diagnostic evaluation of the patient prior to use to determine proper catheter size should be carried out.

Do not use this device for any purpose other than stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

Pigtail straightener, if any, is not intended for use in accessory channel of endoscope.

This device is not intended for use beyond 29 days.

Sphincterotomy is not required for device placement.

It is recommended that drainage collection bag be available for use during procedure.

3-way stopcock in Drainage connecting tube is not intended for introduction of any medication, contrast media or saline to the patient.

It is recommended that drainage catheter be periodically evaluated for patency.

WARNINGS

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Visually inspect the packaging and device. If the packaging is opened or damaged when received, do not use. If a device abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

MR Safety Information

This symbol means the device is MR Conditional.

Non-clinical testing has demonstrated the Nasal Biliary Drainage Sets which contain radiopaque marker bands are MR Conditional. These devices have part numbers with a –RT suffix. A patient with these devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3.0 T or less
- Maximum spatial field gradient of 1,900 gauss/cm (19 T/m)
- Maximum MR system reported, whole body average specific absorption rate (SAR) of 4 W/kg

In non-clinical testing, the image artifact caused by the device extends approximately 4 mm from the Nasal Biliary Drainage Set when imaged with a gradient echo pulse sequence and a 3 T MR system.

This symbol means the device is MR Safe.

The Nasal Biliary Drainage Sets which do not contain radiopaque marker bands are MR Safe.

In non-clinical testing, the image artifact caused by the device extends approximately 4 mm from the Nasal Biliary Drainage Set when imaged with a gradient echo pulse sequence and a 3 T MR system.

INSTRUCTIONS FOR USE

Illustrations

1. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
2. Visually inspect the device confirming it is free from damage prior to use.
3. Flush drainage catheter using a sterile syringe filled with sterile water or saline attached to Luer lock fitting on end of catheter.
4. Carefully remove Tuohy-Borst connector from end of drainage catheter.
5. If tip of drainage catheter is configured with a pigtail and is preloaded with a pigtail straightener, advance pigtail straightener along shaft and to tip of catheter to straighten pigtail curl.

CAUTION: Care must be exercised when straightening pigtail curl in order to avoid kinking or breaking catheter.

6. Introduce drainage catheter and pigtail straightener, if any, onto pre-positioned wire guide and advance in small increments until the catheter enters the accessory channel of endoscope.

CAUTION: Ensure the pigtail straightener, if any, remains positioned at tip of catheter in order to avoid kinking or breaking catheter.

7. Continue to advance drainage catheter in small increments until pigtail portion of the catheter has fully entered the accessory channel of endoscope. Remove pigtail straightener from the catheter.

8. Continue to advance drainage catheter in small increments until catheter is endoscopically visualized exiting endoscope.

NOTE: Elevator of endoscope must be open to allow drainage catheter to exit.

9. Under fluoroscopic monitoring, slowly advance drainage catheter to desired position in biliary duct above obstruction.
10. When drainage catheter is in desired position, slowly remove endoscope while holding drainage catheter and wire guide in place at patient's mouth

CAUTION: Do not remove endoscope without removing the Tuohy-Borst connector from the drainage catheter. Doing so may result in dislocation of drainage catheter from the desired position in the biliary duct.

CAUTION: Do not remove endoscope without holding the drainage catheter and wire guide. Doing so may result in dislocation of drainage catheter from the desired position in the biliary duct.

11. Fluoroscopically confirm that drainage catheter and wire guide remain positioned in biliary duct, then slowly remove wire guide in small increments while maintaining position of drainage catheter.

CAUTION: Do not remove wire guide without maintaining the position of the drainage catheter. Doing so may result in dislocation of drainage catheter from the desired position in the biliary duct.

12. Lubricate the nasal transfer tube with a water-soluble lubricant approved for medical use and advance through nostril until it exits in back of throat.
 13. Grasp connection end of nasal transfer tube in throat and pull it out through mouth (**See Fig. 2**).
 14. If placing an 8.5 French or smaller catheter, lubricate the proximal end of drainage catheter with water soluble lubricant if required. Thread the proximal end of the catheter through the oral end of the nasal transfer tube and advance until it exits through the nasal end or until it cannot be advanced any further. If placing a 10 French or larger tube, securely connect the tip of the tube to the oral end of the nasal transfer tube (**See Fig. 3**).
 15. Hold drainage catheter, to prevent kinking and removal, and slowly pull nasal transfer tube until nasal transfer tube is fully removed from nostril and drainage catheter emerges through nostril (**See Fig. 4**).
- CAUTION:** Do not pull nasal transfer tube without holding the drainage catheter. Doing so may result in kinking of drainage catheter or dislocation of drainage catheter from the desired position in the biliary duct.
16. Disconnect/remove the nasal transfer tube from the drainage catheter.
 17. Wipe any remaining lubricant from the proximal end of the drainage catheter.
 18. Reattach Tuohy-Borst connector onto end of drainage catheter. Contrast may be injected to check placement of catheter within duct.

19. Secure drainage catheter to patient, then connect end of catheter to Luer lock fitting on drainage connecting tube. Attach universal connector to appropriate drainage bag. Ensure the handle of the 3-way stopcock is turned to the correct position to begin drainage.

Removal

20. Carefully remove the drainage catheter after indwell time by gently pulling it through the nasal cavity until it exits through the nose

Upon completion of procedure, dispose of components per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages.

This symbol on the label indicates that this device contains phthalates. Specific phthalates contained in the device are identified beside or below the symbol by the following acronyms

- DEHP: Di(2-ethylhexyl) phthalate

The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

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© COOK 2021
William A. Cook Australia Pty. Ltd.
COOK IRELAND LTD.