

- EN 4 **Pushing Catheter and Guiding Catheter**  
 Instructions for Use
- BG 7 **Избутващ катетър и водещ катетър**  
 Инструкции за употреба
- HR 11 **Potisni kateter i vodeći kateter**  
 Upute za uporabu
- CS 14 **Tlačný katetr a vodič katetr**  
 Návod k použití
- DA 17 **Skubbekateter og styrekateter**  
 Brugsanvisning
- NL 20 **Pushing-katheter en geleidekatheter**  
 Gebruiksaanwijzing
- ET 24 **Suruv kateeter ja juhtkateeter**  
 Kasutusjuhised
- FR 27 **Cathéter pousse-prothèse et cathéter guide**  
 Mode d'emploi
- DE 31 **Platzierungskatheter und Führungskatheter**  
 Gebrauchsanweisung
- EL 35 **Καθετήρας ώθησης και οδηγός καθετήρας**  
 Οδηγίες χρήσης
- HU 39 **Tolókatéter és vezetőkátéter**  
 Használati utasítás
- IT 42 **Catetere spingitore e catetere guida**  
 Istruzioni per l'uso
- LV 46 **Bīdītājkatētrs un vadītājkatētrs**  
 Lietošanas instrukcija
- LT 49 **Stumiamasis kateteris ir kreipiamasis kateteris**  
 Naudojimo nurodymai
- NO 52 **Skyvekateter og ledekateter**  
 Bruksanvisning
- PL 55 **Cewnik popychający i cewnik prowadzący**  
 Instrukcja użycia
- PT 59 **Cateter posicionador e cateter guia**  
 Instruções de utilização
- RO 63 **Cateter de împingere și cateter de ghidaj**  
 Instrucțiuni de utilizare
- SK 67 **Zatlačací katéter a vodiaci katéter**  
 Návod na použitie
- ES 71 **Catéter empujador y catéter guía**  
 Instrucciones de uso
- SV 75 **Påförarkateter och guidingskateter**  
 Bruksanvisning

- Figure 1: Pushing Catheter (PC)**  
**Фигура 1: Избутващ катетър (PC)**  
**Slika 1: Potisni kateter (PC)**  
**Obrázek 1: Tlačný katetr (PC)**  
**Figur 1: Skubbekateter (PC)**  
**Afbeelding 1: Pushing-katheter (PC)**  
**Joonis 1: Suruv kateeter (PC)**  
**Figure 1 : Cathéter pousse-prothèse (PC)**  
**Abbildung 1: Platzierungskatheter (PC)**  
**Εικόνα 1: Καθετήρας ώθησης (PC)**  
**1. ábra: Tolókatéter (PC)**  
**Figura 1 – Catetere spingitore (PC)**  
**1. attēls Biditājkatētrs (PC)**  
**1 pav.: Stumiamasis kateteris (PC)**  
**Figur 1: Skyvekateter (PC)**  
**Rysunek 1: Cewnik popychający (PC)**  
**Figura 1: Cateter posicionador (PC)**  
**Figura 1: Cateter de împingere (PC)**  
**Obrázok 1: Zatláčací katéter (PC)**  
**Figura 1: Catéter empujador (PC)**  
**Figur 1: Páförarkateter (PC)**

Pushing Catheter  
 Избутващ катетър  
 Potisni kateter  
 Tlačný katetr  
 Skubbekateter  
 Pushing-katheter  
 Suruv kateeter  
 Cathéter pousse-prothèse  
 Platzierungskatheter  
 Καθετήρας ώθησης  
 Tolókatéter  
 Catetere spingitore  
 Biditājkatētrs  
 Stumiamasis kateteris  
 Skyvekateter  
 Cewnik popychający  
 Cateter posicionador  
 Cateter de împingere  
 Zatláčací katéter  
 Catéter empujador  
 Páförarkateter

**Figure 2: Guiding Catheter (GC)**  
**Фигура 2: Водещ катетър (GC)**  
**Slika 2: Vodeći kateter (GC)**  
**Obrázek 2: Vodící katetr (GC)**  
**Figur 2: Styrekateter (GC)**  
**Afbeelding 2: Geleidekatheter (GC)**  
**Joonis 2: Juhtkateeter (GC)**  
**Figure 2 : Cathéter guide (GC)**  
**Abbildung 2: Führungskatheter (GC)**  
**Εικόνα 2: Οδηγός καθετήρας (GC)**  
**2. ábra: Vezetőkatéter (GC)**  
**Figura 2 – Catetere guida (GC)**  
**2. attēls Vadītājkatētrs (GC)**  
**2 pav.: Kreipiamasis kateteris (GC)**  
**Figur 2: Ledekateter (GC)**  
**Rysunek 2: Cewnik prowadzący (GC)**  
**Figura 2: Cateter guia (GC)**  
**Figura 2: Cateter de ghidaj (GC)**  
**Obrázok 2: Vodiaci katéter (GC)**  
**Figura 2: Catéter guía (GC)**  
**Figur 2: Guidingkateter (GC)**

**Guiding Catheter Hub**  
**HubВтулка на водещия катетър**  
**Čvoršte vodećeg katetera**  
**Ústí vodicího katetru**  
**Styrekatetermuffe**  
**Aanzetstuk geleidekatheter**  
**Juhtkateetri muhv**  
**Embase du cathéter guide**  
**Ansatz des Führungskatheters**  
**Ομφαλός οδηγού καθετήρα**  
**Vezetőkatéter kónusza**  
**Connettore del catetere guida**  
**Vadītājkatetra galviņa**  
**Kreipiamojo kateterio įvorė**  
**Ledekatetermuffe**  
**Złączka cewnika prowadzącego**  
**Conector do cateter guia**  
**Ambou cateter de ghidaj**  
**Hrdlo vodiaceho katétra**  
**Conector del catéter guía**  
**Guidingkateterns fattning**

**Guiding Catheter**  
**Водещ катетър**  
**Vodeći kateter**  
**Vodící katetr**  
**Styrekateter**  
**Geleidekatheter**  
**Juhtkateeter**  
**Cathéter guide**  
**Führungskatheter**  
**Οδηγός καθετήρας**  
**Vezetőkatéter**  
**Catetere guida**  
**Vadītājkatētrs**  
**Kreipiamasis kateteris**  
**Ledekateter**  
**Cewnik prowadzący**  
**Cateter guia**  
**Cateter de ghidaj**  
**Vodiaci katéter**  
**Catéter guía**  
**Guidingkateter**

**Radiopaque Band**  
**Рентгеноположителна лента**  
**Rendgenski vidljiva traka**  
**Rentgenkontrastní proužek**  
**Röntgenfast marking**  
**Radiopake band**  
**Röntgenkontrastne riba**  
**Bande radio-opaque**  
**Röntgendichte Markierung**  
**Ακτινοσκοπική ταινία**  
**Sugárfogó sáv**  
**Banda radiopaca**  
**Rentgenkontrastėjoša josla**  
**Rentgenokontrastinis žiedas**  
**Radioopakt bånd**  
**Cieniodajny pasek**  
**Banda radiopaca**  
**Bandā radiopacā**  
**Rádioopakny krúžok**  
**Banda radiopaca**  
**Röntgentätt band**

**CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

**Table 1: Compatible Cook Biliary Stents for the Pushing Catheter and Guiding Catheter**

| Device Family           | Cook Biliary Stents |                       |                |                           |                        |         |
|-------------------------|---------------------|-----------------------|----------------|---------------------------|------------------------|---------|
|                         | Cotton-Huibregtse®  | Soehendra Tannenbaum® | Cotton-Leung®  | Cotton-Leung® (Sof-Flex®) | Plastic Biliary Stent* | Zimmon® |
|                         | CHBSO               | TTSO                  | CLSO           | CLSO-SF                   | PBS*                   | ZSO     |
| Pushing Catheter (PC)   | √                   | √                     | √              | √                         | √                      | √       |
| Guiding Catheter (GC-5) | 8.5Fr               | 8.5Fr                 | 8.5Fr          | 10Fr                      | 8.5Fr                  | X       |
| Guiding Catheter (GC-6) | 10Fr<br>11.5Fr      | 10Fr<br>11.5Fr        | 10Fr<br>11.5Fr | X                         | 10Fr<br>11.5Fr         | X       |

\* Note: PBS devices are not available in all markets, including the US.

**Table 2: Compatible Cook Pancreatic Stents for the Pushing Catheter and Guiding Catheter**

| Device Family           | Cook Pancreatic Stents |                      |                                     |
|-------------------------|------------------------|----------------------|-------------------------------------|
|                         | Geenen®                | Geenen® (Sof-Flex®)  | Zimmon®                             |
|                         | GPSO/<br>GPSOS         | GPSO-SF/<br>GPSOS-SF | SPSOF/<br>SPSOS/<br>ZPSOF/<br>ZPSOS |
| Pushing Catheter (PC)   | √                      | √                    | √                                   |
| Guiding Catheter (GC-5) | 8.5Fr                  | X                    | X                                   |
| Guiding Catheter (GC-6) | 10Fr<br>11.5Fr         | X                    | X                                   |

Convention is as follows;

- The symbol "X" means not indicated for use with.
- Checkmark "√" means compatible with corresponding stent sizes.
- Fr refers to the stent diameter in French size.

**Note:** Pushing catheters are compatible with the stents indicated in **Table 1** and **Table 2** e.g. PC-7 is compatible with a 7Fr stent, PC-5 is compatible with a 5Fr stent etc.

#### **INTENDED USE**

This device is used for endoscopic biliary and pancreatic stent placement.

## **DEVICE DESCRIPTION**

This IFU covers Pushing Catheter and Guiding Catheter introduction systems represented by **Figures 1** and **2**. See **Table 1** for the compatible Cook Biliary Stents and **Table 2** for the compatible Cook Pancreatic Stents. For all stent information including orientation, refer to the stent instructions for use.

## **CONTRAINDICATIONS**

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Inability to pass wire guide / guiding catheter / stent through the obstructed area.

## **POTENTIAL ADVERSE EVENTS**

Those associated with ERCP include, but are not limited to: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with biliary stent placement include, but are not limited to: fever • obstruction of the biliary/pancreatic duct • obstruction of the common bile duct • pain/discomfort • stent migration • stent occlusion • trauma to the biliary tract or duodenum • trauma to the pancreatic tract or duodenum.

## **PRECAUTIONS**

A complete diagnostic evaluation must be performed, prior to use, to determine the proper stent size.

Sphincterotomy is not necessary for the device placement.

Refer to the label for the minimum channel size required for this device.

Refer to the label regarding the selection of the appropriate wire guide.

Do not use excessive force to advance the stent. The positioning sleeve/pigtail straightener (where supplied with stent) is not intended for use in the accessory channel of the endoscope.

The device must be used under fluoroscopic monitoring.

Dislodgement of a placed stent is possible when attempting additional procedures.

Do not use this device for any purpose other than stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

## **WARNINGS**

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

## INSTRUCTIONS FOR USE

## Illustrations

### I. FOR NON-PIGTAIL STENTS:

1. Gently ensure the full extension of all the flaps.
2. Load the positioning sleeve onto the duodenal flap end of the stent.

#### For 5Fr to 7Fr stents:

- 3a. Introduce the stent, ductal end first, and the positioning sleeve onto a pre-positioned wire guide.
- 3b. Advance the pushing catheter over the wire guide to advance the stent into the accessory channel. As the front flap of the stent enters the accessory channel, keep the positioning sleeve over the back flap until the stent is completely in the accessory channel.

#### For 8.5 Fr and larger stents:

- 3c. Remove the hub from the end of the guiding catheter, then introduce the guiding catheter into the accessory channel over a pre-positioned wire guide.
- 3d. With the elevator open, advance until a sufficient length of the guiding catheter is above the obstruction.
- 3e. Introduce the stent, ductal end first, and the positioning sleeve onto the guiding catheter and the pre-positioned wire guide.
- 3f. Advance the pushing catheter over the wire guide and guiding catheter to advance the stent into accessory channel. As the front flap of the stent enters the accessory channel, keep the positioning sleeve over the back flap until the stent is completely in the accessory channel.

**Then refer to steps 4-9 below**

### II. FOR PIGTAIL STENTS:

1. Straighten the ductal pigtail. Use the pigtail straightener where provided.
2. Introduce the stent, ductal end first, and the pigtail straightener onto a pre-positioned wire guide until the straightener reaches second curl.
3. Advance the pushing catheter over the wire guide to advance the pigtail stent into the accessory channel.

**Then refer to steps 4-9 below**

4. As the pushing catheter advances the stent completely into the accessory channel, slide the positioning sleeve/pigtail straightener (where supplied) back over the pushing catheter until it reaches the end of catheter, keeping it clear of the accessory channel.
5. Advance the pushing catheter in small increments until the stent is in the desired position.
6. Fluoroscopically and endoscopically confirm the desired stent position. Note: For non-pigtail stents of 8.5 Fr or greater only: If desired, replace the hub on the guiding catheter before flushing with the contrast media.
7. After confirming the stent position, gently remove the wire guide and the guiding catheter (where supplied) from the endoscope while maintaining the position of the stent with the pushing catheter.
8. Gently remove the pushing catheter from the accessory channel.
9. These stents may be removed using standard endoscopic techniques.

**Upon completion of the procedure, remove all the components of the device from the working channel of the endoscope and dispose of the introducer components per institutional guidelines for biohazardous medical waste.**

### HOW SUPPLIED

These devices are supplied Ethylene Oxide (EtO) sterilized in a peel open pouch.

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© COPYRIGHT COOK 2019 ©
**COOK IRELAND LTD.**