
EN 4	Pushing Catheter and Guiding Catheter Instructions for Use
BG 8	Избутващ катетър и водещ катетър Инструкции за употреба
CS 13	Tlačný katetr a vodící katetr Návod k použití
DA 17	Skubbekateter og styrekateter Brugsanvisning
DE 21	Platzierungskatheter und Führungskatheter Gebrauchsanweisung
EL 26	Καθετήρας ώθησης και οδηγός καθετήρας Οδηγίες χρήσης
ES 30	Catéter empujador y catéter guía Instrucciones de uso
ET 35	Suruv kateeter ja juhtkateeter Kasutusjuhend
FR 39	Cathéter pousse-prothèse et cathéter guide Mode d'emploi
HR 43	Potisni kateter i vodeći kateter Upute za uporabu
HU 47	Tolókatéter és vezetőkátéter Használati utasítás
IT 52	Catetere spingitore e catetere guida Istruzioni per l'uso
LT 56	Stumiamasis kateteris ir kreipiamasis kateteris Naudojimo nurodymai
LV 60	Bīdītājkatetrs un vadītājkatetrs Lietošanas instrukcija
NL 64	Pushing-katheter en geleidekatheter Gebruiksaanwijzing
NO 69	Skyvekateter og ledekateter Bruksanvisning
PL 73	Cewnik popychający i cewnik prowadzący Instrukcja użycia
PT 78	Cateter posicionador e cateter guia Instruções de utilização
RO 82	Cateter de împingere și cateter de ghidaj Instrucțiuni de utilizare
SK 86	Zatláčací katéter a vodiaci katéter Návod na použitie
SL 90	Potisni kateter in vodilni kateter Navodila za uporabo
SR 94	Potisni kateter i vodeći kateter Uputstvo za upotrebu
SV 98	Påförarkateter och guidingkateter Bruksanvisning
TR 102	İtici Kateter ve Kılavuz Kateter Kullanma Talimatı

EN Figure 1: Pushing Catheter (PC)

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
PUSHING CATHETER AND GUIDING CATHETER

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

This IFU covers Pushing Catheter and Guiding Catheter introducers.

The guiding catheter is available in 5 Fr and 6 Fr in a length of 320 cm. The pushing catheter is available in 3 Fr, 4 Fr, 5 Fr, 6 Fr, 7 Fr, 8 Fr, 8.5 Fr, 10 Fr and 11 Fr in lengths of 170 cm.

Performance Characteristics

The function and key features of the device are described below.

Guiding Catheter - The function of the guiding catheter is to guide the stent to its intended location.

Guiding catheters contain one or more of the following features:

- **Atraumatic tip** - The end of the guiding catheter is designed to be atraumatic to the anatomy.
- **Radiopacity** - The guiding catheter tubing material is radiopaque. The guiding catheter also contains 2 marker bands to further enhance its fluoroscopic visibility.
- **Tuohy-Borst** - The guiding catheter contains a Tuohy-Borst with a connection to facilitate the attachment of a standard Luer syringe.

Pushing Catheter - The function of the Pushing Catheter is to advance the stent, over a pre-positioned wire guide and a Guiding Catheter, to its intended location, and to maintain the position of the stent as it is being deployed. It contains the following feature:

- **Proximal Identifier** - An identifier at its proximal end. It is used to identify the orientation in which the device must be used.

Device Compatibility

The Pushing Catheter and Guiding Catheter introducers are compatible with the following:

- Endoscopes with 3.2 mm, 3.7 mm and 4.2 mm accessory channel (Refer to label for selection of the appropriate accessory channel)
- 0.018", 0.021" and 0.035" wire guides (Refer to label for selection of the appropriate wire guide size)
- Pigtail Straighteners/Positioning Sleeve (Flap Protector)
- Endoscopic Cap/ Wire Guide Locking Device
- Contrast Media/Saline
- Standard Luer syringe
- Cook Pancreatic & Biliary Stents (**see Table 1**)

Table 1: Compatible Cook Pancreatic & Biliary Stents for the Pushing Catheter and Guiding Catheter

Device Family	Cook Biliary Stents							Cook Pancreatic Stents		
	Cotton-Huibregtse®	Soehendra Tannenbaum®	Cotton-Leung®	Cotton-Leung® (Sof-Flex®)	Plastic Biliary Stent*	Zimmon®	Compass	Geenen®	Geenen® (Sof Flex)	Zimmon®
	CHBSO	TTSO	CLSO	CLSO-SF	PBS*	ZSO	CBBSO	GPSO/ GPSOS	GPSO-SF/ GPSOS-SF	SPSOF/ SPSOS/ ZPSOF/ ZPSOS
Pushing Catheter (PC)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Guiding Catheter (GC-5)	8.5 Fr	8.5 Fr	8.5 Fr	10 Fr	8.5 Fr	X	X	8.5 Fr	X	X
Guiding Catheter (GC-6)	10 Fr 11.5 Fr	10 Fr 11.5 Fr	10 Fr 11.5 Fr	X	10 Fr 11.5 Fr	X	X	10 Fr 11.5 Fr	X	X

***Note:** PBS devices are not available in all markets, including the US.

Convention is as follows;

- The symbol "X" means not indicated for use with.
- Checkmark "✓" means compatible with corresponding stent sizes.
- Fr refers to the stent diameter in French size.

Note: Pushing catheters are compatible with the stents indicated in **Table 1** e.g. PC-7 is compatible with a 7 Fr stent, PC-5 is compatible with a 5 Fr stent etc.

Intended User

The product is intended for use by physicians trained and experienced in ERCP techniques.

Standard techniques for the endoscopic placement of stents should be employed.

Patient Population(s)

Adult patients requiring biliary stenting for obstruction. The nature of the underlying pathology is prevalent in varying patients, therefore the devices are indicated for patients with biliary duct obstruction caused by common bile duct stones, malignant biliary obstruction and benign or malignant strictures. These devices can be used in all patient populations (≥10 kg) requiring pancreatic stenting. The nature of the underlying pathology is prevalent in varying patients, therefore the devices are indicated for patients with pancreatitis, stricture, stones, cancer, anatomic anomalies, disrupted duct, fistula/leak, pancreatic fluid collection and patients in whom prophylactic stenting for prevention of post-ERCP pancreatitis is required.

Contact with Body Tissue

The device is tissue contacting in line with intended use.

Operating Principles

The guiding catheter tracks over a pre-positioned wire guide as it guides the stent to its intended location. It has radiopaque marker bands which enhance the fluoroscopic visibility of the guiding catheter. The pushing catheter advances the stent over the guiding catheter where applicable to its intended location while also maintaining the position of the stent as it is being deployed.

INTENDED USE

This device is used for endoscopic biliary/pancreatic stent placement.

INDICATIONS FOR USE

These devices are indicated for use with biliary and pancreatic stents for the following indications.

For endoscopic biliary stent placement for biliary drainage of obstructed ducts that could be caused by common bile duct stones, malignant biliary obstruction, benign or malignant strictures or other obstructed biliary conditions requiring drainage.

Endoscopic pancreatic stent placement is used for pancreatic duct drainage that could be caused by pancreatitis, stricture, pancreatic cancer, anatomic anomalies of the pancreas, pancreatic fluid collection, pancreatic stones, disrupted duct, fistula / pancreatic leak. Pancreatic stents are also used prophylactically for prevention of post-ERCP pancreatitis.

CLINICAL BENEFITS

The biliary/pancreatic stent introducers have an indirect therapeutic impact as they function to facilitate the successful placement of the stent. The introducers are essential for stent placement.

CONTRAINDICATIONS

Those specific to ERCP.

Inability to pass wire guide / guiding catheter / stent through the obstructed area.

WARNINGS

This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.

Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.

Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

PRECAUTIONS

A complete diagnostic evaluation must be performed, prior to use, to determine the proper stent size.

Sphincterotomy is not necessary for the device use.

Refer to the label for the minimum channel size required for this device.

Refer to the label regarding the selection of the appropriate wire guide.

Care must be exercised when straightening the pigtail curls, if any, in order to avoid kinking or breaking the stent.

Do not use excessive force to advance the stent.

The positioning sleeve/pigtail straightener is not intended for use in the accessory channel of the endoscope.

The device must be used under fluoroscopic monitoring.

Dislodgement of a placed stent is possible when attempting additional procedures.

Do not use this device for any purpose other than stated intended use.

Store the device in a dry location.

Use of this device restricted to a trained healthcare professional.

The product is intended for use by physicians trained and experienced in ERCP techniques.

Standard techniques for the endoscopic placement of stents should be employed.

These introducer sets are only compatible with Cook pancreatic and biliary stents as outlined in **Table 1**.

POTENTIAL ADVERSE EVENTS

Those associated with ERCP: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • liver abscess • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with biliary/pancreatic stent placement: fever • obstruction of the pancreatic duct • obstruction of the common bile duct • pain/discomfort • stent migration • trauma to the biliary duct or duodenum • trauma to the pancreatic duct or duodenum.

HOW SUPPLIED

These devices are supplied Ethylene Oxide (EO) sterilized in a peel open pouch.

DEVICE PREPARATION

1. Prior to introducing any devices ensure working channel of endoscope is lubricated with water or a water based lubricant.
2. Referring to the device package label, ensure that the endoscope has a channel size greater than or equal to the minimum channel size required to operate the device.
3. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
4. Prior to using the device(s) visually inspect it/them for abnormalities that could result in improper function.
5. Attach endoscopic cap or wire guide locking device to endoscope.

INSTRUCTIONS FOR USE

I. FOR NON-PIGTAIL STENTS:

Illustrations

1. Gently ensure the full extension of all the flaps.
2. Load the positioning sleeve onto the duodenal flap end of the stent.

For 5Fr to 7Fr stents:

- 3a. Introduce the stent, ductal end first, and the positioning sleeve onto a pre-positioned wire guide.
- 3b. Advance the pushing catheter over the wire guide to advance the stent into the accessory channel. As the front flap of the stent enters the accessory channel, keep the positioning sleeve over the back flap until the stent is completely in the accessory channel.

For 8.5 Fr and larger stents:

- 3a. Remove the hub from the end of the guiding catheter, then introduce the guiding catheter into the accessory channel over a pre-positioned wire guide.
- 3b. With the elevator open, advance until a sufficient length of the guiding catheter is above the obstruction.
- 3c. Introduce the stent, ductal end first, and the positioning sleeve onto the guiding catheter and the pre-positioned wire guide.
- 3d. Advance the pushing catheter over the wire guide and guiding catheter to advance the stent into accessory channel. As the front flap of the stent enters the accessory channel, keep the positioning sleeve over the back flap until the stent is completely in the accessory channel.

Then refer to steps 4-9 below

II. FOR PIGTAIL STENTS:

1. Straighten the ductal pigtail. Use the pigtail straightener where provided.
2. Introduce the stent, ductal end first, and the pigtail straightener onto a pre-positioned wire guide until the straightener reaches second curl.
3. Advance the pushing catheter over the wire guide to advance the pigtail stent into the accessory channel.

Then refer to steps 4-9 below

4. As the pushing catheter advances the stent completely into the accessory channel, slide the positioning sleeve/pigtail straightener (where supplied) back over the pushing catheter until it reaches the end of catheter, keeping it clear of the accessory channel.
5. Advance the pushing catheter in small increments until the stent is in the desired position.
6. Fluoroscopically and endoscopically confirm the desired stent position. Note: For non-pigtail stents of 8.5 Fr or greater only: If desired, replace the hub on the guiding catheter before flushing with the contrast media.

7. After confirming the stent position, gently remove the wire guide and the guiding catheter (where supplied) from the endoscope while maintaining the position of the stent with the pushing catheter.
8. Gently remove the pushing catheter from the accessory channel.
9. These stents may be removed using standard endoscopic techniques.

DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be coiled for disposal in accordance with institutional guidelines.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2026
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.