
- EN 10** **EchoTip® Insight™ Endoscopic Ultrasound Needle Connecting Tube Female to Male Luer Lock**
 Instructions for Use
- BG 15** **Ендоскопска ултразвукова игла EchoTip® Insight™**
Свързваща тръба, мъжки към женски луер лок
 Инструкции за употреба
- CS 21** **Endoskopická ultrazvuková jehla EchoTip® Insight™**
Spojovací hadička, ženská-mužská spojka Luer lock
 Návod k použití
- DA 27** **EchoTip® Insight™ endoskopisk ultralydssnål**
Forbindelsesslange, hun- til han-Luer lock
 Brugsanvisning
- DE 32** **EchoTip® Insight™ endoskopische Ultraschallnadel**
Verbindungsschlauch, Luer-Lock-Buchse auf -Stecker
 Gebrauchsanweisung
- EL 38** **Ενδοσκοπική βελόνα υπερήχων EchoTip® Insight™**
Συνδετικός σωλήνας, θηλυκή σε αρσενική ασφάλιση Luer
 Οδηγίες χρήσης
- ES 44** **Aguja de ecoendoscopia EchoTip® Insight™**
Tube conector con conexiones Luer Lock de hembra a macho
 Instrucciones de uso
- ET 50** **EchoTip® Insight™-i endoskoopiline ultrahelinõel**
Ühendustoru emase ja isase Luer-lukuga
 Kasutusjuhised
- FR 55** **Aiguille ultrasonore endoscopique EchoTip® Insight™**
Tubulure de connexion, Luer lock femelle raccordé à un Luer lock mâle
 Mode d'emploi
- HU 61** **EchoTip® Insight™ endoszkópos echogén tű**
Összekötő cső hüvelyes-dugós Luer-zárral
 Használati utasítás
- IT 67** **Ago ecoendoscopico EchoTip® Insight™**
Tube connettore, Luer Lock da femmina a maschio
 Istruzioni per l'uso

LT 73	„EchoTip® Insight™“ echoendoskopijos adata Jungiamasis vamzdelis, lizdinė ir kištukinė Luerio jungtys Naudojimo nurodymai
LV 78	EchoTip® Insight™ endoskopiskās ultraskaņas adata Savienojošā caurulīte ar sievišķo un vīrišķo Luer Lock savienotāju Lietošanas instrukcija
NL 83	EchoTip® Insight™ endoscopische ultrasone naald Verbindingsslang, vrouwelijke-naar-mannelijke Luerlock Gebruiksaanwijzing
NO 89	EchoTip® Insight™ endoskopisk ultralydnål Forbindelsesslange, hunn-til-hann-luer-lock Bruksanvisning
PL 95	Endoskopowa igła ultrasonograficzna EchoTip® Insight™ Dren łączący z żeńskim i męskim złączem Luer lock Instrukcja użycia
PT 101	Agulha de ultrassons endoscópica EchoTip® Insight™ Tubo de ligação, Luer-lock fêmea a macho Instruções de utilização
RO 106	Ac de ultrasonografie endoscopică EchoTip® Insight™ Tub conector de la Luer lock „mamă” la Luer lock „tată” Instrucțiuni de utilizare
SK 112	Endoskopická ultrazvuková ihla EchoTip® Insight™ Spojovacia hadička so zásuvkovým a zástrčkovým spojom Luer Lock Návod na použitie
SV 117	EchoTip® Insight™ endoskopisk ultraljudsnål Kopplingsslang med Luer-lås av typ hona till hane Bruksanvisning


Fig. 1

A. **EN Male Luer**

- BG Мъжки лuer
- CS Mužská spojka Luer
- DA Han-Luer
- DE Luer-Stecker
- EL Αρσενικός σύνδεσμος Luer
- ES Conexión Luer macho
- ET Isane Luer
- FR Luer mâle
- HU Dugós Luer
- IT Luer maschio
- LT Kištukinė Luerio jungtis
- LV Vīrišķais Luer savienotājs
- NL Mannelijke Luer
- NO Hann-luer
- PL Męski złącze Luer
- PT Luer macho
- RO Luer „tată”
- SK Zástrčkový spoj Luer
- SV Luer-lås av hantyp

Fig. 2

A. EN Connecting tube female Luer

- BG Свързваща тръба, женски луер
CS Ženská spojka Luer spojovací hadičky
DA Forbindelsesslangens hun-Luer
DE Luer-Buchse des Verbindungsschlauchs
EL Θηλυκός σύνδεσμος Luer συνδετικού σωλήνα
ES Conexión Luer hembra del tubo conector
ET Ühendustoru emane Luer-lukk
FR Luer femelle de la tubulure de connexion
HU Összekötő cső hüvelyes Luer-illesztéke
IT Luer femmina del tubo connettore
LT Jungiamojų vamzdelio lizdinė Luerio jungtis
LV Savienojošās caurulītes sievišķais Luer savienotājs
NL Vrouwelijke Luer van verbindingsslang
NO Forbindelsesslange, hunn-luer
PL Żeńskie złącze Luer drenu łączącego
PT Luer fêmea do tubo de ligação
RO Tub conector cu mecanism de închidere Luer „mamă”
SK Zásuvkový spoj Luer spojovacej hadičky
SV Kopplingsslangens Luer-lås av hontyp

B. EN Connecting tube male Luer

- BG Свързваща тръба, мъжки луер
CS Mužská spojka Luer spojovací hadičky
DA Forbindelsesslangens han-Luer
DE Luer-Stecker des Verbindungsschlauchs
EL Άρσενικός σύνδεσμος Luer συνδετικού σωλήνα
ES Conexión Luer macho del tubo conector
ET Ühendustoru isane Luer-lukk
FR Luer mâle de la tubulure de connexion
HU Összekötő cső dugós Luer-illesztéke
IT Luer maschio del tubo connettore
LT Jungiamojų vamzdelio kištukinė Luerio jungtis
LV Savienojošās caurulītes vīrišķais Luer savienotājs
NL Mannelijke Luer van verbindingsslang
NO Forbindelsesslange, hann-luer
PL Męskie złącze Luer drenu łączącego
PT Luer macho do tubo de ligação
RO Tub conector cu mecanism de închidere Luer „tată”
SK Zástrčkový spoj Luer spojovacej hadičky
SV Kopplingsslangens Luer-lås av hantyp


C. EN Stopcock sideport

- BG Страничен порт на спирателно кранче
CS Postranní port uzavíracího kohoutu
DA Stophanens sidehul
DE Seitenloch des Absperrhahns
EL Πλευρική θύρα στρόφιγγας
ES Orificio lateral de la llave de paso
ET Sulgurventiili küljeava
FR Orifice latéral du robinet
HU Elzárócsap oldalnyílása
IT Raccordo laterale del rubinetto
LT Čiaupo šoninė anga
LV Krāna sānu atvere
NL Zijopening van afsluitkraan
NO Sideport for stoppekran
PL Port boczny kraniaka
PT Orifício lateral da torneira de passagem
RO Orificiu lateral al robinetului
SK Bočný port kohútika
SV Infusionskranens sidoport

D. EN Stopcock male Luer

- BG Мъжки луер на спирателно кранче
CS Mužská spojka Luer uzavíracího kohoutu
DA Stophanens han-Luer
DE Luer-Stecker des Absperrhahns
EL Άρσενικός σύνδεσμος Luer στρόφιγγας
ES Conexión Luer macho de la llave de paso
ET Sulgurventiili isane Luer
FR Luer mâle du robinet
HU Elzárócsap dugós Luer-illesztéke
IT Luer maschio del rubinetto
LT Čiaupo kištukinė Luerio jungtis
LV Krāna vīrišķais Luer savienotājs
NL Mannelijke Luer van afsluitkraan
NO Hann-luer for stoppekran
PL Męskie złącze Luer kraniaka
PT Luer macho da torneira de passagem
RO Luer „tată” al robinetului
SK Zástrčkový spoj Luer kohútika
SV Infusionskranens Luer-lås av hantyp

Fig. 3

**A. EN Sliding sheath adjuster**

BG	Плъзгащ се регулатор на защитната обвивка
CS	Posuvné nastavovací zařízení sheathu
DA	Glidende sheath-justering
DE	Verschiebbarer Schleuseneinsteller
EL	Συρόμενος προσαρμογέας θηκαρίου
ES	Ajustador deslizante de la vaina
ET	Liughülsi regulaator
FR	Ajusteur de gaine coulissant
HU	Elcsúsztható hüvelyhossz-beállító
IT	Regolatore scorrevole della guaina
LT	Slankusis vamzdelio reguliatorius
LV	Pārbidāmais ievadapvalka regulētājs
NL	Schuifregeling van de sheath
NO	Glidende hylsejusteringsmekanisme
PL	Przesuwany regulator koszulki
PT	Regulador deslizante da bainha
RO	Ajustor glisant al tecii
SK	Posuvný nastavovač puzdra
SV	Glidande hylsjusterare

B. EN Needle handle

BG	Дръжка на иглата
CS	Rukojeť jehly
DA	Nålegreb
DE	Nadelgriff
EL	Λαβή βελόνας
ES	Mango de la aguja
ET	Nõela käepide
FR	Poignée de l'aiguille
HU	Tű fogantyúja
IT	Impugnatura dell'ago
LT	Adatos kotelis
LV	Adatas rokturis
NL	Naaldhandvat
NO	Nålhandtak
PL	Uchwyt igły
PT	Punho da agulha
RO	Mânerul acului
SK	Rúčka ihly
SV	Nålhandtag

C. EN Sheath reference mark

BG	Означение за дължината на защитната обвивка
CS	Referenční značka sheathu
DA	Referencemærke på sheath
DE	Schleusenbezugsmarke
EL	Σήμανση αναφοράς θηκαρίου
ES	Marca de referencia de la vaina
ET	Hülsi tähis
FR	Repère de référence sur la gaine
HU	Hüvelyhossz referenciajele
IT	Indicatore di riferimento della guaina
LT	Vamzdelio atskaitos žymė
LV	Ievadapvalka atskaites atzīme
NL	Referentiemarkering van de sheath
NO	Referansemerke på hylsen
PL	Znacznik referencyjny wysunięcia koszulki
PT	Marca de referência da bainha
RO	Marcaj de referință al tecii
SK	Referenčná značka puzdra
SV	Hylsans referensmarkering


D. EN Safety ring

BG	Предпазен пръстен
CS	Pojistný prstenec
DA	Sikkerhedsring
DE	Sicherheitsring
EL	Δακτύλιος ασφαλείας
ES	Anillo de seguridad
ET	Kaitserõngas
FR	Bague de sécurité
HU	Biztonsági gyűrű
IT	Anello di sicurezza
LT	Apsauginis žiedas
LV	Drošības gredzens
NL	Veiligheidsring
NO	Sikkerhetsring
PL	Pierścień bezpieczeństwa
PT	Anel de segurança
RO	Inel de siguranță
SK	Bezpečnostný krúžok
SV	Säkerhetslåsring

E. EN Female Luer

BG	Женски лuer
CS	Ženská spojka Luer
DA	Hun-Luer
DE	Luer-Buchse
EL	Θηλυκός σύνδεσμος Luer
ES	Conexión Luer hembra
ET	Emane Luer
FR	Luer femelle
HU	Hüvelyes Luer-illeszték
IT	Luer femmina
LT	Lizdinė Luerio jungtis
LV	Sievšķais Luer savienotājs
NL	Vrouwelijke Luer
NO	Hunn-luer
PL	Żeńskie złącze Luer
PT	Luer fêmea
RO	Luer „mamă”
SK	Zásuvkový spoj Luer
SV	Luer-lås av hontyp

Fig. 4

**A. EN Syringe**

BG Женски лuer
CS Ženská spojka Luer
DA Hun-Luer
DE Luer-Buchse
EL Θηλυκός σύνδεσμος Luer
ES Conexión Luer hembra
ET Emame Luer
FR Luer femelle
HU Hűvelyes Luer-illeszték
IT Luer femmina
LT Lizdinė Luerio jungtis
LV Sievišķais Luer savienotājs
NL Vrouwelijke Luer
NO Hunn-luer
PL Żeńskie złącze Luer
PT Luer fêmea
RO Luer „mamă”
SK Zásuvkový spoj Luer
SV Luer-lås av hontyp

B. EN Compass pressure transducer

BG Трансдюсер за налягане Compass
CS Snímač tlaku Compass
DA Compass-tryktransducer
DE Compass Druckwandler
EL Μορφοτροπέας πίεσης Compass
ES Transductor de presión Compass
ET Rõhuandur Compass
FR Transducteur de pression Compass
HU Compass nyomásjeladó
IT Trasduttore di pressione Compass
LT „Compass” spaudimo jutiklis
LV Compass spiediena devējs
NL Compass druktransducer
NO Compass-trykktansducer
PL Przetwornik ciśnienia Compass
PT Transdutor de pressão Compass
RO Traductor de presiune Compass
SK Snímač tlaku Compass
SV Compass tryckomvandlare

C. EN Connecting tube

BG Свързваща тръба
CS Spojovací hadička
DA Forbindelsesslange
DE Verbindungsschlauch
EL Συνδετικός σωλήνας
ES Tubo conector
ET Ühendustoru
FR Tubulure de connexion
HU Összekötő cső
IT Tubo connettore
LT Jungiamasis vamzdelis
LV Savienojošā caurulīte
NL Verbindingsslang
NO Forbindelsesslange
PL Dren łączący
PT Tubo de ligação
RO Tub de conectare
SK Spojovacia hadička
SV Kopplingslang


D. EN Three-way stopcock

BG Трйпътнo спйрателнo кранчe
CS Trojcestný uzavírací kohout
DA Trevejsstophane
DE Dreiweg-Absperrhahn
EL Στρόφιγγα τριών οδών
ES Clave de paso de tres vías
ET Kolmepoolne sulgurventiil
FR Robinet à trois voies
HU Háromutas elzárócsap
IT Rubinetto a tre vie
LT Trijų krypčių čiaupas
LV Trīsvirzienu noslēgkrāns
NL Driewegafsluitkraan
NO Treveis stoppekran
PL Kranik trójdrożny
PT Torneira de passagem de três vias
RO Robinet tridirecțional
SK Trojsmerný kohútik
SV Trevägs infusionskran

E. EN EchoTip Insight needle

BG Игла EchoTip Insight
CS Jehla EchoTip Insight
DA EchoTip Insight-nål
DE EchoTip Insight Nadel
EL Βελόνα EchoTip Insight
ES Aguja EchoTip Insight
ET Nõel EchoTip Insight
FR Aiguille EchoTip Insight
HU EchoTip Insight tű
IT Ago EchoTip Insight
LT „EchoTip Insight” adata
LV EchoTip Insight adata
NL EchoTip Insight naald
NO EchoTip Insight-nål
PL Igła EchoTip Insight
PT Agulha EchoTip Insight
RO Ac EchoTip Insight
SK Ihla EchoTip Insight
SV EchoTip Insight-nål

Fig. 5

A. EN Stopcock 'OFF' position while priming connecting tube

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

Fig. 8

ECHOTIP® INSIGHT™ ENDOSCOPIC ULTRASOUND NEEDLE CONNECTING TUBE FEMALE TO MALE LUER LOCK

INTENDED USE

EchoTip Insight Needle: This device is used to access the venous vasculature related to the portal circulation and hepatic outflow through the accessory channel of an ultrasound endoscope to allow for direct pressure measurement.

Connecting Tube: The tube is used for the transfer of liquids between the EchoTip Insight components.

DEVICE DESCRIPTION

The EchoTip Insight Needle and Connecting Tube are packaged together with the Compass pressure transducer in an outer box. Once assembled, these devices are used with an ultrasound endoscope in patients with portal hypertension.

The EchoTip Insight Needle is a device consisting of a handle with adjustable components to allow the user to adjust the extension of the needle and sheath. The device has a 5.2 Fr (1.73 mm) outer sheath for protection of the 25 ga (0.56 mm) needle and is attached to the accessory channel of an endoscope. The device has a female Luer fitting on the needle handle to allow for the connection of the stopcock. The device is supplied with a 10 ml syringe that is attached to a Compass pressure transducer that is connected to both the needle and connecting tube when assembled.

The connecting tube consists of a 90 cm tube, a female Luer fitting, a male Luer fitting and a stopcock. The connecting tube is used to attach the transducer to the needle handle. The stopcock is used to aid priming of the assembled components.

CONTRAINDICATIONS

Those specific to the primary endoscopic procedure to be performed in gaining access to desired target.

This device should only be used to access the hepatic vein or portal vein. Any other use is contraindicated.

POTENTIAL ADVERSE EVENTS

Those associated with an endoscopic ultrasound procedure include: allergic reaction to medication • aspiration • cardiac arrhythmia or arrest • death • fever • hypotension • infection • pain/discomfort • perforation • respiratory depression or arrest.

Those associated with the device include: arteriovenous fistula • cholangitis • damage to blood vessels • embolism • fever • hemorrhage • hypotension • infection • inflammation • nerve damage • pain/discomfort • perforation • pneumoperitoneum • sepsis • septicemia/bacteremia • thrombosis • vessel occlusion • vessel trauma.

WARNINGS

The tip of the needle is sharp and could cause injury to the patient or user if not used with caution.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

This device contains nickel, which may cause allergic reaction in individuals with nickel sensitivity.

Product is sterile if package is unopened or undamaged. Do not use the product if the peel-open packages, or the box they are supplied in, are damaged. Do not use the product if there is doubt as to whether the product is sterile.

Visually inspect the packaging and device. If the packaging is opened or damaged when received, do not use. Visually inspect the device with particular attention to kinks, bends, or breaks. If a device abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Ensure both vessels are visible and accessible prior to initiating procedure. Access to the vasculature should be approached through the liver parenchyma to reduce risk of bleeding.

The risks and potential adverse effects of the procedure in patients with pre-existing conditions (such as those listed below) must be considered relative to the potential benefits of this procedure:

- Patients with active gastrointestinal bleeding which is not controlled by medical intervention.
- Patients with known infection which is not controlled by medical intervention.
- Patients with known portal vein thrombosis.
- When using with heparinized saline in patients with known heparin-induced thrombocytopenia.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Do not use the device supplied for any purpose other than stated intended use.

Refer to the Compass pressure transducer Instructions for Use for technical specifications and detailed operation instructions including warnings, precautions and operating range of transducer.

Use of the devices supplied are restricted to a trained healthcare professional. The EchoTip Insight Needle should only be used by physicians trained in Endoscopic Ultrasound procedures.

The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of device. Failure to retract the needle may result in damage to the endoscope.

The transducer must be positioned at the phlebostatic axis (fourth intercostal space at the mid anterior posterior midpoint of the chest wall) while obtaining a pressure measurement to ensure reading is accurate.

A sterile fluid (i.e. saline or heparinized saline) is used to flush the system once all parts are assembled during system preparation. When using heparinized saline, refer to flushing solution labelling prior to use to ensure all appropriate cautions, warnings and contraindications are observed.

In a well-lit area, ensure no large air bubbles that may occlude the system are present. If air bubbles are observed, they must be removed from the system before proceeding. Air in the system may increase the response time of the transducer to measure pressure within the portal or hepatic vein.

Once assembled, ensure the closed system remains intact and is not opened to atmospheric pressure during the procedure as doing so will affect the measurements.

Device damage such as needle kinking may result in a non-responsive/inaccurate measurement or higher forces when priming fluid. Do not continue to use the product if there is doubt as to whether the product is damaged.

Vessel penetration may lead to intrahepatic bleeding and non-responsive/inaccurate pressure measurement.

SYSTEM PREPARATION

Illustrations

1. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
2. Visually inspect the device components confirming they are free from damage prior to use.
3. Remove all components from their individual packaging.
4. Remove the Luer cap from the transducer.
5. Switch on the transducer as outlined in the Compass pressure transducer Instructions for Use.
6. Fill 10 ml syringe (Fig. 1) with sterile fluid, ensuring that all large air bubbles are removed from the syringe.
7. Attach the prepared syringe to the female Luer fitting of the transducer.
8. Attach the male Luer lock fitting of the transducer to the white female Luer fitting of the connecting tube and secure tightly.
9. Attach the male Luer lock fitting of the connecting tube stopcock (Fig. 2) to the female Luer fitting on the needle handle (Fig. 3) and secure tightly. **Note:** The system is now fully assembled and should correspond with the image shown in Fig. 4.
10. Prime the system up to the stopcock by slowly advancing plunger into syringe ensuring no large air bubbles are present (Fig. 5 & Fig. 6). **Note:** Fluid should be visible exiting the sideport of the stopcock.

11. Adjust the stopcock to close the sideport before proceeding to prime the needle with minimum 1 ml of sterile fluid (Fig. 6 & Fig. 7). **Note:** Fluid should be visible exiting the distal tip of the sheath.
12. Inspect system to ensure all connections are tight, fully engaged and there are no large air bubbles present in the system. **Note:** If air is present in the tubing or Luer connections, repeat Steps 10-12.
13. Identify the patient's phlebostatic axis and level the transducer at this point. **Caution:** Position the transducer at this level while obtaining the pressure measurement. Any movement will change current or subsequent pressure measurements.

If a syringe exchange is required during the procedure, complete the following steps:

- a. If the device is within the endoscope, remove the device by disconnecting the Luer lock fitting from the accessory channel port by rotating the device handle counter clockwise and withdraw the entire device from the endoscope.
- b. Disconnect the syringe from the female Luer fitting of the transducer.
- c. Fill 10 ml syringe (Fig. 1) with sterile fluid, ensuring that all large air bubbles are removed from the syringe and attach the syringe to the female Luer fitting of the transducer.
- d. Repeat Steps 10 through 12 of "System Preparation".

INSTRUCTIONS FOR USE

1. Identify desired vasculature by endoscopic ultrasound where a triplicate pressure measurement will be taken. **Note:** Ensure both vessels are visible and accessible through the liver parenchyma prior to proceeding with the following steps.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the device into the accessory channel of the endoscope. **Caution:** If resistance is encountered on device introduction, reduce angulation of scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
4. Attach the device to the endoscope accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is endoscopically visible, confirming that the sheath has emerged from the working channel of the scope. To adjust sheath length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (Fig. 3). Tighten thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.

6. While maintaining position of the ultrasound endoscope, set the needle to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of safety ring (Fig. 8). Tighten thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the endoscope. Failure to attach the device prior to the needle adjustment or extension may result in damage to the endoscope.
7. Extend the needle by advancing the needle handle to the pre-positioned safety ring into desired vessel under ultrasound guidance, either hepatic or portal vein. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction. **Note:** Following puncture continue to visualize puncture site to ensure no leakage and continuously observe throughout the procedure.
8. Using the prepared syringe, instill up to 0.5 ml of sterile fluid while observing the ultrasound view. Confirm location within the vessel and ensure that the needle tip is not against the vessel wall as this may alter pressure readings.
9. Maintaining position of both the needle tip and the transducer, allow time for pressure reading to equilibrate. Record the pressure reading displayed on the transducer. **Caution:** Ensure sufficient time is given (minimum 60 seconds) to the device to allow it to stabilize before recording the measurement from the attached transducer. **Note:** For subsequent measurements instill up to 0.5 ml of sterile fluid to flush the system for the next measurement.
10. Without retracting the needle, repeat Step 9 to obtain triplicate measurements while maintaining this position in the vessel. **Note:** Measurements in the portal and hepatic vein are completed in triplicate to ensure values obtained are reproducible and to record the average of the three values. If they are not reproducible, additional measurements should be obtained. If the additional measurements are still not reproducible or additional sterile fluid is required, follow Steps a–d of “System Preparation” before obtaining additional measurements.
11. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place. (Fig. 8)
12. Reposition the endoscope and identify the second vessel to be accessed (opposite of the vessel accessed in Step 7). Obtain triplicate measurements in this vessel by repeating Steps 6 through 10 of the “Instructions for Use”.
13. Upon completion of the procedure, retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.

14. To remove the device from the endoscope, disconnect the Luer lock fitting from the accessory channel port by rotating the device handle counter clockwise and withdraw the entire device from the endoscope.
15. Calculate and record the portal pressure gradient (the difference between the average hepatic vein and average portal vein measurements).

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Supplied together in non-sterile box. Individual devices supplied sterilized by ethylene oxide gas in peel-open packages. Only devices supplied should be used for procedure. Store in a dark, dry, cool place. Avoid extended exposure to light.

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© COOK 2024
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.