- EN 11** **EchoTip Insight® Endoscopic Ultrasound Needle Connecting Tube Female to Male Luer Lock**
 Instructions for Use
- CS 16** **Endoskopická ultrazvuková jehla EchoTip Insight®**
Spojovací hadička, ženská-mužská spojka Luer lock
 Návod k použití
- DA 21** **EchoTip Insight® endoskopisk ultralydsnål**
Forbindelsesslange, hun- til han-Luer lock
 Brugsanvisning
- DE 25** **EchoTip Insight® endoskopische Ultraschallnadel**
Verbindungsschlauch, Luer-Lock-Buchse auf -Stecker
 Gebrauchsanweisung
- EL 31** **Ενδοσκοπική βελόνα υπερήχων EchoTip Insight®**
Συνδετικός σωλήνας, θηλυκή σε αρσενική ασφάλιση Luer
 Οδηγίες χρήσης
- ES 36** **Aguja de ecoendoscopia EchoTip Insight®**
Tubo conector con conexiones Luer Lock de hembra a macho
 Instrucciones de uso
- FI 42** **Endoskooppisen EchoTip Insight® -ultraäänineulan**
Liitäntäletku, naaras-urokseen-luer-lukko
 Käyttöohjeet
- FR 46** **Aiguille ultrasonore endoscopique EchoTip Insight®**
Tubulure de connexion, Luer lock femelle raccordé à un Luer lock mâle
 Mode d'emploi
- IT 52** **Ago ecoendoscopico EchoTip Insight®**
Tubo connettore, Luer Lock da femmina a maschio
 Istruzioni per l'uso
- NL 57** **EchoTip Insight® endoscopische ultrasone naald**
Verbindingslang, vrouwelijke-naar-mannelijke Luerlock
 Gebruiksaanwijzing
- NO 63** **EchoTip Insight® endoskopisk ultralydsnål**
Forbindelsesslange, hunn-til-hann-luer-lock
 Bruksanvisning

PL
67

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






Fig. 6

EN Vertical device orientation while priming system

- CS Vertikální orientace prostředku při plnění systému
 DA Produktet holdes lodret under priming af systemet
 DE Vertikale Haltung des Instruments beim Entlüften des Systems
 EL Προσανατολισμός της συσκευής σε κάθετη θέση κατά την αρχική πλήρωση του συστήματος
 ES Orientación vertical del dispositivo durante el cebado del sistema
 FI Laitteen pystysuunta järjestelmän esitäytön aikana
 FR Orientation verticale du dispositif pendant l'amorçage du système
 IT Orientamento verticale del dispositivo durante il priming del sistema
 NL Verticale oriëntatie van hulpmiddel tijdens voorvullen
 NO Vertikal innretting av anordning under priming av system
 PL Pionowe ustawienie urządzenia podczas wypelniania systemu
 PT Orientação vertical do dispositivo ao encher o sistema
 SV Produkten hålls lodrätt när systemet fylls


Fig. 7

A. **EN Stopcock 'OFF' position while priming needle**

- CS Uzavírací kohout v poloze „OFF“ (Vypnuto) při plnění jehly
 DA Stophane i 'OFF'-position under priming af nålen
 DE Stellung „OFF“ (Aus) des Absperrhahns beim Entlüften der Nadel
 EL Θέση «OFF» (Κλειστή) της στρόφιγγας κατά την αρχική πλήρωση της βελόνας
 ES Posición «OFF» (cerrado) de la llave de paso durante el cebado de la aguja
 FI Sulkuhanan "OFF"-asento esitytettäessä neulaa
 FR Position « OFF » (fermée) du robinet pendant l'amorçage de l'aiguille
 IT Rubinetto in posizione "OFF" (chiuso) durante il priming dell'ago
 NL 'OFF'-positie van afsluitkraan tijdens voorvullen van naald
 NO Stoppekranens OFF-posisjon under priming av nål
 PL Kranik w pozycji „OFF” (zamknięty) podczas wypelniania igły
 PT Posição “OFF” (Fechado) da torneira de passagem ao purgar a agulha
 SV Infusionskran läge "OFF" när nålen fylls

Fig. 8

ECHOTIP INSIGHT® ENDOSCOPIC ULTRASOUND NEEDLE CONNECTING TUBE FEMALE TO MALE LUER LOCK

DEVICE DESCRIPTION

The EchoTip Insight needle and connecting tube are packaged together with the Compass Pressure Transducer in an outer box. Once assembled, these devices are used with an ultrasound endoscope. The polyetheretherketone (PEEK) outer sheath protects the stainless steel needle and ultrasound endoscope from damage.

Performance Characteristics

The EchoTip Insight needle is attached to the accessory channel of an ultrasound endoscope. The handle contains adjustable components to allow the user to adjust the extension of the needle and sheath. The needle has a 5.2 Fr (1.73 mm) outer sheath which has an adjustable length of 0-5 cm for protection of the 25 Ga (0.56 mm) needle and ultrasound endoscope. The needle has an extension range of 0-8 cm.

The EchoTip Insight needle has Luer fittings on the handle to attach the connecting tube and to connect the needle to the ultrasound endoscope accessory channel. The needle has an echogenic needle pattern at its distal tip to aid visibility under endoscopic ultrasound (EUS) when puncturing the venous vasculature. The EchoTip Insight needle is supplied with a 10 mL standard Luer syringe. The syringe is attached to the Compass Pressure Transducer.

The 90 cm connecting tube contains a female Luer fitting, a male Luer fitting and a stopcock. The connecting tube is used to attach the transducer to the needle handle. The stopcock aids priming of the assembled components and transfers sterile fluid, i.e. saline or heparinized saline, between the Compass Pressure Transducer and EchoTip Insight needle.

Device Compatibility

The EchoTip Insight needle and connecting tube are compatible with the following:

- Ultrasound endoscope with 2.8 mm accessory channel
- Compass Pressure Transducer (Compass UniversalHg - CUHG1-INT)
- Sterile fluid (i.e. saline or heparinized saline)
- 10 mL standard Luer syringe

Patient Population

Adult patients with portal hypertension.

Intended User

The product is intended for use by physicians trained and experienced in endoscopic ultrasound techniques.

Standard techniques for endoscopic ultrasound procedures should be employed.

Contact with Body Tissue

The device has contact with body tissue and blood in line with its intended use.

Operating Principles

The EchoTip Insight needle and connecting tube provide a lumen through which a column of sterile fluid can be primed. The pressure applied at the distal tip of the needle is measured by the attached Compass Pressure Transducer. Once the EchoTip Insight needle and connecting tube are connected to the Compass Pressure Transducer, the needle is inserted into the accessory channel of an ultrasound endoscope.

The Luer fitting on the base of the needle handle is attached to the ultrasound endoscope accessory channel. The sheath can be adjusted using the sheath adjuster thumbscrew on the needle handle to ensure it exits the accessory channel. The proximal thumbscrew on the needle handle is used to set the desired needle extension length. The needle is extended into the vasculature by advancing the

handle. The echogenic needle pattern is visible under endoscopic ultrasound to confirm the location of the needle tip.

The pressure applied from the blood flow within the vasculature to the column of fluid within the system is transferred to the transducer at the proximal end of the needle handle. Once measurements are recorded, the needle may be retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark.

INTENDED USE

EchoTip Insight Needle: This device is used to access the hepatic and portal veins or their branches, via a transgastric-transhepatic approach or transduodenal-transhepatic approach at the duodenal bulb, through the accessory channel of an ultrasound endoscope to allow for direct pressure measurement.

Connecting Tube: The tube is used for the transfer of liquids between the EchoTip Insight components.

INDICATIONS FOR USE

Adult patients with portal hypertension.

CLINICAL BENEFITS

The EchoTip Insight Endoscopic Ultrasound Needle and Connecting Tube allow for direct portal pressure measurement under endoscopic ultrasound.

CONTRAINDICATIONS

Those specific to the primary endoscopic procedure to be performed in gaining access to desired target.

This device should only be used to access the hepatic or portal veins or their branches. Any other use is contraindicated.

WARNINGS

- The tip of the needle is sharp and could cause injury to the patient or user if not used with caution.
- This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of the device.
- This device contains nickel, which may cause allergic reaction in individuals with nickel sensitivity.
- Product is sterile if package is unopened and undamaged. Do not use the product if the peel-open packages, or the box they are supplied in, are damaged. Do not use the product if there is doubt as to whether the product is sterile.
- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect the device with particular attention to kinks, bends, or breaks. If a device abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for a return authorization.
- Ensure both vessels are visible and accessible prior to initiating procedure. Access to the vasculature should be approached through the liver parenchyma to reduce risk of bleeding.
- The risks and potential adverse effects of the procedure in patients with pre-existing conditions (such as those listed below) must be considered relative to the potential benefits of this procedure:
 - Patients with active gastrointestinal bleeding which is not controlled by medical intervention.
 - Patients with known infection which is not controlled by medical intervention.
 - Patients with known portal vein thrombosis.
 - When using with heparinized saline in patients with known heparin-induced thrombocytopenia.
 - Device not indicated for use in pregnant/nursing women.

PRECAUTIONS

- Refer to package label for minimum channel size required for this device.
- Do not use the device supplied for any purpose other than stated intended use.
- Refer to the Compass Pressure Transducer Instructions for Use for technical specifications and detailed operation instructions including warnings, precautions and operating range of transducer.
- Use of the devices supplied are restricted to a trained healthcare professional. The EchoTip Insight needle should only be used by physicians trained in endoscopic ultrasound procedures.
- The product is intended for use by physicians trained and experienced in endoscopic ultrasound techniques.
- Standard techniques for endoscopic ultrasound procedures should be employed.
- The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of device. Failure to retract the needle may result in damage to the endoscope.
- The transducer must be positioned at the phlebostatic axis (fourth intercostal space at the mid anterior posterior midpoint of the chest wall) while obtaining a pressure measurement to ensure reading is accurate.
- A sterile fluid (i.e. saline or heparinized saline) is used to flush the system once all parts are assembled during system preparation. When using heparinized saline, refer to flushing solution labelling prior to use to ensure all appropriate cautions, warnings and contraindications are observed.
- In a well-lit area, ensure no large air bubbles that may occlude the system are present. If air bubbles are observed, they must be removed from the system before proceeding. Air in the system may increase the response time of the transducer to measure pressure within the portal or hepatic veins or their branches.
- Once assembled, ensure the closed system remains intact and is not opened to atmospheric pressure during the procedure as doing so will affect the measurements.
- Device damage such as needle kinking may result in a non-responsive/inaccurate measurement or higher forces when priming fluid. Do not continue to use the product if there is doubt as to whether the product is damaged.
- Vessel penetration may lead to intrahepatic bleeding and non-responsive/inaccurate pressure measurement.
- Store the device in a dry location and keep out of sunlight.

POTENTIAL ADVERSE EVENTS

Those associated with an endoscopic ultrasound procedure: allergic reaction to medication • aspiration • cardiac arrhythmia or arrest • death • fever • hypotension • infection • pain/discomfort • perforation • respiratory depression or arrest.

Those associated with the device: allergic reaction to nickel • arteriovenous fistula • cholangitis • damage to blood vessels • embolism • fever • hemorrhage • hypotension • infection • inflammation • nerve damage • pain/discomfort • perforation • pneumoperitoneum • sepsis • septicemia/bacteremia • thrombosis • vessel occlusion • vessel trauma.

HOW SUPPLIED

The EchoTip Insight Endoscopic Ultrasound Needle and Connecting Tube are supplied ethylene oxide (EO) sterilized in a sealed tray and a peel-open pouch, respectively, and packaged together in a non-sterile box. Only devices supplied should be used for procedure.

SYSTEM PREPARATION

Illustrations

1. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
2. Visually inspect the device components confirming they are free from damage prior to use.
3. Remove all components from their individual packaging.
4. Remove the Luer cap from the transducer.

5. Switch on the transducer as outlined in the Compass Pressure Transducer Instructions for Use.
6. Fill 10 mL syringe (Fig. 1) with sterile fluid, ensuring that all large air bubbles are removed from the syringe.
7. Attach the prepared syringe to the female Luer fitting of the transducer.
8. Attach the male Luer lock fitting of the transducer to the white female Luer fitting of the connecting tube and secure tightly.
9. Attach the male Luer lock fitting of the connecting tube stopcock (Fig. 2) to the female Luer fitting on the needle handle (Fig. 3) and secure tightly. **Note:** The system is now fully assembled and should correspond with the image shown in Fig. 4.
10. Prime the system up to the stopcock by slowly advancing plunger into syringe ensuring no large air bubbles are present (Fig. 5 & Fig. 6). **Note:** Fluid should be visible exiting the sideport of the stopcock.
11. Adjust the stopcock to close the sideport before proceeding to prime the needle with minimum 1 mL of sterile fluid (Fig. 6 & Fig. 7). **Note:** Fluid should be visible exiting the distal tip of the sheath.
12. Inspect system to ensure all connections are tight, fully engaged and there are no large air bubbles present in the system. **Note:** If air is present in the tubing or Luer connections, repeat Steps 10-12.
13. Identify the patient's phlebostatic axis and level the transducer at this point. **Caution:** Position the transducer at this level while obtaining the pressure measurement. Any movement will change current or subsequent pressure measurements.

If a syringe exchange is required during the procedure, complete the following steps:

- a. If the device is within the endoscope, remove the device by disconnecting the Luer lock fitting from the accessory channel port by rotating the device handle counterclockwise and withdraw the entire device from the endoscope.
- b. Disconnect the syringe from the female Luer fitting of the transducer.
- c. Fill 10 mL syringe (Fig. 1) with sterile fluid, ensuring that all large air bubbles are removed from the syringe and attach the syringe to the female Luer fitting of the transducer.
- d. Repeat Steps 10 through 12 of "System Preparation".

INSTRUCTIONS FOR USE

1. Identify desired vasculature by endoscopic ultrasound where a triplicate pressure measurement will be taken. **Note:** Ensure both vessels are visible and accessible through the liver parenchyma prior to proceeding with the following steps.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the device into the accessory channel of the endoscope. **Caution:** If resistance is encountered on device introduction, reduce angulation of scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
4. Attach the device to the endoscope accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is endoscopically visible, confirming that the sheath has emerged from the working channel of the scope. To adjust sheath length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (Fig. 3). Tighten thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining position of the ultrasound endoscope, set the needle to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of safety ring (Fig. 8). Tighten thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure

the device has been attached to the accessory channel of the endoscope. Failure to attach the device prior to the needle adjustment or extension may result in damage to the endoscope.

7. Extend the needle by advancing the needle handle to the pre-positioned safety ring into desired vessel under ultrasound guidance, either hepatic or portal veins or their branches.
Caution: If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction. **Note:** Following puncture continue to visualize puncture site to ensure no leakage and continuously observe throughout the procedure.
8. Using the prepared syringe, instill up to 0.5 mL of sterile fluid while observing the ultrasound view. Confirm location within the vessel and ensure that the needle tip is not against the vessel wall as this may alter pressure readings.
9. Maintaining position of both the needle tip and the transducer, allow time for pressure reading to equilibrate. Record the pressure reading displayed on the transducer. **Caution:** Ensure sufficient time is given (minimum 60 seconds) to the device to allow it to stabilize before recording the measurement from the attached transducer. **Note:** For subsequent measurements instill up to 0.5 mL of sterile fluid to flush the system for the next measurement.
10. Without retracting the needle, repeat Step 9 to obtain triplicate measurements while maintaining this position in the vessel. **Note:** Measurements in the portal and hepatic veins or their branches are completed in triplicate to ensure values obtained are reproducible and to record the average of the three values. If they are not reproducible, additional measurements should be obtained. If the additional measurements are still not reproducible or additional sterile fluid is required, follow Steps a–d of “System Preparation” before obtaining additional measurements.
11. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place. (Fig. 8)
12. Reposition the endoscope and identify the second vessel to be accessed (opposite of the vessel accessed in Step 7). Obtain triplicate measurements in this vessel by repeating Steps 6 through 10 of the “Instructions for Use”.
13. Upon completion of the procedure, retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
14. To remove the device from the endoscope, disconnect the Luer lock fitting from the accessory channel port by rotating the device handle counterclockwise and withdraw the entire device from the endoscope.
15. Separate the device components and prepare for disposal. **Note:** Adequate care should be taken when disconnecting the device components to prevent contamination or spillage of fluid remaining in the needle and connecting tube.
16. Calculate and record the portal pressure gradient (the difference between the average hepatic vein and average portal vein measurements).

DISPOSAL OF DEVICES

Refer to the Compass Pressure Transducer Instructions for Use for transducer disposal guidelines.

The EchoTip Insight needle and connecting tube may be contaminated with potentially infectious substances of human origin and should be coiled for disposal in accordance with institutional guidelines.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2025
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.