EN
 3
 BG
 5
 HR
 8
 CS
 11
 DA
 14
 NL
 16
 ET
 19
 FR
 22
 DE
 25
 EL
 28
 HU
 31
 IT
 34
 LV
 37
 LT
 39
 NO
 41
 PL
 43
 PT
 46
 RO
 49
 SK
 52
 ES
 55
 SV
 58

Pushing Catheter

Instructions for Use

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Zatlačací katéter

Catéter empujador

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CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Table 1: Compatible Cook Biliary Stents for the Pushing Catheter

Device Family	Cook Biliary Stents				
	Cotton-Huibregtse®	Cotton-Leung®	Cotton-Leung® (Sof-Flex®)	Plastic Biliary Stent*	Zimmon®
	CHBSO	CLSO	CLSO-SF	PBS*	ZSO
Pushing Catheter (PC-E)	√	√	√	√	√

* Note: PBS devices are not available in all markets, including the US.

Convention is as follows:

- Checkmark "√" means compatible with corresponding stent sizes.

Note: Pushing catheters are compatible with the stents indicated in **Table 1** e.g. PC-7E is compatible with a 7Fr stent, PC-5E is compatible with a 5Fr stent etc.

INTENDED USE

This device is used for endoscopic biliary stent placement.

DEVICE DESCRIPTION

This IFU covers Pushing Catheter (PC-E) introducers represented by **Figure 1**. See **Table 1** for compatible Cook Biliary Stents. For all stent information including orientation, refer to stent instructions for use.

CONTRAINDICATIONS

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Inability to pass wire guide / stent through the obstructed area.

POTENTIAL ADVERSE EVENTS

Those associated with ERCP include, but are not limited to: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with biliary stent placement include, but are not limited to: fever • obstruction of the pancreatic duct • pain/discomfort • stent migration • trauma to the biliary duct or duodenum.

PRECAUTIONS

A complete diagnostic evaluation must be performed, prior to use, to determine the proper stent size.

Refer to the label for the minimum channel size required for this device.

Refer to the label regarding the selection of the appropriate wire guide.

Do not use excessive force to advance the stent. The positioning sleeve/pigtail straightener (where supplied with stent) is not intended for use in the accessory channel of the endoscope.

Dislodgement of a placed stent is possible when attempting additional procedures.

Do not use this device for any purpose other than stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

WARNINGS

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

INSTRUCTIONS FOR USE

Illustrations

I. FOR NON-PIGTAIL STENTS:

1. Gently ensure the full extension of all the flaps.
2. Load the positioning sleeve onto the duodenal flap end of the stent.
- 3a. Introduce the stent, ductal end first, and the positioning sleeve onto a pre-positioned wire guide.
- 3b. Advance the pushing catheter over the wire guide to advance the stent into the accessory channel. As the front flap of the stent enters the accessory channel, keep the positioning sleeve over the back flap until the stent is completely in the accessory channel.

Then refer to steps 4-9 below

II. FOR PIGTAIL STENTS:

1. Straighten the ductal pigtail. Use the pigtail straightener where provided.
2. Introduce the stent, ductal end first, and the pigtail straightener onto a pre-positioned wire guide until the straightener reaches second curl.
3. Advance the pushing catheter over the wire guide to advance the pigtail stent into the accessory channel.
4. As the pushing catheter advances the stent completely into the accessory channel, slide the positioning sleeve/pigtail straightener (where supplied) back over the pushing catheter until it reaches the end of catheter, keeping it clear of the accessory channel.
5. Advance the pushing catheter in small increments until the stent is in the desired position.
6. Fluoroscopically and endoscopically confirm the desired stent position.
7. After confirming the stent position, gently remove the wire guide from the endoscope while maintaining the position of the stent with the pushing catheter.
8. Gently remove the pushing catheter from the accessory channel.
9. These stents may be removed using standard endoscopic techniques.

Upon completion of the procedure, remove all the components of the device from the working channel of the endoscope and dispose of the introducer components per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

These devices are supplied Ethylene Oxide (EtO) sterilized in a peel open pouch.

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© COPYRIGHT COOK 2020 ©
COOK IRELAND LTD.