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
OASIS™ AND FUSION® OASIS™ (ONE ACTION STENT INTRODUCTION SYSTEM)

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

This IFU covers the OASIS One Action Stent Introduction Systems and Fusion OASIS Introduction Systems represented in **Figures 1 and 2**. See **Table 1** for compatible Cook Biliary Stents.

The OASIS One Action Stent Introduction System (OA) is an introduction system that is available in 8.5 Fr, 10 Fr and 11 Fr. It consists of a guiding catheter component of length 205 cm and a pushing catheter component of length 170 cm. The guiding catheter has 4 radiopaque marker bands.

The OASIS One Action Stent Introduction System (OA-E) is an introduction system that is available in 8.5 Fr and 10 Fr. It consists of a guiding catheter component of lengths 318 cm and 320 cm and a pushing catheter component of lengths 293cm and 295 cm. The guiding catheter has 4 radiopaque marker bands.

The Fusion OASIS One Action Stent Introduction System (FS-OA) is an introduction system that is available in 8.5 Fr and 9 Fr. It consists of a guiding catheter component of length 195 cm and a pushing catheter component of length 170 cm. The guiding catheter has 1 radiopaque marker band.

Performance Characteristics

The function and key features of the device are described below.

Guiding Catheter - The function of the guiding catheter is to guide the stent to its intended location.

Guiding catheters contain one or more of the following features:

- Atraumatic tip - The end of the guiding catheter is designed to be atraumatic to the anatomy.
- Radiopacity - The guiding catheter tubing material is radiopaque. The guiding catheter also contains marker bands to further enhance its fluoroscopic visibility.
- Proximal Wire Port (Hub) - Used to facilitate the attachment and detachment of the stylet and allows the wire guide to exit in the long wire configuration.
- IDE Port - The IDE port allows for use with shorter wire guide configurations (FS-OA only).

Pushing Catheter - The function of the Pushing Catheter is to advance the Stent over a pre-positioned wire guide and a Guiding Catheter, to its intended location, and to maintain the position of the Stent as it is being deployed. It contains the following feature:

- Proximal Identifier - An identifier at its proximal end. It is used to identify the orientation in which the device must be used.

Stylet - Enhances the pushability of the device when used in a 'short-wire' configuration (FS-OA only).

Device Compatibility

These devices are compatible with the following:

- Endoscope with 3.2 mm, 3.7 mm and 4.2 mm accessory channel (refer to label for selection of the appropriate accessory channel)
- 0.035" Wire Guide (refer to label for selection of the appropriate wire guide size)
- Positioning Sleeve (Flap Protector)
- Endoscopic Cap/Wire Guide Locking Device
- Cook Biliary Stents (see **Table 1**)

Note: When using the Fusion OASIS Introduction System (FS-OA) with the CLSO and CHBSO stent families; they are not compatible with stent lengths greater than 14 cm.

Note: Introduction Systems are compatible with the stents indicated in **Table 1**, e.g. FS-OA-10 is compatible with a 10 Fr stent, FS-OA-8.5 is compatible with a 8.5 Fr stent.

Table 1: Compatible Cook Biliary Stents for the OASIS™ and Fusion® OASIS™ One Action Stent Introduction Systems

Device Family	Cook Biliary Stents				
	Cotton-Huibregetse®	Soehendra Tannenbaum®	Cotton-Leung®	Cotton- Leung® Sof-Flex®	Plastic Biliary Stent*
	CHBSO	TTSO	CLS0	CLS0-SF	PBS*
OASIS™ Introduction System (OA)	8.5 Fr 10 Fr 11.5 Fr	8.5 Fr 10 Fr 11.5 Fr	8.5 Fr 10 Fr 11.5 Fr	X	8.5 Fr 10 Fr 11.5 Fr
OASIS™ Introduction System (OA-E)	8.5 Fr 10 Fr	8.5 Fr 10 Fr	8.5 Fr 10 Fr	X	8.5 Fr 10 Fr
Fusion® OASIS™ Introduction System (FS-OA)	8.5 Fr 10 Fr	8.5 Fr 10 Fr	8.5 Fr 10 Fr	10 Fr	8.5 Fr 10 Fr

*Note: PBS devices are not available in all markets, including the US.

Convention is as follows:

- The symbol "X" means not indicated for use with.
- Fr refers to stent diameter in French size.

Patient Population(s)

These devices can be used in adult patients requiring placement of a biliary stent. The nature of the underlying pathology is prevalent in varying patients, therefore the devices are indicated for patients with biliary duct obstruction caused by common bile duct stones, malignant biliary obstruction and benign or malignant strictures.

Intended User

The product is intended for use by physicians trained and experienced in ERCP techniques. Standard techniques for the endoscopic placement of stents should be employed.

Contact with body tissue

The device is tissue contacting in line with intended use.

Operating Principles

The function of the Pushing Catheter is to advance the stent, over a pre-positioned wire guide to its intended location, and to maintain the position of the stent as it is being deployed.

The function of the guiding catheter is to guide the stent to its intended location.

The FS-OA may be used for either 'long-wire' or 'short-wire' ERCP techniques.

Long-wire: Upon removal of the stylet the FS-OA may be used in an over-the-wire configuration, whereby the full device is loaded onto the wire guide until the wire guide exits at the proximal wire port (hub).

Short-wire: The FS-OA contains an IDE port. By allowing the wire guide to exit its intended lumen through this port, only the distal section of the device is loaded onto the wire guide during use, thereby facilitating the 'short-wire' ERCP technique. The stylet enhances pushability when placing the stent.

INTENDED USE

This device is used for endoscopic biliary stent placement.

INDICATIONS FOR USE

For endoscopic biliary stent placement for biliary drainage of obstructed ducts that could be caused by common bile duct stones, malignant biliary obstruction, benign or malignant strictures or other obstructed biliary conditions requiring drainage.

CLINICAL BENEFITS

The biliary stent introducers have an indirect therapeutic impact as they function to facilitate the successful placement of a biliary stent. The introducers are essential for stent placement.

CONTRAINDICATIONS

- Those specific to ERCP.
- Inability to pass a wire guide/ guiding catheter/ stent through the obstructed area.

WARNINGS

- Fusion OASIS Introduction System contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.
- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

PRECAUTIONS

- A complete diagnostic evaluation must be performed prior to use to determine the proper stent size.
- Refer to the label for the minimum channel size required for this device.
- Refer to the label for the appropriate wire guide size required for this device.
- Compatible only with Cook non-pigtail biliary stents as outlined in **Table 1**.
- This device must be used under fluoroscopic monitoring.
- Do not use excessive force to advance the stent.
- The positioning sleeve (supplied with the stent) is not intended for use in the accessory channel of the endoscope.
- Sphincterotomy is not necessary for device use.
- Dislodgement of a placed stent is possible when attempting additional procedures.
- Do not use this device for any purpose other than stated intended use.
- Store the device in a dry location.
- Use of this device is restricted to a trained healthcare professional.
- The product is intended for use by physicians trained and experienced in ERCP techniques.
- Standard techniques for the endoscopic placement of stents should be employed.

POTENTIAL ADVERSE EVENTS

Those associated with ERCP: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • liver abscess • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with biliary stent placement: allergic reaction to nickel (FS-OA only) • fever • obstruction of the pancreatic duct • pain/discomfort • stent migration • trauma to the biliary duct or duodenum.

HOW SUPPLIED

These devices are supplied Ethylene Oxide (EO) sterilized in a peel open pouch.

DEVICE PREPARATION

1. Prior to introducing any devices ensure working channel of endoscope is lubricated with water or a water based lubricant.
2. Referring to the device package label, ensure that the endoscope has a channel size greater than or equal to the minimum channel size required to operate the device.
3. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
4. Prior to using the device(s) visually inspect it/them for abnormalities that could result in improper function.
5. Attach endoscopic cap or wire guide locking device to endoscope.
6. Load positioning sleeve onto duodenal end of the stent in order to collapse the flaps. The stent should fit into the flap protector with the entirety of the duodenal flaps collapsing and being covered.
7. Lubricate the wire guide locking device/endoscopic cap with a water-based lubricant.
8. Pre-load the desired stent and positioning sleeve onto the tip of the introducer. Refer to the stent IFU for the intended placed orientation of the stent in the biliary duct.

I. IF USING INTRA DUCTAL EXCHANGE (IDE) PORT AND A PRE-POSITIONED SHORT WIRE.

1. Unlock the short wire from the Wire Guide Locking Device and advance the introducer with the preloaded stent over a pre-positioned wire guide ensuring that the wire guide exits the guiding catheter at the IDE port.
2. Introduce the device into the endoscope accessory channel and relock the wire guide.
3. With the elevator open, advance the device in short increments until it is endoscopically visualized exiting the scope. **Note:** As the front flap of stent enters the accessory channel, keep the positioning sleeve over the back flap until it is completely in the accessory channel. Slide the positioning sleeve over the pushing catheter, keeping it clear of the accessory channel.
4. Advance the device into the appropriate duct.
5. Fluoroscopically monitor the radiopaque band(s) on the guiding catheter. When a sufficient length of the guiding catheter is above the obstruction, disconnect the guiding and pushing catheter hubs on the introducer. **Note:** A sufficient length of the guiding catheter must be above the obstruction prior to the stent positioning.
6. Maintain the position of the guiding catheter and advance the stent into the desired position using the pushing catheter.
7. Fluoroscopically and endoscopically confirm the desired stent placement.
8. Unlock the wire guide from the Wire Guide Locking Device.
9. Retract the wire guide until it exits the guiding catheter at the IDE port. **Note:** Fluoroscopically visualize the radiopaque band at the IDE port. When the radiopaque distal tip of the wire guide passes the band, a disengagement from the wire guide lumen will occur.
10. Advance the disengaged wire guide to maintain ductal access.
11. Relock the wire guide into the Wire Guide Locking Device.
12. Retract the guiding catheter into the endoscope while maintaining position of the stent using the pushing catheter.
13. Fluoroscopically and endoscopically confirm the desired stent placement.
14. Remove the pushing catheter.

II. IF USING PROXIMAL WIRE PORT (PWP) AND A PRE-POSITIONED LONG WIRE GUIDE. (See Figures 1 and 2):

1. Remove the stylet, if applicable.
2. Advance the introducer with the preloaded stent over a pre-positioned wire guide and into the accessory channel of the endoscope until the proximal end of the wire guide exits PWP. **Note:** If applicable, ensure that the wire guide does not exit at the IDE port.

REFER TO STEPS 3-7 IN "SECTION I", THEN RESUME WITH STEP 3 BELOW:

3. Retract and remove both the guiding catheter and the wire guide while maintaining the position of the stent with the pushing catheter.
4. Confirm the position of the stent once the guiding catheter and the wire guide are removed.
5. Remove the pushing catheter.

These stents may be removed using standard endoscopic techniques.

DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be coiled for disposal in accordance with institutional guidelines.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2024
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.