EN	EchoTip® ClearCore™ EUS Biopsy Needle
6	Instructions for Use
CS	Bioptická jehla EchoTip® ClearCore™ EUS
10	Návod k použití
DA	EchoTip® ClearCore™ EUS-biopsinål
14	Brugsanvisning
DE	EchoTip® ClearCore™ EUS-Biopsienadel
18	Gebrauchsanweisung
EL	Βελόνα βιοψίας EUS EchoTip® ClearCore™
23	Οδηγίες χρήσης
ES	Aguja de biopsia para ecoendoscopia EchoTip® ClearCore™
28	Instrucciones de uso
FI	EchoTip® ClearCore™ -EUS-biopsianeula
32	Käyttöohjeet
FR	Aiguille de biopsie écho-endoscopique EchoTip® ClearCore™
36	Mode d'emploi
HR	Igla za biopsiju EchoTip® ClearCore™ EUS
41	Upute za uporabu
IS	EchoTip® ClearCore™ sýnatökunál fyrir holsjáróm skoðun
45	Leiðbeiningar um notkun
IT	Ago per biopsia con ecoendoscopia EchoTip® ClearCore™
49	Istruzioni per l'uso
NL	EchoTip® ClearCore™ EUS biopsienaald
53	Gebruiksaanwijzing
NO	EchoTip® ClearCore™ EUS-biopsinål
58	Bruksanvisning
PL	Igła biopsyjna EchoTip® ClearCore™ EUS
62	Instrukcja użycia
PT	Agulha de Biopsia EUS EchoTip® ClearCore™
66	Instruções de utilização
RO	Ac de biopsie EUS EchoTip® ClearCore™
71	Instrucţiuni de utilizare
RU	Биопсийная игла EchoTip® ClearCore™ для ультразвукового эндоскопа
75	Инструкция по применению

SL

80

Biopsijska igla EUS EchoTip® ClearCore™

Navodila za uporabo

SV

84

EchoTip® ClearCore™ EUS biopsinål

Bruksanvisning

TR

88

EchoTip® ClearCore™ EUS Biyopsi İğnesi

Kullanma Talimatı

1.


EN Plunger

CS Pist
DA Stempel
DE Kolben
EL Έμβολο
ES Émbolo
FI Mäntä
FR Poussoir
HR Klip
IS Stimpill
IT Stantuffo
NL Plunjer
NO Stempel
PL Tłoczek
PT Êmbolo
RO Piston
RU Поршень
SL Bat
SV Kolv
TR Piston

EN Stopcock "open"

CS „Otevřený“ uzavírací kohout
DA Stophane "åben"
DE Absperrhahn „geöffnet“
EL Στρόφιγγα «ανοικτή»
ES Llave de paso «abierta»
FI Sulkuhana "auki"
FR Robinet ouvert
HR Sigurnosni ventil otvoren
IS Krani „opinn“
IT Rubinetto "aperto"
NL Afsluitkraan 'open'
NO Stoppekran åpen
PL Kranik w pozycji „otwartej”
PT Torneira de passagem "aberta"
RO Robinet „deschis”
RU Запорный кран «открыт»
SL Zapiralni ventil „odprt”
SV Kran "öppen"
TR Stopkok "açık"

2.

EN Sliding Sheath Adjuster

CS Posuvné nastavovací zařízení sheathu
 DA Glidende sheath-justering
 DE Verschiebbarer Schleuseneinsteller
 EL Συρόμενος προσαρμογέας θηκαριού
 ES Ajustador deslizante de la vaina
 FI Holkin liukusäädin
 FR Ajusteur de gaine coulissant
 HR Klizni podešivač uvodnice
 IS Rennistilling fyrir slíður
 IT Regolatore scorrevole della guaina
 NL Schuifregeling van de sheath
 NO Glidende hylsejusteringsmekanisme
 PL Przesuwany regulator koszulki
 PT Regulador deslizante da bainha
 RO Ajustor glisant al tecii
 RU Передвижной механизм регулировки гильзы
 SL Drsni nastavljalnik tulca
 SV Glidande hylsjusterare
 TR Kayar Kılıf Ayarlayıcı

EN Sheath Reference Mark

CS Referenční značka sheathu
 DA Referencemærke på sheath
 DE Schleusenbezugsmarke
 EL Σήμανση αναφοράς θηκαριού
 ES Marca de referencia de la vaina
 FI Holkin viitemerkintä
 FR Repère de référence sur la gaine
 HR Referentna oznaka uvodnice
 IS Viðmiðunarkerki á slíðri
 IT Indicatore di riferimento della guaina
 NL Referentiemarkering van de sheath
 NO Referansemerke på hylsen
 PL Znacznik referencyjny wysunięcia koszulki
 PT Marca de referência da bainha
 RO Marcă de referință al tecii
 RU Установочная метка на гильзе
 SL Referenčna oznaka tulca
 SV Hylsans referensmarkering
 TR Kılıf Referans İşareti

EN Needle Handle

CS Rukojeť jehly
 DA Nålegreb
 DE Nadelgriff
 EL Λαβή βελόνας
 ES Mango de la aguja
 FI Neulan kahva
 FR Poignée de l'aiguille
 HR Drška igle
 IS Nálarhandfang
 IT Impugnatura dell'ago
 NL Naaldhandvat
 NO Nålhåndtak
 PL Uchwyt igły
 PT Punho da agulha
 RO Mâner al acului
 RU Рукоятка иглы
 SL Ročaj igle
 SV Nålhåndtag
 TR İğne Sapı


EN Safety Ring

CS Pojistný prstenec
 DA Sikkerhedsring
 DE Sicherheitsring
 EL Δακτύλιος ασφαλείας
 ES Anillo de seguridad
 FI Turvarengas
 FR Bague de sécurité
 HR Sigurnosni prsten
 IS Öryggishringur
 IT Anello di sicurezza
 NL Veiligheidsring
 NO Sikkerhetsring
 PL Pierścień bezpieczeństwa
 PT Anel de segurança
 RO Inel de siguranță
 RU Запорное кольцо
 SL Varnostni obroč
 SV Säkerhetslåsring
 TR Güvenlik Halkası

EN Stylet Hub

CS Ústí styletu
 DA Stiletmuffe
 DE Mandrinsatz
 EL Ομφαλός του στείλεου
 ES Conector del estilete
 FI Mandriinin kanta
 FR Embase du stylet
 HR Čvorište sonde
 IS Nafkjarni
 IT Pomello del mandrino
 NL Aanzetstuk van het stilet
 NO Stilettfeste
 PL Złączka mandrynu
 PT Conector do estilete
 RO Acoperitoarea stiletului
 RU Разъем стилета
 SL Pesto stileta
 SV Mandrängfattning
 TR Stile Göbeği

3.

4.

ECHOTIP® CLEARCORE™ EUS BIOPSY NEEDLE

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The EchoTip® ClearCore™ EUS Biopsy Needle is available in two needle sizes (19 Ga and 22 Ga) and contains the following components:

- Needle
- Syringe

Performance Characteristics

The function and key features of these components are described below:

EchoTip® ClearCore™ EUS Biopsy Needle:

- **Handle** – The handle contains adjustable components to allow the user to adjust the advancement or extension of the needle and the sheath. The handle safety ring slides and locks the needle at the desired needle length advancement (0–8 cm). The '0 cm' reference mark ensures complete needle retraction within the sheath. The sliding sheath adjuster allows the user to set the required sheath length extension (0–5 cm) for their ultrasound endoscope (EUS). The '0 cm' reference mark denotes the minimum sheath length available. The base of the sliding sheath adjuster has a Luer lock fitting, which facilitates the attachment of the device to the scope accessory channel.
- **Sheath** – The sheath covers the needle during introduction, advancement, and withdrawal of the device from the scope accessory channel. It protects the needle and the scope accessory channel from damage by the needle tip. ECHO-BX-22 has a flexible coiled sheath (PTFE-coated stainless steel coil spring). ECHO-BX-19 has a flexible PEEK sheath. Table 1 outlines the sheath size (Fr) for the device:

Table 1. Sheath size (Fr) for the device

Device	Sheath Size (Fr)
ECHO-BX-22	5.3
ECHO-BX-19	5.4

- **Needle** – The cobalt chromium needle has a dimpling pattern to provide visualization of the needle tip under endoscopic ultrasound. The needle cannula has a Franseen tip design to puncture and facilitate sample collection from the target site.
- **Stylet** – The stylet is inserted through the length of the needle. The stylet provides support to the needle during advancement into the target site, prevents sample collection outside of the target site, and may also be used for sample retrieval following collection, if desired. The stylet is straight (ECHO-BX-19) or coiled (ECHO-BX-22) upon removal. The stylet can be removed prior to or following puncture or when aspiration is required at the target site.

Syringe:

- The device is supplied with a syringe to utilize the standard vacuum technique for sample collection, if desired. The Luer lock fitting at the tip of the syringe can be attached to the Luer fitting on the handle. The syringe has a plunger and a stopcock to create a vacuum for sample collection. It is supplied with clips attached to the plunger to allow the user to select the preferred volume for sample aspiration.

Device Compatibility

The EchoTip® ClearCore™ EUS Biopsy Needles are compatible with:

- ultrasound endoscope with 2.8 mm minimum accessory channel
- standard Luer syringe

Patient Population

Adult patients requiring fine needle biopsy (FNB) of submucosal and extramural lesions, mediastinal masses, lymph nodes, and intraperitoneal masses within or adjacent to the gastrointestinal tract.

Intended User

The product is intended for use by physicians trained and experienced in ultrasound endoscopy techniques. Standard techniques for ultrasound endoscopy should be employed.

Contact with Body Tissue

This device is tissue contacting in line with intended use.

Operating Principles

The device is inserted into the ultrasound endoscope accessory channel and attached when the fittings meet. The sheath length can be adjusted depending on the length of the EUS scope using the sliding sheath adjuster. The needle length advancement can be adjusted as required. The device is used to puncture the target site to obtain a sample. If using a standard vacuum syringe technique, the stylet is then removed, and the syringe is attached to the Luer fitting on the needle handle. Turning the syringe stopcock to the "open" position allows negative pressure to aspirate the sample. The handle is gently moved in small increments back and forth within the target site to obtain a sample. If using the syringe, the stopcock is returned to the "closed" position. Upon completion of the procedure the needle is retracted into the sheath and the device is disconnected from the scope accessory channel.

INTENDED USE

This device is used with an ultrasound endoscope for fine needle biopsy (FNB), of submucosal and extramural lesions, mediastinal masses, lymph nodes, and intraperitoneal masses within or adjacent to the gastrointestinal tract.

INDICATIONS FOR USE

The device is indicated for ultrasound guided tissue sampling of submucosal and extramural lesions, mediastinal masses, lymph nodes, and intraperitoneal masses within or adjacent to the gastrointestinal tract, which is used for subsequent pathological examination to aid in disease diagnosis and patient management.

CLINICAL BENEFITS

This device facilitates successful sampling of submucosal and extramural lesions, mediastinal masses, lymph nodes, and intraperitoneal masses within or adjacent to the gastrointestinal tract to aid in disease diagnosis and patient management.

CONTRAINDICATIONS

- Those specific to primary endoscopic procedure to be performed to gain access to the target site.
- Coagulopathy.

WARNINGS

- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect the device with particular attention to kinks, bends, and breaks. If an abnormality is detected that would prohibit proper working condition do not use. Please notify Cook Medical for return authorization.
- Not for use in the heart or vascular system.
- The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.
- This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.
- The stylet contains nickel titanium alloy (nitinol), which may cause an allergic reaction in individuals with nickel sensitivity.

PRECAUTIONS

- When targeting multiple sites, replace device for each site.
- Needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement, or withdrawal of the device. Failure to retract needle may result in damage to the endoscope.
- Device not indicated for use in pregnant/nursing women.
- Do not use this device for any purpose other than the stated intended use.
- Use of this device is restricted to a trained healthcare professional.
- The product is intended for use by physicians trained and experienced in ultrasound endoscopy techniques.
- Standard techniques for ultrasound endoscopy should be employed.
- Store the device in a dry location.
- This symbol on the label indicates that the needle contains cobalt (Co) at a level above 0.1% w/w,

which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the needle contains cobalt as a cobalt alloy which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

POTENTIAL ADVERSE EVENTS

Those associated with gastrointestinal endoscopic ultrasound procedure: allergic reaction to medication • aspiration • cardiac arrhythmia or arrest • damage to blood vessels • fever • hemorrhage • hypotension • infection • pain/discomfort • perforation • respiratory depression or arrest.

Those associated with the device: acute pancreatitis • allergic reaction to nickel • fever • hemorrhage • infection • pain/discomfort • perforation • peritonitis • portal vein gas and thrombosis • pneumoperitoneum • tumor seeding (through the needle tract).

HOW SUPPLIED

These devices are supplied ethylene oxide (EO) sterilized in a tray with a Tyvek lid.

SYSTEM PREPARATION

1. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
2. Support the handle at both ends while removing it from the tray. Once handle is released from the tray, the sheath can be released from the secure clips.
3. Visually inspect the device components confirming they are free from damage prior to use.
4. Examine syringe. It has two plunger locks that must be depressed to advance the plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when the stopcock is in the "open" position, aligned with the syringe. (See Fig. 1.)
5. Prepare syringe as follows:
 - a. With the stopcock in the "open" position, depress the plunger locks and fully advance the plunger into the syringe.
 - b. Turn the stopcock 90° to the "closed" position.
 - c. Pull the plunger back until it is locked into place at the desired setting, creating a vacuum.
 - d. Set the prepared syringe aside until biopsy is required.

INSTRUCTIONS FOR USE

1. Identify the desired target site by endoscopic ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope. **Caution:** If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting on the accessory channel port. (See Fig. 3.)
4. Attach the device to the scope accessory channel port by rotating the device handle clockwise until the fittings are connected.

5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust the length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window. (See Fig. 2.) Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining the position of the scope, set the needle advancement to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring. (See Fig. 4.) Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the advancement of the needle in centimeters. **Caution:** During needle adjustment or advancement, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or advancement may result in damage to the scope.
7. Advance the needle into the target site by advancing the handle of the device to the pre-positioned safety ring. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
8. Standard vacuum syringe techniques may be applied for biopsy (Steps 7, 9–11) or, if desired, other techniques that may or may not incorporate the use of the stylet may be used.
9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional biopsy is required later.
10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer fitting on the needle handle.
11. Turn the stopcock to the “open” position aligned with the syringe, allowing negative pressure in the syringe to facilitate biopsy. Move the needle handle back and forth in small increments within the target site. **Note:** Do not remove the needle from the target site during biopsy.
12. Upon completion of biopsy, if using a syringe, turn stopcock to the “closed” position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
13. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
14. Unlock the thumbscrew on the safety ring and advance the needle.
15. Use a new syringe or the stylet to expel the specimen, then prepare per institutional guidelines.
16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.
17. For additional biopsy from the same target site, gently reinsert the stylet into the Luer fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.
18. Additional samples may be obtained by repeating Step 2 through 16 of the “Instructions for Use”.
19. Upon completion of procedure disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.

DEVICE DISPOSAL

This device may be contaminated with potentially infectious substances of human origin and should be coiled for disposal in accordance with institutional guidelines.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2024
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.