EN**Single Use Sterile Surgical Laser Fiber****IT****Fibre Laser Chirurgiche Sterili Monouso****FR****Fibre Laser Chirurgicale Stérile à Usage Unique****DE****Sterile Chirurgische Einweg-Laserfaser****ES****Fibra láser quirúrgica estéril desechable**

Instruction for use

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DGM001332.02

Table 1 - Fiber Specification / Specifiche della fibra / Spécification de la fibre / Faser -spezifikationen / Especificación de fibra

Fiber / Fibra / Fibre / Faser / Fibra	Ball-Tip	200 μm	272 μm	365 μm	550 μm	800 μm	1000 μm
Core diameter (μm) Diámetro del núcleo (μm) Diamètre du noyau (μm) Kerndurchmesser (μm) Diámetro núcleo (μm)	272 $\pm 2\%$	200 $\pm 2\%$	272 $\pm 2\%$	365 $\pm 2\%$	550 $\pm 2\%$	800 $\pm 2\%$	1000 $\pm 2\%$
Outer diameter (μm) Diámetro externo (μm) Diamètre externe (μm) Außendurchmesser (μm) Diámetro exterior (μm)	420 $\pm 5\%$ (jacket)*	375 $\pm 5\%$	420 $\pm 5\%$	570 $\pm 5\%$	750 $\pm 5\%$	1200 $\pm 5\%$	1450 $\pm 5\%$
Short-term Min bending radius (mm) Raggio di curvatura minimo a breve termine (mm) Rayon de courbure min. à court terme (mm) Kurzfristiger Mindestbiegeradius (mm) Radio de curvado mínimo a corto plazo (mm)	≥ 17	≥ 13	≥ 17	≥ 21	-	-	-

*The Max outer diameter for the Ball-Tip fiber is given by the distal ball diameter, which is 0.7 ± 0.1 mm

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
EN Single-use Sterile Surgical Laser Fiber

Scope of validity

These instructions detail the care and the use of the Quanta System Single-Use Sterile Surgical Laser Fibers. They are not recommendations for the medical or the surgical application of this device. Any physician using this device should be thoroughly familiar with the surgical procedures performed before the use of the device.

⚠ WARNING: This product shall be used by trained and qualified users only.

⚠ CAUTION: Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

In order to ensure safe handling of the laser probe, you must first fully read and understand this User Manual. Refer to the laser system User Manual for specific instructions concerning warnings, cautions, contraindications and clinical use of the laser.

Product description and specification

The Sterile Surgical Laser Fibers are sterile packed instruments for short-time, invasive and non-invasive use. To avoid risks due to infection and injury, the fiber should be used only if the sterile packaging has not been damaged and not expired. The Sterile Surgical Laser Fibers are sterilized with Ethylene Oxide. Refer to the Fiber Specification Table 1.

Indications for Use

Quanta System Surgical Laser Fibers are:

- intended to be used in conjunction with any cleared surgical laser manufactured by Quanta System equipped with SMA 905 or SMA 906 or compatible connector.
- indicated for use in general surgical applications for: incision, excision, vaporization, ablation, hemostasis or coagulation of soft tissue in a contact or non-contact mode (with a compatible laser marketed for use in the desired application).
- indicated for use in general surgery, urology, gastroenterology, gynecology, dermatology, vascular surgery, neurosurgery, plastic surgery, ENT/otolaryngology, endovenous occlusion of the greater saphenous vein in the patient with superficial vein reflux and laser assisted lipolysis with an approved compatible laser marketed for use in the desired application.
- also indicated for use in open or closed endoscopic applications where incision, excision, tissue dissection, excision of external tumors and lesions, complete or partial resection of internal organs, tumors or lesions, tissue vaporization, hemostasis and or coagulation may be indicated with an approved compatible laser marketed for use in the desired application.
- also intended as an aid for otologic procedures, for use in incision, excision, coagulation and vaporization of soft and fibrous tissue including osseous tissue with an approved compatible laser marketed for use in the desired application.
- also indicated for use in lithotripsy with an approved compatible laser marketed for use in the desired application.
- indicated for use with laser devices emitting radiation from 532 nm to 2100 nm, with pulsed and continuous wave (CW) emission mode.
- indicated, but not limited, for use with Diode laser, Argon, KTP/532, Ho:YAG, Nd:YAG, Tm:YAG pulsed and continuous wave CW laser devices.

Quanta System Surgical Laser Fibers may be used in surgical specialties or procedures for which compatible lasers have received regulatory clearance: for a complete information about applications, contraindications, precautions and warnings when using optical fiber it is necessary to refer to the applicable laser device User Manual.

Application and Safety Notes

Refer to the specific laser User Manual for the maximum input power values as well for warnings, cautions, contraindications and clinical use of the fiber in conjunction with the laser itself. Please do not hesitate to contact Quanta System for technical assistance. In order to ensure safe handling of the laser probe, you must first fully read and understand this User Manual. Please take care to read and understand the following Safety Notes:

- To ensure the proper and safe handling of the Sterile Surgical Laser Fiber, they should be used only by physicians familiar with the use of medical laser devices.
- Apply the general guidelines and instructions for working with laser radiation. The precautionary measures are indicated on the labels and in the User Manual of the laser device.
WARNING: All present people must wear safety goggles while the laser system is in use. The requirements for the protective equipment depend on the particular application. Refer to the User Manual of the laser device used.
- **WARNING: the Sterile Surgical Laser Fiber may be used only under surgical-room procedures.** The sterile packaging must be checked for any defects. Fibers from opened or damaged packaging are **not-sterile and must not be used.** Please check the expiration date.
- Perform a visual inspection after the fiber has been removed from the packaging. Ensure that the distal end is intact and that the SMA plug is clean. Do not use fiber with a damaged distal end or SMA plug.
- The minimum bending radius of the laser fiber must be maintained at all times throughout the entire time of use. The smallest permissible short-time bending radius corresponds to 100 times the radius of the fiber cladding. For the use of the fiber with flexible endoscopes, it is recommended considering the short-term minimum bending radius as given in Table 1.
- Before and after removing the laser fiber from its packaging, check the fiber for visible signs of damage, in particular, breakage. If any damage is found, do not use the product!
- **WARNING: Do not use a defective laser fiber,** nor should you use the laser fiber in an incorrect manner because this could cause severe eye or tissue injury, unintended exposure of patients or operating staff to laser radiation, or fire in the treatment area. Observe the detailed safety information in the User Manual of the respective laser device and the instructions for protection against laser radiation.
- Always prevent the distal end of the laser fiber from coming into contact with reflecting surfaces of other instruments or products used in order to prevent any uncontrolled leakage radiation or damage to the fiber or tissue.
- The laser may be activated only after a precise targeting and inspection of the tissue to be treated and with a correctly positioned fiber tip.
- Remnants of blood and tissues or defects can lead to overheating and changes in the fiber tip, thus directly influencing the power density of the laser radiation and the resulting therapeutic tissue effects. This can lead to a significant risk for the patient and the treating physician. If such an effect is observed the distal end must be cleaned with the laser in **Standby state.** If the full power of the laser cannot be restored, the fiber must be replaced.
- During the treatment the laser fiber can be cleaned with a sterile pad soaked with hydrogen-peroxide or sterilized water.
WARNING: during cleaning the laser must be in Standby state and the distal end of the fiber has been cooled down.
- In contact surgery treatments the distal end of the laser fiber should be moved over or placed on the tissue surface without pressure. Avoid any lateral pressure.

Operating instruction

Refer to the laser system User Manual for specific instructions concerning warnings, cautions, contraindications and clinical use of the laser.

- Check first the expiration date. The sterile packaging must be checked for any defects. Fibers from opened or damaged packaging are **not-sterile and must not be used.**

- Perform a visual inspection after the fiber has been removed from the packaging. Ensure that the distal end is intact and that the SMA plug is clean. Do not use fibers with a damaged distal end or SMA plug.
- Check the faces of the connector and of the distal fiber end for soiling or foreign matter. Damaged or soiled faces can damage or destroy the product and/or the laser system used. **NOTE: only the use of a suitable microscope objective (or magnifying glass at a minimum 40x or 20x magnifications) can assure an accurate inspection of the connector face.**
- Insert carefully the fiber connector in straight line into the fiber coupler of the laser device and tighten that by hand until you feel a resistance in doing that. Do not use tools. Please note that the fiber coupler of the laser system may be damaged, if the fiber connector is handled improperly or tightened with too much force.
- Switch on the laser device as shown in the User Manual and activate the aiming beam.
- Check the laser fiber once again for buckling, breakage and other defects. Pay special attention to the radiation leakage of the aiming beam outside the distal end surface. Direct the distal end at a non-reflecting surface: the target beam should generate a sharply defined, non-frayed, circular light spot (Fig. A2). Any irregular output profile may result in lower emission quality (Fig. A1).

⚠ WARNING: If any damage is found, do not use the laser fiber!
If used in conjunction with a surgical instrument, make sure that the fiber tip extends out of the instrument.

⚠ WARNING: Please, consider that the tip of Ball-Tip fibers is typically degraded by use in a short time during laser treatment. If necessary, the tip can be refurbished (intraoperative cleaving: see the Tip Refurbishing Section) and the fiber can be used as a standard, frontal-emission, fiber.

Tip Refurbishing

Please refer to the following instructions to refurbish the fiber tip during the treatment, when the degraded fiber tip results in a bad quality of laser emission (intraoperative cleaving)

1. Use a pair of scissors, or ceramic cutter, in order to cut off around 30 mm.
2. Use a pair of stripping pliers suitable for the fiber diameter in order to strip the fiber jacket by around 40 mm.
3. Manually remove the coating material (silicone).
4. Use a ceramic or equivalent tool in order to score a crack around 5 mm distant from the outer plastic jacket.
5. Pull the fiber tip off in the direction of the longitudinal fiber axis; do not apply a twisting or bending action (Fig. B).
6. Check the beam pattern.

Properly dispose of separated quartz tip, then inspect the cleaved tip to assess cleave quality. If a better cleave is required or the obtained one is not satisfactory, repeat the cleavage procedure from the beginning till a satisfactory cleave quality is achieved. A good cleave quality consists in homogeneous and circular profile of the aiming beam (Fig. A2), whereas a bad one is characterized by a blurred profile with irregular regions (Fig. A1).

⚠ WARNING: When performing intraoperative trimming, only sterile or sterilized devices can be used. If properly sterilized, the same device can be used intraoperatively more than once exclusively for trimming the same fiber.

⚠ WARNING: Once a Ball-Tip fiber has been refurbished, it actually behaves as a standard, frontal-emission, fiber. In this case, take care during the re-insertion of the cleaved fiber inside flexible endoscopes.

Disposal

⚠ WARNING: Do not reuse, reprocess or re-sterilize.

Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to a device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness and death of the patient.

⚠ WARNING: Dispose of the laser fiber as a hospital waste following local laws.

Possible risks and/or complication

In order to ensure safe handling of the laser probe, you must first fully read and understand this User Manual. Refer to the laser system User Manual for specific instructions concerning warnings, cautions, contraindications and clinical use of the laser. Complications that could occur during laser treatments include local and/or systemic infection, thermal changes to the surrounding structures, local hematoma, dissection and perforation, tissue adhesion, distal tip detachment, and discomfort during and/or after (laser) energy application. In the unlikely event of a detached tip, it may be visually located through an appropriate scope and removed using forceps. Irrigate the area thoroughly to remove any traces of the tip.

Reporting of serious incident

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Quanta System S.p.A.