IFU

Quanta System
LASER IN OUR DNA

EN

Reusable Sterile Surgical Laser Fiber

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
Fig. C *

Fig. D *

*Pictures and data from STERRAD® 100NX™ User's Guide Ref 99970, 08-54119-0-001, June 2007.

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Sterile Surgical Laser Fiber

Scope of validity

These instructions detail the care and the use of the Quanta System Reusable Sterile Surgical Laser Fibers. They are not recommendations for the medical or the surgical application of this device. Any physician using this device should be thoroughly familiar with the surgical procedures prior to using the device.

⚠ WARNING: This product shall be used by trained and qualified users only.

⚠ CAUTION: Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

In order to ensure safe handling of the laser probe, you must first fully read and understand this User Manual. Refer to the laser system User Manual for specific instructions concerning warnings, cautions, contraindications and clinical use of the laser.

Product description and specification

The Sterile Surgical Laser Fibers are sterile packed instruments for short-time, invasive and non-invasive use. To avoid risks due to infection and injury the fiber should be used only if the sterile packaging has not been damaged and not expired. The Sterile Surgical Laser Fibers are sterilized with Ethylene Oxide. The re-usable fibers have been designed for steam re-sterilization or hydrogen peroxide gas-plasma. Cleaning and re-sterilization is required before each further use. Re-sterilization of the product is permissible only in compliance with the reprocessing guidelines. See Table 1.

Indications for Use

Quanta System Surgical Laser Fibers are:

- intended to be used in conjunction with any cleared surgical laser manufactured by Quanta System equipped with SMA 905 or SMA 906 or compatible connector.
- indicated for use in general surgical applications for: incision, excision, vaporization, ablation, hemostasis or coagulation of soft tissue in a contact or non-contact mode (with a compatible laser marketed for use in the desired application).
- (surgical fiber optic laser delivery devices) are indicated for use in general surgery, urology, gastroenterology, gynecology, dermatology, vascular surgery, neurosurgery, plastic surgery, ENT/otolaryngology, endovenous occlusion of the greater saphenous vein in the patient with superficial vein reflux and laser assisted lipolysis with an approved compatible laser marketed for use in the desired application.
- also indicated for use in open or closed endoscopic applications where incision, excision, tissue dissection, excision of external tumors and lesions, complete or partial resection of internal organs, tumors or lesions, tissue vaporization, hemostasis and/or coagulation may be indicated with an approved compatible laser marketed for use in the desired application.
- also intended as an aid for otologic procedures, for use in incision, excision, coagulation and vaporization of soft and fibrous tissue including osseous tissue with an approved compatible laser marketed for use in the desired application.
- also indicated for use in lithotripsy with an approved compatible laser marketed for use in the desired application.
- indicated for use with laser devices emitting radiation from 532 nm to 2100 nm, with pulsed and continuous wave (CW) emission mode.
- indicated, but not limited, for use with Diode laser, Argon, KTP/532, Ho:YAG, Nd:YAG, Tm:YAG pulsed and continuous wave CW laser devices.

Quanta System Surgical Laser Fibers may be used in surgical specialties or procedures for which compatible lasers have received regulatory clearance: for a complete information about applications, contraindications, precautions and warnings when using optical fibers it is necessary to refer to the applicable laser device User Manual.

Application and Safety Notes

Refer to the specific laser User Manual for the maximum input power values as well for warnings, cautions, contraindications and clinical use of the fiber in conjunction with the laser itself. Please do not hesitate to contact Quanta System for technical assistance. In order to ensure safe handling of the laser probe, you must first read fully and understand this User Manual.

Please take care read and understand the following Safety Notes:

- To ensure the proper and safe handling of the Sterile Surgical Laser Fiber, they should be used only by physicians familiar with the use of medical laser devices.
- Apply the general guidelines and instructions for working with laser radiation. The precautionary measures are indicated on the labels and in the User manual of the laser device. **WARNING: All present people must wear safety goggles** while the laser system is in use. The requirements for the protective equipment depend on the particular application. Refer to the user manual of the laser device used.
- **WARNING: the Sterile Surgical Laser Fiber may be used only under surgical-room procedures.** The sterile packaging must be checked for any defects. Fibers from opened or damaged packaging are **not-sterile** and **must not be used**. Please check the expiration date.
- Perform a visual inspection after the fiber has been removed from the packaging. Ensure that the distal end is intact and that the SMA plug is clean. Do not use fiber with a damaged distal end or SMA plug.
- The minimum bending radius of the laser fiber must be maintained at all times throughout the entire time of use. The smallest permissible short-time bending radius corresponds to 100 times the radius of the fiber cladding. For the use of the fiber with flexible endoscopes, it is recommended considering the short-term minimum bending

radius as given in Table 1.

- Before and after removing the laser fiber from its packaging, check the fiber for visible signs of damage, in particular, breakage. If any damage is found, do not use the product!
- **WARNING: Do not use a defective laser fiber**, nor should you use a laser fiber in an incorrect manner because this could cause severe eye or tissue injury, unintended exposure of patients or operating staff to laser radiation, or fire in the treatment area. Observe the detailed safety information in the User Manual of the respective laser device and the instructions for protection against laser radiation.
- Always prevent the distal end of the laser fiber from coming into contact with reflecting surfaces of other instruments or products used in order to prevent any uncontrolled leakage radiation or damage to the fiber or tissue.
- The laser may be activated only after a precise targeting and inspection of the tissue to be treated and with a correctly positioned fiber tip.
- Remnants of blood and tissues or defects can lead to overheating and changes in the fiber tip, thus directly influencing the power density of the laser radiation and the resulting therapeutic tissue effects. This can lead to a significant risk for the patient and the treating physician. If such an effect is observed the distal end must be cleaned with the laser in the **Standby mode**. If the full power of the laser cannot be restored, the fiber must be replaced.
- During the treatment the laser fiber can be cleaned with a sterile pad soaked with hydrogen-peroxide or sterilized water. **WARNING: during cleaning the laser must be in Standby mode and the distal end of the fiber has been cooled down.**
- In contact surgery treatments the distal end of the laser fiber should be moved over or placed on the tissue surface without pressure. Avoid any lateral pressure.

Operating instruction

Refer to the laser system User Manual for specific instructions concerning warnings, cautions, contraindications and clinical use of the laser.

- For initial use check first the expiration date. The sterile packaging must be checked for any defects. Fibers from opened or damaged packaging are not-sterile and must not be used.
- Perform a visual inspection after the fiber has been removed from the packaging. Ensure that the distal end is intact and that the SMA plug is clean. Do not use fibers with a damaged distal end or SMA plug.
- Check the faces of the connector and of the distal fiber end for soiling or foreign matter. Damaged or soiled faces can damage or destroy the product and/or the laser system used. **NOTE: only the use of a suitable microscope objective (or magnifying glass at a minimum 40x or 20x magnifications) can assure an accurate inspection of the connector face.**
- Insert carefully the fiber connector in straight line into the fiber coupler of the laser device and tighten that by hand until you feel a resistance in doing that. Do not use tools. Please note that the fiber coupler of the laser system may be damaged, if the fiber connector is handled improperly or tightened with too much force.
- Switch on the laser device as shown in the User Manual and activate the aiming beam.
- Check the laser fiber once again for buckling, breakage and other defects. Pay special attention to the radiation leakage of the aiming beam outside the distal end surface. Direct the distal end at a non-reflecting surface: the target beam should generate a sharply defined, non-frayed, circular light spot (**Fig. A2**). Any irregular output profile may result in lower emission quality (**Fig. A1**).

⚠ WARNING: If any damage is found, do not use the laser fiber!

If used in conjunction with a surgical instrument, make sure that the fiber tip extends out of the instrument.

Tip Refurbishing

Please refer to the following instructions to refurbish the fiber tip at the end of the treatment, before reprocessing the fiber (*postoperative cleaving*), or during the treatment when the degraded fiber tip results in a bad quality of laser emission (*intraoperative cleaving*):

1. Use a pair of scissors, or blade, or ceramic cutter, in order to cut off around 30 mm.
2. Use a pair of stripping pliers suitable for the fiber diameter in order to strip the fiber jacket by around 40 mm.
3. Manually remove the coating material (silicone).
4. Use a ceramic or equivalent tool in order to score a crack around 5 mm distant from the outer plastic jacket.
5. Pull the fiber tip off in the direction of the longitudinal fiber axis; do not apply a twisting or bending action (**Fig. B**).
6. Check the beam pattern.

Properly dispose of separated quartz tip, then inspect the cleaved tip to assess cleave quality. If a better cleave is required or the obtained one is not satisfactory, repeat the cleavage procedure from the beginning till a satisfactory cleave quality is achieved. A good cleave quality consists in homogeneous and circular profile of the aiming beam (**Fig. A2**), whereas a bad one is characterized by a blurred profile with irregular regions (**Fig. A1**).

⚠ WARNING: When performing intraoperative trimming, only sterile or sterilized devices can be used. If properly sterilized, the same device can be used intraoperatively more than once exclusively for trimming the same fiber.

⚠ WARNING: When performing postoperative trimming, also non-sterile devices (cleaver and fiber strippers) can be used (only prior to the reprocessing and sterilization of a multiple use fiber).

Reprocessing

Reprocessing of the surgical laser fiber means the mechanical re-processing of the distal laser fiber tip as well as the cleaning and sterilization of the laser fiber.

⚠ WARNING: The fiber must be cleaned, prepared, checked, packaged and sterilized before each re-use in accordance with the instruction of this document.

⚠ WARNING: If the fiber has been used ten (10) times, a further re-sterilization is not permitted as its safety and effectiveness may not be ensured. The number of the re-sterilization cycles already carried out must be checked and recorded by the user.

⚠ WARNING: The instructions provided below have been validated by the Manufacturer as being CAPABLE of preparing the optical fibers for re-use. The customer is responsible to follow the reprocessing instructions using validated processes.

⚠ WARNING: On patients with confirmed or suspected Transmissible spongiform encephalopathies (TSEs), also known as prion disease, use only Single-use Sterile Surgical fibers.
If a Reusable Sterile Surgical fiber has been accidentally used on patients with confirmed or suspected Transmissible spongiform encephalopathies (TSEs), the fiber must be wasted according to local laws.
The reprocessing parameters provided below have not been validated by the Manufacturer as being CAPABLE of sterilizing the optical fibers in case of iatrogenic and nosocomial exposure to TSE agents.

Reprocessing – Pre treatment

In order to avoid any contamination of the connector during cleaning phases put the protection cap on the fiber connector immediately after it is removed from the laser.

Immediately after use (within 1 hour), remove coarse dirt from the instruments using fully demineralized water (thoroughly rinse the connector back in particular for at least 1 minute) and place the fibers individually into an active non-protein-fixing solution, following the concentration, temperature and contact time defined in instruction for use of the manufacturer of the solution.

Before the check and preparation for cleaning, please follow these steps:

1. Remove the instrument from the bath and rinse it thoroughly with water at least five times.
2. Dry with compressed and filtered air, in order to blow the instruments dry inside and outside.
3. Do the next step of check and preparation as soon as possible.

After wiping clean, coil up the fiber to form a ring (see below for permissible bending radii) and wrap the connector end two to three times around the coiled-up coil in order to secure the connector end in place and prevent the coil from unwinding.

Please note that this phase does not replace the steps to be carried out later.

Reprocessing – Checking and preparing the optical fiber

Use a suitable microscope or a magnifying lens in order to check the input face of the optical fiber (from the connector side). The face must be free from scratches, discoloration, deformation or foreign particles. Continued use in the presence of such defects can damage the optical fiber and/or laser device and is hence not permitted. Do not use the optical fiber any further in such a case.

Use a suitable microscope or a magnifying lens in order to check the distal fiber end. The distal fiber end must be free from scratches, discoloration, deformation, flaking or foreign particles. If necessary, use a test laser in order to control the quality of the beam pattern. The requirements defined in the Operating instruction section must be taken into consideration. If necessary, prepare the distal end of the optical fiber following the instruction of the Tip Refurbishing section.

Reprocessing - Cleaning

Cleaning can be carried out manually or by using automatic washing disinfectors. Inappropriate automated cleaning may affect the longevity or damage the fiber. When using automated cleaning follow the special handling instructions.

Before manual or automated cleaning, in order to avoid any contamination of the connector during cleaning phases, put the protection cap on the fiber connector immediately after it is removed from the laser.

Manual Cleaning

To prevent water spots, use only fully demineralized water, which is microbiologically at a minimum of drinking water quality. The relevant water treatment systems must be serviced in accordance to the manufacturer's instructions. The instruments must be dried thoroughly immediately after rinsing. Compressed air drying is preferred to other drying methods (such as drying with a cloth), because it is not only very gentle but also a highly effective technique. For manual cleaning, we recommend using soft, lint-free cloths or towels, plastic brushes or cleaning guns for cleaning; active non-protein-fixing process chemicals with or without anti-microbial effect and/or enzymes should be used.

When selecting the cleaning agents to be used, make sure that

- the agents are suitable for cleaning metal and plastic made instruments
- the agents are suitable for ultrasonic cleaning (no foaming)
- the cleaning agent has a tested efficiency and it is compatible with the cleaning agent used before
- the agents are compatible with the instruments (refer to the "Material strength" section).

The cleaning solution should be aldehyde-free (in order to prevent adhesion of blood contamination), of proven efficiency (such as VAH/DGFM or FDA approval and/or CE marking), suitable for cleaning the instruments, and compatible with the instruments (refer to the "Material strength" section). Always use a dedicated, clean, soft piece of cloth, never metal brushes or steel wool, in order to remove dirt manually.

The main reasons for mechanical damage during manual reprocessing include:

- Use of metal brushes
- Use of coarse scouring agents
- Use of too much force.

The cleaning solution Neodisher MediClean has been validated by the Manufacturer, according to the manufacturer's instructions concerning concentration, temperature and contact time. This solution should always be used strictly following the manufacturer's instructions for use, under all conditions. The cleaning solution should be freshly prepared on a daily basis; where contaminations levels are high, it is advisable to prepare fresh solutions at even

shorter intervals. The cleaning solution should be sterile or almost sterile (maximum of 10 germs per ml) as well as low-endotoxin water (maximum of 0.25 endotoxin units per ml - purified water / highly purified water) and/or filtered air for drying.

In order to clean the optical fibers, these steps should be followed:

1. Place the wound-up fiber individually into a sufficiently large wire basket. Make sure that the fiber end is not damaged by the wire.
2. Place the wire basket with the instrument into the cleaning bath (with ultrasonic support, if necessary) and ensure the fibers are immersed for 5 minutes in 0,5 % Neodisher MediClean (Dr. Weigert, Hamburg) cleaning solution as validated solution by the manufacturer.
3. Wipe gently the surface of the fibers with a soft cloth soaked in the 0,5% cleaning solution until no visible soil remains.
4. Brush gently in cold fully demineralized water with a soft nylon bristle brush until no visible soil remains.
5. Subsequently remove the wire basket containing the instrument from the cleaning bath and thoroughly rinse the fiber at least three times with fully demineralized water.
6. Dry with compressed and filtered air, in order to blow the instruments dry inside and outside.
7. Check the instruments (see Section Reprocessing – Checking and preparing optical fiber).

Reprocessing – Automated cleaning

Cleaning can be carried out by using automatic washer-disinfectors, meeting the requirements of the ISO 15883 series. Inappropriate automated cleaning may affect the longevity or damage the fiber.

When using automatic washer-disinfectors, use pH mild alkaline Neodisher Mediclean (Dr. Weigert, Hamburg) detergents. The concentration and exposure times specified by the manufacturers of the cleaning agents must be observed under all conditions.

In order to clean the optical fiber, the following steps shall be applied:

1. Place the coiled-up instrument individually into a sufficiently large cleaning basket observing the minimum bending radius. Make sure that the fiber tip does not become damaged by touching the cleaning basket.
2. Do no overlap fibers in the same cleaning basket.
3. Place the cleaning basket with the fiber on the washer disinfecter trays.
4. The load must be homogeneous by type of material.
5. Make sure that the fibers are loaded according to the homogeneous load configuration, defined in your validated cleaning process; between a tray and the other must remain the necessary space to allow the passage of an hand.
6. The load must be distributed so as to allow the sterilizing agent to freely circulate and penetrate trays placed inside the washer disinfecter.
7. Make sure that the fiber could not become damaged by other instruments or medical devices and the fiber must not touch the walls or the door of the washer disinfecter.
8. Set an automated washing cycle, according to the following washing parameters, validated by the Manufacturer*:
 - a. 2 minutes pre-cleaning with cold tap water
 - b. Draining
 - c. 5 minutes cleaning with 55 °C fully demineralized water and 0.5% Neodisher Mediclean (Dr. Weigert, Hamburg) cleaning solution
 - d. Draining
 - e. 3 minutes neutralization with cold fully demineralized water
 - f. Draining
 - g. 2 minutes rinsing with fully demineralized water
 - h. Draining
 - i. Drying for 10 minutes
9. Do the next step as soon as possible (see Section Reprocessing – Checking and preparing optical fiber).

*The validation testing of automated cleaning process has been conducted with a Miele G7836 CD washer, with artificially contaminated fibers and left to dry at ambient conditions for 1 hour. Hospital loads should not exceed the maximum drying time of contaminated fibers, validated by this testing.

Reprocessing – Final check and packaging

Before packaging, check the optical fiber for corrosion of the connector, surface damage, flaking and soiling. Do not use further damaged optical fiber. Soiled ones must be cleaned one more time.

Final check must verify that the distal fiber end has not been damaged during cleaning. Whenever the distal end is damaged, it is necessary to repeat all steps from Reprocessing – Checking and preparing the optical fiber. The cleaning processes must be also repeated.

Lubrication is not required – Do not use any oil.

As soon as the final check has been passed, pack the optical fiber (still in a coiled-up condition) in sufficiently large disposable sterilization packaging unit (single or double packaging). The packaging of the medical device shall be compatible with steam and/or gas-plasma sterilization and with the optical fiber. The packaging materials and process shall be validated, according to EN ISO 11607 series.

Reprocessing - Sterilization

The reusable optical fiber has been tested for the steam sterilization process or hydrogen peroxide gas-plasma. Any other sterilization process is not allowed.

⚠ WARNING: The instantaneous sterilization process is NOT ALLOWED.

⚠️ WARNING: Hot-air sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization processes are NOT ALLOWED.

⚠️ WARNING: The parameters provided below have been validated by the Manufacturer as being CAPABLE of sterilizing the optical fibers for re-use. The customer is responsible to follow the reprocessing instructions using validated processes.

⚠️ WARNING: This manual contains a number of references to accessories and solutions for the Surgical Laser Fibers use. The Surgical Laser Fibers have been tested and validated for use with the accessories and disposables listed in this manual. The Manufacturer has not validated the use of accessories or disposables other than those specified in this manual. The Manufacturer does not assume responsibility or liability for use of accessories or disposable other than those specified in this manual. It is the responsibility of the user to validate that other accessories or disposables provide safe and effective performance.

Recommended parameters for Steam sterilization

Steam Sterilization Load:

- Use a steam sterilizer cleared by the US FDA or according to EN 13060 and/or EN 285, with validated sterilization process, according to ANSI/AAMI/ISO 17665.
- The fibers must be cleaned, perfectly dried and packaged before being placed in the sterilizer.
- The load must be distributed so as to allow the sterilizing agent to freely circulate and penetrate each package placed inside the sterilization chamber.
- The material must not touch the walls or the door of the sterilization chamber.
- Packs should be placed in non-compressed grilled baskets and placed vertically, to facilitate the removal of air and drying.
- A necessary space to allow the passage of a hand must remain between a pouch and the other.
- In the case of paper and polypropylene packs, arrange the packages by matching the polypropylene part with polypropylene of the other package in such a way that the contents do not overlap.
- The load must be homogeneous by type of material, because the different type of material can produce room heating and cooling at different times.
- Sterilization load cannot exceed the validated quantities per tray, as defined in your steam sterilization validation data.
- Sterilize without any additional load.
- Maximum exposure temperature of 138 °C (280 °F) according to ANSI/AAMI/ISO 17665 or EN ISO 17665.
- Use steam sterilization process parameters according to your validated methods, according to EN ISO 17665 requirements.

The following steam sterilization parameters have been validated by the manufacturer.

Pre-vacuum mode:

- 3 pre-vacuum phases
- 134 °C, holding time 3 minutes (commonly used in Europe) or 135 (275 °F) holding time 3 minutes and 132 °C (270 °F) holding time 4 minutes (commonly used in US)
- 20 minutes drying time

Or

Gravity mode:

- 15 minutes minimum exposure time for gravity-displacement steam sterilization at 134 °C (commonly used in Europe), 135 °C (275 °F) / 132 °C (270 °F) (commonly used in US)
- 20 minutes drying time

Recommended parameters for hydrogen peroxide sterilization (Gas-plasma sterilization)

Low temperature hydrogen peroxide, using the gas-plasma sterilizer Sterrad® 100NX™, Sterrad® Flex cycle, with set of validated parameters as defined in your validated gas-plasma sterilization process, according to EN ISO 14937 and the sterilizer and Hydrogen peroxide solution manufacturer process specifications.

Gas-plasma sterilizing solution Hydrogen peroxide shall be stored and used at the temperature defined by the manufacturer's instructions concerning concentration, temperature and exposure time. This solution should always be used strictly following the manufacturer's instructions for use, under all conditions.

The following gas-plasma sterilization parameters have been validated for the optical fibers by Quanta System, using the gas-plasma sterilizer Sterrad® 100NX™, manufacturer Advanced Sterilization Products, a Johnson&Johnson company, Division of Eticon, Inc. Irvine, CA 92618-9824, Sterrad® Flex cycle:

- Cycle time 42 minutes exposure – Flex Cycle
- Hydrogen peroxide (H2O2) concentration: 58 % - 59,5 %, at temperature and the relative humidity required for the Sterrad® Flex cycle
- Cycle temperature: 47 °C - 56 °C (116.6 °F – 132.8 °F)
- Other process parameters as required for the Sterrad® Flex cycle.

Gas-plasma Sterilization Load:

- Use a steam sterilizer cleared by the US FDA or European Notified Body, with validated sterilization process, according to ANSI/AAMI/ 14937 or EN ISO 14937.
- The fibers must be clean, perfectly dry and packaged before being placed in the sterilizer.
- The load must be distributed so as to allow the sterilizing agent to freely circulate and penetrate each package placed inside the sterilization chamber.
- When placing the load on the shelves, make certain that you do not block the ultraviolet lamp beam in the front right (your right) side of the chamber. Make sure the load is centered on the shelves and that the shelves are centered in the chamber (see Fig. C - Do not block the UV lamp).

- Do not allow any part of the load to touch the electrode, the back wall of the chamber, or the inside of the door.
- Leave at least 1 inch (25 mm) of free space between the load and the electrode to allow hydrogen peroxide to diffuse around the load (see Fig. D - Do not stack trays).
- The material must not touch the walls or the door of the sterilization chamber.
- Packs should be placed in non-compressed grilled baskets and placed vertically, to facilitate the removal of air and drying.
- A necessary space to allow the passage of a hand must remain between a pouch and the other.
- In the case of paper packs and polypropylene arrange the packages by matching the polypropylene part with polypropylene of the other package in such a way that the contents do not overlap.
- The load must be homogeneous by type of material, because the different type of material can produce room heating and cooling at different times.
- Sterilization load cannot exceed the validated quantities per tray, as defined in your Gas-plasma sterilization validation data.
- Sterilize without any additional load.
 - Maximum sterilization temperature of 47 °C - 56 °C.

Use gas-plasma sterilization process parameters according to your validated methods, according to ANSI/AAMI/ 14937 or EN ISO 14937 requirements.

Post-sterilization techniques/activities

After sterilization, store the fiber in the sterilization packaging in a dry and dust-free condition.

Reprocessing - Material strength

When selecting cleaning agents, please ensure that the following ingredients are not present:

- organic, mineral and oxidizing acids
- lye (neutral/enzymatic cleaning agent recommended)
- organic solvents (such as alcohols, ethers, ketones, petroleum ethers)
- oxidants (such as hydrogen peroxides)
- halogens (chlorine, iodine, bromine)
- aromatic/halogenated hydrocarbons
- oils.

Never use metal brushes or steel wool to clean any of the instruments.

⚠️ WARNING: Do not expose the optical fiber to temperatures above 138 °C (280 °F).

⚠️ WARNING: Do not reprocess the optical fiber more than 10 times.

Reuse and disposal

If a suitable care and an undamaged and clean condition are observed, the laser fiber can be reused up to the maximum value stated on the label (10x or 10 times). No liability whatsoever is assumed in the case of non-compliance.

⚠️ WARNING: Any further reuse and/or the use of damaged and/or soiled instruments is forbidden.

⚠️ WARNING: Dispose of the laser fiber as a hospital waste following local laws.

⚠️ WARNING: The instructions provided above have been validated by the Manufacturer as being CAPABLE of preparing a medical device for re-use. It remains the responsibility of the processor to follow our instructions through its validated processes. This requires validation and routine monitoring of the process. Likewise any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

Possible risks and/or complication

In order to ensure safe handling of the laser probe, you must first fully read and understand this User Manual. Refer to the laser system User Manual for specific instructions concerning warnings, cautions, contraindications and clinical use of the laser. Complications that could occur during laser treatments include local and/or systemic infection, thermal changes to the surrounding structures, local hematoma, dissection and perforation, tissue adhesion, distal tip detachment, and discomfort during and/or after (laser) energy application. In the unlikely event of a detached tip, it may be visually located through an appropriate scope and removed using forceps. Irrigate the area thoroughly to remove any traces of the tip.

Reporting of serious incident

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Quanta System S.p.A.