EN

**Single Use Sterile Surgical Laser
Lateral Fiber**

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EN

Single-use Sterile Surgical Laser Lateral Fiber

In order to ensure the safe use of the lateral fiber, you must first fully read and understand this manual. Refer to the laser system User Manual for specific instructions concerning warnings, cautions, contraindications and clinical use of the laser. The use of this device is restricted to medical personnel already qualified and skilled with the use of laser devices.

Indications for Use

Surgical Laser fibers are intended to be used to deliver the laser radiation to the target tissue when used with any cleared/certified surgical laser with operational wavelengths between 532 nm – 2200 nm equipped with SMA 905 or SMA 906 or compatible connector, as per the indications of the laser device used with.

Restrictions for use

Their use is restricted to lasers:

- capable of generating a beam in the fiber input face that is smaller than the input face and is well aligned to the center;
- with peak power up to 15KW and peak power density up to 50MW/cm²;
- with numeric aperture not higher than the specific fiber numeric aperture.

Intended users and Intended patient population are the ones of the laser device intended to be used with the fibers.

Description and Technical data

The lateral laser fiber makes the positioning towards the tissue to be treated more precise and comfortable through the lateral emission characteristics of the laser energy.

Single use / reusable	Single use
Sterility	Fibers are supplied sterile (Ethylene Oxide)
Shelf life	5 years
Recognition connector	Code plug or RFID (if present)
Connector	SMA 905 or SMA 906 or compatible connector
Angle of laser emission	80° (relative to the axis of the fiber)
Numerical Aperture	0.22 ± 0.02
Length	3 m ± 0,1 m
Internal Diameter	600 µm ± 10 µm
Outer diameter	2000 µm ± 60 µm
Short-term Min bending radius	≥ 38.3 mm
Transport and storage temperature range	-29° C to 60° C

Warnings

- The use with a laser with characteristics not compliant with the intended use is **strictly forbidden** as it may result in risks for the patient and for bystanders.
- The fiber is single use only and it can only be used on a single patient during one treatment. The reuse is **strictly forbidden**.
- Apply the general guidelines and instructions for working with laser radiation. The precautionary measures are indicated on the labels and in the User manual of the laser device. Refer to the user manual of the laser device used.
- The incorrect use of the fiber may result in severe eye or tissue injury, unintended exposure of patients or operating staff to laser radiation, or fire in the treatment area. Observe the detailed safety information in the user manual of the respective laser device and the instructions for protection against laser radiation.
- Always prevent the distal end of the laser fiber from coming into contact with reflecting surfaces of other instruments or products used in order to prevent any uncontrolled leakage radiation or damage to the fiber or tissue.
- In contact surgery treatments the distal end of the optical fiber should be moved over or placed on the tissue surface without pressure. **Avoid any lateral pressure.**
Note: in contact mode, the fiber could overheat due to the burn-in tissue on the distal end. Pay attention to ensure a sufficiently high flow rate of the rinsing to cool the tissue and the fiber.
- If the fiber is used in conjunction with a surgical instrument, make sure that the fiber tip extends out of the instrument.
- In the unlikely event of a detached tip, it may be visually located through an appropriate scope and removed using forceps. Irrigate the area thoroughly to remove any traces of the tip.
- in case of intraoperative refurbishing of the distal tip, the laser emission will not be lateral anymore, but would become unpredictable. The use of the optical fiber in this condition could cause serious injuries and it is not allowed. Please change the fiber with a new one in case of damage of the distal tip.

Pre-use Inspection

- To avoid the risk of infection the fiber should be used only if, before using, the user makes sure that the sterile packaging has not been damaged and that the fiber is not expired. Fibers from opened or damaged packaging are **not-sterile and must not be used**, as well as expired fibers.
- Before and after removing the optical fiber from its packaging, check the fiber for visible signs of damage, in particular, breakage. If any damage is found, do not use the product.
- Perform a visual inspection after the fiber has been removed from the packaging. Make sure that the distal end is intact and that the SMA input face is clean. **Do not use fiber** with a damaged distal end or SMA input face, as damaged or soiled faces may seriously damage the product and/or the laser device.

Operating instruction

- For the fibers with Code Plug, connect the recognizing connector to the laser device in order to recognize the fiber (Fig. A).
- For the fibers with RFID, the laser automatically recognizes the fiber (Fig. B).
- Carefully insert the fiber connector in straight line into the fiber coupler of the laser device and tighten that by hand until you feel a resistance in doing that. Do not use tools. The fiber coupler of the laser system may get damaged, if the fiber connector is handled improperly or tightened with too much force.
- Check the laser fiber once again for bucking, breakage and other defects. Pay special attention to the radiation leakage of the aiming beam outside the distal end surface. Direct the distal end at a non-reflecting surface: the target beam should generate a sharply defined, non-frayed, circular light spot. Any irregular output profile may result in lower emission quality.
- During the treatment the optical fiber can be cleaned with a sterile pad soaked with hydrogen-peroxide or sterilized water; during cleaning the laser device must be placed in Standby mode and the distal end of the fiber must have been cooled down.
- If during a procedure tissue debris stick on optical fiber tip, place the laser in Standby mode, remove the optical fiber from the endoscopic instrument and carefully clean the tip with a moist sterile gauze.
The correct method to remove debris from the optical fiber is to gently clean it from the tip towards the handpiece.
Remnants of blood and tissues or defects can lead to overheating and changes in the fiber tip, thus directly influencing the power density of the laser radiation and the resulting therapeutic tissue effects. If the full power of the laser cannot be restored after the cleaning, the fiber must be replaced.
- In case the fiber tip gets damaged during the treatment, or a bad quality of laser emission is observed change the fiber with a new one.

Disposal

- Do not reuse, reprocess or re-sterilize.

Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to a device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness and death of the patient.

- Dispose the Optical Fiber as a hospital waste following local laws.

Reporting of serious incidents

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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