Document No.	05SDU00825	Secret Level	<input type="checkbox"/> General Lv1 <input checked="" type="checkbox"/> Secrecy Lv2 <input type="checkbox"/> Confidential Lv3
Version	A-0	Doc. Type	DHF
Scope	Ascend™ Single-Use Flexible Ureteroscope (FDA)	Code Model	RP002E00-OD1/RP002P00-OD1 ASUS-86-900R-DC/ASUS-86-300R-DC/ ASUS-86-300S-DC/ASUS-86-900S-DC/ ASUS-86-900R/ASUS-86-300R/ ASUS-86-300S/ASUS-86-900S

Instruction for Use

Ascend™ Single-Use Flexible Ureteroscope

Drafted by		Date	
Reviewed by		Date	
Reviewed by		Date	
Approved by		Date	
Date of Publication		Seal	

--	--	--

D i s t r i b u t i o n	<input type="checkbox"/> Quality Assurance Department	<input type="checkbox"/> Quality Inspection Department
	<input type="checkbox"/> Laboratory	<input type="checkbox"/> Production Department
	<input type="checkbox"/> R&D Department	<input type="checkbox"/> NPI Department
	<input type="checkbox"/> Engineering Department	<input type="checkbox"/> PE Technical Department
	<input type="checkbox"/> Planning Department	<input type="checkbox"/> Procurement Department
	<input type="checkbox"/> Regulatory Affairs Department	<input type="checkbox"/> After-Sales Department
	<input type="checkbox"/> Domestic Sales Department	<input type="checkbox"/> Domestic Marketing Department
	<input type="checkbox"/> Overseas Sales Department	<input type="checkbox"/> Overseas Marketing Department
	<input type="checkbox"/> Clinical Medical Product Department	<input type="checkbox"/> Administration & HR Department

Change History

Rev.	Date	Revised contents	Prepared by	Change No.
A-0	2023.12.07	Initial release.	Xueping Zeng	20231207001

Instructions for Use (IFU)

Ascend™ Single-Use Flexible Ureteroscope

Model:

ASUS-86-900R-DC, ASUS-86-300R-DC, ASUS-86-300S-DC, ASUS-86-900S-DC,
ASUS-86-900R, ASUS-86-300R, ASUS-86-300S, ASUS-86-900S

Contents

Page

English (Instructions for Use).....2
Français (Manuel d'utilisation).....21

Table of Contents

1. Precautions, Warnings and Tips.....	3
2. Product Overview.....	5
3. Other products which are intended to be used in combination with the device.....	9
4. Product Performance	9
5. Installation and Operation	10
6. Troubleshooting.....	14
7. Other Items.....	16
8. Electromagnetic Compatibility Information.....	17
9. Disposal.....	20

1. Precautions, Warnings and Tips

1.1 Warnings

1	Prior to the use of this product, please read this IFU carefully, including all precautions, warnings and tips. Failure to do so may result in injury to the patient.
2	No modification of this product is allowed without the permission of the manufacturer.
3	Damage to the device due to improper use and incorrect maintenance by the user, as well as their consequences and responsibilities, shall be taken by the user.
4	The Ascend™ Single-Use Flexible Ureteroscope (hereinafter referred to as "endoscope") shall be disposed of in accordance with the hospital waste disposal procedures after using, and shall not be reused. Improper handling or reuse may result in cross-contamination or infection of patients.
5	Personal protective equipment shall be worn during operation, so as to prevent the operators from exposure to potentially infectious substances. Personal protective equipment should prevent from skin exposure with appropriate size and enough length.
6	Do not hit, collide or drop the endoscope, and do not bend, pull or twist the endoscope, as this may cause damage to the endoscope or dislodge parts into the patient, resulting in injury, bleeding or perforation.
7	Do not pull, twist or rotate the bending section which is already bent at an angle, as this would cause injury, bleeding or perforation of the patient and would prevent the bending section from being straightened during the examination.
8	Do not insert/withdraw the insertion portion suddenly or with excessive force as this would result in injury, bleeding or perforation to patient.
9	Do not press the bending section hard, as the outer skin of the bending section would be stretched or cracked, then may result in water leakage.
10	Do not insert the video cable connector if it becomes damp or contaminated, as this would result in electric shock and serious damage to the endoscope, and also jeopardize the safety of the patients or operators.
11	Do not pull the video cable or knock the connector during examination, as it would result in abnormal images.
12	Do not excessively coil the video cable, as it would result in abnormal images.
13	Please plug/unplug and connect the video cable with the endoscopic video image processor off, otherwise it would result in damage to the device.
14	If there is no live image shown, the camera module may have been damaged, please stop using the device immediately. Do not use the endoscopic diagnosis and treatment accessories, inserting, withdrawing the insertion portion of the endoscope, or performing rotational movement, as this would result in injury, bleeding or perforation to patient.
15	When the endoscopic live image is frozen, do not use the endoscopic accessories for diagnosis and treatment, inserting, withdrawing the insertion portion of the endoscope, or performing rotational movement, as this would result in injury, bleeding or perforation to patient.
16	If the live image disappears unexpectedly during the examination, or if the frozen images cannot be recovered, do not use the endoscopic accessories for diagnosis and treatment, insert or withdraw the insertion portion of the endoscope, or perform angle control, and perform troubleshooting as described in section "6.Troubleshooting" of this IFU, as it would result in patient injury, bleeding or perforation.
17	LED illumination is required for endoscopic observation, and improper use would result in viable tissues alterations such as protein denaturation and tissue perforation of viable tissues, please follow the precautions as follows: (1) The maximum temperature limit of the tip is 43°C, if the surface temperature exceeds 41 °C, as it would result in mucosal burns. Do maintain the proper distance required for adequate observation and do not constantly observe tissue at close range or allow the distal tip of the endoscope to remain in

	unnecessary contact with the mucosa for prolonged durations, as this would result in an elevated temperature at the tip (possibly exceeding 41°C). Always power off the endoscope when it is not used for a long time to avoid unnecessary illumination of the endoscope; (2) Do not expose the light directly to the eyes to avoid damage to the eyes.
18	Electromagnetic interference may be generated when using this product near equipment marked with the following symbols or other portable mobile RF (radio frequency) communication devices (such as cellular phones). In case of electromagnetic interference, take appropriate measures, such as adjusting the direction, location or shielding the using site of this device. 
19	To eliminate electromagnetic interference from other devices (any device other than this product or the components constituting this system), observe whether this system is operating properly in the environment in which it is used.
20	The influence from strong static electricity or X-ray would result in abnormal endoscopic image display.
21	The camera module and LED at the distal tip of this device are low-voltage sensitive electronic devices. Therefore, non-conventional contact shall be avoided. Prior to the use of this endoscope, it is recommended to touch the ground with the hands to ensure that the static electricity of the touching body is released, and to wear anti-static bracelets or other anti-static devices to avoid electrostatic breakdown of sensitive electronic devices.
22	This product is not an MR safe product and cannot be used safely in an MRI environment. Simultaneous use of nuclear magnetic resonance imaging (NMRI), endoscopes and accessories would result in hazards and generate artificial artifacts; please follow the instructions and safety precautions specified by the manufacturer.
23	If the product is used under extremely poor conditions, LED may be damaged.
24	Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.
25	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
26	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
27	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
28	Please contact with manufacturer and local competent authority if have any serious incident related with the device.
29	Electromagnetic interference* (EMI) may occur on this instrument near equipment marked with the following symbol  . Electromagnetic interference may occur on this instrument near other portable and mobile radio frequency (RF) communications equipment or specific RF emitters for example, cellular phones, diathermy, electrocautery, and RFID or other specific potential electromagnetic disturbance sources. To check for EMI, verify the system's operation in which it will be used. Should EMI occur, employ mitigation measures like reorienting or repositioning the instrument, or shielding its location. Placing this instrument near other medical electrical equipment or mobile RF communications equipment may result in EMI, which may degrade the video image.

1.2  Cautions and Tips

1	In many therapeutic areas, this product may be used in conjunction with laser surgery, pneumatic or electro-hydraulic endoscopic accessories. Please follow the instructions and safety precautions for the products and accessories used.
2	Please perform operations by following the precautions below, as the endoscopic image may disappear unexpectedly during the examination or the frozen image cannot be recovered. (1) Insert the video cable into the endoscopic video image processor and pull it gently to verify that the connection is firm; (2) Do not bend, tap, pull or twist the endoscope, as this would damage the endoscope and cause water leakage or damage to other parts; (3) Prior to connecting the video cable to the endoscopic video image processor, verify that the connector is completely dry, otherwise the device may malfunction.
3	Darkening of the endoscope images during use may indicate that blood or mucous membranes are adhering to the camera module at the distal tip of the endoscope. Carefully withdraw the endoscope from the patient, wipe away the blood or mucous membranes, and restore optimal illumination to ensure a safe examination.
4	Do not exclusively rely on endoscopic observation as the initial diagnostic basis for a lesion to make any decisions about potential diagnostic or device interventions.
5	This endoscope contains no memory chip, and cannot store endoscopic information. If you need to access the information photographed and recorded by the endoscope, please access it via the endoscopic video image processor attached.
6	Spare device: Be sure to prepare a spare endoscope to avoid interrupting the inspection in case of device failure or abnormal function.
7	The operation must be carried out in a visual state and placed in a place that meets the requirements.
8	The operator must be physicians or healthcare providers under supervision of physicians and they must have received sufficient training on clinical endoscope technology.
9	The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
10	The input voltage of the endoscope is 3.3VDC, and the input current is 150mA.

2. Product Overview

2.1 Product name

Ascend™ Single-Use Flexible Ureteroscope

2.2 Model

ASUS-86-900R-DC, ASUS-86-300R-DC, ASUS-86-300S-DC, ASUS-86-900S-DC,

ASUS-86-900R, ASUS-86-300R, ASUS-86-300S, ASUS-86-900S

2.3 Scope of application

The device has been developed to be used with endoscopic accessories such as a biopsy forceps and Endoscopic Video Image Processor for endoscopy and endoscopic surgery within urinary tract and interior of the ureter and kidney.

2.4 Intended patient population

Adult


2.5 Requirement of operator

- a) Education:
 - A physician or a medical staff under the supervision of a physician.
- b) Knowledge:
 - Must have had enough clinical endoscopic technical training.
- c) Language understanding:
 - Can understand English or local language in the User Manual and displayed on touchscreen console.
- d) Experience:
 - Trained doctors
 - Professional doctors or medical staff who have studied and fully understood the IFU.
- e) Permissible impairments:
 - Mild reading vision impairment or vision corrected to log MAR 0, 2 (6/10 or 20/32)






2.6 Contraindications

- (1) Uncontrolled systemic hemorrhagic disease.
- (2) Weak constitution, cannot tolerate anesthesia.
- (3) Benign prostatic hyperplasia affects ureteral ureteroscopy insertion.
- (4) Lower urinary tract tumor.
- (5) Bladder contracture.
- (6) Uncontrolled hemorrhagic disease of the whole body.
- (7) Uncontrolled urinary tract infection.
- (8) Severe cardiopulmonary insufficiency, unable to tolerate surgery.
- (9) Severe urethral stricture and ureteral stricture.
- (10) Severe hip deformity with difficult lithotomy position.


2.7 Structure composition





Structural differences

Model	Decorative cover	Distal tip part
ASUS-86-900R-DC		
ASUS-86-300R-DC		
ASUS-86-900S-DC		
ASUS-86-300S-DC		

ASUS-86-900R

Structural differences

Model	Decorative cover	Distal tip part
ASUS-86-900R		
ASUS-86-300R		
ASUS-86-900S		
ASUS-86-300S		

No.	Name	Functions
1	Rocker	When the rocker is rotated in the “U” direction, the bending section is rotated upwards; When the rocker is rotated in the “D” direction, the bending section is rotated downwards.
2	Snap lock	Rotate towards the direction shown by the “LOCK→” arrow to lock the bending section, which locks the desired position, and rotate in the opposite direction to unlock it.
3	Handle	Used to operate the endoscope, e.g. control the angle.
4	Instrument channel port	Waterproof cap can be installed, endoscopic diagnosis and treatment accessories can be inserted through this channel.
5	Instrument channel	The endoscopic accessories for diagnosis and treatment and perfusate enter the patient's body cavity through the instrument channel.
6	Insertion portion	Including insertion tube and distal tip, for insertion into the body cavity of the patient.

7	Bending section	For adjusting and moving the distal tip of the endoscope.
8	Distal tip	Export of camera modules, LED carriers and accessories for diagnosis and treatment.
9	Video cable interface	Insert video cable, connect to medical endoscopic video image processor to display endoscope image.

2.8 Sterile level

The product is sterilized with ethylene oxide, it is a single use sterile product, and SAL is 10^{-6} .

3. Other products which are intended to be used in combination with the device

Other Product	Product Name
Image processor in combination with the device (ASUS-IPD-V100)	Endoscopic Video Image Processor
The recommended endoscopic accessories (Endoscopic accessories cleared or approved by FDA and compatible with the endoscope are required)	Nitinol Tipless Stone Extractor
	Medical laser fiber

4. Product Performance

Item		Parameter	
Model		ASUS-86-900R-DC, ASUS-86-300R-DC, ASUS-86-300S-DC, ASUS-86-900S-DC	ASUS-86-900R, ASUS-86-300R, ASUS-86-300S, ASUS-86-900S
Basic Size	Working length	670mm	670mm
	Maximum width of insertion section	3.15mm (9.5Fr)	3.18mm (9.5Fr)
	Maximum OD of flexible insertion portion	Φ2.85mm(8.6Fr)	Φ2.85mm(8.6Fr)
	Minimum width of instrument channel	1.2mm(3.6Fr)	1.2mm(3.6Fr)
Optical Properties	Field of view	120°	120°
	Viewing angle	0°	0°
	Depth of field	5mm~50mm	5mm~50mm
Angle range of bending section		275° up - 275° down	275° up - 275° down

Operating Environment	Ambient temperature	10°C ~ 40°C	10°C ~ 40°C
	Relative humidity	30% ~ 85%	30% ~ 85%
	Air pressure	70kPa ~ 106kPa	70kPa ~ 106kPa
		Prohibited to use in MRI environment	Prohibited to use in MRI environment
Transportation Environment	Ambient temperature	-29°C ~ 55°C	-29°C ~ 55°C
	Relative humidity	10% ~ 95%	10% ~ 95%
	Air pressure	70kPa ~ 106kPa	70kPa ~ 106kPa
Storage Environment	Ambient temperature	-29°C ~ 55°C	-29°C ~ 55°C
	Relative humidity	10% ~ 95%	10% ~ 95%
	Air pressure	70kPa ~ 106kPa	70kPa ~ 106kPa

5. Installation and Operation

5.1 Preparation and check

Please prepare and check the product in accordance with the following instructions before using. In case of any abnormality, please refer to “6.Troubleshooting” section of this IFU for troubleshooting. If there is still a problem with the device, please do not use this product and contact Cook Medical.

	The product is sterilized with ethylene oxide and it is a single use sterile product. Ensure that the packaging is intact, and do not use in case of any package damages before using.
	Do not use it if the product is beyond its shelf life of two years.

5.1.1 Endoscopy check

- (1) Check the endoscope handle section for any deformation, loose parts or other abnormalities;
- (2) Check the camera module at the distal tip of the endoscope for any scratches, cracks, stains, gaps or other abnormalities;
- (3) Check the entire insertion portion of the endoscope for any abnormalities such as damages, stains and cracks.

5.1.2 Bending function check

- (1) Straighten the bending section;

- (2) Verify that the snap lock is in the unlocked position as shown in “LOCK→”;
- (3) Slowly adjust the rocker in both directions to the maximum angle to verify that the bending section bends smoothly and that the maximum angle can be reached;
- (4) Slowly adjust the rocker to the middle position to verify that the bending section can return to the straightened state smoothly.

	In case of abnormal bending functions, do not use the endoscope, as the bending section may not be straightened, resulting in patient injury, bleeding or perforation.
-----------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------

5.1.3 Connection of the power supply with Endoscopic Video Image Processor

- (1) Connect the endoscopic video image processor to a DC power supply;
- (2) Verify that the video cable connectors are connected the corresponding interface according to the pair markings;
- (3) Connect one end connector of the video cable to the endoscope and the other to the endoscopic video image processor;
- (4) Turn on the endoscopic video image processor and verify that the endoscope image is displayed properly.

	Don't move the endoscopic video image processor , video cable, adapter, power line, etc, as it would result in poor connection, which affects the signal interruption.
------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------

5.1.4 Endoscopic live image check

- (1) Verify that the endoscope LED illumination is functioning properly;
- (2) Properly adjust the screen brightness;
- (3) Verify that the images are free from noisy point, blurring, fogging or other abnormalities;
- (4) Slowly adjust the maximum angle of rocker up and down respectively, and observe whether the image is abnormal;
- (5) Verify that the image shows no transient disappearance or other abnormalities.

	Do not use the endoscopic video image processor if real-time endoscopic images cannot be observed, as it may result in patient injury.
-------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------

5.2 Operation

	Slowly turn the rocker during using, as it may result in damage to the bending function and injury, bleeding or perforation to patient.
	Do not insert or withdraw the endoscope under the following conditions, otherwise, it would result in patient injury, bleeding or perforation. -When the endoscopic accessories for diagnosis and treatment protrude from the distal tip of the endoscope; -When the bending section is at a bending angle.

	<p>In case of any of the following conditions during the examination, the examination shall be suspended immediately and the endoscope shall be withdrawn from the patient as described in Section “6.Troubleshooting,” the continued use of which may result in injury, bleeding, or perforation to the patient.</p> <ul style="list-style-type: none"> -Any malfunction of endoscope; -Unexpected disappearance, freezing, abnormal color or incomplete display of live image; -Noisy point, blurring or fogging on the live image; -Rocker cannot be adjusted.
	<p>If the endoscope is dropped or the distal tip is hit by a hard object, stop using the endoscope even if the camera module is not damaged by visual check.</p>

5.2.1 Grip

- (1) Grip the handle of the endoscope with one hand firstly;
- (2) Adjust the rocker with the thumb;
- (3) Hold the insertion portion with the other hand.

5.2.2 Perfusion

Infuse perfusate through the forceps/perfusion plug and enter the patient's body cavity through the instrument channel.

	<p>In the case of perfusion, when the endoscopic accessories for diagnosis and treatment pass through the instrument channel, if there is no sufficiently clear image shown, stop the examination and withdraw the endoscopic accessories for diagnosis and treatment, and then perform perfusion.</p>
-------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

5.2.3 Insert endoscope

While observing the endoscope image, slowly insert the insertion portion into the patient's body cavity.

If necessary, slowly adjust the rocker to adjust the position of the tip to continue the insertion.

	<p>Do not apply products such as olive oil or products containing petroleum-based lubricants (e.g. Vaseline®) to the endoscope, as these products may strain and damage the outer skin of the bending section.</p>
	<p>Do not over bend the insertion portion near the handle connection, as the insertion portion may be damaged.</p>
	<p>If the accessories for diagnosis and treatment are used with the bending section in a bent position, the rocker shall be locked with a snap lock to maintain the angle of the distal tip.</p>
	<p>Do not adjust the rocker with force when the snap lock is locked, as it would result in damage to the endoscope.</p>

5.2.4 Use of endoscopic accessories for diagnosis and treatment

- (1) Verify that the endoscopic accessories for diagnosis and treatment are closed or included in the sheath;
- (2) Slowly and straightly insert the endoscopic accessories for diagnosis and treatment and straightly into the instrument channel port;
- (3) Hold the accessories for diagnosis and treatment by hand approximately 4 cm from the entrance of the instrument channel, and slowly advance it while observing the endoscopic image and inserting it flat into the

instrument channel;

(4) When withdrawing the accessories for diagnosis and treatment from the instrument channel, the accessories for diagnosis and treatment shall be closed or retracted into the sheath, and then withdrawn from the instrument channel slowly and straightly while observing the live image.

	When applying the accessories for diagnosis and treatment, please ensure that the distance between the camera module and the mucosa is greater than the minimum visual distance of the endoscope to guarantee that the accessories for diagnosis and treatment always appear in the image, as it would result in patient injury or device damage.
	If it is difficult to insert or withdraw the accessories for diagnosis and treatment, the bending section shall be straightened as much as possible without losing the image. Forceful insertion or withdrawal of the accessories for diagnosis and treatment would result in damage to the instrument channel or the accessories for diagnosis and treatment, or cause patient injury, bleeding or perforation.
	If the accessories for diagnosis and treatment are not visible in the endoscopic image, discontinue the use of the accessories for diagnosis and treatment, as it would result in patient injury, bleeding or perforation or device damage.
	Do not insert the accessories for diagnosis and treatment by force or abruptly, as it would result in accidental protrusion of the accessories for diagnosis and treatment from the distal tip of the endoscope, causing damage to the device and patient injury, bleeding or perforation.
	Liquid may be ejected when withdrawing the accessories for diagnosis and treatment from the instrument channel. To prevent this situation, cover the accessories and the channel surroundings with a piece of gauze during the withdrawal.
	If it is not possible to withdraw the endoscopic accessories for diagnosis and treatment from the endoscope, close the accessories for diagnosis and treatment, and then withdraw the endoscope and the accessories for diagnosis and treatment carefully while observing the image, so as not to cause soft tissue trauma.
	Do not fire the laser until the laser probe is put in the correct position displayed in the endoscopic image, keep a proper distance between the target and the tip of the endoscope, and use the lowest possible power output in case of patient injury, burn, bleeding, perforation, or device damage.

5.2.5 Withdrawal of endoscope

- (1) Rotate the snap lock to the unlocked position as indicated by the “LOCK→” arrow;
- (2) Adjust the rocker to return the bending section to its straightened state;
- (3) Observe the endoscopic image while carefully withdrawing the endoscope;
- (4) Please dispose of this product in accordance with hospital waste disposal procedures after surgery, and do not reuse it.

	If the endoscope is unable to be withdrawn smoothly from the patient, do not attempt to withdraw it with force. Forced withdrawal would result in injury, bleeding or perforation of the patient.
	Do not allow the body fluid of the patient sticking to the endoscope to come into contact with the bed or floor while withdrawing the endoscope, otherwise, the patient’s body fluids would pose a risk of infection to the patient or medical staff.



6. Troubleshooting

Item	Failure	Possible cause	Solution
Image quality or brightness	No image display	Failure to open the Endoscopic Video Image Processor	Check the power connection and turn on the endoscopic video image processor.
		The video cable is not connected firmly	Verify that the side with the pair markings is facing upwards, insert the video cable firmly into the endoscopic video image processor until hearing the “click” sound.
	Image is too dark or too bright	Improper brightness setting	Adjust the brightness setting of the Endoscopic Video Image Processor as described in the IFU.
Infuse perfusate	No inflow of perfusate/leakage of perfusate	Improper installation of the clamp/perfusion plug	Rotate the fastening ring to align the instrument channel port with the fastening ring.
		Instrument channel leakage	Stop the use of this endoscope.
		Instrument channel blockage	Stop the use of this endoscope.
Endoscopic accessories for diagnosis and treatment	The accessories for diagnosis and treatment cannot pass through the instrument channel smoothly	Instrument channel blockage	Withdraw the accessories for diagnosis and treatment, try to insert carefully again, if it still cannot pass smoothly, please stop the use of the endoscope.
		Use of incompatible accessories for diagnosis and treatment	Select compatible accessories for diagnosis and treatment.
Bending angle	Rocker adjustment obstruction	Snap lock is locked	Rotate the snap lock to the unlocked position as indicated by the “LOCK→” arrow.
	Bending angle cannot be adjusted	Endoscope damage	Stop the use of the endoscope , if it has been inserted into the body cavity of the patient, carefully withdraw the endoscope while observing the image.
Withdraw abnormal endoscope	Endoscope is withdrawn when the image is displayed normally	/	<p>(1) Switch the endoscopic video image processor to the real-time image display interface, and if the image magnification function is applied, please disable the function first;</p> <p>(2) If accessories for diagnosis and treatment are used, first close the accessories for diagnosis and treatment and retract them into the sheath, then adjust the rocker to return the bending section to the straightened state, and then carefully and slowly withdraw the accessories for diagnosis and treatment;</p> <p>(3) While keeping the bending section straightened, slowly and carefully withdraw the endoscope while observing the image.</p>

<p>Withdraw abnormal endoscope</p>	<p>The endoscope is withdrawn when the image cannot be displayed/the image cannot be recovered after freezing</p>	<p>/</p>	<p>(1) Turn off the endoscopic video image processor and then turn it on again, if the real-time image display function is restored, it can continue to be used normally; (2) If the real-time image still cannot be displayed, first close the accessories for diagnosis and treatment and retract them into the sheath, then adjust the rocker to the middle to return the bending section to the straightened state, and then carefully and slowly withdraw the accessories for diagnosis and treatment; (3) While keeping the bending section straightened, withdraw the endoscope slowly; (4) Turn off endoscopic video image processor.</p>
------------------------------------	-------------------------------------------------------------------------------------------------------------------	----------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

7. Other Items

7.1 Description of markings

Symbol	Description	Symbol	Description
	Refer to instruction manual/booklet		Do not use if package is damaged and consult instructions for use
	Endoscopic equipment		The device should be sent to the special agencies according to local regulations
	Type BF applied part		Use-by date
	Do not re-use		Temperature limit
	Do not re-sterilize		Humidity limitation
	Sterilized using ethylene oxide		Atmospheric pressure limitation
	Unsafe under nuclear magnetic resonance examination environment		Medical device
	Batch code		Unique Device Identifier
	Serial number		Single sterile barrier system
	Manufacturer		Keep dry
	Date of manufacture		Warning
	Fragile, handle with care		Caution
	Upward		U.S. Federal law restricts this device to sale by or on the order of a licensed practitioner

8. Electromagnetic Compatibility Information

	This product shall be used only in the following specified electromagnetic environment. Users and medical staff shall ensure that the device is used in such an environment.
	Leakage current multiplication - patient risk: In case that the product is used in combination with medical electrical devices and/or motorized endoscopic accessories, the leakage current may be multiplied.
	If this product or endoscopic accessories are used with products of different manufacturers, as well as with medical electrical products, it shall be ensured that the conditions of BF (insulated, ungrounded part of the application) are satisfied, and products with the following symbols indicate compliance with BF regulations: 
	Using accessories, transducers and cables beyond the requirements of manufacturer can increase the emission of the equipment or system or reduce the immunity except they are sold by the manufacturer of the equipment or system as spare parts of internal components. There are no other accessories and cables in the endoscope itself. It must be used together with the endoscopic video image processor (Model:ASUS-IPD-V100) distributed by Cook Medical to ensure that it meets and complies with the interference radiation and immunity standards. See the manual of endoscopic video image processor (Model:ASUS-IPD-V100) for the requirements of relevant cable.
	Devices or systems should not be used in close proximity or stacked with other devices, and if they must be used in close proximity with other devices or stacked, they shall be observed to verify that they can operate properly in the configuration in which they are used.
	Active medical devices are subject to special EMC precautions and therefore shall be installed and used in accordance with these guidelines.
	Portable and mobile communication RF devices may affect the normal use of medical electrical devices.

Table 1 Guidance and manufacturer's declaration - electromagnetic emissions

Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class A
Harmonic emissions IEC 61000-3-2	Not Applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable

Table 2 Guidance and manufacturer's declaration - electromagnetic Immunity

Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV power supply lines ±1 kV signal input/output 100 kHz repetition frequency	±2 kV power supply lines ±1 kV signal input/output 100 kHz repetition frequency
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U _T ; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % U _T ; 1 cycle and 70 % U _T ; 25/30 cycles; Single phase: at 0°. 0 % U _T ; 250/300 cycle	0 % U _T ; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % U _T ; 1 cycle and 70 % U _T ; 25/30 cycles; Single phase: at 0°. 0 % U _T ; 250/300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF IEC61000-4-3	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
NOTE U _T is the a.c. mains voltage prior to application of the test level.		

Table 3 Guidance and manufacturer's declaration - electromagnetic Immunity

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380–390	TETRA 400	Pulse modulation 18 Hz	27	27
	450	430–470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	28	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	9	9
	745					
	780					
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28	28
	870					
	930					
	1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	28	28
	1845					
	1970					
	2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28	28
	5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9	9
5500						
5785						

Table 4 Guidance and manufacturer's declaration - electromagnetic Immunity

Radiated RF IEC61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	Test Frequency	Modulation	IEC 60601-1-2 Test Level(A/m)	Compliance level (A/m)
	30 kHz	CW	8	8
	134,2 kHz	Pulse modulation 2,1 kHz	65	65
	13,56 MHz	Pulse modulation 50 kHz	7,5	7,5

9. Disposal

The requirement not to dispose of WEEE as unsorted municipal waste and to collect such WEEE separately.

In accordance with Directive 2012/19/EU and national disposal regulations regarding old electrical and electronic devices, please be advised that such items must be disposed of in a special way within the European Union (EU). These regulations require the environmentally friendly recycling/disposal of old electrical and electronic devices. Such items must not be disposed of as domestic refuse. This has been expressed using the icon of the “crossed out trash can”.

The Ascend™ Single-Use Flexible Ureteroscope (Model: ASUS-86-900R-DC, ASUS-86-300R-DC, ASUS-86-300S-DC, ASUS-86-900S-DC, ASUS-86-900R, ASUS-86-300R, ASUS-86-300S, ASUS-86-900S) contains recyclable electronic waste. To protect the environment, do not dispose of it in household waste as unsorted municipal waste, but collect separately and take it to appropriate local collection center.

The Ascend™ Single-Use Flexible Ureteroscope (Model: ASUS-86-900R-DC, ASUS-86-300R-DC, ASUS-86-300S-DC, ASUS-86-900S-DC, ASUS-86-900R, ASUS-86-300R, ASUS-86-300S, ASUS-86-900S) is classified to small equipment (5th category) according to ANNEX IV of Directive 2012/19/EU.

Other countries

For country-specific information on disposal, contact your local dealers or importer.

Manufactured By:

Guangzhou Red Pine Medical Instrument Co.,Ltd.

12 F, No.87 Luoxuan Avenue, Guangzhou International Bioisland, Huangpu District,
510000 Guangzhou City, Guangdong Province, PEOPLE'S REPUBLIC OF CHINA

Distributed By:

Cook Medical LLC
400 Daniels Way
Bloomington, IN 47404, USA
www.cookmedical.com

Not all RPNs are available in specific markets. Contact your Cook representative for information.

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