Instructions for Use (IFU)

Ascend™ Single-Use Flexible Ureteroscope

Model: ASUS-75-900R, ASUS-75-300R, ASUS-75-300S, ASUS-75-900S

Contents

Page

English (Instructions for Use).....2
Français (Manuel d'utilisation).....22

TABLE OF CONTENTS

1. Indications for Use	3
2. Product Overview	3
3. Contraindications	6
4. Warnings	7
5. Preparation and Operation	12
6. Disposal Instructions	14
7. Storage and Transportation Conditions	14
8. Troubleshooting	15
9. Description of Symbols	17
10. Electromagnetic Compatibility Information	18

1 Indications for Use

The device has been developed to be used with endoscopic accessories and Endoscopic Video Image Processor for endoscopy and endoscopic surgery within urinary tract and interior of the ureter and kidney.

2 Product Overview

2.1 Product name

Ascend™ Single-Use Flexible Ureteroscope

2.2 Model

ASUS-75-900R, ASUS-75-300R, ASUS-75-300S, ASUS-75-900S

2.3 Intended patient population

Adult

2.4 Requirement of operator

(1) Background:

- A physician or a medical staff under the supervision of a physician, and he or she must receive sufficient training on clinical endoscope technology.

(2) Language understanding:


- Can understand English or local language in the IFU and displayed on touchscreen console.

(3) Experience:

- Trained doctors.

















- Professional doctors or medical staff who have studied and fully understood the IFU.

2.5 Structure composition

ASUS-75-900R

Structural differences

Model	Mark	Deflection	Distal tip	Position of endoscopic accessories in the image
ASUS-75-900R		Reverse deflection 		
ASUS-75-300R		Reverse deflection 		
ASUS-75-900S		Standard deflection 		
ASUS-75-300S		Standard deflection 		

No.	Name	Functions
1	Rocker	<p>When the rocker is rotated in the "U" direction, the bending section is deflected up;</p> <p>When the rocker is rotated in the "D" direction, the bending section is deflected down;</p> <p>When the rocker is rotated to the "-", the bending section returns to the straightened state.</p>

2	Handle	The handle enables the user to grip and operate the endoscope.
3	Instrument channel	This channel is used for the insertion of endoscopic accessories and for feeding fluid.
4	T-valve	For injecting fluid and delivering the fluid to the distal tip.
5	Instrument channel port	For inserting endoscopic accessories.
6	Insertion portion	For inserting into the patient's body cavity.
7	Insertion tube	For connecting the handle and the bending section and protecting instrument channel, module cables, etc. from contacting the patient's body.
8	Bending section	For adjusting and moving the distal tip of the endoscope.
9	Distal tip	Contains the camera module, LED carriers and outlet of endoscopic accessories.
10	Video cable interface	Connects the endoscope to the Endoscopic Video Image Processor through the video cable.
11	Video cable	Transmits the electrical signal to the Endoscopic Video Image Processor through the video cable to display images captured by the endoscope.

2.6 Product performance

Item		Parameter
Model		ASUS-75-900R/ASUS-75-300R ASUS-75-300S/ASUS-75-900S
Basic Size	Working length	670 mm
	Maximum insertion portion width	2.75 mm (8.25 Fr)
	OD of flexible insertion portion	Φ2.50 mm (7.5 Fr)
	ID of instrument channel	Φ1.20 mm (3.6 Fr)
	Minimum instrument channel width	1.10 mm (3.3 Fr)
Optical Properties	Field of view	120°
	Direction of view	0°
	Depth of field	3 mm~50 mm
Angle range of bending section		285° up~285° down
Operating Environment	Ambient temperature	10°C~40°C
	Relative humidity	30%~85%
	Air pressure	70 kPa~106 kPa

		Prohibited to use in MRI environment
--	---	--------------------------------------

2.7 Product compatibility

The endoscope is intended to be used in combination with the image processor as follows.

Product Name	Model
Endoscopic Video Image Processor HD	ASUS-IPD-V200

The endoscope is intended to be used in combination with the recommended endoscopic accessories with proper parameters as follows.

Endoscopic Accessories	Parameters
Stone Extraction Basket	Maximum OD \leq 2.2 Fr
Laser Fiber	Maximum OD \leq 276 μ m

2.8 Sterile level

The product is sterilized with ethylene oxide and is a single-use sterile product with a SAL of 10^{-6} .

2.9 Shelf-life

3 years.

3 Contraindications

- (1) Uncontrolled systemic hemorrhagic disease.
- (2) Weak constitution, cannot tolerate anesthesia.
- (3) Benign prostatic hyperplasia affects ureteral ureteroscopy insertion.
- (4) Lower urinary tract tumor.
- (5) Bladder contracture.
- (6) Uncontrolled hemorrhagic disease of the whole body.
- (7) Uncontrolled urinary tract infection.
- (8) Severe cardiopulmonary insufficiency, unable to tolerate surgery.
- (9) Severe urethral stricture and ureteral stricture.
- (10) Severe hip deformity with difficult lithotomy position.
- (11) Pregnancy.

4 Warnings

4.1 Preoperative warnings

4.1.1 General preoperative warnings

- (1) Prior to the use of this Ascend™ Single-Use Flexible Ureteroscope (hereinafter referred to as "endoscope"), please read this IFU carefully, including all precautions, warnings and tips. Failure to do so may result in injury to the patient.
- (2) No modification of this product is allowed without the permission of the manufacturer.
- (3) This endoscope does not contain any user-serviceable parts. Do not disassemble, modify or attempt to repair it, patient or user injury or equipment damage may happen.
- (4) Do not attempt to clean and reuse the device, as it is a single use device. Reuse of the product can cause contamination leading to infection.
- (5) In many therapeutic areas, this endoscope may be used in conjunction with endoscope accessories. For specific operation of endoscope accessories and related safety information, refer to the IFU of the corresponding accessories.
- (6) Always make sure to select the endoscopic accessory with the correct size. There is no guarantee that instruments selected solely using the working channel size will be compatible in combination. The compatibility of selected instruments should be tested before the procedure.

4.1.2 Preoperative warnings regarding usage environment

- (1) This endoscope shall be used only in the electromagnetic environment specified in "10. Electromagnetic Compatibility Information". Users and medical staff shall ensure that the device is used in such an environment.
- (2) Electromagnetic interference may be generated when using this product near equipment marked with the following symbols or other portable mobile RF (radio frequency) communication devices (such as cellular phones). In case of electromagnetic interference, take appropriate measures, such as adjusting the direction, location or shielding the using site of this device.

- (3) To eliminate electromagnetic interference from other devices (any device other than this product or the components constituting this system), observe whether this system is operating properly in the environment in which it is used.
- (4) If the endoscope is affected by strong static electricity or X-ray, it would result in abnormal endoscopic image display.
- (5) This endoscope is not an MR safe product and cannot be used safely in an MRI environment. Simultaneous use of nuclear magnetic resonance imaging (NMRI), endoscopes and accessories would result in hazards and generate artificial artifacts; please follow the instructions and safety precautions specified by the manufacturer.
- (6) Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the endoscope, including cables specified by the manufacturer. Otherwise, it would result in degradation of the product performance.

- (7) This endoscope or systems should not be used in close proximity or stacked with other devices, and if they must be used in close proximity with other devices or stacked, they shall be observed to verify that they can operate properly in the configuration in which they are used.
- (8) The operation must be carried out in a visual state and placed in a place that meets the requirements.
- (9) Leakage current multiplication - patient risk: In case that the product is used in combination with medical electrical devices and/or motorized endoscopic accessories, the leakage current may be multiplied.
- (10) If this endoscope or endoscopic accessories are used with products of different manufacturers, as well as with medical electrical products, it shall be ensured that the conditions of BF (insulated, ungrounded part of the application) are satisfied, and products with the following symbols indicate compliance with BF regulations:

- (11) Using accessories, transducers and cables beyond the requirements of manufacturer can increase the emission of the equipment or system or reduce the immunity except they are sold by the manufacturer of the equipment or system as spare parts of internal components.

There are no other accessories and cables in the endoscope itself. It must be used together with the Endoscopic Video Image Processor HD (Model: ASUS-IPD-V200) distributed by Cook Medical to ensure that it meets the interference radiation and immunity standards. See the IFU of Endoscopic Video Image Processor HD (Model: ASUS-IPD-V200) for the requirements of relevant cables.

- (12) Active medical devices are subject to special EMC precautions and therefore shall be installed and used in accordance with these guidelines.
- (13) The emissions characteristics of this endoscope make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this endoscope might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

4.1.3 Preoperative warnings regarding the operator

- (1) The operator must be physicians or healthcare providers under supervision of physicians and they must have received sufficient training on clinical endoscope technology.
- (2) Personal protective equipment shall be worn during operation, so as to prevent the operators from exposure to potentially infectious substances. Personal protective equipment should prevent skin exposure with appropriate size and enough length.

4.1.4 Preoperative warnings regarding the product preparation

Have backup system ready in case of malfunction.

4.1.5 Preoperative warnings regarding product inspection

- (1) Prior to use always perform an inspection and functionality check according to the “5.1. Preparation and check”. Inspect the diagnostic and treatment accessories in conjunction with this device according to their respective instruction manuals. Do not use the device if the device or the diagnostic and treatment accessories of the device is damaged in any way or the functionality check still fails after having been corrected referring to “8.Troubleshooting”. Otherwise, it would result in injury or infection on the patient or user.
- (2) This endoscope is a single-use device. Before initial use, make sure the packaging is intact. If the packaging is damaged, do not use the product.
- (3) Do not use it if the product is beyond its shelf life of three years.
- (4) Do not use the endoscope in case of abnormal bending functions or the bending section cannot be straightened. Otherwise, it could result in patient injury, bleeding or perforation.
- (5) Do not insert the video cable if it becomes damp or contaminated, as this could result in electric shock and serious damage to the endoscope, and also jeopardize the safety of the patients or the operators.

4.2 Intraoperative warnings

4.2.1 Intraoperative warnings regarding operation tips

- (1) To prevent potential thermal injury, begin at the lowest brightness setting.
- (2) Do not exert excessive force when operating.
- (3) Please perform operations by following the precautions below, or the endoscopic image may disappear unexpectedly during the examination or the frozen image cannot be recovered.
 - ①Insert the video cable into the Endoscopic Video Image Processor and pull it gently to verify that the connection is firm;
 - ②Do not bend, tap, pull or twist the endoscope, as this would damage the endoscope and cause water leakage or damage to other parts;
 - ③Prior to connecting the video cable to the Endoscopic Video Image Processor, verify that the video cable is completely dry, otherwise it would result in device malfunctions.
- (4) In case of any of the following conditions during the examination, do not use endoscopic accessories, insert or withdraw the insertion portion of the endoscope, or perform angle control, as this could result in injury, bleeding or perforation to patient. Perform troubleshooting as described in section “8. Troubleshooting” of this IFU.
 - ①No real-time image shown;
 - ②Real-time image is frozen.
- (5) LED illumination is required for endoscopic observation, and improper use would result in viable tissues alterations such as protein denaturation and tissue perforation of viable tissues, please follow the precautions as follows:
 - ①When the distal tip contacts body, there will be no normal live image shown. The maximum temperature

of the distal tip may reach up to 60 °C at the 40°C environment, which would result in mucosal burns. Do maintain the proper distance required for adequate observation and do not constantly observe tissue at close range or allow the distal tip of the endoscope to remain in unnecessary contact with the mucosa for prolonged duration, as this would result in an elevated temperature at the tip (possibly exceeding 41°C);

②Be sure to disconnect the endoscope when it is not used for a long time to avoid unnecessary illumination;

③Do not direct light into the eyes to avoid light damage to the eyes.

- (6) If the endoscope is used under extremely poor conditions, LED may be damaged.
- (7) The endoscope is connected to the power supply network through the isolation measures of the Endoscopic Video Image Processor HD (Model: ASUS-IPD-V200).
- (8) Slowly turn the rocker during using, or it may result in damage to the bending section or cause injury, bleeding or perforation to patient.
- (9) Please plug/unplug and connect the video cable with the Endoscopic Video Image Processor off, otherwise it would result in damage to the device.
- (10) Do not excessively twist the video cable, as it would result in abnormal images.
- (11) Before each use or after a change of viewing modes/settings, the operator should check to ensure the view observed through the device provides a live image (rather than a stored one) and has the correct image orientation. Otherwise, it could result in patient injury or equipment damage.
- (12) Before using endoscopic accessories, the outer surface of the accessories which are intended to be inserted into the patient should be inspected for any rough surfaces, sharp edges, protrusions, or other irregularities that may cause injury to the patient.
- (13) When using the endoscopic accessory, ensure that the distance between the camera module and the mucosa is greater than the minimum visual distance of the endoscope to guarantee that the accessory for diagnosis and treatment always appears in the image, or it could result in patient injury or device damage.
- (14) If it is difficult to insert or withdraw the accessory for diagnosis and treatment, the bending section shall be straightened as much as possible without losing the image. Forceful insertion or withdrawal of the accessory for diagnosis and treatment would result in damage to the instrument channel or the accessory for diagnosis and treatment, or cause patient injury, bleeding or perforation.
- (15) Fluid may be ejected during the withdraw of the endoscopic accessory from the instrument channel. To prevent this situation, cover the accessory and the channel surroundings with a piece of sterile gauze during the withdrawal.
- (16) When using compatible laser equipment:
 - ①User shall take steps to protect eye and skin, including but not limited to using personal protective equipment;
 - ②Prepare, inspect and connect the laser unit and laser probe as described in their IFU;
 - ③Do not fire the laser until the laser probe is put in the correct position displayed in the endoscopic image, keep a proper distance between the target and the distal tip, and use the lowest possible power output to prevent patient injury, burn, bleeding, perforation, or device damage;
 - ④Do not use a damaged laser probe. A laser probe with a damaged sheath or distal tip may cause patient injury and/or equipment damage.

4.2.2 Intraoperative warnings regarding operation prohibition

- (1) Do not damage the insertion portion during use, as it may leave parts of the product inside the patient or expose sharp surfaces that may cause damage to mucosa.
- (2) Do not hit, collide or drop the endoscope, and do not bend, pull or twist the endoscope, as this may cause damage to the endoscope or dislodge parts into the patient, resulting in possible injury, bleeding or perforation.
- (3) Do not pull, twist or rotate the bending section which is already bent at an angle, as this may cause injury, bleeding or perforation of the patient.
- (4) Do not press the bending section or adjust the bending angle with excessive force, or the endoscope may be damaged and affect the bending function.
- (5) Do not over bend the insertion portion near the handle connection, as the insertion portion may be damaged.
- (6) Do not insert or withdraw the insertion portion or endoscopic accessory suddenly or with excessive force as this could result in injury, bleeding or perforation to patient.
- (7) Do not insert or withdraw the endoscope under the following conditions, otherwise, it could result in patient injury, bleeding or perforation.
 - When the endoscopic accessory protrudes from the distal tip of the endoscope;
 - When the bending section is at a bending angle.
- (8) Do not pull the video cable during examination, as it would result in abnormal images.
- (9) Do not move the Endoscopic Video Image Processor, video cable, adapter, power line etc. during use, as it could result in abnormal image.
- (10) To avoid interference and injury, do not use diathermy device in conjunction with the endoscope.

4.2.3 Intraoperative warnings regarding stopping operation

- (1) In case of any of the following conditions during the examination, suspend the examination immediately and withdraw the endoscope from the patient as described in Section “8. Troubleshooting”, the continued use of which may result in injury, bleeding, or perforation to the patient.
 - Any malfunction of endoscope;
 - Unexpected disappearance, freezing, abnormal color or incomplete display of live image;
 - Noisy point, blurring or fogging on the live image;
 - Rocker cannot be adjusted.
- (2) Stop using the endoscope if the endoscope is dropped or the distal tip is hit by a hard object.
- (3) The amount of fluid will be reduced when an endoscopic accessory is inserted in the instrument channel. Confirm the amount of fluid in advance.

4.2.4 Intraoperative warnings regarding troubleshooting

- (1) Darkening of the endoscope images during use may indicate that blood or mucous membranes are adhering to the camera module at the distal tip of the endoscope. Carefully withdraw the endoscope from the patient, wipe away the blood or mucous membranes, and restore optimal illumination to ensure a safe examination.
- (2) Many problems (e.g. voltage dips/interruptions testing) can be solved by restarting the Endoscopic Video Image Processor manually during procedure. Restarting the image processor should be tried in the following order:
 - a. Stop the operation and ensure that the power is restored and the interference disappears.
 - b. Restart the Ureteroscope system to display image.If the problem cannot be resolved by the actions, contact Cook Medical.

4.3 Postoperative warnings

- (1) Do not exclusively rely on endoscopic observation as the initial diagnostic basis for a lesion to make any decisions about potential diagnostic or device interventions.
- (2) This endoscope is single-use medical device. Please dispose it according to the regulations of the local hospital after use.
- (3) In the event of any serious incident related to the device, please contact Cook Medical.

5 Preparation and Operation

5.1 Preparation

Please prepare and check the product in accordance with the following instructions before use. In case of any abnormality, please refer to “8. Troubleshooting” of this IFU for troubleshooting. If there is still a problem with the device, please do not use this product and contact Cook Medical.

5.1.1 Check the endoscope

- (1) Check the handle for any deformation, loose parts or other abnormalities;
- (2) Check the camera module at the distal tip of the endoscope for any scratches, cracks, stains, gaps or other abnormalities;
- (3) Check the entire insertion portion of the endoscope for any abnormalities such as damages, stains and cracks;
- (4) Check whether the functional components of the endoscope are complete.

5.1.2 Check the bending function

- (1) Straighten the bending section;
- (2) Slowly adjust the rocker in both directions to the maximum angle to verify that the bending section bends smoothly and that the maximum angle can be reached;

(3) Slowly adjust the rocker to the middle position to verify that the bending section can return to the straightened state smoothly.

5.1.3 Connect the endoscope to the Endoscopic Video Image Processor

- (1) Connect the Endoscopic Video Image Processor to the power supply;
- (2) Connect the video cable interface to the Endoscopic Video Image Processor according to the pair markings;
- (3) Turn on the Endoscopic Video Image Processor and verify that the endoscope image is displayed properly.

5.1.4 Check the endoscopic live image

- (1) Verify that the endoscope LED illumination is functioning properly;
- (2) Properly adjust the screen brightness;
- (3) Verify that the images are free from noisy point, blurring, fogging or other abnormalities;
- (4) Slowly adjust the maximum angle of rocker up and down respectively, and observe whether the image is abnormal;
- (5) Verify that the image shows no transient disappearance or other abnormalities;
- (6) Verify that real-time endoscopic image appears on the screen by pointing the distal tip of the endoscope at an object, such as the palm.

5.1.5 Check the instrument channel

- (1) Loosen the turncap of the T-valve and insert the endoscopic accessory through the instrument channel port. Verify that the endoscopic accessory extends smoothly from the instrument channel outlet and that no foreign matters are present at the distal tip.
- (2) Verify that the endoscopic accessory can be smoothly withdrawn from the instrument channel port.

5.2 Operation

5.2.1 Grip

- (1) Grip the handle of the endoscope with one hand firstly;
- (2) Hold the insertion tube with the other hand;
- (3) Adjust the rocker with the thumb to control the bending section.

5.2.2 Infusion

Infuse fluid through the T-valve.

5.2.3 Insert the endoscope

While observing the endoscope image, slowly insert the insertion portion into the patient's body cavity. If necessary, slowly adjust the rocker to adjust the position of the distal tip to continue the insertion.

5.2.4 Adjust the bending angle

Operate the rocker as needed to guide the insertion and observation of the tip.

5.2.5 Use the endoscopic accessory

- (1) Verify that the endoscopic accessory is closed, turned off or retracted into its sheath;
- (2) Slowly and straightly insert the endoscopic accessory into the instrument channel port;
- (3) Hold the endoscopic accessory and slowly advance it while observing the endoscopic image and inserting it flat into the instrument channel;
- (4) Ensure that the endoscopic accessory is closed, turned off or retracted into its sheath first. Then withdraw it from the instrument channel slowly and straightly while observing the live image.

5.2.6 Withdraw the endoscope

- (1) Adjust the rocker to return the bending section to its straightened state;
- (2) Observe the endoscopic image while carefully withdrawing the endoscope.

6 Disposal Instructions

When disposing of the endoscope or the endoscopic accessories, follow hospital waste disposal procedures and all applicable national and local laws and guidelines, and do not reuse it.

7 Storage and Transportation Conditions

Store the endoscope in a clean and dry condition as shown in the table below.

Storage Environment	Ambient temperature	3°C ~ 30°C
	Relative humidity	25% ~ 85%
	Air pressure	70 kPa ~ 106 kPa

Transportation environment condition should be in accordance with the table below.


Transportation Environment	Ambient temperature	-29°C ~ 55°C
	Relative humidity	10% ~ 95%
	Air pressure	70 kPa ~ 106 kPa

8 Troubleshooting

Item	Situation	Possible cause	Solution
Image quality or brightness	No image display	Failure to turn on the Endoscopic Video Image Processor	Check the power connection and turn on the Endoscopic Video Image Processor.
		The video cable is not connected firmly	Verify that the side with the pair markings is facing upwards, unplug both ends of the video cable and plug them in back.
	Image is too dark or too bright	Improper brightness setting	Adjust the brightness setting of the Endoscopic Video Image Processor as described in the IFU.
	Image is not clear	The camera module is dirty	Withdraw the endoscope and wipe the camera module with sterile gauze.
Infusion	Leakage of fluid/no inflow of fluid	Instrument channel port is not closed	Tighten the turncap of T-valve.
		Instrument channel leakage	Stop the use of this endoscope.
		T-valve is not opened.	Open the T-valve.
		Instrument channel blockage	Stop the use of this endoscope.
Endoscopic accessories for diagnosis and treatment	Accessories for diagnosis and treatment cannot pass through the instrument channel smoothly	Instrument channel blockage	Withdraw the accessory for diagnosis and treatment, and try to insert it carefully again. If it still cannot pass through smoothly, please stop the use of the endoscope.
		Use of incompatible accessories for diagnosis and treatment	Select compatible accessories for diagnosis and treatment.
Bending angle	Bending angle cannot be adjusted	Endoscope damage	Stop the use of the endoscope. If it has been inserted into the body cavity of the patient, carefully withdraw it while observing the image.
Withdrawal of abnormal endoscope	Endoscope is withdrawn when the image is displayed normally	/	<p>(1) Switch the Endoscopic Video Image Processor to the real-time image display interface, and if the image magnification function is applied, please disable the function first;</p> <p>(2) If an accessory for diagnosis and treatment is used, first close the accessory for diagnosis and treatment, then adjust the rocker to return the bending section to the straightened state, and then carefully and slowly withdraw the accessory for</p>

			<p>diagnosis and treatment;</p> <p>(3) While keeping the bending section straightened, slowly and carefully withdraw the endoscope while observing the image.</p>
	<p>Endoscope is withdrawn when the image cannot be displayed/the image cannot be recovered after freezing</p>	<p>/</p>	<p>(1) Turn off the Endoscopic Video Image Processor and then turn it on again. If the real-time image display function is restored, it can continue to be used normally;</p> <p>(2) If the real-time image still cannot be displayed, first close the accessory for diagnosis and treatment, then adjust the rocker to the middle to return the bending section to the straightened state, and then carefully and slowly withdraw the accessory for diagnosis and treatment;</p> <p>(3) While keeping the bending section straightened, withdraw the endoscope slowly;</p> <p>(4) Turn off Endoscopic Video Image Processor.</p>

9 Description of Symbols

Symbol	Description	Symbol	Description
	Refer to instruction manual/booklet		Do not use if package is damaged and consult instructions for use
	Endoscopic equipment		The device should be sent to the special agencies according to local regulations
	Type BF applied part		Use-by date
	Do not re-use		Temperature limit
	Do not re-sterilize		Humidity limitation
	Sterilized using ethylene oxide		Atmospheric pressure limitation
	Unsafe under nuclear magnetic resonance examination environment		Medical device
	Batch code		Unique Device Identifier
	Serial number		Single sterile barrier system
	Manufacturer		Keep dry
	Date of manufacture		Keep away from sunlight
	Fragile, handle with care		Caution
	Upward		Not made with Phthalate (DEHP)
	Stacking limit by number		U.S. Federal law restricts this device to sale by or on the order of a licensed practitioner

10 Electromagnetic Compatibility Information

Essential performance: 1) Endoscopic real-time image can be transmitted normally. In case of interference, the image can be restored after the interference is removed. 2) Illumination can be provided normally. In case of interference, the lighting can be restored after the interference is removed. 3) Be able to display image with correct orientation.

	This endoscope shall be used only in the following specified electromagnetic environment. Users and medical staff shall ensure that the device is used in such an environment.
	Leakage current multiplication - patient risk: In case that the product is used in combination with medical electrical devices and/or motorized endoscopic accessories, the leakage current may be multiplied.
 	If this endoscope or endoscopic accessories are used with products of different manufacturers, as well as with medical electrical products, it shall be ensured that the conditions of BF (insulated, ungrounded part of the application) are satisfied, and products with the following symbols indicate compliance with BF regulations:
	Using accessories, transducers and cables beyond the requirements of manufacturer can increase the emission of the equipment or system or reduce the immunity except they are sold by the manufacturer of the equipment or system as spare parts of internal components. There are no other accessories and cables in the endoscope itself. It must be used together with the Endoscopic Video Image Processor HD (Model: ASUS-IPD-V200) distributed by Cook Medical to ensure that it meets and complies with the interference radiation and immunity standards. See the IFU of Endoscopic Video Image Processor HD (Model: ASUS-IPD-V200) for the requirements of relevant cables.
	Devices or systems should not be used in close proximity or stacked with other devices, and if they must be used in close proximity with other devices or stacked, they shall be observed to verify that they can operate properly in the configuration in which they are used.
	Active medical devices are subject to special EMC precautions and therefore shall be installed and used in accordance with these guidelines.
	The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
	Many problems(e.g. voltage dips/interruptions testing) can be solved by restarting the Endoscopic Video Image Processor manually during procedure. Restarting the image processor should be tried in the following order: a. Stop the operation and ensure that the power is restored and the interference disappears. b. Restart the Ureteroscopy system to display image. If the problem cannot be resolved by the actions, contact Cook Medical .

Table 1 Guidance and manufacturer's declaration - electromagnetic emissions

Emissions test	Compliance	Compliance level		
RF conducted emissions CISPR 11	Group 1 Class A ^{a*}	Frequency range	(Test voltage:120V-230V) ^{b*}	
		MHz	Quasi-peak	Average
		0.15-0.50	79 dB(μV)	66 dB(μV)

		0.50-30	73 dB(μ V)	60 dB(μ V)
RF radiated emissions CISPR 11	Group 1 Class A ^{a*}	Frequency range MHz	Quasi-peak (Test voltage:120V-230V) ^{b*}	
			10 m measuring distance	3 m measuring distance ^{c*}
		30-230	40 dB(μ V/m)	50 dB(μ V/m)
		230-1000	47 dB(μ V/m)	57 dB(μ V/m)
Harmonic emissions IEC 61000-3-2	Not Applicable			
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable			
<p>a* The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.</p> <p>b* The endoscope system is rated power of ≤ 20 k VA, therefore, selection of the appropriate set of limits shall be based on the rated a.c. power stated by the manufacturer.</p> <p>c* The 3 m separation distance applies only to small size equipment meeting the size criterion defined in 3.17 of CISPR 11 Standard.</p>				

Table 2 Guidance and manufacturer's declaration - electromagnetic immunity

Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4	± 2 kV power supply lines ± 1 kV signal input/output 100 kHz repetition frequency	± 2 kV power supply lines ± 1 kV signal input/output 100 kHz repetition frequency
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV differential mode ± 0.5 kV, ± 1 kV, ± 2 kV common mode	Not Applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U_T ; 0.5 cycle. At 0° , 45° , 90° , 135° , 180° , 225° , 270° and 315° . 0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycles; Single phase: at 0° . 0 % U_T ; 250/300 cycle	Not Applicable
Power frequency magnetic field	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz

IEC 61000-4-8		
Conducted RF IEC61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF IEC61000-4-3	3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz
NOTE U _T is the a.c. mains voltage prior to application of the test level.		

Table 3 Guidance and manufacturer’s declaration - electromagnetic immunity

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380–390	TETRA 400	Pulse modulation 18 Hz	27	27
	450	430–470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	28	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	9	9
	745					
	780					
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28	28
	870					
	930					
	1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT;	Pulse modulation 217 Hz	28	28
	1845					
1970						

			LTE Band 1, 3, 4, 25; UMTS			
	2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28	28
	5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9	9
	5500					
	5785					

Table 4 Guidance and manufacturer’s declaration - electromagnetic immunity

Radiated RF IEC61000-4-39	Test Frequency	Modulation	IEC 60601-1-2 Test Level(A/m)	Compliance level (A/m)
(Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	30 kHz	CW	8	8
	134.2 kHz	Pulse modulation 2.1 kHz	65	65
	13.56 MHz	Pulse modulation 50 kHz	7.5	7.5

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Guangzhou Red Pine Medical Instrument Co.,Ltd.
Cook Medical LLC
William A. Cook Australia Pty. Ltd.