Instructions for Use (IFU)

Ascend™ Single-Use Flexible Cysto-Nephroscope HD

Model: ASUS-HD-15R-CYS, ASUS-HD-15S-CYS

Contents	Page
English (Instructions for Use).....	2
Français (Manuel d'utilisation).....	26

TABLE OF CONTENTS

1. Indications for use	3
2. Product Overview	3
3. Contraindications	6
4. Warnings	6
5. Installation and Operation	13
6. Storage conditions	18
7. Disposal instructions	18
8. Troubleshooting	19
9. Description of symbols	21
10. Electromagnetic Compatibility Information	22

1. Indications for use

Intended purpose: This product shall be used in conjunction with the Endoscopic Video Image Processor HD and other peripheral equipment for observation, diagnosis, photography and treatment within the bladder, urethra, and kidney.

The device is suitable for professional healthcare facility environments such as hospitals and clinics.

This product shall not be used for other purposes.

2. Product Overview

2.1 Product name

Ascend™ Single-Use Flexible Cysto-Nephroscope HD

2.2 Model

ASUS-HD-15R-CYS, ASUS-HD-15S-CYS

2.3 Intended patient population

Adult

2.4 Requirement of operator

(1) Background:

- A physician or a medical staff under the supervision of a physician, and he or she must receive sufficient training on clinical endoscope technology.

(2) Knowledge:

- Must have had enough clinical endoscopic technical training.

(3) Language understanding:

- Can understand English or local language in the IFU and displayed on touchscreen console.

(4) Experience:


- Trained doctors
- Professional doctors or medical staff who have studied and fully understood the IFU.

(5) Permissible impairments:





- Mild reading vision impairment or vision corrected to log MAR 0, 2 (6/10 or 20/32)

2.5 Structure composition

ASUS-HD-15R-CYS/ASUS-HD-15S-CYS

Structural differences

Model	Decorative cover	Deflection
ASUS-HD-15R-CYS		Reverse deflection 
ASUS-HD-15S-CYS		Standard deflection 

No.	Name	Functions
1	Rocker	When the rocker is rotated in the “U” direction, the bending section moves up; When the rocker is rotated in the “D” direction, the bending section moves down. When the rocker is rotated to the “-”, the bending section is in a straight position.
2	Handle	The handle enables the user to grip the endoscope.
3	Instrument channel	The endoscopic accessories for diagnosis and treatment and fluid enter the patient's body cavity through the instrument channel.
4	T-valve	For injecting liquid through the T-value and delivering to the distal tip.
5	Instrument channel port	For inserting endoscopic diagnosis and treatment accessories into the instrument channel.
6	Insertion portion	For inserting into the patient's body cavity.
7	Insertion tube	For connecting the handle and the bending section and protecting instrument channel, module cables, etc. from contacting the patient's body.
8	Bending section	For adjusting and moving the distal tip of the endoscope.
9	Distal tip	The distal tip is the export of camera module, LED carriers and accessories for diagnosis and treatment.
10	Video cable interface	Connects the endoscope to the Endoscopic Video Image Processor HD to display image captured.
11	Connector section	Transmits the electrical signal to the Endoscopic Video Image Processor HD to display images captured by the endoscope.

2.6 Product performance

Item		Parameter
Type		ASUS-HD-15R-CYS/ASUS-HD-15S-CYS
Basic Size	Working length	380 mm
	Maximum insertion portion width	Φ5.8 mm(17.4Fr)
	Insertion tube width	Φ5.0 mm(15.0 Fr)
	Outer diameter of the distal tip	Φ5.4 mm(16.2 Fr)
	Minimum instrument channel width	Φ2.12mm(6.36Fr)
	Instrument channel width	Φ2.2 mm(6.6 Fr)
Optical Properties	Field of view	120°
	Direction of view	0°
	Depth of field	3 mm~100 mm
Angle range of bending section		Up 225° Down 225°
Operating	Ambient temperature	10℃~40℃

Environment	Relative humidity	30%~85%
	Air pressure	70 kPa~106 kPa
		Prohibited to use in MRI environment

2.7 Product compatibility

Confirm that the endoscopic accessory for diagnosis and treatment with proper parameter can be inserted into the instrument channel(such as biopsy forceps).

Ho:YAG laser treatment accessories

The Endoscope is compatible with high-frequency electrosurgical units that have a Maximum Rated Repetitive Peak Voltage: Electrocision: 1400Vp (2800Vp-p); Electrocoagulation: 2500Vp (5000Vp-p).

2.8 Sterile level

The product is sterilized with ethylene oxide and is a single-use sterile product with a SAL of 10^{-6} .

2.9 Shelf-life

3 years.

3. Contraindications

- (1) Cases not covered by the "Indications for use" in the IFU;
- (2) Cases that are contraindications prescribed by the hospital's administration or other official institutions such as academic societies on endoscopy, or those not covered by the indications specified by these authorities;
- (3) Cases in which any equipment that is incompatible with this device is used;
- (4) Severe urethral stricture;
- (5) Patients with acute infection (acute urethritis, acute prostatitis, acute epididymitis);
- (6) Patients on anticoagulants or with coagulopathies.

Note: For details about the packing list, please refer to the IFU for each product.

4. Warnings

4.1 Preoperative warnings

4.1.1 General preoperative warnings

- (1) Prior to the use of the Ascend™ Single-Use Flexible Cysto-Nephroscope HD (hereinafter referred to as "endoscope"), please read this IFU carefully, including all precautions, warnings and tips. Failure to do

so may result in injury to the patient.

- (2) No modification of this product is allowed without the permission of the manufacture.
- (3) This product does not contain any user-serviceable parts. Do not disassemble, modify or attempt to repair it; patient or user injury or equipment damage can result. Damage to the device due to improper use and incorrect maintenance by the user, as well as their consequences and responsibilities, shall be taken by the user.
- (4) The endoscope shall be disposed of in accordance with the hospital waste disposal procedures after using, and shall not be reused. Improper handling or reuse may result in cross-contamination or infection of patients.
- (5) In the event of any serious incident related to the device, please contact Cook Medical and local competent authority.
- (6) In many therapeutic areas, this product may be used in conjunction with endoscope therapy accessories. Please follow the instructions and safety precautions for the products and accessories used.
- (7) It is important to select the correct size of endoscope attachment. For more information on the combination of endoscope and accessory for diagnosis and treatment, please refer to the IFU for the accessory. For instructions on the operation of the accessory, refer to respective IFU.

4.1.2 Preoperative warnings regarding usage environment

- (1) This product shall be used only in the specified electromagnetic environment according to "10. Electromagnetic Compatibility Information". Users and medical staff shall ensure that the device is used in such an environment.
- (2) Electromagnetic interference may be generated when using this product near equipment marked with the following symbols or other portable mobile RF (radio frequency) communication devices (such as cellular phones). In case of electromagnetic interference, take appropriate measures, such as adjusting the direction, location or shielding the using site of this device.

- (3) To eliminate electromagnetic interference from other devices (any device other than this product or the components constituting this system), observe whether this system is operating properly in the environment in which it is used.
- (4) If the endoscope is affected by strong static electricity or X-ray, it would result in abnormal endoscopic image display.
- (5) This product is not an MR safe product and cannot be used safely in an MRI environment. Simultaneous use of nuclear magnetic resonance imaging (NMRI), endoscopes and accessories would result in hazards and generate artificial artifacts; please follow the instructions and safety precautions specified by the manufacturer.
- (6) Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the endoscope, including cables specified by the manufacturer. Otherwise, it could result in degradation of the product performance.
- (7) This endoscope should not be used in close proximity or stacked with other devices, and if they must be used in close proximity with other devices or stacked, they shall be observed to verify that they can operate properly in the configuration in which they are used.

- (8) Before endoscopy and endoscopic treatment, remove all metal objects (watches, glasses, necklaces, etc.) from the patient's body. During high-frequency cauterization, if the patient is wearing a metal object, the tissue around the metal object may burn.
- (9) The operation must be carried out in a visual state and placed in a place that meets the requirements.
- (10) Leakage current multiplication - patient risk: In case that the product is used in combination with medical electrical devices and/or motorized endoscopic accessories, the leakage current may be multiplied.
- (11) Before the use of endoscope with the high-frequency accessories, it is necessary to install and connect the equipment according to the IFU, and make sure that noise will not affect the observation and surgery. The surface of the high-frequency accessories entering the patient should be checked before use to ensure that there are no rough surfaces, sharp edges or protrusions that will cause harm. If not, may result in injury to the patient.
- (12) If this product or endoscopic accessories are used with products of different manufacturers, as well as with medical electrical products, it shall be ensured that the conditions of BF (insulated, ungrounded part of the application) are satisfied, and products with the following symbols indicate compliance with BF regulations:

- (13) Using accessories, transducers and cables beyond the requirements of manufacturer may increase the emission of the equipment or system or reduce the immunity to emissions. Spare parts of internal components produced by the manufacturer of the equipment/system may be an exception.

There are no other accessories and cables in the endoscope itself. It must be used together with the Endoscopic Video Image Processor HD (Model: ASUS-IPD-V200) distributed by Cook Medical to ensure that it meets and complies with the interference radiation and immunity standards. See the IFU of Endoscopic Video Image Processor HD (Model: ASUS-IPD-V200) for the requirements of relevant cables.

- (14) Active medical devices are subject to special EMC precautions and therefore shall be installed and used in accordance with these guidelines.
- (15) The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

4.1.3 Preoperative warnings regarding the operator

- (1) The operator must be physicians or healthcare providers under supervision of physicians and they must receive sufficient training on clinical endoscope technology.
- (2) Personal protective equipment shall be worn during operation, so as to prevent the operators from exposure to potentially infectious substances. Personal protective equipment should prevent skin exposure with appropriate size and enough length.
- (3) Using high-frequency current may cause accidental burns. Please wear insulation gloves.

4.1.4 Preoperative warnings regarding the product preparation

- (1) Use of accessories, transducers and cables other than those specified or provided by the manufacturer of the endoscope could result in increased electromagnetic emissions or decreased electromagnetic immunity of the endoscope and result in improper operation.
- (2) This endoscope contains no memory chip, and cannot store endoscopic information. If you need to access the information photographed and recorded by the endoscope, please access it via the Endoscopic Video Image Processor HD.
- (3) Spare device: Be sure to prepare a spare endoscope to avoid interrupting the inspection in case of device failure or abnormal function.

4.1.5 Preoperative warnings regarding product inspection

- (1) Prior to use always perform an inspection and functionality check according to the “5.1. Preparation and check”. Inspect the diagnostic and treatment accessories in conjunction with this device according to their respective instruction manuals. Do not use the device if the device or the diagnostic and treatment accessories of the device is damaged in any way or the functionality check still fails after having been corrected referring to the chapter “8.Troubleshooting”. Otherwise, it would result in injury or infection on the patient or user.
- (2) The product is sterilized with ethylene oxide and it is a single use sterile product. Before initial use, make sure the packaging is intact. If the packaging is damaged, do not use the product.
- (3) The shelf life of this product is three years. Do not use it if the product is beyond its shelf life of three years.
- (4) In case of abnormal bending functions, do not use the endoscope, as the bending section may not be straightened, resulting in patient injury, bleeding or perforation.
- (5) Do not insert the video cable if it becomes damp or contaminated, as this would result in electric shock and serious damage to the endoscope, and also jeopardize the safety of the patients or operators.

4.2 Intraoperative warnings

4.2.1 Intraoperative warnings regarding operation tips

- (1) Please perform operations by following the precautions below, or the endoscopic image may disappear unexpectedly during the examination or the frozen image cannot be recovered.
 - ①Insert the video cable into the Endoscopic Video Image Processor HD and pull it gently to verify that the connection is firm;
 - ②Do not bend, tap, pull or twist the endoscope, as this would damage the endoscope and cause water leakage or damage to other parts;

- ③Prior to connecting the video cable to the Endoscopic Video Image Processor HD, verify that the video cable is completely dry, otherwise it could result in device malfunctions.
- (2) In case of any of the following conditions during the examination, do not use endoscopic diagnosis and treatment accessories, insert or withdraw the insertion portion of the endoscope, or perform angle control, as this could result in injury, bleeding or perforation to patient. Perform troubleshooting as described in section “8. Troubleshooting” of this IFU.
 - ①No real-time image shown
 - ②Real-time image is frozen
- (3) The patient should be closely observed throughout the use of the product.
- (4) LED illumination is required for endoscopic observation, and improper use could result in viable tissue alterations such as protein denaturation and tissue perforation of viable tissues, please follow the precautions as follows:
 - ①If the distal tip contacts the body, there will be no normal live image shown. The maximum temperature of the distal tip may reach to 60 °C at the 40°C environment, which would result in mucosal burns. Do maintain the proper distance required for adequate observation and do not constantly observe tissue at close range or allow the distal tip of the endoscope to remain in unnecessary contact with the mucosa for prolonged duration, as this would result in an elevated temperature at the tip (possibly exceeding 41°C).
 - ②Be sure to power off the endoscope when it is not used for a long time to avoid unnecessary illumination.
 - ③Do not direct light into the eyes to avoid light damage to the eyes
- (5) If the product is used under extremely poor conditions, LED may be damaged.
- (6) The product is connected to the power supply network through the isolation measures of the Endoscopic Video Image Processor HD (Model: ASUS-IPD-V200).
- (7) Slowly turn the rocker during using, or it may result in damage to the bending function and injury, bleeding or perforation to patient.
- (8) Please plug/unplug and connect the video cable with the Endoscopic Video Image Processor HD off, otherwise it may result in damage to the device.
- (9) Do not excessively coil the video cable, as it would result in abnormal images.
- (10) Ensure that before each use or after a change of viewing modes/settings, the OPERATOR should check to ensure the view observed through the device provides a live image (rather than a stored one) and has the correct image orientation. Otherwise, it could result in patient injury or equipment damage.
- (11) When applying the accessory for diagnosis and treatment, ensure that the distance between the camera module and the mucosa is greater than the minimum visual distance of the endoscope to guarantee that the accessory for diagnosis and treatment always appears in the image, or it could result in patient injury or device damage.
- (12) If it is difficult to insert or withdraw the accessory for diagnosis and treatment, the bending section should be straightened as much as possible without losing the image. Forceful insertion or withdrawal of the accessory for diagnosis and treatment could result in damage to the instrument channel or the accessory for diagnosis and treatment, or cause patient injury, bleeding or perforation.
- (13) If the accessory for diagnosis and treatment is not visible in the endoscopic image, discontinue the use of

the accessory for diagnosis and treatment, or it could result in patient injury, bleeding or perforation or device damage.

- (14) Liquid may be ejected during the withdraw of the accessory for diagnosis and treatment from the instrument channel. To prevent this situation, cover the accessory and the channel surroundings with a piece of gauze during the withdrawal.
- (15) Set the electrosurgical equipment to the minimum necessary output level. If the output level is too high, the insulation of the endoscope or accessory may be damaged and cause operator or patient burns.
- (16) Do not activate the high-frequency current until the distal tip of the high-frequency accessory is within the scope of the endoscope. In addition, make sure that the distance between the electrodes and the mucosa near the target area and the distal tip of the endoscope is appropriate. Emitting the high-frequency current when the distal tip of the accessory is not visible or too close to the distal tip of the endoscope may cause injury to the patient, bleeding or perforation, and damage to the device.
- (17) The use of high-frequency currents may interfere with the normal endoscopic image. This interference is not indicative of a malfunction of the endoscope.
- (18) Apply only electrical conductive jelly to the insertion tube and into the urethra when operating the high-frequency cauterization.
- (19) Use nonelectrolytic solution for infusion. Once urine or blood is mixed into the solution, the nonelectrolytic solution will be diluted, and electrosurgical faculty may be degraded. In this case, do not raise the power of the electrosurgical equipment. Instead, replace the diluted solution with new solution. If the power setting is raised, patient injury can result.
- (20) When using laser equipment, the laser probe should be placed in the correct position shown in the endoscope image. The target should be kept at the appropriate distance from the distal tip, and the output power should be used as low as possible.
- (21) When using compatible laser equipment, the user shall take steps to protect eye and skin, including, but not limited to wearing goggles.

4.2.2 Intraoperative warnings regarding operation prohibition

- (1) Do not hit, collide or drop the endoscope. Do not bend, pull or twist the endoscope. These actions may cause damage to the endoscope or dislodge parts into the patient; which may result in injury, bleeding or perforation.
- (2) Do not pull, twist or rotate the bending section which is already bent at an angle, as this would cause injury, bleeding or perforation of the patient. It would also prevent the bending section from being straightened during the examination.
- (3) Do not press the bending section hard, as the outer skin of the bending section would be stretched or cracked, which may result in water leakage.
- (4) Do not adjust the bending angle with excessive force. Otherwise, the endoscope may be damaged and affect the bending function.
- (5) Do not over bend the insertion portion near the handle connection, as the insertion portion may be damaged.
- (6) Do not insert or withdraw the insertion portion suddenly or with excessive force as this could result in

injury, bleeding or perforation to patient.

- (7) Do not insert or withdraw the endoscope under the following conditions. It could result in patient injury, bleeding or perforation.
 - When the endoscopic accessory for diagnosis and treatment protrudes from the distal tip of the endoscope;
 - When the bending section is at a bending angle.
- (8) Do not pull the video cable during examination, as it could result in abnormal images.
- (9) Do not move the Endoscopic Video Image Processor HD, video cable, adapter, power line etc. during use, as it would result in poor connection which would cause signal interruption.
- (10) Do not use the Endoscopic Video Image Processor HD if real-time endoscopic images cannot be observed, as it may result in patient injury.
- (11) Endoscopes should not be used when delivering oxygen or highly flammable anesthetic gases to patients. This can cause harm to patients.
- (12) Do not perform electrosurgery or laser cauterization while supplying oxygen. This may result in combustion during cauterization.
- (13) Do not insert the accessory for diagnosis and treatment by force or abruptly, or it may result in accidental protrusion of the accessory for diagnosis and treatment from the distal tip of the endoscope, potentially causing damage to the device and patient injury, bleeding or perforation.
- (14) When performing electrosurgery, do not allow the external surface of the control section and its surroundings to become too wet. Operator and patient injury can result.
- (15) Do not use the spray solidification mode of the electrosurgical equipment when performing high-frequency cauterization. It may result in damage to the endoscope and could cause operator or patient burns.
- (16) Never use a nonisolated forceps for high-frequency treatment.
- (17) To avoid interference and/or injury, do not use diathermy device in conjunction with the endoscope.
- (18) Do not apply products such as olive oil or products containing petroleum-based lubricants (e.g. Vaseline®) to the endoscope, as these products may strain and damage the outer skin of the bending section.
- (19) If the endoscope is unable to be withdrawn smoothly from the patient, do not attempt to withdraw it with force. Forced withdrawal could result in injury, bleeding or perforation of the patient.
- (20) While withdrawing the endoscope, be cautious of patient's body fluid adhered. Avoid contacting bed or floor as patient's body fluids could pose a risk of infection to the patient or medical staff.

4.2.3 Intraoperative warnings regarding stopping operation

- (1) In the case of any of the following conditions during the examination, the examination shall be suspended immediately and the endoscope shall be withdrawn from the patient as described in Section 8. "Troubleshooting". The continued use of which may result in injury, bleeding, or perforation to the patient.
 - Any malfunction of endoscope;
 - Unexpected disappearance, freezing, abnormal color or incomplete display of live image;

- Noisy point, blurring or fogging on the live image;
- Rocker cannot be adjusted.

- (2) If the endoscope is dropped or the distal tip is hit by a hard object, stop using the endoscope even if visual check reveals no damage to the camera module.
- (3) The amount of fluid will be low when the endoscopic accessory for diagnosis and treatment is inserted into the instrument channel. Confirm the amount of fluid with that condition in advance.

4.2.4 Intraoperative warnings regarding troubleshooting

- (1) Darkening of the endoscope images during use may indicate that blood or mucous membranes are adhering to the camera module at the distal tip of the endoscope. Carefully withdraw the endoscope from the patient, wipe away the blood or mucous membranes, and restore optimal illumination to ensure a safe examination.
- (2) If it is not possible to withdraw the endoscopic accessory for diagnosis and treatment from the endoscope, close the accessory for diagnosis and treatment, and then withdraw the endoscope and the accessory for diagnosis and treatment carefully while observing the image, so as not to cause soft tissue trauma.
- (3) Some issues encountered (e.g. voltage dips/interruptions testing) can be solved by manually restarting the Endoscopic Video Image Processor HD. Restarting the processor should be performed in the following step:
 - a. Stop the operation and ensure that the power is restored and the interference disappears.
 - b. Restart the Cysto-Nephroscope system to display image.

If the problem cannot be resolved by the actions, contact Cook Medical.

4.3 Postoperative warnings

- (1) Do not exclusively rely on endoscopic observation as the initial diagnostic basis for a lesion to make any decisions about potential diagnostic or device interventions.
- (2) This endoscope is single-use medical device. Please dispose it according to the regulations of the local hospital after use.

5. Installation and Operation

5.1 Preparation and check

Please prepare and check the product in accordance with the following instructions before use. In case of any abnormality, please refer to “8. Troubleshooting” of this IFU for troubleshooting. If there is still a problem with the device, please do not use this product and contact Cook Medical.

	The product is sterilized with ethylene oxide and it is a single use sterile product. Before initial use, make sure the packaging is intact. If the packaging is damaged, do not use the product.
	The shelf life of this product is three years. Do not use it if the product is beyond its shelf life of three years.

5.1.1 Check the endoscope

- (1) Check the handle for any deformation, loose parts or other abnormalities;
- (2) Check the camera module at the distal tip of the endoscope for any scratches, cracks, stains, gaps or other abnormalities;
- (3) Check the entire insertion portion of the endoscope for any abnormalities such as damages, stains and cracks.
- (4) Check whether the functional components of the endoscope are complete.

5.1.2 Check the bending function

- (1) Straighten the bending section;
- (2) Slowly adjust the rocker in both directions to the maximum angle to verify that the bending section bends smoothly and that the maximum angle can be reached;
- (3) Slowly adjust the rocker to the middle position to verify that the bending section can return to the straightened state smoothly.

	In case of abnormal bending functions, do not use the endoscope, as the bending section may not be straightened. This may result in patient injury, bleeding or perforation.
---	--

5.1.3 Connect the endoscope to the Endoscopic Video Image Processor HD

- (1) Connect the Endoscopic Video Image Processor HD to the power supply;
- (2) Connect the video cable to the Endoscopic Video Image Processor HD according to the pair markings;
- (3) Turn on the Endoscopic Video Image Processor HD and verify that the endoscope image is displayed properly.

	Do not move the Endoscopic Video Image Processor HD, video cable, adapter, power line etc. during use, as it would result in poor connection which could cause signal interruption.
---	---

5.1.4 Check the endoscopic live image

- (1) Verify that the endoscope LED illumination is functioning properly;
- (2) Properly adjust the screen brightness;
- (3) Verify that the images are free from noisy point, blurring, fogging or other abnormalities;
- (4) Slowly adjust the maximum angle of rocker up and down respectively, and observe whether the image is abnormal;

- (5) Verify that the image shows no transient disappearance or other abnormalities.
- (6) Verify that real-time endoscopic image appears on the screen by pointing the distal tip of the endoscope at an object, such as the palm of the hand.

	Do not use the Endoscopic Video Image Processor HD if real-time endoscopic images cannot be observed, as it may result in patient injury.
---	---

5.1.5 Check the instrument channel

- (1) Loosen the turncap of the T-valve and insert the endoscopic accessory through the instrument channel port. Verify that the endoscopic diagnosis and treatment accessories extends smoothly from the exit of the instrument channel and no foreign matter are present at the distal tip.
- (2) Verify that the endoscopic diagnosis and treatment accessories can be smoothly withdrawn from the entrance of the instrument channel.

5.2 Operation

	Slowly turn the rocker during use, or it may result in damage to the bending function and injury, bleeding or perforation to patient.
	Do not insert or withdraw the endoscope under the following conditions; otherwise, it may result in patient injury, bleeding or perforation. -When the endoscopic accessory for diagnosis and treatment protrudes from the distal tip of the endoscope; -When the bending section is at a bending angle.
	In case of any of the following conditions during the examination, the examination shall be suspended immediately and the endoscope shall be withdrawn from the patient as described in Section "8. Troubleshooting", the continued use of which may result in injury, bleeding, or perforation to the patient. -Any malfunction of endoscope; -Unexpected disappearance, freezing, abnormal color or incomplete display of live image; -Noisy point, blurring or fogging on the live image; -Rocker cannot be adjusted.
	If the endoscope is dropped or the distal tip is hit by a hard object, stop using the endoscope even if visual check reveals no damage to the camera module.
	When performing electrosurgery, do not allow the external surface of the control section and its surroundings to become too wet. Operator and patient injury can result.

5.2.1 Grip

- (1) Grip the handle of the endoscope with one hand firstly;
- (2) Adjust the rocker with the thumb. The bending section of the endoscope can be controlled with the rocker. When advanced in the "U" direction, the bending section moves up; When advanced in the "D" direction, the bending section moves down; When turned in the "-" direction, the bending section is in a straight state. The

insert tube should always be held as straight as possible to ensure the best angle of bend at the distal tip.

(3) Hold the insertion portion with the other hand.

5.2.2 Infusion

Infuse liquid through the T-valve.

	The amount of fluid will be low when the endoscopic accessory for diagnosis and treatment is inserted into the instrument channel. Confirm the amount of fluid with that condition in advance
---	---

5.2.3 Insert the endoscope

While observing the endoscope image, slowly insert the insertion portion into the patient's body cavity. If

necessary, slowly adjust the rocker to adjust the position of the distal tip to continue the insertion.

	Do not apply products such as olive oil or products containing petroleum-based lubricants (e.g. Vaseline®) to the endoscope, as these products may strain and damage the outer skin of the bending section.
	Do not over bend the insertion portion near the handle connection, as the insertion portion may be damaged.

5.2.4 Adjust bending angle

Operate the rocker as needed to guide the insertion and observation of the tip.

	Do not adjust the bending angle with excessive force. Otherwise, the endoscope may be damaged and affect the bending function.
---	--

5.2.5 Use the endoscopic accessory

- (1) Verify that the endoscopic accessory for diagnosis and treatment is closed;
- (2) Slowly and straightly insert the endoscopic accessory for diagnosis and treatment into the instrument channel port;
- (3) Hold the endoscopic accessory for diagnosis and treatment by hand approximately 4 cm from the instrument channel port, and slowly advance it while observing the endoscopic image and inserting it flat into the instrument channel;
- (4) Ensure that the endoscopic accessory is closed first. Then withdraw it from the instrument channel slowly and straightly while observing the live image.

	When applying the accessory for diagnosis and treatment, ensure that the distance between the camera module and the mucosa is greater than the minimum visual distance of the endoscope to guarantee that the accessory for diagnosis and treatment always appears in the image, or it may result in patient injury or device damage.
	If it is difficult to insert or withdraw the accessory for diagnosis and treatment, the bending section shall be straightened as much as possible without losing the image. Forceful insertion or withdrawal of the accessory for diagnosis and treatment would result in damage to the instrument channel or the accessory for diagnosis and treatment, or cause patient injury, bleeding or perforation.

	If the accessory for diagnosis and treatment is not visible in the endoscopic image, discontinue the use of the accessory for diagnosis and treatment, or it could result in patient injury, bleeding or perforation or device damage.
	Do not insert the accessory for diagnosis and treatment by force or abruptly, or it could result in accidental protrusion of the accessory for diagnosis and treatment from the distal tip of the endoscope, causing damage to the device and patient injury, bleeding or perforation.
	Liquid may be ejected during the withdraw of the accessory for diagnosis and treatment from the instrument channel. To prevent this situation, cover the accessory and the channel surroundings with a piece of gauze during the withdrawal.
	If it is not possible to withdraw the endoscopic accessory for diagnosis and treatment from the endoscope, close the accessory for diagnosis and treatment, and then withdraw the endoscope and the accessory for diagnosis and treatment carefully while observing the image, so as not to cause soft tissue trauma.
	It is important to select the correct size of endoscope attachment. For more information on the combination of endoscope and accessory for diagnosis and treatment, please refer to the IFU for the accessory. For instructions on the operation of the accessory, refer to respective IFU.
	Before the use of endoscope with the high-frequency accessories, it is necessary to install and connect the equipment according to the IFU, and make sure that noise will not affect the observation and surgery. The surface of the high-frequency accessories entering the patient should be checked before use to ensure that there are no rough surfaces, sharp edges or protrusions that will cause harm. If not, may result in injury to the patient.
	Using high-frequency current may cause accidental burns. Please wear insulation gloves.
	Do not activate the high-frequency current until the distal tip of the high-frequency accessory is within the scope of the endoscope. In addition, make sure that the distance between the electrodes and the mucosa near the target area and the distal tip of the endoscope is appropriate. Emitting the high-frequency current when the distal tip of the accessory is not visible or too close to the distal tip of the endoscope may cause injury to the patient, bleeding or perforation, and damage to the device.
	The use of high-frequency currents may interfere with the normal endoscopic image. This interference is not indicative of a malfunction of the endoscope.
	Do not perform electrosurgery or laser cauterization while supplying oxygen. This may result in combustion during cauterization.
	Do not use the spray solidification mode of the electrosurgical equipment when performing high-frequency cauterization. It may result in damage to the endoscope and could cause operator or patient burns.
	Set the electrosurgical equipment to the minimum necessary output level. If the output level is too high, the insulation of the endoscope or accessory may be damaged and cause operator or patient burns.
	Never use a nonisolated forceps for high-frequency treatment.
	Apply only electrical conductive jelly to the insertion tube and into the urethra when operating the high-frequency cauterization.
	Use nonelectrolytic solution for infusion. Once urine or blood is mixed into the solution, the nonelectrolytic solution will be diluted, and electrosurgical faculty may be degraded. In this case, do not raise the power of the electrosurgical equipment. Instead, replace the diluted solution with new solution. If the power setting is raised, patient injury can result.

	When using laser equipment, the laser probe should be placed in the correct position shown in the endoscope image. The target should be kept at the appropriate distance from the distal tip, and the output power should be used as low as possible.
	When using compatible laser equipment, user shall take steps to protect eye and skin, including, but not limited to, wear goggles.

5.2.6 Withdraw the endoscope

- (1) Adjust the rocker to return the bending section to its straightened state;
- (2) Observe the endoscopic image while carefully withdrawing the endoscope;
- (3) Dispose of this product in accordance with hospital waste disposal procedures after surgery, and do not reuse it.

	If the endoscope is unable to be withdrawn smoothly from the patient, do not attempt to withdraw it with force. Forced withdrawal could result in injury, bleeding or perforation of the patient.
	Do not allow the body fluid of the patient sticking to the endoscope to come into contact with the bed or floor while withdrawing the endoscope, otherwise, the patient's body fluids would pose a risk of infection to the patient or medical staff.

6. Storage conditions

Store the endoscope in a clean and dry condition as below table.

Storage Environment	Ambient temperature	3°C ~ 30°C
	Relative humidity	25% ~ 85%
	Air pressure	70 kPa ~ 106 kPa

Transportation environment condition should be accordance with the below table.

Transportation Environment	Ambient temperature	-29°C ~ 55°C
	Relative humidity	10% ~ 95%
	Air pressure	70 kPa ~ 106 kPa

7. Disposal instructions

When disposing of the endoscope or the instrument for biopsy, follow all applicable national and local laws and guidelines.


	The endoscope shall be disposed of in accordance with the hospital waste disposal procedures after using, and shall not be reused. Improper handling or reuse may result in cross-contamination or infection of patients.
---	---

8. Troubleshooting

Item	Situation	Possible cause	Solution
Image quality or brightness	No image display	Failure to turn on the Endoscopic Video Image Processor HD	Check the power connection and turn on the Endoscopic Video Image Processor HD.
		The video cable is not connected firmly	Verify that the side with the pair markings is facing upwards, plug the video cable to the Endoscopic Video Image Processor HD firmly.
	Image is too dark or too bright	Improper brightness setting	Adjust the brightness setting of the Endoscopic Video Image Processor HD as described in the IFU.
	Image is not clear	The distal tip of the endoscope is dirty	Withdraw the endoscope and wipe the lens with sterile gauze.
Infusion	No inflow of fluid/leakage of fluid	Instrument channel port is not closed	Rotate the fastening ring of T-valve.
		Instrument channel leakage	Stop the use of this endoscope.
		Instrument channel blockage	Stop the use of this endoscope.
Endoscopic accessories for diagnosis and treatment	Accessories for diagnosis and treatment cannot pass through the instrument channel smoothly	Instrument channel blockage	Withdraw the accessory for diagnosis and treatment, and try to insert it carefully again. If it still cannot pass through smoothly, please stop the use of the endoscope.
		Use of incompatible accessories for diagnosis and treatment	Select compatible accessories for diagnosis and treatment.
Bending angle	Bending angle cannot be adjusted	Endoscope damage	Stop the use of the endoscope. If it has been inserted into the body cavity of the patient, carefully withdraw it while observing the image.
Withdrawal of abnormal endoscope	Endoscope is withdrawn when the image is displayed normally	/	(1) Switch the Endoscopic Video Image Processor HD to the real-time image display interface, and if the image magnification function is applied, please disable the function first; (2) If an accessory for diagnosis and treatment is used, first close the accessory for diagnosis and treatment, then adjust the rocker to return the bending section to the straightened state, and then carefully and slowly withdraw the accessory for diagnosis and treatment;

			(3) While keeping the bending section straightened, slowly and carefully withdraw the endoscope while observing the image.
	Endoscope is withdrawn when the image cannot be displayed/the image cannot be recovered after freezing	/	(1) Turn off the Endoscopic Video Image Processor HD and then turn it on again. If the real-time image display function is restored, it can continue to be used normally; (2) If the real-time image still cannot be displayed, first close the accessory for diagnosis and treatment, then adjust the rocker to the middle to return the bending section to the straightened state, and then carefully and slowly withdraw the accessory for diagnosis and treatment; (3) While keeping the bending section straightened, withdraw the endoscope slowly; (4) Turn off Endoscopic Video Image Processor HD.

9. Description of symbols

Symbol	Description	Symbol	Description
	Refer to instruction manual/booklet		Do not use if package is damaged and consult instructions for use
	Endoscopic equipment		The device should be sent to the special agencies according to local regulations
	Type BF applied part		Use-by date
	Do not re-use		Temperature limit
	Do not re-sterilize		Humidity limitation
	Sterilized using ethylene oxide		Atmospheric pressure limitation
	Unsafe under nuclear magnetic resonance examination environment		Medical device
	Batch code		Unique Device Identifier
	Serial number		Single sterile barrier system
	Manufacturer		Keep dry
	Date of manufacture		Keep away from sunlight
	Fragile, handle with care		Caution
	Upward		Product is phthalate(DEHP) free
	Stacking limit by number		U.S. Federal law restricts this device to be sold by or on the order of a licensed practitioner.

10. Electromagnetic Compatibility Information

Essential performance: 1) Endoscopic real-time image can be transmitted normally. In case of interference, the image can be restored after the interference is removed. 2) Illumination can be provided normally. In case of interference, the lighting can be restored after the interference is removed. 3) Be able to display image with correct orientation.

	This product shall be used only in the following specified electromagnetic environment. Users and medical staff shall ensure that the device is used in such an environment.
	Leakage current multiplication - patient risk: In case that the product is used in combination with medical electrical devices and/or motorized endoscopic accessories, the leakage current may be multiplied.
	If this product or endoscopic accessories are used with products of different manufacturers, as well as with medical electrical products, it shall be ensured that the conditions of BF (insulated, ungrounded part of the application) are satisfied, and products with the following symbols indicate compliance with BF regulations: 
	Using accessories, transducers and cables beyond the requirements of manufacturer can increase the emission of the equipment or system or reduce the immunity except they are sold by the manufacturer of the equipment or system as spare parts of internal components. There are no other accessories and cables in the endoscope itself. It must be used together with the Endoscopic Video Image Processor HD (Model: ASUS-IPD-V200) distributed by Cook Medical to ensure that it meets and complies with the interference radiation and immunity standards. See the IFU of Endoscopic Video Image Processor HD (Model: ASUS-IPD-V200) for the requirements of relevant cables.
	Devices or systems should not be used in close proximity or stacked with other devices, and if they must be used in close proximity with other devices or stacked, they shall be observed to verify that they can operate properly in the configuration in which they are used.
	Active medical devices are subject to special EMC precautions and therefore shall be installed and used in accordance with these guidelines.
	The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
	Many problems(e.g. voltage dips/interruptions testing) can be solved by restarting the device by manual during procedure. Restarting the device should be tried in the following order: - Stop the operation and ensure that the power is restored and the interference disappears. - Restart the Cysto-Nephroscope system to display image. If the problem cannot be resolved by the actions, contact Cook Medical.

Table 1 Guidance and manufacturer's declaration - electromagnetic emissions

Emissions test	Compliance	Compliance level		
RF conducted emissions CISPR 11	Group 1 Class A ^a *	Frequency range	(Test voltage:120V-230V) ^b *	
		MHz	Quasi-peak	Average
		0.15-0.50	79 dB(μV)	66 dB(μV)

		0.50-30	73 dB(μV)	60 dB(μV)
RF radiated emissions CISPR 11	Group 1 Class A ^{a*}	Frequency range MHz	Quasi-peak (Test voltage:120V-230V) ^{b*}	
			10 m measuring distance	3 m measuring distance ^{c*}
		30-230	40 dB(μV/m)	50 dB(μV/m)
		230-1000	47 dB(μV/m)	57 dB(μV/m)
Harmonic emissions IEC 61000-3-2	Not Applicable			
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable			
^{a*} The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.				
^{b*} The endoscope system is rated power of ≤20 k VA, therefore, selection of the appropriate set of limits shall be based on the rated a.c. power stated by the manufacturer.				
^{c*} The 3 m separation distance applies only to small size equipment meeting the size criterion defined in 3.17 of CISPR 11 Standard.				

Table 2 Guidance and manufacturer's declaration - electromagnetic Immunity

Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV power supply lines ±1 kV signal input/output 100 kHz repetition frequency	±2 kV power supply lines ±1 kV signal input/output 100 kHz repetition frequency
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode	Not Applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U _T ; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % U _T ; 1 cycle and 70 % U _T ; 25/30 cycles; Single phase: at 0°. 0 % U _T ; 250/300 cycle	Not Applicable
Power frequency magnetic field	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz

IEC 61000-4-8		
Conducted RF IEC61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF IEC61000-4-3	3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz
NOTE U _T is the a.c. mains voltage prior to application of the test level.		

Table 3 Guidance and manufacturer's declaration - electromagnetic Immunity

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380–390	TETRA 400	Pulse modulation 18 Hz	27	27
	450	430–470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	28	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	9	9
	745					
	780					
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28	28
	870					
	930					
	1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT;	Pulse modulation 217 Hz	28	28
	1845					
	1970					

			LTE Band 1, 3, 4, 25; UMTS			
	2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28	28
	5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9	9
	5500					
	5785					

Table 4 Guidance and manufacturer's declaration - electromagnetic Immunity

Radiated RF IEC61000-4-39	Test Frequency	Modulation	IEC 60601-1-2 Test Level(A/m)	Compliance level (A/m)
(Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	30 kHz	CW	8	Not Applicable
	134.2 kHz	Pulse modulation 2.1 kHz	65	65
	13.56 MHz	Pulse modulation 50 kHz	7.5	7.5

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Guangzhou Red Pine Medical Instrument Co.,Ltd.
Cook Medical LLC
William A. Cook Australia Pty. Ltd.