Harrison Fetal Bladder Stent Set	EN 3
Souprava stentu Harrison pro fetální močový měchýř	CS 7
Harrison føtalt blærestentsæt	DA 11
Harrison Fetus-Blasenstentset	DE 15
Σετ ενδοπρόσθεσης εμβρυικής ουροδόχου κύστης Harrison	EL 20
Equipo de stent vesical fetal Harrison	ES 25
Set d'endoprothèse de Harrison pour vessie fœtale	FR 29
Harrison magzati hólyagsztent-készlet	HU 34
Set standard con stent per vescica fetale Harrison	IT 38
Harrison stentset voor foetale blaas	NL 42
Føtalt Harrison-blærestentsett	NO 47
Zestaw Harrison do stentowania pęcherza płodu	PL 51
Conjunto de stent para bexiga fetal Harrison	PT 55
Harrison urinblåsestentset för foster	SV 60
Harrison胎儿膀胱支架套装	ZH 64

HARRISON FETAL BLADDER STENT SET

PROFESSIONAL LABELING AND INSTRUCTIONS FOR USE

Humanitarian Device: Authorized by Federal (U.S.A.) law for use in the treatment of fetal post-vesicular obstructive uropathy. The effectiveness of this device for this use has not been demonstrated.

Caution: Federal law restricts this device to sale, distribution, and use by or on the order of a physician with appropriate training and experience.

Drs. Michael R. Harrison, Mitchell S. Golbus, and Roy A. Filley, Fetal Treatment Program at the University of California Medical Center in San Francisco, were the first to successfully develop this technique for suprapubic catheterization of the fetal bladder. Dr. Curtis Lowery later requested this modification which permits a 5.0 French multi-length stent to be placed through a 13 gage needle.

DEVICE DESCRIPTION

The Harrison Fetal Bladder Stent is a double pigtail stent with an outer diameter of 5.0 French and a length ranging from 1.5 cm to 3.5 cm. The proximal pigtail is a double coil that is oriented perpendicular to the stent to allow the pigtail to lie flat along the fetal abdomen. The distal pigtail is a single coil, and it is oriented perpendicular to the stent in order to enhance retention. Both the proximal and distal coils have 4 sideports.

The Harrison Fetal Bladder Stent Set includes the following:

- Multi-length double pigtail stent with a usable length between the pigtails of 1.5 cm to 3.5 cm
- .038 inch (0.97 mm) diameter TFE coated stainless steel wire guide 40 cm long
- 5.0 French positioner 24 cm long
- 13 gage Echotip® trocar needle 18 cm long

The trocar tip is mechanically processed for enhanced visualization during ultrasound.

INDICATIONS FOR USE

The Harrison Fetal Bladder Stent Set is intended for use in fetal urinary tract decompression following the diagnosis of fetal post-vesicular obstructive uropathy in fetuses of 18 to 32 weeks gestational age.

CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is broken.

CONTRAINDICATIONS

The Harrison Fetal Bladder Stent should not be applied if any of the following conditions are present:

- Severe congenital abnormalities that jeopardize neonatal survival
- Abnormal karyotype
- Renal cortical cysts or evidence of renal failure

WARNINGS

1. Implantation of the Harrison Fetal Bladder Stent may result in leakage of amniotic fluid and/or complete rupture of the membranes.
2. Implantation of the Harrison Fetal Bladder Stent carries the possible risk of infection and the development of chorioamnionitis. This could lead to the need for intervention, including termination of pregnancy, and in rare cases, loss of the uterus.
3. The procedure to implant the Harrison Fetal Bladder Stent carries the risk of premature labor which could lead to premature delivery, and in rare cases, damage to the uterus.
4. Once the stent is implanted, there is the possibility that it may become obstructed or dislodged, resulting in the need for repeated stent placements.
5. The Harrison Fetal Bladder Stent is designed and intended for decompression of the fetal urinary tract. There are no data regarding this stent's safety and effectiveness in decompressing, draining, or treating other fetal cavities or conditions.

PRECAUTIONS

1. Patient Evaluation

A complete medical history should be obtained to determine conditions that might influence the selection of the procedure or to identify conditions that are absolute or relative contraindications to surgery.

Physicians should carefully evaluate each case, and only use the stent where there is the potential for significant renal or pulmonary impairment if no intervention is made.

Evidence of residual renal function must be established according to established protocol prior to stent placement.

The following procedures should be conducted to determine if a fetus is a candidate for the Harrison Fetal Bladder Stent:

- a. ultrasound evaluation demonstrating a lower obstructive uropathy (e.g., bilateral hydronephrosis, ureterostasis, megacystitis or oligohydramnios), and ruling out the presence of other congenital anomalies;
- b. fetal karyotype to rule out chromosomal anomalies;
- c. serial vesicocentesis to evaluate the fetal renal function through fetal urinary biochemical parameters. The parameters and their respective cut-off values are shown below:
 - Na⁺ < 100 mg/dL
 - Ca⁺⁺ < 8 mg/dL
 - Osmolarity < 200 mOsm/L
 - β-2-μglobulin < 4 mg/L
 - protein < 20 mg/L

2. Patient Counseling

Prior to implantation of the stent, the patient should be fully informed about alternative procedures and possible side effects and complications. Once patient assessment is complete and a patient is identified as a potential candidate, it is strongly recommended that the patient be given counseling by the physician and other allied health personnel, including genetic counselors, perinatal nurses, perinatal nurse practitioners, or social workers. The patient should be encouraged to discuss openly and fully any questions she may have concerning the Harrison Fetal Bladder Stent. The following may serve as an outline of the issues for discussing and advising the patient:

- a. Explain the patient's individual condition in depth.
 - b. Due to the puncture of the uterus and membranes, this procedure carries the risk of leakage of amniotic fluid and/or complete rupture of membranes.
 - c. Because this is an invasive procedure, this procedure carries the risk of infection and the development of chorioamnionitis. This could lead to the need for intervention, including termination of the pregnancy, and, in rare cases, loss of the uterus.
 - d. Because there is instrumentation of the uterus, this procedure carries the risk of premature labor which could lead to premature delivery, and, in rare cases, damage to the uterus.
 - e. In spite of efforts to assess their baby's renal function prior to placement of the stent, the potential exists for additional renal failure or additional compromise that was either not recognized or that develops following stent placement.
 - f. Even if the procedure is successful and the baby survives birth, there might already be some renal damage which could lead to renal transplant at some point in the baby's life.
 - g. Stents may become obstructed or dislodged resulting in the need for repeated stent placements, either at the time of the original stent placement or later in the pregnancy depending on the time when the original stent becomes obstructed or dislodged.
 - h. Provide the patient with a copy of the Patient Information Booklet.
- ## 3. Clinical Usage
- a. For each procedure, use only stents enclosed in the sealed, sterile package.
 - b. Physicians should be thoroughly familiar with the use of the Harrison Fetal Bladder Stent prior to performing any procedure.

Fetal bladder stents are intended for use by physicians trained and experienced in techniques for fetal bladder stent placement.

ADVERSE EVENTS

The following adverse events have been reported with the use of the Harrison Fetal Bladder Stent:

- Maternal chorioamnionitis

The following side effects have been reported with the use of the Harrison Fetal Bladder Stent:

- Stent migration and stent blockage

Potential adverse events and side effects associated with the use of the stent could include:

- Abdominal wall herniation and its sequelae (including gastroschisis)
- Ascites
- Maternal sepsis
- Amniotic fluid leak
- Direct trauma to the fetus, such as fetal intestinal perforation
- Uterine injury or bleeding, placental bleeding
- Preterm labor
- Spontaneous abortion

CLINICAL STUDIES

At the time of the HDE submission, 148 Harrison Fetal Bladder Stent Sets had been placed in 75 patients. Outcome data from 20 fetuses who were candidates for the stent is included in the HDE. ***It is important to note that the data presented in the HDE are based on a nonrandom group of patients and physicians, of which only a small proportion have responded to the sponsor's retrospective survey. The rates and proportions computed from these data are subject to bias, and the ability to generalize these rates to a larger patient population is limited.***

Table 1. Outcome Data (n=20)

6 fetuses did not receive the stent

- 1 survived to delivery
- 5 did not survive to delivery

14 fetuses received the stent

- 11 survived to delivery
 - 6 survived to discharge
 - 5 did not survive to discharge
- 1 did not survive to delivery
- 2 were electively terminated

The clinical data presented in the HDE were obtained from retrospective case reports. No information is available about the cause of death for the one fetus who died prior to delivery or the five fetuses who died prior to discharge following stent placement. However, all of these fetuses are known to have had significant renal disease, as evidenced by biochemical parameters and/or anatomical abnormalities. These conditions are associated with an extremely high fetal mortality rate.

The submission also includes information on the frequency of need for replacement stents. Six replacements were needed due to stent migration, and one was needed due to partial stent obstruction.

Table 2. Replacement Stents

11 fetuses received a single stent

3 fetuses required a replacement stent

2 required a total of two stents

1 required a total of three stents

The length of time that the stent was in place was available for 9 stent placements in 7 fetuses. The data are shown below:

Table 3. Duration of Stent Placement (days)

0	2-3	3	5	8	14	21	42	70
---	-----	---	---	---	----	----	----	----

As previously noted, due to the limited nature of the clinical data provided in the submission, it is not possible to extrapolate these results to the population in general. However, the limited clinical data suggests that the device will not expose patients to an unreasonable or significant risk of illness or injury, and that the probable benefit to health from using the device outweighs the risk of injury or illness.

INSTRUCTIONS FOR USE

PATIENT PREPARATION

1. Perform final ultrasound to ascertain exact position of fetus.

Manipulation of the fetus *in utero* may be necessary to attain a more ideal fetal position for the performance of the procedure.

2. Sedate the mother if deemed advisable.

Sedation of the fetus is usually not necessary. However, if fetus will require significant manipulation or there is excessive fetal movement, sedation of the fetus *in utero* may be advisable.

3. Ascertain during the ultrasound if there is an adequate window of fluid within the amniotic space for the proximal stent coil to form. Often, due to the presence of oligohydramnios, there is an inadequate window of fluid for the performance of the procedure. If fluid window is insufficient, amnio-infusion may be required.

4. For amnio-infusion, instill between 500 ml and 1000 ml of fluid within the amniotic space. Warmed saline is suggested.

CAUTION: Do not microwave IV bags. Warm under running water or in an approved warming device.

Actual amount of fluid required to achieve an adequate fluid window is dependent on individual patient condition.

5. At the time of amnio-infusion, it is advisable to administer intra-amniotic antibiotics. Due to the threat of chorioamnionitis, a broad spectrum antibiotic such as Nafcillin (500 mg-1000 mg) or a cephalosporin (1 gr-2 gr) that covers penicillin-resistant staph is recommended.

- **Administration of specific antibiotics and dosages is dependent on individual patient conditions, and should be determined by the physician.**
- **The use of tocolytic agents during and after the procedure may be advisable. Close observation of the fetal condition and any signs of pre-term labor are mandatory.**

ASSEMBLY INSTRUCTIONS

1. Just prior to the placement procedure, load the positioner onto the wire guide.
2. Load the stent onto the wire guide by introducing the wire guide into the multi-length coil of the stent until 3 mm-4 mm of the wire guide extends beyond the distal coil of the stent. Place non-flared end of the positioner against the end of the multi-length (proximal) stent coil (Figure A).

Figure A

NOTE: The single coil of the stent is the distal coil which is placed in the fetal bladder. Be sure stent is properly loaded on wire guide before placement procedure begins.

CAUTION: To assure stent coils maintain optimum memory, do not load the stent onto the wire guide until immediately prior to placement through the needle.

IMPLANTING HARRISON FETAL BLADDER STENT

1. Prepare skin puncture site with an antibiotic skin prep. A small wheal of local anesthetic should be applied to the puncture site.
2. Using a scalpel with a No. 11 blade, make a small nick in the skin just large enough to allow passage of the 13 gage trocar needle.
3. Using ultrasonic guidance, advance the trocar percutaneously through the mother's abdominal and uterine walls and into the fetal bladder.

If amnio-infusion was not necessary to produce adequate fluid window, antibiotics may be administered through the needle after amniotic placement but prior to fetal bladder puncture.

Make sure that a 2 cm window of fluid is available between the fetal bladder and the uterine wall.

4. Trocar tip should be advanced 5mm-10mm into the fetal bladder.

CAUTION: If the trocar is advanced too far into an over-distended bladder, proximal stent coil may not have room to form and improper placement may result.

NOTE: Once proper trocar and needle positioning within the fetal bladder has been verified, withdraw the trocar and immediately place thumb over the top of the needle hub to prevent premature decompression. A full fetal bladder is desirable to enhance visualization of coil formation and positioning.

5. While stabilizing the needle, advance the stent assembly into the needle.
6. After the positioner has entered the hub of the needle and is within the needle cannula, stabilize the positioner and remove the wire guide.

The entire stent must be advanced through needle hub and into the needle before removing wire guide.

7. The positioner can now be utilized to advance the stent coil into the fetal bladder. The first ink mark on the positioner indicates that the single pigtail has completely exited the needle when the mark is in line with the needle hub (Figure B).


Figure B

8. Once the straightened pigtail has coiled in the fetal bladder, the needle is very slowly withdrawn over the positioner to the second ink mark, while stabilizing the positioner with one hand (Figure C).


Figure C

CAUTION: This action is important to prevent withdrawal of the stent from the fetal bladder and ensure coiling of the remaining pigtail coil in the amniotic fluid.

CAUTION: If the needle tip is out of the fetus and at the edge of the uterus (fluid window) before reaching the second ink mark, the stent may be advanced into the fluid window with the positioner.

CAUTION: Exact care must be exercised to avoid leaving any segment of the stent in the uterine wall. The extra-fetal segment of the stent will coil to prevent entry of the entire stent into the fetal bladder. The intra-vesical coil will prevent expulsion of the stent from the fetus.

9. Once proper positioning of the stent is assured, the positioner and the needle are removed together.
10. Document proper stent positioning with ultrasound evaluation.
11. Observe fetal bladder to assure stent is functioning properly. Proper positioning is normally evidenced by rapid decompression of the fetal bladder.
12. Monitor the fetus to assure there is no fetal distress or premature labor.

FOLLOW-UP

1. Follow-up is performed by serial ultrasound examinations.
2. Initial follow-up examination should be performed within 48 - 72 hours following the placement procedure.
3. Continuing follow-up examinations should be performed on a weekly basis for the remainder of the pregnancy.

REMOVAL OF THE STENT

The stent should remain in place until appropriate neonatal and pediatric urologic evaluation has taken place. Once the time for removal of the stent has been determined, it should be removed using standard aseptic techniques.

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**ZH
44**

T _ F B S S _ P _ R E V 1

RELIEVING PRESSURE IN MY BABY'S OBSTRUCTED URINARY TRACT USING THE HARRISON FETAL BLADDER STENT SET

Humanitarian Device: Authorized by Federal (U.S.A.) law for use in the treatment of fetal post-vesicular obstructive uropathy. The effectiveness of this device for use has not been demonstrated.

PATIENT INFORMATION

DEFINITIONS

- **Amniocentesis:** A procedure in which a sample of the amniotic fluid is taken and studied.
- **Amniotic Sac:** The space around your baby which is filled with fluid.
- **Chorioamnionitis:** An inflammation of the membrane surrounding your baby.
- **Chorionic Villus Sampling (CVS):** A procedure in which a sample of the tissue in the placenta, which is the membrane surrounding your baby, is taken and studied.
- **Fetal Karyotype:** A study of your baby's chromosomes.
- **Gestational Age:** The time that has passed since your baby's conception.
- **Lower Obstructive Uropathy:** A blockage in the urinary tract below the bladder.
- **Maternal Sepsis:** An infection in the mother.
- **Open Fetal Surgery:** The partial removal of a fetus from the uterus so surgery can be performed to correct a defect.
- **Urinary Ascites:** The leaking of urine into your baby's abdomen.
- **Urinary Tract:** The system which removes urine from the body.

WHY IS THERE PRESSURE IN MY BABY'S URINARY TRACT?

Normally, your baby's urine drains from the kidneys into the bladder. Then it drains through the urinary tract and into the amniotic sac. This is not happening with your baby. Your baby has a condition known as lower obstructive uropathy. This means that there is probably a blockage in your baby's urinary tract. Your baby's urine can't flow into your amniotic sac. Pressure builds up behind this blockage. If left untreated, this pressure could cause damage to your baby's lungs and kidneys. This damage could lead to stillbirth. It could cause your baby to die shortly after being born because the lungs or the kidneys aren't working. It could also cause severe physical deformities. A Harrison Fetal Bladder Stent can decompress, or relieve the pressure in, your baby's urinary tract. This happens because the urine drains through the stent from your baby's bladder into your amniotic sac. It doesn't have to go through your baby's urinary tract so there is no pressure build-up.

ARE MY BABY AND MYSELF CANDIDATES FOR THIS PROCEDURE?

Your doctor can tell you if you and your baby could benefit from this procedure.

You and your baby might be candidates for this procedure if:

- Your baby is 18 to 32 weeks gestational age and has a blocked urinary tract.
- Your baby shows no physical deformities under a detailed ultrasound.
- A study of your baby's chromosomes, also called a fetal karyotype, shows no other serious defects.

WHEN IS THIS PROCEDURE USUALLY DONE?

The doctor can perform this procedure if a baby has reached at least 18 weeks but no older than 32 weeks gestational age. Since each case is different, your doctor will be able to give you more information about you and your baby.

WHAT WILL MY DOCTOR DO TO RELIEVE THE PRESSURE IN MY BABY'S URINARY TRACT?

You probably found out your baby might have an obstructed urinary tract during a routine ultrasound examination. You will need to have another detailed ultrasound examination. This is to make sure your baby has no physical problems.

Then, the doctor will perform a study of your baby's chromosomes. In order to do this, the doctor will have to obtain a sample of your baby's tissue to examine. The doctor does this in one of two ways.

The first is chorionic villus sampling or CVS. To perform CVS, the doctor inserts a catheter into your uterus and takes a tissue sample. The second way is amniocentesis. To perform amniocentesis, the doctor inserts a needle through your abdomen and takes a sample of your amniotic fluid. This tissue or fluid is studied to make sure your baby has no other serious problems.

Finally, the doctor may sample your baby's urine several times over a few days. The doctor inserts a needle through your abdomen and into your baby's bladder. The doctor studies the urine to find out how well your baby's kidneys are working.


The doctor will give you antibiotics to help protect you against infection. You will receive a local anesthesia. The doctor may give you an additional sedative through an I.V. if it is necessary. Your doctor may also give your baby a sedative. The doctor does that if your baby is moving a lot or needs to be moved a great deal to find the best position.

They may also perform an amnio-infusion. To do so, the doctor will infuse, or insert, fluid if there is not enough in the space between your baby and the wall of your uterus. If you need an amnio-infusion, you will probably receive your antibiotics at the same time.

The doctor will place the Harrison Fetal Bladder Stent under ultrasound. The doctor will insert a needle through your abdomen and into the baby's bladder. This is the same as the urine sampling procedure.


If you did not need amnio-infusion, you will receive your antibiotics through this needle. The doctor will inject the antibiotics while the tip of the needle is in the amniotic space, just before it enters your baby's bladder.

The Harrison Fetal Bladder Stent then goes through the needle over a wire guide until the end of it curls up inside your baby's bladder.

Coil Forms in Baby's Bladder


The doctor then pulls the needle back until the other end of the stent is free. This end will curl up just outside your baby's abdomen in the amniotic sac. The doctor then completely removes the needle and leaves the stent in place to drain.

Coil forms in amniotic space as needle is removed

The whole process should only take about 15-30 minutes. The doctor will want you to stay at the hospital for observation long enough to make sure you and your baby are all right. During this time, the doctor will monitor your uterus to make sure there aren't any contractions or irritability. If this happens, the doctor will give you medicine that stops contractions. They will also monitor your baby's heart rate. They do this with an external fetal monitor. Once the doctor is sure everything is fine, you can go home.

You will need to have your first ultrasound examination within 48-72 hours after the procedure. You will need follow-up ultrasound examinations every week until your baby is born. The doctor will be examining the stent to make sure that it is working properly and hasn't moved. Sometimes, through the baby's movement or by grasping it, the stent may come out of your baby's bladder or become kinked. The stent could also become blocked. If any of these things happen, the baby's urine will no longer be able to flow into the amniotic sac. Then the doctor may need to repeat the procedure and replace the stent.

Parts of your body and your baby's body affected by the procedure

WHAT ARE THE RISKS ASSOCIATED WITH RELIEVING THE PRESSURE IN MY BABY'S URINARY TRACT USING THE HARRISON FETAL BLADDER STENT?

There are some possible risks involved with the placement of the Harrison Fetal Bladder Stent.

- You might have chorioamnionitis, an inflammation of the fetal membrane. This may cause you to lose the fluid around your baby. It may cause infection in your baby or cause your baby to be stillborn. This inflammation can happen after any procedure, including placing a fetal bladder stent, in which a doctor places an instrument into your uterus during pregnancy.

- You might have urinary ascites which is the leakage of urine into your baby's abdomen. This happens when the baby's bladder stretches too far. It can result from fetal stent placement. It usually goes away by itself following the placement of a fetal bladder stent.
- Placing the fetal bladder stent could cause you to go into preterm labor. It can happen after any surgery which goes into the uterus during pregnancy.
- The doctor places a needle through your abdomen and uterus and into your baby's bladder. This makes a path for the stent. There might be some minor bleeding from your uterus or the placenta and minor injury to your uterus caused by the passing of the needle. The bleeding will usually stop after a short time.
- You might have maternal sepsis. This is an infection in you, not your baby. It can happen any time a doctor places an instrument into your uterus while you're pregnant. This is why the doctor will give you antibiotics before and after placing the stent.
- Your amniotic fluid may leak from the space between your baby and the wall of your uterus. This is another problem that can happen any time a doctor places an instrument into your uterus while you're pregnant. It may go away by itself. However, there is no other way your doctor can treat the amniotic fluid leak.
- The needle could perforate, or make a hole in, your baby's intestine. It could cause other damage if the doctor doesn't place it accurately. The doctor watches the procedure on ultrasound to make sure this doesn't happen.
- Placing a fetal bladder stent could cause you to abort or miscarry. This could happen if there is an inflammation of the fetal membrane.
- After the stent placement, your baby might move and displace, or pull the stent out of the bladder. Your baby might grab the stent and pull it out. The stent might become kinked or blocked. If this happens, your doctor may need to place another stent in your baby's bladder.

The discomfort you will have when your doctor places this stent is about the same as you had when they took a fluid sample or drained urine from your baby's bladder.

ARE THERE ANY OTHER WAYS TO RELIEVE THE PRESSURE IN MY BABY'S BLADDER?

Your doctor can insert a needle into your baby's bladder above the blockage and drain the urine to relieve the pressure. The doctor will need to do this as often as the pressure builds up to prevent damage to your baby's kidneys and lungs.

Some doctors can perform an open fetal surgery. The doctor removes your baby from your uterus far enough to repair the blockage. After the repair, the doctor would place your baby back in your uterus and your pregnancy would continue.

DECIDING WHETHER OR NOT TO RELIEVE THE PRESSURE IN MY BABY'S URINARY TRACT

Deciding whether or not to relieve the pressure in your baby's urinary tract is an important decision. The decision to have the placement of the Harrison Fetal Bladder Stent or any other procedure to relieve the pressure is up to you and your partner. You do not have to have this procedure performed. You have read about the potential risks of the procedure. You must understand these risks. You should also know the potential risks of leaving this condition untreated. The pressure build-up could damage your baby's lungs and/or kidneys. This could lead to physical deformity or to your baby's death.

If you are thinking about having this procedure performed, you should discuss it with your doctor as soon as possible. Your doctor will be able to explain this procedure to you in more detail. If you decide to have the procedure, your doctor will help you schedule an appointment.

Before you make this decision, you must understand that relieving the pressure in your baby's urinary tract *will not* correct the original defect. It will not get rid of the blockage in your baby's urinary tract. It will allow urine to drain from your baby's bladder without being backed up by the blockage. Your baby's lungs and kidneys can continue to grow without pressure. After your baby is born, your doctor will still need to correct the original defect.

HOW CAN I FIND OUT MORE ABOUT RELIEVING THE PRESSURE IN MY BABY'S URINARY TRACT?

If you have any other questions, please ask your doctor.

WARNING

Following decompression of the fetal urinary tract, **IMMEDIATELY REPORT** any pain, bleeding or fluid loss to your doctor. These abnormal conditions should be closely monitored.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd